



THE CITY LIFE.

With sleek good looks and high caliber performance, the C300 Corpus offers the perfect blend of stability and agility. Features include: adjustable suspension; front and rear battery access; and a seat-to-floor height of 460 mm (18 inches) with up to 45 degrees of tilt. No other product in its class can even compete. Small, nimble and quick, the C300 Corpus is designed for navigating busy streets and getting in and out of tight spaces. Plus, the Corpus is the world's most popular powered seating system with an ergonomic back featuring built-in, adjustable lumbar and trunk supports. The modular design of the Corpus also allows for growth and added seat functions. An ideal combination, the C300 Corpus is ready to take you in any direction you choose. Lead the way.

C300
CORPUS EXPORT

permobil

The Art of Mobility™

POWER FUNCTIONS



Power adjustable tilt.



Power adjustable seat height.



Power adjustable legrests.



Power adjustable recline.



Rnet (Standard) Alternative Drive Controls available (Option).



UniTrack: Complete and flexible support system (Option).



Flip-up armrests for easy transfer (Standard).

SHROUD COLORS



Crystal Blue



Vintage Wine



Metallic Green



Platinum Silver



Black Onyx



Chili Pepper Red



Popstar Pink



Titan Blue

C300 CORPUS™ SPECIFICATIONS

Weight incl. Batteries	approx. 145 kg (320 lbs)	Seat elevator (option)	200 mm max. (8")	Max Speed	7 km/h (4.5 mph)
Base length	1080 mm (42,5")	Power tilt angle (option)	45° max. (EU=30°)	Driving speed tilted/reclined	½ speed at 135° backrest angle
Base width	620 mm (24")	Power recline (option)	145° max	Range	30–35 km (18–22 miles)
Turning radius	625 mm (24,5")	Seat Widths	430/480 mm (17"/19")	Battery Capacity	2 x 60 Ah
Seat to floor height	460 mm (18")	Seat Depths	410–560 mm (16"–22")	Max. user weight	120 kg / 265 lbs
		Backrest heights	535/605 mm (21"/24")		



The Art of Mobility™

Permobil AB, Box 120, SE 861 23 Timrå, Sweden. Phone: +46 60 59 59 00.

www.permobil.com