

C300
PS

EXPORT



The Art of Mobility™

C300 PS™

The C300 is a compact front wheel drive power wheelchair with sleek good looks. Designed for getting in and out of tight spaces, it is easy to manoeuvre indoors and perfect for getting out and about in the city.

The new PS seat is simple and smart. The telescopic frame and expanding seat pan allow for a full range of growth. The backrest upholstery incorporates adjustable slingbacks to provide fine adjustment of seating comfort.


POWER FUNCTIONS (OPTION)




Adjustable slingbacks to allow fine adjustment of seating comfort.




SHROUD COLORS




Crystal Blue




Vintage Wine



Metallic Green



Platinum Silver



Black Onyx



Chili Pepper Red



Popstar Pink




Titan Blue



VR2 electronics.



4-wheel adjustable independent suspension for comfort and grip.



Tri-spoke split rims with color matching hubcaps. Solid or air filled.

C300 PS™ SPECIFICATIONS

Weight incl. batteries	approx. 130 kg (290 lbs)	Seat elevator (option)	200 mm max. (8")	Max Speed VR2	7 km/h (4.5 mph)
Base length	1200 mm (47.25")	Power tilt angle (option)	45° max. (EU 30°)	Driving speed tilted	½ speed at >20° tilt
Base width	620 mm (24.4")	Seat Widths	400/450 mm (16"/18")	Range	30–35 km (18–22 miles)
Turning radius	637 mm (25")	Seat Depths	410–535 mm (16"–21")	Battery Capacity	2 × 60 Ah
Seat to floor height	460 mm (18")	Backrest heights	360–605 mm (14"–24")	Max. user weight	120 kg (265 lbs)



The Art of Mobility™

Permobil AB, Box 120, SE 861 23 Timrå, Sweden. Phone: +46 60 59 59 00.

www.permobil.com