



**C350**

EXPORT

# C350 CORPUS

Rear-wheel driven and fast.



The Art of Mobility™

# C350 Corpus™

**If you prefer a rear wheel drive chair and demand top class seating comfort...**

The C350 offers you a fast, easy to drive, powerbase with an attractive modern design.

With a Corpus seating system on your C350 you have the possibility to customise your seating to fit you perfectly. The almost limitless adjustment of depth, width, height and angle of the Corpus' multiple seat components contributes to an ergonomically pleasing package with a series of power seat functions at your fingertips.

POWER FUNCTIONS (OPTION)



SHROUD COLORS



Power tilt and recline.



Power seat elevator.



Power legrest.



R-Net option (10 km/h).  
Alternative Drive Controls available. (Option)



4-wheel adjustable independent suspension for comfort and grip.



Tri-spoke split rims with color matching hubcap. Solid or air filled.



UniTrack – complete and flexible support system. (Option)

## C350 CORPUS™ SPECIFICATIONS

Weight incl. batteries	142 kg (313 lbs)	Seat elevator (option)	200 mm max. (8")	Max Speed	10 km/h (6.5 mph)
Base length	1065 mm (42")	Power tilt angle (option)	45° max. (EU=30°)	Driving tilted/reclined ½ speed at >135° backrest angle	
Base width	625 mm (24.6")	Power recline (option)	150° max.	Range	30–35 km (18–22 miles)
Turning radius	945 mm (37")	Seat Widths	430/480 mm (17"/19")	Battery Capacity	2 x 60 Ah
Seat to floor height	460 mm (18")	Seat Depths	410–560 mm (16, 18, 20, 22")	Max. user weight	136 kg (300 lbs)
		Backrest heights	535mm/605 mm (21"/ 24")		



**The Art of Mobility™**

Permobil AB. Box 120, SE 861 23 Timrå, Sweden. Phone: +46 60 59 59 00.

[www.permobil.com](http://www.permobil.com)