



C350

EXPORT

C350 PS

Rear-wheel driven and fast.



The Art of Mobility™

C350 PS™

If you prefer a rear wheel drive chair...


The C350 offers you a fast, easy to drive, powerbase with an attractive modern design.

The new PS seat is simple and smart. The telescopic frame and expanding seat pan allow for a full range of growth. The backrest upholstery incorporates adjustable slingbacks to allow fine adjustment of seating comfort.

POWER FUNCTIONS (OPTION)




SHROUD COLORS



VR2 electronics. (Standard)



R-Net option (10 km/h).
Alternative Drive Controls available. (Option)



4-wheel adjustable independent suspension for comfort and grip.



Tri-spoke split rims with color matching hubcaps. Solid or air filled.

C350 PS™ SPECIFICATIONS

Weight incl. Batteries	approx. 140 kg (310 lbs)	Seat elevator (option)	200 mm max. (8")	Max Speed VR2	8 km/h (5 mph)
Base length	1170 mm (46")	Power tilt angle (option)	45° max. (EU=30°)	Max Speed Rnet	10 km/h (6.5 mph)
Base width	620 mm (24.4")	Seat Widths	400/450 mm (16"/18")	Driving speed tilted	½ speed at >20° tilt
Turning radius	890 mm (37")	Seat Depths	410–535 mm (16"–21")	Range	30–35 km (18–22 miles)
Seat to floor height	460 mm (18")	Backrest heights	360–605 mm (14"–24")	Battery Capacity	2 x 60 Ah
				Max. user weight	136 kg (300 lbs)



The Art of Mobility™

Permobil AB, Box 120, SE 861 23 Timrå, Sweden. Phone: +46 60 59 59 00.

www.permobil.com