



## FLEXIBILITY IS BOUNDLESS.

The streamlined design and robust power that you expect from Permobil's high-grade power wheelchairs is once again presented in the new rear wheel drive C350 TS. No other wheelchair in this category offers the perfect combination of stability, seating comfort and versatility that this model does at the best overall price. Thanks to Permobil's passion for engineering you now have an even broader choice of seating systems.

The TS (Tension Seat) system's strongpoint is its adjustability. The built-in, individually tensionable, straps in the backrest and under the seating surface allow the seating to be shaped exactly to the user's needs. Integrated UniTrack profiles provide fast and seamless adjustment of seat depth, armrest height and positioning of optional supports to provide optimum seating comfort. Up to 4 powered seat functions: seat elevator, tilt, recline and legrest, are available as options.

# C350 TS

## POWERED SEAT FUNCTIONS



Power tilt and recline available as option.

SHROUD COLOR



Black



Seamless adjustment of seat depth and arm rest height.



The built-in, individually tensionable, straps in the backrest and under the seating surface allow the seating to be shaped exactly to the user's needs.



Split-rim wheels with color matched hub caps.



Powered seat elevator.  
(Available as option)



Power leg rest angle adjustment.  
(Available as option)



VR2 electronics with options for 4 powered seat functions.

## C350 TS™ SPECIFICATIONS

Weight incl. Batteries	151 kg (333 lbs)	Seat Widths	380–530 mm (15"–21")
Base length	1065 mm (42")	Seat Depths	380–550 mm (15"–22")
Base width	625 mm (24.6")	Backrest heights	450–600 mm (18"–23.5")
Turning radius	890 mm (35")	Max Speed	8 km/h (5 mph)
Seat to floor height	460–510 mm (18"–20")	Range	30–35km (18–22 miles)
Seat elevator (option)	200 mm max. (8")	Battery Capacity	2 × 60 Ah
Power tilt angle (option)	30° (0–30°, –5–25°)	Max. user weight	136 kg (300 lbs)
Power recline (option)	125° max		



The Art of Mobility™

Permobil AB.  
Box 120, SE-861 23 Timrå, Sweden.  
Phone: +46 60 59 59 00.

[www.permobil.com](http://www.permobil.com)