



ON THE MOVE.

Designed for very active people, the C500 Corpus™ is perfect for both indoor and outdoor use. Powerful motors give it enough torque to navigate through rough terrain, while the adjustable independent suspension provides a smooth ride. Users can also maneuver in tight spaces with a 715 mm (28") turning radius, perfect for crowded clubs or restaurants. The Corpus seating system of the C500 is ergonomically pleasing with built-in, adjustable lumbar and trunk supports. Power functions include: tilt, recline, seat elevator, legrest elevation and many other custom accessories for the person on-the-go.



C500
CORPUS EXPORT



THE POWER OF MOBILITY




POWER FUNCTIONS



Power adjustable tilt, seat height, legrests and recline.

SHROUD COLORS



Co-pilot: Easy to use, power assisted attendant control. (Option)




Magic Drive: Opens up a world of alternative control possibilities. (Option)



ESP (Enhanced Steering Performance) virtually eliminates the need for constant manual correction when driving on uneven ground, soft surfaces, angled and sloped areas or at high speeds.



UniTrack: Complete and flexible support system. (Option)



Permlock C: Lock your Permobil into position for safer travel by car. (Option)



R-Net (Standard). Alternative Drive Controls available. (Option)

C500 CORPUS™ SPECIFICATIONS

Weight incl. Batteries	158 kg (348 lbs)
Base length	1150 mm (45.2")
Base width	675 mm (26.6")
Turning radius	715 mm (28")
Seat to floor height*	520 mm (20.5")

Seat elevator (option)**	200 mm max. (8")
Power tilt angle (option)***	45° max. (EU=30°)
Power recline (option)	150° max.
Seat Widths	430/480 mm (17"/19")
Seat Depths	360-560 mm (14-22")
Backrest heights	535/605 mm (21"/24")

Max Speed Standard	7.5 km/h (4.7 mph)
Max Speed Option S model	10 km/h (6.5 mph)
Driving speed	
tilted	½ speed at 135° backrest angle
Range	30-40 km (18-25 miles)
Battery Capacity	2 x 73 Ah
Max. user weight	136 kg (300 lbs) (EU 120 kg / 265 lbs)

* Lowrider 460mm (18")
 **Lowrider 150mm (6")
 ***Lowrider 20°



THE POWER OF MOBILITY

Permobil AB, Box 120, SE 861 23 Timrå, Sweden. Phone: +46 60 59 59 00.

www.permobil.com