



The Power of Mobility

A PRACTICALLY INVISIBLE CO-PILOT

It is important for a wheelchair to be easy for an attendant to control. Permobil's new wheelchair Co-pilot, which is used by attendants, is very different to earlier versions of attendant controls. The main advantage of the Co-pilot is that it is power assisted and the wheelchair is highly intuitive to control; it's as easy as steering a pushchair. The control system responds immediately to any change of speed or direction. The Co-pilot also has a unique function that allows the attendant to control the seat functions and speed.



Gently grip the control handle and the wheelchair will move in the desired direction. Push lightly to the right and the wheelchair will move to the right without any difficulty at all. If you start to walk a little faster, the wheelchair will also accelerate. Thanks to the proportional control the wheelchair is easier and safer to operate for both the wheelchair user and their attendant.



The ergonomically designed handle can be folded away when it is not in use.



The height and angle of the handle is adjustable, making it even more user friendly for the attendant.

For more information and prices please contact your Permobil sales representative.

Call 01484722888 or visit our website at www.permobil.com

COPILLOT122009UK