



THE CITY LIFE

With sleek good looks and high calibre performance, the C300 Corpus offers the perfect blend of stability and agility. Features include: adjustable suspension; front and rear battery access; and a seat-to-floor height of 46 cm with a 45 degree tilt. No other product in its class can even compete. Small, nimble and quick, the C300 Corpus is designed for navigating busy streets and getting in and out of tight spaces. In addition to which, the Corpus is the world's most popular powered seating system with an ergonomic back featuring built-in, adjustable lumbar and trunk supports. The modular design of the Corpus also allows for growth and added seat functions. A perfect match, the C300 Corpus is ready to take you wherever you want to go. Lead the way.

C300
CORPUS

POWER FUNCTIONS



Power adjustable tilt



R-Net electronics with ICS



Power adjustable recline



Power seat height adjustment



Power adjustable legrests

SPECIFICATIONS

Length	108 cm	Legrest angle - adjustable	80° - 150°
Width	62 cm	Battery capacity	2 x 60 AH
Min. transport length	88,5 cm	Range	25 - 35 km
Seat height	51 cm	Max. speed	7 km/hr
Seat height with electric seat lift	51 - 71 cm	Turning radius (ISO 7176-5)	62,5 cm
Weight incl. batteries	112 kg	Obstacle capability	7 cm
Max. user weight	120 kg	Suspension	Yes
Seat width	43 - 48 cm	Electronics	R-Net
Seat depth	41 - 46 - 51 - 56 cm	Electric seat lift	Yes
Backrest width	34 - 39 - 44 cm	Electric tilt adjustment	Yes
Backrest height	53,5 - 60,5 cm	Electric legrest adjustment	Yes
Distance between armrests	35 - 40 - 45 cm	Electric backrest adjustment	Yes
Backrest angle - adjustable	90° - 165°	Warranty	1 Year



Dark blue metallic



Red metallic



Dark green metallic



Silver metallic



Black metallic



Pink



Light blue metallic

FRAME COLOURS

ISO 7176-19 / ISO 10542-3
EN12182/EN12184



The Power of Mobility

Permobil Ltd.
Oceanic, Ground Floor
Navigation Business Park
Waters Meeting Road
Bolton, BL1 8SW
info@permobil.co.uk

www.permobil.com