



AN ORGANIC EXPERIENCE

The C400 Corpus combines innovative design, styling and advanced engineering concepts. The C400 maintains a small profile but combines better indoor manoeuvrability with superior outdoor performance. With a rugged, advanced suspension system, the C400 has advanced obstacle climbing ability. The Corpus seating system of the C400 is modular in design, allowing for growth and additional seat functions. Ergonomically pleasing with power tilt and recline, customization of the Corpus is practically unlimited... an important feature when you are becoming one with nature.

C400

CORPUS

POWER FUNCTIONS



Power adjustable tilt

FRAME COLOURS



Dark blue metallic



Red metallic



Green metallic



Black metallic



R-Net electronics with ICS



Power adjustable recline



Power seat height adjustment



Power adjustable legrests

ISO 7176-19 / ISO 10542-3
EN12182/EN12184



SPECIFICATIONS

Length	108,5 cm	Legrest angle - adjustable	80° - 150°
Width	62 cm	Battery capacity	2 x 73 AH
Min. transport length	89 cm	Range	30 - 35 km
Seat height	56 cm	Max. speed	7 km/hr
Seat height with electric seat lift *	56 - 76 cm	Turning radius (ISO 7176-5)	66,5 cm
Weight incl. batteries	140 kg	Obstacle capability	6 cm
Max. user weight	120 kg	Suspension	Yes
Seat width	43 - 48 cm	Electronics	R-Net
Seat depth	36 - 41 - 46 - 51 - 56 cm	Electric seat lift	Yes
Backrest width	34 - 39 - 44 cm	Electric tilt adjustment	Yes
Backrest height	53,5 - 60,5 cm	Electric legrest adjustment	Yes
Distance between armrests	35 - 40 - 45 cm	Electric backrest adjustment	Yes
Backrest angle - adjustable	80° - 165°	Warranty	1 Year



The Power of Mobility

Permobil Ltd.
Oceanic, Ground Floor
Navigation Business Park
Waters Meeting Road
Bolton, BL1 8SW
info@permobil.co.uk

www.permobil.com

* Lowrider 47,5 - 62,5