



**standard**  
**ESP** Enhanced Steering Performance

# ON THE MOVE

Designed for very active people, the C500 Corpus is perfect for both indoor and outdoor use. Powerful motors give it enough torque to navigate rough terrain, while the adjustable independent suspension provides a smooth ride. Users can also manoeuvre in tight spaces with a 71,5 cm turning radius, perfect for crowded clubs or restaurants. The Corpus seating system of the C500 is ergonomically pleasing with built-in, adjustable lumbar and trunk supports. Power functions include: tilt, recline, seat lift, legrest elevation and many other custom accessories for people on-the-go.

## C500 CORPUS

**POWER FUNCTIONS**



Power adjustable tilt



R-Net electronics with ICS



Power adjustable recline



Power seat height adjustment



Power adjustable legrests

FRAME COLOURS



Dark blue metallic



Red metallic



Green metallic



Black metallic

**SPECIFICATIONS**

Length	115 cm	Legrest angle - adjustable	80° - 150°
Width	65 cm	Battery capacity	2 x 73 AH
Min. transport length	90,5 cm	Range	30 - 40 km
Seat height	56 cm	Max. speed	7,5 - 10 km/hr
Seat height with electric seat lift *	56 - 76 cm	Turning radius (ISO 7176-5)	71,5 cm
Weight incl. batteries	143 kg	Obstacle capability	6 cm
Max. user weight	120 kg	Suspension	Yes
Seat width	43 - 48 cm	Electronics	R-Net
Seat depth	36 - 41 - 46 - 51 - 56 cm	Electric seat lift	Yes
Backrest width	34 - 39 - 44 cm	Electric tilt adjustment	Yes
Backrest height	53,5 - 60,5 cm	Electric legrest adjustment	Yes
Distance between armrests	35 - 40 - 45 cm	Electric backrest adjustment	Yes
Backrest angle - adjustable	80° - 165°	Warranty	1 Year

ISO 7176-19 / ISO 10542-3  
EN12182/EN12184



**The Power of Mobility**

Permobil Ltd.  
Oceanic, Ground Floor  
Navigation Business Park  
Waters Meeting Road  
Bolton, BL1 8SW  
info@permobil.co.uk

www.permobil.com

\* Lowrider 47,5 - 62,5