



# CHAIRMAN HD CORPUS



**MAXIMUM SPEED**

6,5 KM/H



**RANGE**

30-40 KM



**MOTOR POWER**

2 X 300 WATTS



Twin castors provide improved stability and bearing pressure distribution



The small support wheels make a big difference to safety



Available with 20 cm seat lift as an option

*The Power of Mobility*

# CHAIRMAN HD CORPUS

The Chairman Heavy Duty is highly robust and built to enable people up to 200 kg to move around as free as a bird. What is striking about this wheelchair is that it hardly any bigger than other units. The chassis is wider and stronger, the backrest and armrests have been adapted to cope with heavier loads and the seating system is deeper, wider and higher. Twin castors provide improved stability and bearing pressure distribution.

The Corpus seating system's modular design provides enhanced functionality through adaptable seat sizing and electric adjustments. The Corpus is designed by one of the leading ergonomic experts in the world resulting in a seating system that is fully adapted to the natural contours of the body. The adaptability of the Corpus seating system is limitless. Depth, width, height, and the angles of the various seat components can be adjusted or adapted. The ergonomically shaped seating and optional electric adjustments make the Corpus the world's most popular seating system.



Optional tilt adjustment enables the seating system to tilt up to 25 degrees



High supports are available in fabric or imitation leather as an option



Optional seat lift enables the seating system to lift up to 20 cm



Calf supports are available in fabric or imitation leather as an option


## TECHNICAL SPECIFICATIONS

LENGTH, MM	1120	LEGREST ANGLE – ADJUSTABLE	80°-135°
WIDTH, MM	700	BATTERY CAPACITY, AH	2x73
MIN. TRANSPORT LENGTH, MM	920	CHARGING TIME, HOURS	8
SEAT HEIGHT, MM	550	RANGE, KM	30-40
SEAT HEIGHT WITH ELECTRIC SEAT LIFT, MM	550-750	MAX. SPEED, KM/H	6,5
WEIGHT INCL. BATTERIES, KG	165	BRAKING DISTANCE, MM	1200
MAX. USER WEIGHT, KG	200	TURNING RADIUS PER ISO 7176-5, MM	750
SEAT WIDTH, MM	480/530/580	OBSTACLE CAPABILITY, MM	60
SEAT DEPTH, MM	435-635	SUSPENSION	No
BACKREST HEIGHT, MM	720	ELECTRIC SEAT LIFT	Yes
ARMREST HEIGHT, MM	190-300	ELECTRIC TILT ADJUSTMENT	Yes
DISTANCE BETWEEN ARMRESTS, MM	470/520/570	ELECTRIC LEGREST ADJUSTMENT	Yes
DISTANCE BETWEEN SEAT CUSHION AND FOOTPLATE, MM	400-500	ELECTRIC BACKREST ADJUSTMENT	Yes
BACKREST ANGLE – ADJUSTABLE	90°-135°	CONTROL SYSTEMS	Yes
		WARRANTY	1 Year

## ELECTRIC ADJUSTMENTS



## COLOURS



BLACK

ISO 7176-19 / ISO 10542-3  
EN12182/EN12184



Permobil Ltd.  
Oceanic, Ground Floor  
Navigation Business Park  
Waters Meeting Road  
Bolton, BL1 8SW  
info@permobil.co.uk

www.permobil.com