

# The C350 makes the difference!



The Power of Mobility

Permobil Ltd.  
Oceanic, Ground Floor • Navigation Business Park • Waters Meeting Road  
Bolton, BL1 8SW • info@permobil.co.uk

www.permobil.com

C300C350TS122010UK



**C300 and C350**



The Power of Mobility

**TS**



**C300 TS**



**C350 TS**

# FLEXIBILITY IS BOUNDLESS

The streamlined design and robust power that you expect from Permobil's high-grade electric wheelchairs is once again present in the brand new C300 TS and C350 TS. No other electric wheelchair in this category offers the perfect combination of stability, seating comfort and versatility that these two do at the best overall price. Thanks to Permobil's passion for engineering you now have an even broader choice of seating systems. The TS seating system is noted for being multifaceted and adjustable. Thanks to the integrated unitrack profile in the durable tubular material you can quickly adjust the seat depth, armrest height and the position of the lateral supports. The standard low, 51-cm seat height ensures optimum functionality. Options include 4 electric power functions: seat lift, and adjustments to tilt, back and leg angles. Furthermore, the TS seat can be adjusted to 3 different heights. The potential is unequalled – 2 types of leg rests are available, the height and angle of armrests can be adjusted and various positioning accessories such as thigh supports, body supports, calf supports and headrest are also available to provide the user with optimum support.



The Power of Mobility


www.permobil.com



VR2 electronics with extension for 4 joystick-operated power functions.



Three spoke rims with hubcaps in the same colour as the capping.



### POWER FUNCTIONS



### TECHNICAL SPECIFICATIONS C300 TS C350 TS

	C300 TS	C350 TS	C300 TS	C350 TS
Length	1080 mm	1065 mm	Charging time	8 hour
Width	620 mm	625 mm	Range	25 - 35 Km
Minimum transport length	885 mm	900 mm	Maximum speed	7 Km/h
Seat height	510 mm	510 mm	Braking distance	1200 mm
Seat height with electric seat lift	510/710 mm	510/710 mm	Turning circle	625 mm
Maximum user's weight	136 kg	136 kg	Obstacle climbing ability	70 mm
Weight (including batteries)	128 kg	130 kg	Suspension	yes
Seat width	380 - 530 mm	380 - 530 mm	Seat lift	yes
Seat depth	380 - 550 mm	380 - 550 mm	Seat tilt	yes
Backrest height	450/600 mm	450/600 mm	Electric leg rest adjustment	yes
Armrest height	190 - 390 mm	190 - 390 mm	Electric backrest adjustment	yes
Distance between armrests	380 - 530 mm	380 - 530 mm	Control systems	VR-2
Backrest angle – adjustment range	90 - 125°	90 - 125°	Guarantee	1 Year
Leg rest angle – adjustment range	90 - 150°	90 - 150°	ISO-standard crash test	yes
Battery capacity	2x60 Ah	2x60 Ah	4-point locking system	yes

ISO 7176-19 / ISO 10542-3  
EN12182/EN12184



Electric back and/or seat tilt



Electric seat lift



Electric leg rest adjustment

