

# K450 MX



**\* Indicates a Required Field**

## Dealer Information



\* **Contact:** \_\_\_\_\_

Dealer Code: \_\_\_\_\_

\* **Dealer Name:** \_\_\_\_\_

Address: \_\_\_\_\_

\* **City:** \_\_\_\_\_

\* **State/Zip:** \_\_\_\_\_

\* **Phone#:** \_\_\_\_\_

\* **Fax#:** \_\_\_\_\_

PO#: \_\_\_\_\_

Email Address: \_\_\_\_\_

## Client Information

*Permobil recommends that the client is evaluated by a certified rehab specialist.*

\* **First Name:** \_\_\_\_\_

\* **Last Name:** \_\_\_\_\_

Diagnosis: \_\_\_\_\_

Funding Source: \_\_\_\_\_

Client Age: \_\_\_\_\_

## Client Measurements

\* **Weight:** \_\_\_\_\_

\* **Height:** \_\_\_\_\_

A) Top of Shoulders: \_\_\_\_\_

B) Chest Depth: \_\_\_\_\_

C) Chest Width: \_\_\_\_\_

D) Seat Depth: \_\_\_\_\_


E) Top of Head: \_\_\_\_\_

F) Elbow to Hand: \_\_\_\_\_

G) Seat Pan to Elbow: \_\_\_\_\_

H) Hip Width: \_\_\_\_\_

\* **I) Knee to Foot:** \_\_\_\_\_



Three column pricing on order/quote:

Display HCPCS Codes on order/quote:

Please send order/quote to fax#: **(800) 231-3256**

Email to: [sales@permobilus.com](mailto:sales@permobilus.com)

**Permobil Inc.**  
6961 Eastgate Blvd.  
Lebanon, TN 37090  
Tel: (800) 736-0925  
Fax: (800) 231-3256  
[www.permobil.com](http://www.permobil.com)

Prices effective **March 8, 2010.**

## Chair Model and Colors

## Select Options:

Part Number	Description	Price	HCPCS Code
<input checked="" type="checkbox"/> I104524-99-0	K450 Base w/ R-net 5.2mph - K450 PS2 <i>Price Includes: K450 Power Base Set Up For Multiple Power Seat Functions, Tie-Down Hardware for Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below.</i>	11,845.00	K0891

## Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104525-99-0	Flat Free Drive Tires - K450	No Charge	
<input type="checkbox"/> I10607	Pneumatic Drive Tires & Flat Free Casters - K450	No Charge	

## Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104534-99-0	Shroud for K450 - Chili Pepper Red	No Charge	
<input type="checkbox"/> I104535-99-0	Shroud for K450 - Onyx Black	No Charge	
<input type="checkbox"/> I104536-99-0	Shroud for K450 - Platinum Silver	No Charge	
<input type="checkbox"/> I104537-99-0	Shroud for K450 - Pop Star Pink	No Charge	
<input type="checkbox"/> I104538-99-0	Shroud for K450 - Titan Blue	No Charge	
<input type="checkbox"/> I104539-99-0	Shroud for K450 - Extreme Yellow	No Charge	

## Base Options

## Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input checked="" type="checkbox"/> I10510	Power Adjustable Seat Height - 8" Travel <i>This item must be selected.</i>	2,495.00	E2300

## Select Options:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> IM34 SLDG	Batteries - Group 34, 12V 60 Amp Hours, Installed	762.00	K0108
<input type="checkbox"/> I104550-99-0	Anti-Tippers for K450	395.00	
<input type="checkbox"/> I104540-99-0	Light Kit - K450 R-net	985.00	
<input type="checkbox"/> I104553-99-0	WC 19 Transit Kit - K450 <i>This kit includes additional anchor points for tie down straps, seat stabilizing bars and hardware, and a WC-19 Approved integrated lap and shoulder belt. This kit has been approved according to ANSI/RESNA WC19. Cannot be ordered with a 10"W Seat Size (I104317-99-0). An ICS Switchbox must be installed on the chair in order to access WC19 mode.</i>	1,585.00	
<input type="checkbox"/> I104554-99-0	Fixed Ventilator Tray - K450 <i>Tray is 9.75"D x 13"W, and increases the overall length of a K450 by approximately 11". The ventilator tray is designed to accommodate an LP10 style ventilator and one U1 Battery. A custom charge will need to be applied if ordering vent tray with any of the following items: WC-19 Transit Kit (I104553-99-0), Attendant Control (ID50882), Push Handles (I104576-99-0), or Back Pack Clips (I10598).</i>	525.00	E1029
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	435.00	

## Permolock

## Select Options:

Part Number	Description	Price	HCPCS Code
-------------	-------------	-------	------------

## Select Options:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104680-99-0	Permolock C3 Locking Base C300, C350, K300, K450 <i>Kit Includes: Permolock C3 Locking Base, Owner's Manual, Manual Release Handle, and Vehicle Mounted Release Button. Does NOT include Dual Locking Pin Kit (I10608). Locking Base must be installed by a NMEDA Dealer. Please provide the NMEDA Dealer and Address where the locking Base will ship to: NAME: _____ ADDR: _____</i>	1,775.00	
<input type="checkbox"/> I104652-99-0	Permolock C3 Dual Locking Pin Kit - K450 <i>This kit includes two pins that are mounted underneath the base of the wheelchair. The pins are fixed and do not push into the wheelchair base.</i>	595.00	

## Joystick Options

## Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard.	No Charge	
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	94.00	E2324
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	139.00	E2323
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	94.00	E2323
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	127.00	E2323
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	127.00	E2323
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	127.00	E2323
<input type="checkbox"/> IPC101B	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	127.00	E2323
<input type="checkbox"/> IPC102B	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	127.00	E2323
<input type="checkbox"/> IPC107B	Bodypoint J/S Handle Dome Shaped, Rubber	127.00	E2323

## Select Options:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10412	Joystick Mount Right - Fixed, R-net Only	No Charge	
<input type="checkbox"/> I10413	Joystick Mount Left - Fixed, R-net Only	No Charge	
<input type="checkbox"/> I102537-99-0	Joystick Mount Right - Swing Away, R-net & VR2	292.00	E1028
<input type="checkbox"/> I102538-99-0	Joystick Mount Left - Swing Away, R-net & VR2	292.00	E1028
<input type="checkbox"/> ID51108	R-net LED Remote Joystick Module <i>Standard Joystick Module. This Joystick does NOT have an LCD Display and does NOT contain mono jacks. Limited to five profiles with this Joystick Module. Must order expandable controller (PN I10411) with this part Please note if ordering this item, the ICS Switchbox must be mounted onto either the armrest pad or the backrest. Cannot order PN: I1822845 or I103469-99-0 with this style of joystick.</i>	No Charge	
<input type="checkbox"/> ID51036	R-net Remote Joystick w/ Mono Jacks <i>Expandable/remote joystick used with multiple power option bases or as an upgrade for no power option or single power option bases. Must order expandable controller (PN: I10411) with this part. Comes standard w/ mono jacks.</i>	530.00	
<input type="checkbox"/> I103814-99-0	Compact Advanced Joystick for R-net <i>Compact sized joystick that does not contain an LCD Display. This is an expandable/remote joystick used with multiple power option bases or as an upgrade for no power option or single power option bases. Must order expandable controller (PN: I10411) with this part. Contains two mono jacks for mode &amp; on/off. Does NOT contain a wrist support pad. If you want a wrist support pad, please select I10621, below.</i>	605.00	E2373
<input type="checkbox"/> I10621	Add Wrist Pad to Compact Joystick Advanced	195.00	
<input type="checkbox"/> ID51040	R-net Remote Joystick w/ Toggle Switches <i>Expandable/remote joystick with toggle switches. Used with multiple power option bases or as an upgrade for no power option or single power option bases. Must order expandable controller (PN: I10411) with this part.</i>	605.00	
<input type="checkbox"/> ID51042	R-net Remote Joystick w/ Speed Pot & Single Toggle <i>Expandable/remote joystick with a single toggle switch for modes/profiles &amp; a speed pot. Used with multiple power option bases or as an upgrade for no power option or single power option bases. Must order expandable controller (PN: I10411) with this part.</i>	730.00	

## Select Options:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> ID51153	R-net Remote Color Screen Joystick w/ Mono Jacks <i>Must order expandable controller (PN: I10411) with this part. Comes standard w/ mono jacks.</i>	835.00	
<input type="checkbox"/> ID51149	R-net Remote Color Screen JS w/ Toggle Switches <i>Must order expandable controller (PN: I10411) with this part.</i>	875.00	
<input type="checkbox"/> ID51190	R-net Remote Color JS w/ Speed Pot & Single Toggle <i>Must order expandable controller (PN: I10411) with this part.</i>	875.00	
<input checked="" type="checkbox"/> I10411	Expandable Controller - R-net <i>Must be ordered with the expandable/remote joysticks (PN: I10410, I103544-99-0, ID51040, ID51108, ID51042, and I103814-99-0).</i>	525.00	E2377
<input checked="" type="checkbox"/> I10513	Harness for Expandable Controller <i>Must be ordered with any Expandable Controller (I10411 or I10450)</i>	445.00	E2313
<input type="checkbox"/> ID50882-A	Attendant Control for R-net Control System Only <i>On the K450, the mounting location for the Attendant Control is off of the Push Handle Bracket. If ordering an attendant control, you must also order I104576-99-0 (Push Handles - MX Seat).</i>	1,084.00	E2331
<input type="checkbox"/> I104571-99-0	Co-Pilot Attendant Control - K450	2,675.00	
<input type="checkbox"/> ID51111	R-net BT Mouse Module <i>The R-net BT Mouse Module allows for mouse control through a standard R-net Joystick or other type of input device (ex: Mini Joystick w/ Omni). It is only usable as a mouse on a computer. It CANNOT be paired to work with other Bluetooth devices such as cell phones, wireless headsets, etc. This kit includes: an R-net BT Mouse Module, a 1.5m R-net Bus Cable, and a USB Style Bluetooth Adapter for the user's computer. Note: The Bluetooth Adapter may not work with all computers.</i>	795.00	

## Seating System

## Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input checked="" type="checkbox"/> I104324-99-0	MX Power Seating System <i>Price Includes: MX Seating System, Adjustable Seat Frame, Power Seat to Floor System, Adjustable Backrest Angle, and One Piece Footplate. Weight Limit 125lbs.</i>	6,455.00	K0108

## Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104529-99-0	No Tilt <i>Weight Limit 125lbs.</i>	No Charge	
<input type="checkbox"/> I104551-99-0	45° Power Tilt <i>Weight Limit 125lbs.</i>	5,309.00	E1002

## Seating System Options

## Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input checked="" type="checkbox"/> I10585	Multiple Seat Function Control Kit for R-net <i>MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive control and/or ICS switch box. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the Seat Function Control Kit.</i>	2,875.00	E2311

## Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104136-99-0	Toggles on Seat Function Control Kit	225.00	
<input type="checkbox"/> I104135-99-0	Push Buttons on Seat Function Control Kit	No Charge	

## Select Options:

Part Number	Description	Price	HCPCS Code
-------------	-------------	-------	------------

Select Options:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I1822845	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge	
<input type="checkbox"/> I1822846	Mount ICS Switchbox onto End of Armrest Pad <i>ICS Switchbox will be mounted to the end of the armrest pad. If this item is selected, the ICS Switchbox will NOT swing away with the Joystick Module. Please indicate below if you would like it mounted onto the left or right arm pad.</i> Mount ICS Switchbox on Left Armpad _____ Mount ICS Switchbox on Right Armpad _____	No Charge	
<input type="checkbox"/> I103470-99-0	Mount ICS Switchbox on Backrest	125.00	
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Not available with VR2 Electronics. Seat function interface box which accommodates up to 8 momentary switches. Provides an additional method to accessing the seat functions. You must still order the Multiple Seat Function Control Kit if ordering this part. 1/8" mini-mono switches not included.</i>	985.00	K0108
<input type="checkbox"/> I104746-99-0	Remote Stop System - R-net, Pilot+ <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Each transmitter and receiver is a serialized matched pair, meaning they must be used together. Custom charge may be required if ordering with a low backrest height or a pediatric wheelchair. If ordering Remote Stop System with a K450, push handles must also be selected.</i>	985.00	
<input type="checkbox"/> I10515	Memory Seat Program for K450 MX <i>If this is selected, the chair will be programmed to allow the user to store and save up to 3 different memory seat positions through the ICS Switchbox. On the K450 MX, ALL seat functions are capable of recalling a saved seat position. Please refer to the K450 Owner's Manual for more information on saving and recalling a memory seat position.</i>	No Charge	

**Backrest Options**

Choose One of the Following

Part Number	Description	Price	HCPCS Code
	<b>MX STYLE BACKRESTS:</b> <i>Use this section to select MX Style Backrests. Back heights measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Please be sure to pick both a backrest shell size and a backrest cushion style (Leatherette or Black Mesh Fabric). Note: A 16"W MX Backrest is NOT available.</i>		
<input type="checkbox"/> I104321-99-0	MX Backrest - 10"W x 14"H <i>Backrest can be adjusted to either 14" or 16" high. Backrest will be mounted in the 16" high position if ordering a headrest.</i>	No Charge	
<input type="checkbox"/> I104322-99-0	MX Backrest - 12"W x 16"H <i>Backrest can be adjusted to either 16" or 18" high.</i>	No Charge	
<input type="checkbox"/> I104323-99-0	MX Backrest - 14"W x 20"H <i>Backrest can be adjusted to either 20" or 22" high.</i>	No Charge	

Choose One of the Following

Part Number	Description	Price	HCPCS Code
	Please Select Cushion Style Below:		
<input type="checkbox"/> I104309-99-0	MX Back Cushion 270x325, 10"W x 14"T - Leatherette	No Charge	
<input type="checkbox"/> I104311-99-0	MX Back Cushion 320x375, 12"W x 16"T - Leatherette	No Charge	
<input type="checkbox"/> I104313-99-0	MX Back Cushion 370x475, 14"W x 20"T - Leatherette	No Charge	
<input type="checkbox"/> I104308-99-0	MX Back Cushion 270x325, 10"W x 14"T - Black Mesh	No Charge	
<input type="checkbox"/> I104310-99-0	MX Back Cushion 320x375, 12"W x 16"T - Black Mesh	No Charge	
<input type="checkbox"/> I104312-99-0	MX Back Cushion 370x475, 14"W x 20"T - Black Mesh	No Charge	

Choose One of the Following

Part Number	Description	Price	HCPCS Code
-------------	-------------	-------	------------

## Choose One of the Following

Part Number	Description	Price	HCPCS Code
	<b>PS JR TUBULAR STYLE BACKRESTS</b>		
	<i>Use this section to select PS Jr Tubular Style Backrest Canes. Note: A 10"W Tubular Backrest is NOT available.</i>		
<input type="checkbox"/> I103717-99-0	PS Jr Back Canes, Height & Angle Adjustable <i>Must also select PS Jr Style Armrests if ordering this backrest. Choose PS Jr Style Backrest Type (Solid Back, Sling Back, or Omit Back) below. Backrest width should match seat width. Back height is adjustable from 16.5"T - 22.5"T, measured from the seat pan to the top of the back.</i>	No Charge	

## Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103719-99-0	PS Jr Solid Backrest - 12"W	No Charge	
<input type="checkbox"/> I103721-99-0	PS Jr Solid Backrest - 14"W	No Charge	
<input type="checkbox"/> I103723-99-0	PS Jr Solid Backrest - 16"W	No Charge	
Standard headrests will not mount to sling backs.			
<input type="checkbox"/> I103725-99-0	PS Jr Sling Backrest - 12"W	No Charge	
<input type="checkbox"/> I103727-99-0	PS Jr Sling Backrest - 14"W	No Charge	
<input type="checkbox"/> I103729-99-0	PS Jr Sling Backrest - 16"W	No Charge	
<input type="checkbox"/> I10449	PS Jr - Omit Backrest	No Charge	

## Seat Sizing

## Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104317-99-0	MX Seat Size 10"W x 10"D <i>Seat depth can be adjusted to 10", 11", or 12". If ordering this seat size, the lowest allowed seat position when using the seat to floor function will be approximately 3.5" from the floor to the top of the seat pan. This seat size is not compatible with the WC 19 Kit (I104553-99-0).</i>	No Charge	
<input type="checkbox"/> I104318-99-0	MX Seat Size 12"W x 12"D <i>Seat depth can be adjusted to 12", 13", or 14". If ordering this seat size, the lowest allowed seat position when using the seat to floor function will be approximately 3.5" from the floor to the top of the seat pan.</i>	No Charge	
<input type="checkbox"/> I104319-99-0	MX Seat Size 14"W x 14"D <i>Seat depth can be adjusted to 14", 15", or 16". If ordering this seat size, with MX Arms, the lowest allowed seat position when using the seat to floor function will be approximately 3.5" from the floor to the top of the seat pan. If ordering PS Jr Arms, the lowest allowed seat position will be 10".</i>	No Charge	
<input type="checkbox"/> I104320-99-0	MX Seat Size 16"W x 16"D <i>Seat depth can be adjusted to 16", 17", or 18". If ordering this seat size with MX Arms, the lowest allowed seat position when using the seat to floor function will be approximately 10" from the floor to the top of the seat pan. If ordering PS Jr Arms, the lowest allowed seat position will be 16".</i>	No Charge	

## Select Options:

Part Number	Description	Price	HCPCS Code
	<b>Seat Cushions Below Are Approximately 2" Thick:</b>		
<input type="checkbox"/> I306384-99-90-0	Leatherette Seat Cush MX - 10"W x 10"D (250x250)	475.00	
<input type="checkbox"/> I306385-99-90-0	Leatherette Seat Cush MX - 12"W x 12"D (300x300)	475.00	
<input type="checkbox"/> I306386-99-90-0	Leatherette Seat Cush MX - 14"W x 14"D (350x350)	475.00	
<input type="checkbox"/> I101818-99-0	Leatherette Seat Cush MX - 16"W x 16"D (400x400)	475.00	
<input type="checkbox"/> I104305-99-0	Black Mesh Seat Cush MX - 10"W x 10"D (250x250)	485.00	
<input type="checkbox"/> I104306-99-0	Black Mesh Seat Cush MX - 12"W x 12"D (300x300)	485.00	
<input type="checkbox"/> I104307-99-0	Black Mesh Seat Cush MX - 14"W x 14"D (350x350)	485.00	
<input type="checkbox"/> I104702-99-0	Black Mesh Seat Cush MX - 16"W x 16"D (400x400)	485.00	

### Armrest Options

Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10512	Height Adjustable Armrest Assembly MX Seat <i>If you are not selecting Tilt, you MUST order this item.</i>	284.00	E0973
<input type="checkbox"/> I10514	Height Adjustable Armrest Assembly MX Seat - NC <i>No Charge if ordering Power Tilt (I104551-99-0)</i>	No Charge	

Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104303-99-0	MX Armrest Pad Short Leatherette 8"L x 4"W	No Charge	
<input type="checkbox"/> I104304-99-0	MX Armrest Pad Long Leatherette 12.5"Lx3.5"W	No Charge	

Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104648-99-0	Height Adj Arm Assembly - PS Jr w/ UT Arm Pads <i>Newly Redesigned! Includes: 40° of angle adjustability, unlimited rotational adjustment, and 2" of offset adjustment medially and laterally. Height adjustment in 1" increments from 7" - 11" from the seat pan. Includes "C" Sized UT Armrest Pad (8"L x 3.5"W). Height Adjustable Armrest Assembly is no charge if ordering Power Tilt. Other sized armpads are be available upon request.</i>	284.00	E0973
<input type="checkbox"/> I10519	Height Adj Arm Assembly - PS Jr w/ UT Arm Pads - NC <i>No Charge when ordering Powered Tilt.</i>	No Charge	

### Legrest Options

Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104314-99-0	MX Footplate - Small (11.5" W x 8"D)	No Charge	
<input type="checkbox"/> I104315-99-0	MX Footplate - Large (14" W x 9.5"D)	No Charge	

Select Options:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103427-99-0	Footstrap Kit	165.00	
<input type="checkbox"/> I103104-99-0	Footstrap Kit - Padded	185.00	
<input type="checkbox"/> I104338-99-0	Small Anti-Skid Overlay for Footplate - K450	125.00	
<input type="checkbox"/> I104337-99-0	Large Anti-Skid Overlay for Footplate - K450	125.00	
<input type="checkbox"/> IFT202SR	BP Ankle Huggers, Med Side Release, Adj (Pair) <i>This item will be sent as a part with the wheelchair. It will not be installed.</i>	225.00	K0108

### Accessories

Select Options:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10400	Bodypoint Non-Padded Hip Belt, UniTrack - Medium	No Charge	
<input type="checkbox"/> I10404	Bodypoint Padded Hip Belt UniTrack (C300) - Medium	113.00	K0108
<input type="checkbox"/> I10623	Bodypoint Padded Hip Belt UniTrack - Extra Small	113.00	
<input type="checkbox"/> I104339-99-0	MF Headrest for MX Seat - Leatherette <i>Must also select PN I10552. Also includes a UT Bracket for lateral adjustment of the headrest. Please indicate below which size you need.</i>	263.00	E0955
<input type="checkbox"/> I104344-99-0	UT Bracket for Lateral Head Adj (10"W MX Back) <i>Allows for Lateral Adjustment of the headrest.</i>	No Charge	
<input type="checkbox"/> I104345-99-0	UT Brack for Lateral Head Adj (12-14"W MX Back) <i>Allows for Lateral Adjustment of the headrest.</i>	No Charge	
<input type="checkbox"/> IHEADADAPTER	Universal Headrest Adapter <i>For aftermarket headrests.</i>	197.00	K0108

**K450 MX****MSRP Price List***Select Options:*

<b>Part Number</b>	<b>Description</b>	<b>Price</b>	<b>HCPCS Code</b>
<input type="checkbox"/> I10552	Adj Removable - Headrest Hardware (1 Each)	212.00	E1028
<input type="checkbox"/> I103077-99-0	UT Hip Suppt (B) w/ Mesh Cover - Low (Pair) <i>4.5"W x 3.5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	225.00	E0956
<input type="checkbox"/> I103079-99-0	UT Hip Suppt (A) w/ Mesh Cover - Low (Pair) <i>6"W x 5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	225.00	E0956
<input type="checkbox"/> I10550	Adj Removable - Thigh Support Hardware (Pair)	424.00	E1028
<input type="checkbox"/> I104575-99-0	UT Trunk Suppt (B) for MX Seat - 3.5"T x 4.5"D <i>Must also select the Swing Away Trunk Support Hardware (I10551).</i>	285.00	E0956
<input type="checkbox"/> I104574-99-0	UT Trunk Suppt (F) for MX Seat - 3.5"T x 8.5"D <i>Must also select the Swing Away Trunk Support Hardware (I10551).</i>	285.00	E0956
<input type="checkbox"/> I10551	Adj Removable - Trunk Support Hardware (Pair)	424.00	E1028
<input type="checkbox"/> I104576-99-0	Push Handles - K450 <i>Mounted off the chassis. Height and angle adjustable. Must be selected if ordering bag clips or an attendant control.</i>	345.00	
<input type="checkbox"/> I10598	Medical Necessity Bag Clips - Carabiner Style <i>Push Handles (I104576-99-0) must also be selected in order to mount bag clips.</i>	75.00	
<input type="checkbox"/> I100598-99-0	Chest Harness <i>Chest harness kit includes a standard positioning belt.</i>	371.00	E0960
<input type="checkbox"/> ISH212	BP Shoulder Harness Medium, Front Pull <i>This item will be sent as a part with the wheelchair. It will not be installed.</i>	196.00	E0960
<input type="checkbox"/> I100578-99-0	Upper Extremity Support Small, Angle Adjustable <i>Must also select PN I10556.</i>	238.00	E0950
<input type="checkbox"/> I103426-99-0	Upper Extremity Support Angle Adjust - UniTrack <i>Side mounted tray that attaches to the UniTrack Accessory Rail. Dimensions are: 19.5"W x 13"D. Can only be ordered with chairs that offer the UniTrack System. Must also select PN I10556.</i>	296.00	E0950
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each)	212.00	E1028

**Alternative Drive Systems***Choose One of the Following*

<b>Part Number</b>	<b>Description</b>	<b>Price</b>	<b>HCPCS Code</b>
	<b>R-net Alternative Drive Input Devices</b> <i>R-net Alternative Drive Input Devices take the place of, or may be used in addition to, the standard R-net Joystick. They are only usable on wheelchairs that have R-net electronics. You MUST select one of each of the following items: (1) Input Device, (2) Expandable Controller, (3) Adjustable/Removable Alt Drive Hardware, AND (4) Omni or Input Module.</i>		
	<b>(1) INPUT DEVICES FOR ALTERNATIVE DRIVE CONTROL</b>		
<input type="checkbox"/> I10558	Gimble Joystick set up for Armrest Mount	1,388.00	E2373
<input type="checkbox"/> I10559	Gimble Joystick set up for Midline Mount	1,388.00	E2373
<input type="checkbox"/> I10560	Gimble Joystick set up for Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	1,388.00	E2373
<input type="checkbox"/> I10561	Gimble Joystick w/ Haims Harness - Small <i>For neck width up to 5".</i>	1,555.00	E2373
<input type="checkbox"/> I10562	Gimble Joystick w/ Haims Harness - Large <i>For neck width up to 7".</i>	1,555.00	E2373
<input type="checkbox"/> I10564	Mini Joystick set up for Armrest Mount	2,035.00	E2312
<input type="checkbox"/> I10565	Mini Joystick set up for Midline Mount	2,035.00	E2312
<input type="checkbox"/> I10567	Mini Joystick w/ Haims Harness - Small <i>For neck width up to 5".</i>	2,135.00	E2312
<input type="checkbox"/> I10568	Mini Joystick w/ Haims Harness - Large <i>For neck width up to 7".</i>	2,135.00	E2312

**K450 MX****MSRP Price List**

Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10569	Mini Joystick w/ Power Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	5,174.00	E2399
<input type="checkbox"/> I1821750	Tablet Control Pad <i>The Tablet Control is a unique drive input which allows a user to control a power wheelchair in a very similar manner to a mouse being controlled by a touchpad on a laptop computer. The dimensions of the Tablet Control Pad are: 12"W x 9.5"D x 1.5"H. Kit also includes: Buddy Button Switch, and Tablet Control Mounting Bracket.</i>	3,614.00	E2328
<input type="checkbox"/> I10570	Head Array - Adult 5 Switch HDR-A5	3,614.00	E2330
<input type="checkbox"/> I10571	Head Array - Adult 3 Switch HDR-A4 <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	3,614.00	E2330
<input type="checkbox"/> I10572	Head Array - Pediatric 5 Switch HDR-C5-RR	3,614.00	E2330
<input type="checkbox"/> I10573	Sip & Puff System <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	1,409.00	E2325
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	1,929.00	E2325,E2326

Select Options:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10076	Right Side Mount For Joystick	No Charge	
<input type="checkbox"/> I10077	Left Side Mount For Joystick	No Charge	
<input type="checkbox"/> I10074	Right Side Mount For Omni	No Charge	
<input type="checkbox"/> I10075	Left Side Mount For Omni	No Charge	
<input type="checkbox"/> I10078	Mount Seat Control Panel on Right	No Charge	
<input type="checkbox"/> I10079	Mount Seat Control Panel on Left	No Charge	
<input type="checkbox"/> I10080	Mount Seat Control Panel on Backrest	No Charge	

**Alt Drive Systems Setup**

Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10576	(2) EXPANDABLE CONTROLLER (REQUIRED) Expandable Controller - Alternative Drive <i>Must also select Harness for Expandable Controller (I10513), located in "Joystick Options".</i>	525.00	E2377

Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10577	(3) ADJUSTABLE REMOVABLE HARDWARE (REQUIRED) Adjustable Removable Hardware - Alt Drive	212.00	E1028

Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> ID51154IR	(4) OMNI OR INPUT MODULE (REQUIRED) <i>Cannot select Input Module with Sip-n-Puff Input Devices (PN I10573 or I10574). Please specify which type of alternative drive input you will be using if not ordering the input device from Permobil. Input Device: _____</i>	1,470.00	

**Alt Drive Accessories**

Select Options:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> IMICROSWITCH	Microlite Switch	148.00	
<input type="checkbox"/> IBUDDYBUTTON	Buddy Button	111.00	K0108

*Select Options:*

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I5808	Mini Cup Switch	108.00	
<input type="checkbox"/> IZTS-P	Proximity Switch	515.00	
<input type="checkbox"/> IMINIJOY5SW	Mini Joystick 5 Switch OC	593.00	
<input type="checkbox"/> I5018	Wafer Board w/5 Switches	497.00	

Notes & Additional Instructions:

**Chair Order Policy:**

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, Inc. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

**Order Acknowledgement:**

I, \_\_\_\_\_, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: \_\_\_\_\_

Printed Name: \_\_\_\_\_

Title: \_\_\_\_\_

Date: \_\_\_\_\_

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. MID PROFILE® Cushion is a registered trademark of the The ROHO Group. *The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.*