

C300 PS



*** Indicates a Required Field**

Dealer Information



* **Contact:** _____

Dealer Code: _____

* **Dealer Name:** _____

Address: _____

* **City:** _____

* **State/Zip:** _____

* **Phone#:** _____

* **Fax#:** _____

PO#: _____

Email Address: _____

Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

* **First Name:** _____

* **Last Name:** _____

Diagnosis: _____

Funding Source: _____

Client Age: _____

Client Measurements

* **Weight:** _____

* **Height:** _____

A) Top of Shoulders: _____

B) Chest Depth: _____

C) Chest Width: _____

D) Seat Depth: _____


E) Top of Head: _____

F) Elbow to Hand: _____

G) Seat Pan to Elbow: _____

H) Hip Width: _____

* **I) Knee to Foot:** _____



Three column pricing on order/quote:

Display HCPCS Codes on order/quote:

Please send order/quote to fax#: **(800) 231-3256**

Email to: sales@permobilus.com

Permobil Inc.
6961 Eastgate Blvd.
Lebanon, TN 37090
Tel: (800) 736-0925
Fax: (800) 231-3256
www.permobil.com

Prices effective **August 9, 2010.**

C300 PS

Chair Model and Colors

CHOOSE BASE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10594	C300 Base PS - C300 PS0 <i>Price Includes: FWD C300 Power Base Set Up For No Power Seat Functions, Shock Absorbing Suspension System, VR2 Controller, Tie-Down Hardware for Strap Systems, Anti-Tippers, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below.</i>	6,645.00	K0848
<input type="checkbox"/> I103599-99-0	C300 Base PS - C300 PS1 <i>Price Includes: FWD C300 Power Base Set Up For One Power Seat Function, Shock Absorbing Suspension System, VR2 Controller, Tie-Down Hardware for Strap Systems, Anti-Tippers, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below.</i>	6,645.00	K0856
<input type="checkbox"/> I10637	6.5mph High Speed Upgrade Package - C300 <i>Price Includes: ESP Module, High Speed 6.5mph Motor Package, and 120A Power Module. Available with R-net Electronics only. Available for shipment after September 7, 2010.</i>	1,895.00	

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103499-99-0	Tri-Spoke Split Rim w/ Flat Free Drive Tire 3x8"	No Charge	
<input type="checkbox"/> I103498-99-0	Tri-Spoke Split Rim w/ Pneumatic Tires - 3x8"	No Charge	

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103576-99-0	C300/K300 Shroud - Onyx Black	No Charge	
<input type="checkbox"/> I101381-99-0	C300/K300 Shroud - Vintage Wine	No Charge	
<input type="checkbox"/> I101383-99-0	C300/K300 Shroud - Crystal Blue	No Charge	
<input type="checkbox"/> I101384-99-0	C300/K300 Shroud - Metallic Charcoal Gray	No Charge	
<input type="checkbox"/> I103577-99-0	C300/K300 Shroud - Platinum Silver	No Charge	
<input type="checkbox"/> I101385-99-0	C300/K300 Shroud - Green Metallic	No Charge	
<input type="checkbox"/> I103575-99-0	C300/K300 Shroud - Titan Blue	No Charge	
<input type="checkbox"/> I103578-99-0	C300/K300 Shroud - Pop Star Pink	No Charge	
<input type="checkbox"/> I103574-99-0	C300/K300 Shroud - Chili Pepper Red	No Charge	

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10451	Match Hub Cap to Shroud Color <i>If this is selected, the Tri-Spoke Hub Cap will be matched to the selected shroud color.</i>	No Charge	
<input type="checkbox"/> I103537-99-0	Hub Cap Platinum Silver <i>If this is selected, the Tri-Spoke Hub Cap will be Platinum Silver.</i>	No Charge	

Base Options

MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I101994-99-0	Seat Tube Fixed	No Charge	
<input type="checkbox"/> I103779-99-0	Power Adjustable Seat Height - 8" Travel	2,445.00	E2300

SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> IM34 SLDG	Batteries - Group 34, 12V 60 Amp Hours, Installed	785.00	K0108
<input type="checkbox"/> I602307-99-0	P+ PP1A Programmer - VSI, VR2 or Pilot+ Controllers	605.00	
<input type="checkbox"/> I102768-99-0	Fixed Vent Shelf for C300	402.00	E1029
<input type="checkbox"/> I0036-7002A	Enhanced Steering Perform Unit (ESP), R-net <i>If selecting this part, you MUST also pick: PN I0036-7002A - ESP Unit. Fits an LP10 Sized Vent & a U1 Battery. Must be ordered with an expandable controller (ex. R-net Electronics). MUST be selected if ordering any ventilator tray (I102769-99-0, I103241-99-0, or I102768-99-0). Recommended item with switched drive inputs for improved tracking.</i>	1,170.00	K0108

C300 PS

Base Options

SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103564-99-0	Light Kit for C300 w/ R-net Electronics <i>Shrouds with lights cannot be ordered with VR2 Electronics.</i>	972.00	
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	443.00	

Permolock

SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104680-99-0	Permolock C3 Locking Base C300, C350, K300, K450 <i>Kit Includes: Permolock C3 Locking Base, Owner's Manual, Manual Release Handle, and Vehicle Mounted Release Button. Does NOT include Dual Locking Pin Kit (I10608). Locking Base must be installed by a NMEDA Dealer. Please provide the NMEDA Dealer and Address where the locking Base will ship to: NAME: _____ ADDR: _____</i>	1,806.00	
<input type="checkbox"/> I10608	Permolock C3 Dual Locking Pin Kit - C300/K300 <i>This kit includes two spring loaded pins that will be mounted underneath the base of the wheelchair. These pins are spring loaded so they push into the wheelchair base when an obstacle is encountered. The kit also includes spacers that may be used if the Permolock C3 Base is being mounted in a vehicle with a sloped floor. The spacers are only needed for a C300 or K300.</i>	605.00	

Joystick Options

CHOOSE JOYSTICK TYPE:

Part Number	Description	Price	HCPCS Code
	PLEASE NOTE: <i>For all remote joystick modules shown below, you must also select an expandable controller. For additional input devices, refer to Alternative Drive Systems Section.</i>		
<input type="checkbox"/> I10423	VR2 Basic Joystick - Standard	No Charge	
<input type="checkbox"/> I103544-99-0	R-net Remote Joystick w/ Mono Jacks	539.00	
<input type="checkbox"/> ID51040	R-net Remote Joystick w/ Toggle Switches	616.00	
<input type="checkbox"/> ID51042	R-net Remote Joystick w/ Speed Pot & Single Toggle	743.00	
<input type="checkbox"/> ID51153	R-net Remote Color Screen Joystick w/ Mono Jacks	850.00	
<input type="checkbox"/> ID51149	R-net Remote Color Screen JS w/ Toggle Switches	890.00	
<input type="checkbox"/> ID51190	R-net Remote Color JS w/ Speed Pot & Single Toggle	890.00	

CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard.	No Charge	
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	96.00	E2324
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	141.00	E2323
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	96.00	E2323
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	129.00	E2323
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	129.00	E2323
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	129.00	E2323
<input type="checkbox"/> IPC101B	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	129.00	E2323
<input type="checkbox"/> IPC102B	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	129.00	E2323
<input type="checkbox"/> IPC107B	Bodypoint J/S Handle Dome Shaped, Rubber	129.00	E2323

CHOOSE JOYSTICK MOUNT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10424	Fixed Joystick Mount - Right - VR2	No Charge	
<input type="checkbox"/> I10425	Fixed Joystick Mount - Left VR2	No Charge	
<input type="checkbox"/> I10412	Joystick Mount Right - Fixed, R-net Only	No Charge	
<input type="checkbox"/> I10413	Joystick Mount Left - Fixed, R-net Only	No Charge	

C300 PS

Joystick Options

CHOOSE JOYSTICK MOUNT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I102537-99-0	Joystick Mount Right - Swing Away, R-net & VR2	297.00	E1028
<input type="checkbox"/> I102538-99-0	Joystick Mount Left - Swing Away, R-net & VR2	297.00	E1028
<input type="checkbox"/> I103820-99-0	Retractable Joystick Mnt Right - R-net, VR2 Only	495.00	E1028
<input type="checkbox"/> I103819-99-0	Retractable Joystick Mnt Left - R-net, VR2 Only	495.00	E1028
<input type="checkbox"/> I101260-99-0	Midline Control Mount - Pilot+ & R-net Only	268.00	E1028

Please indicate which armrest to mount midline below.

Mount Midline Joystick on Left Armrest: _____

Mount Midline Joystick on Right Armrest: _____

SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10411	Expandable Controller - R-net <i>Must be ordered with the expandable/remote joysticks (PN: I10410, I103544-99-0, ID51040, ID51108, ID51042, and I103814-99-0) OR when ordering Alternative Drive Packages.</i>	545.00	E2377
<input type="checkbox"/> I10513	Harness for Expandable Controller <i>Must be ordered with any Expandable Controller (I10411 or I10450)</i>	475.00	E2313
<input type="checkbox"/> ID50882	Attendant Control for R-net Control System Only	1,103.00	E2331

Not available with VR2 Controllers.

<input type="checkbox"/> ID51111	R-net BT Mouse Module <i>The R-net BT Mouse Module allows for mouse control through a standard R-net Joystick or other type of input device (ex: Mini Joystick w/ Omni). It is only usable as a mouse on a computer. It CANNOT be paired to work with other Bluetooth devices such as cell phones, wireless headsets, etc. This kit includes: an R-net BT Mouse Module, a 1.5m R-net Bus Cable, and a USB Style Bluetooth Adapter for the user's computer. Note: The Bluetooth Adapter may not work with all computers.</i>	809.00	
----------------------------------	--	--------	--

Seating System

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10281	PS Seating System Only - C300/C350 <i>Includes: Seat Frame, Fixed Back Frame, Fixed Legrests and Footplates, Fixed Arm Supports and Standard Armrest Pads. Weight Limit 300 lbs.</i>	No Charge	
<input type="checkbox"/> I10142	PS Seating System With 45 Degree CG Power Tilt <i>Price Includes: Seat Frame, All Seat Sizes, Angle Adjustable Back Canes, Height Adjustable Armrests, Standard Armrest Pads, Swing Away Legrests and Footplates. Weight Limit 300lbs.</i>	4,865.00	E1002

Seating System Options

CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10289	Single Seat Function Control Switch - LiberT <i>This device is a single toggle switch used to control only one seat function. If control of the seat function through the drive control is needed, then the drive electronics must be upgraded to Pilot+ or R-net Electronics, and a Single Seat Function Control Kit (Part Number I10129), must be selected.</i>	No Charge	
<input type="checkbox"/> I10419	Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control and/or ICS Switchbox. Order with R-net Controller or remote joystick options. Please indicate below if you want push buttons or toggles on the Seat Function Control Kit.</i>	1,495.00	E2310
<input type="checkbox"/> I10585	Multiple Seat Function Control Kit for R-net	2,925.00	E2311

MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive control and/or ICS switch box. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the Seat Function Control Kit.

C300 PS

Seating System Options

CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I102798-99-0	Push Buttons on Seat Function Control Kit	No Charge	
<input type="checkbox"/> I102799-99-0	Toggle Switches on Seat Function Control Kit	249.00	

CHOOSE MOUNTING POSITION OF SEAT FUNCTION CONTROL KIT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge	
<input type="checkbox"/> I103470-99-0	Mount ICS Switchbox on Backrest	127.00	

SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Not available with VR2 Electronics. Seat function interface box which accommodates up to 8 momentary switches. Provides an additional method to accessing the seat functions. You must still order the Multiple Seat Function Control Kit if ordering this part. 1/8" mini-mono switches not included.</i>	1,002.00	K0108
<input type="checkbox"/> I104746-99-0	Remote Stop System - R-net, Pilot+ <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Each transmitter and receiver is a serialized matched pair, meaning they must be used together. Custom charge may be required if ordering with a low backrest height or a pediatric wheelchair. If ordering Remote Stop System with a K450, push handles must also be selected.</i>	1,002.00	
<input type="checkbox"/> I104769-99-0	Remote Stop System VSI, VR2, Safegate <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Each transmitter and receiver is a serialized matched pair, meaning they must be used together. Custom charge may be required if ordering with a low backrest height or a pediatric wheelchair. If ordering Remote Stop System with a K450, push handles must also be selected.</i>	1,002.00	

Backrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10144	Sling Upholstery Back Support - Standard <i>Width will be the same as seat width.</i>	No Charge	
<input type="checkbox"/> I10145	Curved Metal Back Support 22" High	417.00	
<input type="checkbox"/> I10146	Curved Metal Back Support 22" High - No Charge <i>No charge if ordering tilt.</i>	No Charge	
<input type="checkbox"/> I10454	Omit Backrest	No Charge	

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10148	Angle Adjustable Back Canes <i>Angle adjustable from 85 to 115 degrees.</i>	129.00	
<input type="checkbox"/> I10149	Angle Adjustable Back Canes - No Charge <i>No charge if ordering tilt.</i>	No Charge	

Seat Sizing

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10150A	Seat Width 16"	No Charge	
<input type="checkbox"/> I10151A	Seat Width 18"	No Charge	
<input type="checkbox"/> I10219A	Seat Width 20"	No Charge	
<input type="checkbox"/> I10220A	Seat Width 22"	No Charge	

C300 PS

Seat Sizing

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10156	Seat Depth 16"	No Charge	
<input type="checkbox"/> I10157	Seat Depth 17"	No Charge	
<input type="checkbox"/> I10158	Seat Depth 18"	No Charge	
<input type="checkbox"/> I10159	Seat Depth 19"	No Charge	
<input type="checkbox"/> I10221	Seat Depth 20"	No Charge	
<input type="checkbox"/> I10222	Seat Depth 21"	No Charge	
<input type="checkbox"/> I10223	Seat Depth 22"	No Charge	

Armrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10163	Height Adjustable Armrest Assembly PS <i>Adjustable in height, with a flip-back design.</i>	289.00	E0973
<input type="checkbox"/> I10230	Height Adjustable Armrest Assembly PS - No Charge <i>Adjustable in height, with a flip-back design. No charge if ordering power tilt.</i>	No Charge	

SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10290	PS Armrest Pads - 10.5"L x 2"W (Pair)	No Charge	
<input type="checkbox"/> I10291	PS Armrest Pads - 14"L x 2"W (Pair)	No Charge	
<input type="checkbox"/> I10292	Waterfall Armrest Pad for PS - 13.5"L x 3"W (Pair)	80.00	
<input type="checkbox"/> I10168	Omit Armrest Pads	No Charge	
<input type="checkbox"/> I10166	Corpus Armrest Pad Adapter Brackets <i>Must be picked if ordering Corpus armrest pads.</i>	51.00	
<input type="checkbox"/> I10208	Corpus Armrest Pad 4"x10" Left Leatherette	94.00	
<input type="checkbox"/> I10209	Corpus Armrest Pad 4"x12" Left Leatherette	94.00	
<input type="checkbox"/> I10210	Corpus Armrest Pad 4"x14" Left Leatherette	94.00	
<input type="checkbox"/> I10211	Corpus Armrest Pad 4"x16" Left Leatherette	94.00	
<input type="checkbox"/> I10212	Corpus Armrest Pad 4"x10" Right Leatherette	94.00	
<input type="checkbox"/> I10213	Corpus Armrest Pad 4"x12" Right Leatherette	94.00	
<input type="checkbox"/> I10214	Corpus Armrest Pad 4"x14" Right Leatherette	94.00	
<input type="checkbox"/> I10215	Corpus Armrest Pad 4"x16" Right Leatherette	94.00	
<input type="checkbox"/> I10169	Armrest Side Support Right Inner <i>Will limit the armrest from flipping back. For Corpus Armrest Pads only.</i>	142.00	
<input type="checkbox"/> I10170	Armrest Side Support Left Inner <i>Will limit the armrests from flipping back. For Corpus Armrest Pads only.</i>	142.00	
<input type="checkbox"/> I101261-99-0	Armrest Side Support Right Outer <i>For Corpus Armrest Pads only.</i>	144.00	
<input type="checkbox"/> I101262-99-0	Armrest Side Support Left Outer <i>For Corpus Armrest Pads only.</i>	144.00	

Legrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10172	70 Degree Swing Away Legrests (Pair)	202.00	
<input type="checkbox"/> I10231	70 Degree Swing Away Legrests (Pair) <i>No charge if ordering Power Tilt.</i>	No Charge	

C300 PS

Legrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10316	Aeron Legrest Assembly for PS Seat <i>90° Legrest Hanger with Manual Adjustment Legrest Elevation. Price includes Corpus Style, Angle Adjustable, Two Piece Footplates. Range of adjustment from the seat pan to the top of the footplate is: 12 - 17" (with calf pads) and 6.5 - 17" (without calf pads).</i>	321.00	

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10173	Footplates (Pair), Standard Small	No Charge	
<input type="checkbox"/> I10174	Footplates (Pair), Standard Large <i>Not available with 16" seat width.</i>	No Charge	
<input type="checkbox"/> I10175	Small Footplates, Multi - Axis Adjustable, Aluminum	198.00	K0108
<input type="checkbox"/> I10176	Large Footplates, Multi - Axis Adjustable, Aluminum <i>Cannot be ordered with 16" or 18" seat widths.</i>	198.00	K0108

SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104684-99-0	UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad. Recommended pad to be used with lower leg lengths of 12" or less. Bracket designed for use with Aeron Center Mount Legrest Assembly. Cannot be ordered with Swing Away Legs (I10586 & I10587).</i>	No Charge	
<input type="checkbox"/> I104686-99-0	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad. Bracket designed for use with Aeron Center Mount Legrest Assembly. Cannot be ordered with Swing Away Legs (I10586 & I10587).</i>	No Charge	
<input type="checkbox"/> I10177	Heel Loops <i>Not available with Aeron Legrest Assembly. (I10316).</i>	92.00	E0951
<input type="checkbox"/> I10197	Calf Supports Leatherette <i>Not available with lower leg lengths of 10" or less.</i>	No Charge	

Accessories

SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10401	Bodypoint Non-Padded Hip Belt, PS - Medium	No Charge	
<input type="checkbox"/> I10406	Bodypoint Padded Hip Belt PS (C300) - Medium	126.00	K0108
<input type="checkbox"/> I10403	Bodypoint Non-Padded Hip Belt, PS - Large	95.00	K0108
<input type="checkbox"/> I10409	Bodypoint Padded Hip Belt PS - Large	158.00	K0108
<input type="checkbox"/> I10199	Headrest-Removable Leatherette <i>Must also select PN I10552.</i>	268.00	E0955
<input type="checkbox"/> I103460-99-0	UT Ergo Headrest w/ Black Fabric Cover <i>Must also select PN I10552.</i>	297.00	E0955
<input type="checkbox"/> I10180	Universal Headrest Adapter <i>Will be sent as part if ordering sling back upholstery.</i>	205.00	K0108
<input type="checkbox"/> I10552	Adj Removable - Headrest Hardware (1 Each)	216.00	E1028
<input type="checkbox"/> I102059-99-0	Medical Necessities Bag <i>Includes Medical Necessity Bag Clips (Carabiner Style).</i>	229.00	
<input type="checkbox"/> I10179	Medical Necessity Bag Clips - Carabiner Style <i>Will be sent as part if ordering sling back upholstery.</i>	76.00	
<input type="checkbox"/> I105144-99-0	Bodypoint Monoflex Large Chest Support Belt Kit <i>21.5"L x 4.5"W Pad Size. Item will be shipped as parts.</i>	145.00	E0960
<input type="checkbox"/> I105143-99-0	Bodypoint Monoflex Medium Chest Support Belt Kit <i>18.25"L x 4"W Pad Size. Item will be shipped as parts.</i>	145.00	E0960
<input type="checkbox"/> I105145-99-0	Bodypoint Monoflex X-Lrge Chest Support Belt Kit <i>24.5"L x 5"W Pad Size. Item will be shipped as parts.</i>	145.00	E0960

C300 PS

Accessories

SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I100578-99-0	Upper Extremity Support Small, Angle Adjustable <i>Must also select PN I10556.</i>	242.00	E0950
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each)	216.00	E1028
<input type="checkbox"/> I101964-99-0	Oxygen Holder for PS or Corpus Seat <i>Cannot be ordered with a Corpus Tubular Backrest. (I10111 or I10112). Custom charge may be applied if ordering in combination with accessories such as medical necessities bags, attendant controls, push handles, back pack clips, etc. This item is mounted to the Ergonomic Backrest. Customization may be required if ordering on a Street.</i>	228.00	E2208
<input type="checkbox"/> I101956-99-0	Crutch Holder for PS or Corpus Seat <i>Cannot be ordered with a Corpus Tubular Backrest. (I10111 or I10112). Accessories such as medical necessities bags and push handles may need to be sent as parts if ordering a crutch holder. This item is mounted to the ergonomic backrest. Customization may be required if ordering on a Street.</i>	228.00	E2207

Alternative Drive Systems

(1) CHOOSE ALTERNATIVE DRIVE INPUT:

Part Number	Description	Price	HCPCS Code
	R-net Alternative Drive Input Devices <i>R-net Alternative Drive Input Devices take the place of, or may be used in addition to, the standard R-net Joystick. They are only usable on wheelchairs that have R-net electronics. You MUST select one of each of the following items: (1) Input Device, (2) Adjustable/Removable Alt Drive Hardware, AND (3) Omni or Input Module.</i>		
<input type="checkbox"/> I10558	Gimble Joystick set up for Armrest Mount	1,412.00	E2373
<input type="checkbox"/> I10559	Gimble Joystick set up for Midline Mount	1,412.00	E2373
<input type="checkbox"/> I10560	Gimble Joystick set up for Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	1,412.00	E2373
<input type="checkbox"/> I10561	Gimble Joystick w/ Haims Harness - Small <i>For neck width up to 5".</i>	1,582.00	E2373
<input type="checkbox"/> I10562	Gimble Joystick w/ Haims Harness - Large <i>For neck width up to 7".</i>	1,582.00	E2373
<input type="checkbox"/> I10563	Gimble Joystick w/ Power Swing Away Chin Boom <i>Custom charge may be applied if ordering smaller backrest sizes and heights.</i>	5,265.00	E2373
<input type="checkbox"/> I10564	Mini Joystick set up for Armrest Mount	2,685.00	E2312
<input type="checkbox"/> I10565	Mini Joystick set up for Midline Mount	2,685.00	E2312
<input type="checkbox"/> I10566	Mini Joystick set up for Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	2,685.00	E2312
<input type="checkbox"/> I10567	Mini Joystick w/ Haims Harness - Small <i>For neck width up to 5".</i>	2,745.00	E2312
<input type="checkbox"/> I10568	Mini Joystick w/ Haims Harness - Large <i>For neck width up to 7".</i>	2,745.00	E2312
<input type="checkbox"/> I10569	Mini Joystick w/ Power Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	5,265.00	E2312
<input type="checkbox"/> I1821750	Tablet Control Pad <i>The Tablet Control is a unique drive input which allows a user to control a power wheelchair in a very similar manner to a mouse being controlled by a touchpad on a laptop computer. The dimensions of the Tablet Control Pad are: 12"W x 9.5"D x 1.5"H. Kit also includes: Buddy Button Switch, and Tablet Control Mounting Bracket.</i>	4,375.00	E2328
<input type="checkbox"/> I10570	Head Array - Adult 5 Switch HDR-A5	4,375.00	E2330
<input type="checkbox"/> I10571	Head Array - Adult 3 Switch HDR-A4 <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	4,175.00	E2330
<input type="checkbox"/> I10572	Head Array - Pediatric 5 Switch HDR-C5-RR	4,375.00	E2330
<input type="checkbox"/> I10573	Sip & Puff System <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	1,475.00	E2325

C300 PS

Alternative Drive Systems

(1) CHOOSE ALTERNATIVE DRIVE INPUT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	1,963.00	E2325,E2326

SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10076	Right Side Mount For Joystick	No Charge	
<input type="checkbox"/> I10077	Left Side Mount For Joystick	No Charge	
<input type="checkbox"/> I10074	Right Side Mount For Omni	No Charge	
<input type="checkbox"/> I10075	Left Side Mount For Omni	No Charge	
<input type="checkbox"/> I10078	Mount Seat Control Panel on Right	No Charge	
<input type="checkbox"/> I10079	Mount Seat Control Panel on Left	No Charge	
<input type="checkbox"/> I10080	Mount Seat Control Panel on Backrest	No Charge	

Alt Drive Systems Setup

(2) CHOOSE ADJ ALT DRIVE HARDWARE (REQUIRED):

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10577	Adjustable Removable Hardware - Alt Drive	216.00	E1028

(3) CHOOSE OMNI OR INPUT MODULE (REQUIRED):

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> ID51154IR	Omni Control Module for R-net w/ Infrared <i>Please specify which type of alternative drive input you will be using if not ordering the input device from Permobil. Input Device: _____</i>	1,496.00	
<input type="checkbox"/> I10581	Input Module (IOM) for R-net Control System <i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head Array or RIM control.</i>	1,012.00	E2351

Alt Drive Accessories

SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> IMICROSWITCH	Microlite Switch Qty: _____	151.00	
<input type="checkbox"/> IBUDDYBUTTON	Buddy Button Qty: _____	113.00	K0108
<input type="checkbox"/> I5808	Mini Cup Switch Qty: _____	110.00	
<input type="checkbox"/> IZTS-P	Proximity Switch Qty: _____	524.00	
<input type="checkbox"/> IMINIJOY5SW	Mini Joystick 5 Switch OC	603.00	
<input type="checkbox"/> I5018	Wafer Board w/5 Switches	506.00	

Notes & Additional Instructions:

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, Inc. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. MID PROFILE® Cushion is a registered trademark of the The ROHO Group. *The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.*