

# C300 Corpus



**\* Indicates a Required Field**

## Dealer Information

\* **Contact:** \_\_\_\_\_

Dealer Code: \_\_\_\_\_

\* **Dealer Name:** \_\_\_\_\_

Address: \_\_\_\_\_

\* **City:** \_\_\_\_\_

\* **State/Zip:** \_\_\_\_\_

\* **Phone#:** \_\_\_\_\_

\* **Fax#:** \_\_\_\_\_

PO#: \_\_\_\_\_

Email Address: \_\_\_\_\_

## Client Information

*Permobil recommends that the client is evaluated by a certified rehab specialist.*

\* **First Name:** \_\_\_\_\_

\* **Last Name:** \_\_\_\_\_

Diagnosis: \_\_\_\_\_

Funding Source: \_\_\_\_\_

Client Age: \_\_\_\_\_

## Client Measurements

\* **Weight:** \_\_\_\_\_

\* **Height:** \_\_\_\_\_

A) Top of Shoulders: \_\_\_\_\_

B) Chest Depth: \_\_\_\_\_

C) Chest Width: \_\_\_\_\_

D) Seat Depth: \_\_\_\_\_


E) Top of Head: \_\_\_\_\_

F) Elbow to Hand: \_\_\_\_\_

G) Seat Pan to Elbow: \_\_\_\_\_

H) Hip Width: \_\_\_\_\_

\* **I) Knee to Foot:** \_\_\_\_\_



Three column pricing on order/quote:

Display HCPCS Codes on order/quote:

Please send order/quote to fax#: **(800) 231-3256**

Email to: [sales@permobilus.com](mailto:sales@permobilus.com)

**Permobil Inc.**  
6961 Eastgate Blvd.  
Lebanon, TN 37090  
Tel: (800) 736-0925  
Fax: (800) 231-3256  
[www.permobil.com](http://www.permobil.com)

Prices effective **August 9, 2010.**

# C300 Corpus

## Chair Model and Colors

### CHOOSE BASE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10592	C300 Base Corpus R-net - C300 PS0 <i>Price Includes: FWD C300 Power Base Set Up For No Power Seat Functions, Shock Absorbing Suspension System, R-net Controller, Tie-Down Hardware for Strap Systems, Anti-Tippers, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below.</i>	6,645.00	K0848
<input type="checkbox"/> I10593	C300 Corpus R-net - C300 PS1 <i>Price Includes: FWD C300 Power Base Set Up For One Power Seat Function, Shock Absorbing Suspension System, R-net Controller, Tie Down Hardware for Strap Systems, Anti-Tippers, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below.</i>	6,645.00	K0856
<input type="checkbox"/> I103260-99-0	C300 Base Corpus R-net - C300 PS2 <i>Price Includes: FWD C300 Power Base Set Up For Multiple Power Seat Functions, Shock Absorbing Suspension System, Tie-Down Hardware for Strap Systems, Anti-Tippers, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below.</i>	6,645.00	K0861
<input type="checkbox"/> I10637	6.5mph High Speed Upgrade Package - C300 <i>Price Includes: ESP Module, High Speed 6.5mph Motor Package, and 120A Power Module. Available with R-net Electronics only. Available for shipment after September 7, 2010.</i>	1,895.00	

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103499-99-0	Tri-Spoke Split Rim w/ Flat Free Drive Tire 3x8"	No Charge	
<input type="checkbox"/> I103498-99-0	Tri-Spoke Split Rim w/ Pneumatic Tires - 3x8"	No Charge	

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103576-99-0	C300/K300 Shroud - Onyx Black	No Charge	
<input type="checkbox"/> I101381-99-0	C300/K300 Shroud - Vintage Wine	No Charge	
<input type="checkbox"/> I101383-99-0	C300/K300 Shroud - Crystal Blue	No Charge	
<input type="checkbox"/> I101384-99-0	C300/K300 Shroud - Metallic Charcoal Gray	No Charge	
<input type="checkbox"/> I103577-99-0	C300/K300 Shroud - Platinum Silver	No Charge	
<input type="checkbox"/> I101385-99-0	C300/K300 Shroud - Green Metallic	No Charge	
<input type="checkbox"/> I103575-99-0	C300/K300 Shroud - Titan Blue	No Charge	
<input type="checkbox"/> I103578-99-0	C300/K300 Shroud - Pop Star Pink	No Charge	
<input type="checkbox"/> I103574-99-0	C300/K300 Shroud - Chili Pepper Red	No Charge	

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10451	Match Hub Cap to Shroud Color <i>If this is selected, the Tri-Spoke Hub Cap will be matched to the selected shroud color.</i>	No Charge	
<input type="checkbox"/> I103537-99-0	Hub Cap Platinum Silver <i>If this is selected, the Tri-Spoke Hub Cap will be Platinum Silver.</i>	No Charge	

## Base Options

### MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I101994-99-0	Seat Tube Fixed	No Charge	
<input type="checkbox"/> I103779-99-0	Power Adjustable Seat Height - 8" Travel	2,445.00	E2300

### SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> IM34 SLDG	Batteries - Group 34, 12V 60 Amp Hours, Installed	785.00	K0108
<input type="checkbox"/> I102768-99-0	Fixed Vent Shelf for C300 <i>If selecting this part, you MUST also pick: PN I0036-7002A - ESP Unit. Fits an LP10 Sized Vent &amp; a U1 Battery. Must be ordered with an expandable controller (ex. R-net Electronics).</i>	402.00	E1029

# C300 Corpus

## Base Options

### SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I102769-99-0	Articulating Vent Tray for C300 <i>If selecting this part, you MUST also pick: PN I0036-7002A - ESP Unit. Fits an LP10 Sized Vent and a Group 34 Battery. Cannot be ordered with any backrest height less than 26.5", a tubular back, or a 13" seat depth. Customization may be required if ordered with any of the following items Co-Pilot, any standard attendant control, push handles, or an oxygen holder. Only available with expandable controllers (R-net Electronics.)</i>	1,318.00	E1030
<input type="checkbox"/> I0036-7002A	Enhanced Steering Perform Unit (ESP), R-net <i>MUST be selected if ordering any ventilator tray (I102769-99-0, I103241-99-0, or I102768-99-0). Recommended item with switched drive inputs for improved tracking.</i>	1,170.00	K0108
<input type="checkbox"/> I103564-99-0	Light Kit for C300 w/ R-net Electronics <i>Shrouds with lights cannot be ordered with VR2 Electronics.</i>	972.00	
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	443.00	

## Permolock

### SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104680-99-0	Permolock C3 Locking Base C300, C350, K300, K450 <i>Kit Includes: Permolock C3 Locking Base, Owner's Manual, Manual Release Handle, and Vehicle Mounted Release Button. Does NOT include Dual Locking Pin Kit (I10608). Locking Base must be installed by a NMEDA Dealer. Please provide the NMEDA Dealer and Address where the locking Base will ship to: NAME: _____ ADDR: _____</i>	1,806.00	
<input type="checkbox"/> I10608	Permolock C3 Dual Locking Pin Kit - C300/K300 <i>This kit includes two spring loaded pins that will be mounted underneath the base of the wheelchair. These pins are spring loaded so they push into the wheelchair base when an obstacle is encountered. The kit also includes spacers that may be used if the Permolock C3 Base is being mounted in a vehicle with a sloped floor. The spacers are only needed for a C300 or K300.</i>	605.00	

## Joystick Options

### CHOOSE JOYSTICK TYPE:

Part Number	Description	Price	HCPCS Code
	PLEASE NOTE: <i>For all remote joystick modules shown below, you must also select an expandable controller. For additional input devices, refer to Alternative Drive Systems Section.</i>		
<input type="checkbox"/> I10410	R-net Remote Joystick w/ Mono Jacks - No Charge	No Charge	
<input type="checkbox"/> I103544-99-0	R-net Remote Joystick w/ Mono Jacks	539.00	
<input type="checkbox"/> ID51040	R-net Remote Joystick w/ Toggle Switches	616.00	
<input type="checkbox"/> ID51042	R-net Remote Joystick w/ Speed Pot & Single Toggle	743.00	
<input type="checkbox"/> ID51153	R-net Remote Color Screen Joystick w/ Mono Jacks	850.00	
<input type="checkbox"/> ID51149	R-net Remote Color Screen JS w/ Toggle Switches	890.00	
<input type="checkbox"/> ID51190	R-net Remote Color JS w/ Speed Pot & Single Toggle	890.00	
<input type="checkbox"/> ID51108	R-net LED Remote Joystick Module	No Charge	
	<i>Standard Joystick Module. This Joystick does NOT have an LCD Display and does NOT contain mono jacks. Limited to five profiles with this Joystick Module. Please note if ordering this item, the ICS Switchbox must be mounted onto either the armrest pad or the backrest. Cannot order PN: I1822845 or I103469-99-0 with this style of joystick.</i>		

### CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard.	No Charge	
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	96.00	E2324
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	141.00	E2323
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	96.00	E2323

# C300 Corpus

## Joystick Options

### CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	129.00	E2323
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	129.00	E2323
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	129.00	E2323
<input type="checkbox"/> IPC101B	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	129.00	E2323
<input type="checkbox"/> IPC102B	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	129.00	E2323
<input type="checkbox"/> IPC107B	Bodypoint J/S Handle Dome Shaped, Rubber	129.00	E2323

### CHOOSE JOYSTICK MOUNT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10412	Joystick Mount Right - Fixed, R-net Only	No Charge	
<input type="checkbox"/> I10413	Joystick Mount Left - Fixed, R-net Only	No Charge	
<input type="checkbox"/> I102537-99-0	Joystick Mount Right - Swing Away, R-net & VR2	297.00	E1028
<input type="checkbox"/> I102538-99-0	Joystick Mount Left - Swing Away, R-net & VR2	297.00	E1028
<input type="checkbox"/> I103820-99-0	Retractable Joystick Mnt Right - R-net, VR2 Only	495.00	E1028
<input type="checkbox"/> I103819-99-0	Retractable Joystick Mnt Left - R-net, VR2 Only	495.00	E1028
<input type="checkbox"/> I101260-99-0	Midline Control Mount - Pilot+ & R-net Only	268.00	E1028

*Please indicate which armrest to mount midline below.*

Mount Midline Joystick on Left Armrest: \_\_\_\_\_

Mount Midline Joystick on Right Armrest: \_\_\_\_\_

### SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input checked="" type="checkbox"/> I10411	Expandable Controller - R-net <i>Must be ordered with the expandable/remote joysticks (PN: I10410, I103544-99-0, ID51040, ID51108, ID51042, and I103814-99-0) OR when ordering Alternative Drive Packages.</i>	545.00	E2377
<input checked="" type="checkbox"/> I10513	Harness for Expandable Controller <i>Must be ordered with any Expandable Controller (I10411 or I10450)</i>	475.00	E2313
<input type="checkbox"/> I104003-99-0	Co-Pilot Attendant Control, Corpus R-net <i>Custom charge may apply if ordering a Corpus Tubular Backrest OR a backrest height of 21" or Lower. Not available with VR2 Controllers.</i>	4,248.00	E2331
<input type="checkbox"/> ID50882	Attendant Control for R-net Control System Only <i>Not available with VR2 Controllers.</i>	1,103.00	E2331
<input type="checkbox"/> ID51111	R-net BT Mouse Module <i>The R-net BT Mouse Module allows for mouse control through a standard R-net Joystick or other type of input device (ex: Mini Joystick w/ Omni). It is only usable as a mouse on a computer. It CANNOT be paired to work with other Bluetooth devices such as cell phones, wireless headsets, etc. This kit includes: an R-net BT Mouse Module, a 1.5m R-net Bus Cable, and a USB Style Bluetooth Adapter for the user's computer. Note: The Bluetooth Adapter may not work with all computers.</i>	809.00	
<input type="checkbox"/> I10582	Output Module (IOM) for R-net Control System <i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). One output module is required for each separate external device to be controlled. Up to three output modules can be installed on one wheelchair.</i>	1,012.00	E2351

## Seating System

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
The Seating Systems Listed Below Include: <i>Corpus Seat Frame, Solid Backrest (Standard Size) with Manual Adjustable Recline Angle, Height Adjustable Armrest Assembly, Solid Seat Pan (17"W x 18"D), Legrest Assembly with Manual Adjustable Elevation. Weight Limit 265 lbs.</i>			
<input type="checkbox"/> I10239	Corpus Seat without Power Tilt or Power Recline	5,995.00	K0108
<input type="checkbox"/> I103525-99-0	Corpus Seat w/ 45° CG Pwr Tilt & Manual Recline	6,295.00	E1002

# C300 Corpus

## Seating System

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10622	Corpus Seat w/ 150° Power Recline	6,850.00	E1004
<input type="checkbox"/> I10417	Corpus Seat w/ 45° CG Pwr Tilt & 150° Pwr Recline	10,395.00	E1007

## Seating System Options

CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10418	Single Seat Function Control Switch for R-net <i>Please choose either Push Buttons or Toggle Switches below.</i>	No Charge	
<input type="checkbox"/> I10419	Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control and/or ICS Switchbox. Order with R-net Controller or remote joystick options. Please indicate below if you want push buttons or toggles on the Seat Function Control Kit.</i>	1,495.00	E2310
<input type="checkbox"/> I10585	Multiple Seat Function Control Kit for R-net	2,925.00	E2311

*MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive control and/or ICS switch box. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the Seat Function Control Kit.*

CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I102798-99-0	Push Buttons on Seat Function Control Kit	No Charge	
<input type="checkbox"/> I102799-99-0	Toggle Switches on Seat Function Control Kit	249.00	

CHOOSE MOUNTING POSITION OF SEAT FUNCTION CONTROL KIT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge	
<input type="checkbox"/> I103471-99-0	Mount ICS Switchbox to End of Arm Pad - Corpus <i>ICS Switchbox will be mounted to the end of the armrest pad. If this item is selected, the ICS Switchbox will NOT swing away with the Joystick Module. Please indicate below if you would like it mounted onto the left or right arm pad.</i> Mount ICS Switchbox on Left Armpad _____ Mount ICS Switchbox on Right Armpad _____	No Charge	
<input type="checkbox"/> I103470-99-0	Mount ICS Switchbox on Backrest	127.00	

SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Not available with VR2 Electronics. Seat function interface box which accommodates up to 8 momentary switches. Provides an additional method to accessing the seat functions. You must still order the Multiple Seat Function Control Kit if ordering this part. 1/8" mini-mono switches not included.</i>	1,002.00	K0108
<input type="checkbox"/> I104746-99-0	Remote Stop System - R-net, Pilot+ <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Each transmitter and receiver is a serialized matched pair, meaning they must be used together. Custom charge may be required if ordering with a low backrest height or a pediatric wheelchair. If ordering Remote Stop System with a K450, push handles must also be selected.</i>	1,002.00	

## Backrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I101853-99-0	BACKREST HEIGHT NOTE: <i>Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°.</i> Ergo Back Fixed 14" W x 21" H - Leatherette	852.00	E2620

# C300 Corpus

## Backrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I101863-99-0	Ergo Back Fixed 14"W x 26.5" High - Leatherette	852.00	E2620
<input type="checkbox"/> I101854-99-0	Ergo Back Sliding 16"W x 21"H - Leatherette	852.00	E2620
<input type="checkbox"/> I101861-99-0	Ergo Back Sliding 16"W x 26.5" H - Leatherette	852.00	E2620
<input type="checkbox"/> I101855-99-0	Ergo Back Sliding 18"W x 21"H - Leatherette	852.00	E2620
<input type="checkbox"/> I101862-99-0	Ergo Back Sliding 18"W x 26.5"H - Leatherette	852.00	E2620
<input type="checkbox"/> I305510-00-90-0	Ergo Back Fixed Tall 16"W x 29.5"H - Leatherette	852.00	E2620
<input type="checkbox"/> I10455	Ergo Back Fixed 14" W x 21" H - Black Mesh	874.00	E2620
<input type="checkbox"/> I10456	Ergo Back Fixed 14"W x 26.5" High - Black Mesh	874.00	E2620
<input type="checkbox"/> I10457	Ergo Back Sliding 16"W x 21"H - Black Mesh	874.00	E2620
<input type="checkbox"/> I10458	Ergo Back Sliding 16"W x 26.5" H - Black Mesh	874.00	E2620
<input type="checkbox"/> I10459	Ergo Back Sliding 18"W x 21"H - Black Mesh	874.00	E2620
<input type="checkbox"/> I10460	Ergo Back Sliding 18"W x 26.5"H - Black Mesh	874.00	E2620
<input type="checkbox"/> I10004	Rectangular Backrest 16"W - Plastic, Fixed	1,139.00	
	<i>Choose backrest height from seat pan (19"-29"): _____ "</i>		
<input type="checkbox"/> I10005	Rectangular Back 18"W - Metal Fixed	1,330.00	
	<i>Choose backrest height from seat pan (19"-29"): _____ "</i>		
<input type="checkbox"/> I10240	Solid Back 16"W X 21"H	No Charge	
	<i>Standard Option with Tilt and/or Recline Systems. Includes Leatherette Back Cushion ONLY.</i>		
<input type="checkbox"/> I10130	Solid Back 16"W x 26.5"H	No Charge	
	<i>Standard Option with Tilt and/or Recline Systems. Includes Leatherette Back Cushion ONLY.</i>		
<input type="checkbox"/> I10241	Solid Back 18"W x 21"H	No Charge	
	<i>Standard Option with Tilt and/or Recline Systems. Includes Leatherette Back Cushion ONLY.</i>		
<input type="checkbox"/> I10131	Solid Back 18"W x 26.5"H	No Charge	
	<i>Standard Option with Tilt and/or Recline Systems. Includes Leatherette Back Cushion ONLY.</i>		
<input type="checkbox"/> I10111	Tubular Backrest Frame Corpus, 16" Wide	363.00	
	<i>Backrest canes to mount aftermarket backrest systems. For tilt only systems. Measurement between the armrests is 17.5". Seat pan to top of canes is 22".</i>		
<input type="checkbox"/> I10112	Tubular Backrest Frame Corpus, 18" Wide	363.00	
	<i>Backrest canes to mount aftermarket backrest systems. For tilt only systems. Measurement between the armrests is 19.5". Seat pan to top of canes is 22"</i>		

SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10014	Custom Backrest Height _____ "	636.00	
	<i>Backrest height measured from seat pan and cannot be lower than 19". Price in addition to above backrest.</i>		
<input type="checkbox"/> I100972-99-0	Split Backrest	544.00	
	<i>Split Backrest 8" from top to allow for transfers and/or transportation. This modification cannot be made on back heights lower than 26.5". Must select either PN: I10429 or I10430 below.</i>		
<input type="checkbox"/> I10429	Split Backrest Cushion	No Charge	
<input type="checkbox"/> I10430	Do NOT Split Backrest Cushion	No Charge	
<input type="checkbox"/> I10284	Flat Backrest Cushion for Ergo Back - Leatherette	91.00	
	<i>This style of cushion is ONLY available in Leatherette. Recommended cushion when ordering lateral supports. Lateral wedges will be omitted if ordering this part.</i>		
<input type="checkbox"/> I10318	ROHO Backrest Cushion For Ergo Back	636.00	
	<i>Price Includes: ROHO MID PROFILE Backrest Cushion which will be matched to fit the Ergonomic Backrest that you select. Approximate height of MID PROFILE cells is 3" in their inflated state. Also includes a Mesh Cover for the ROHO Backrest Cushion, Hand Pump, and Patch Kit. This will NOT fit these part numbers: (I10004, I10005, I10111, I10112, I305510-00-90-0). Please indicate below if you would like to OMIT the Ergo Back Cushion.</i>		
<input type="checkbox"/> I10319	ROHO Backrest Cushion For Ergo Back - Notched	636.00	
	<i>Recommended ROHO cushion when ordering Lateral Trunk Supports. Please indicate below if you would like to OMIT the Ergo Back Cushion.</i>		

# C300 Corpus

## Backrest Options

### SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10320	Omit Ergonomic Backrest Cushion	No Charge	

## Seat Sizing

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I100684-99-0	Corpus Seat 17"W x 16"D	No Charge	
<input type="checkbox"/> I100689-99-0	Corpus Seat 17"W x 18"D	No Charge	
<input type="checkbox"/> I100695-99-0	Corpus Seat 17"W x 20"D	No Charge	
<input type="checkbox"/> I100720-99-0	Corpus Seat 17"W x 22"D	No Charge	
<input type="checkbox"/> I100690-99-0	Corpus Seat 19"W x 18"D	No Charge	
<input type="checkbox"/> I100696-99-0	Corpus Seat 19"W x 20"D	No Charge	
<input type="checkbox"/> I100721-99-0	Corpus Seat 19"W x 22"D	No Charge	
<input type="checkbox"/> I10007	Special Seat Size 21"W x 18"D	No Charge	
<input type="checkbox"/> I10008	Special Seat Size 21"W x 20"D	No Charge	
<input type="checkbox"/> I10009	Special Seat Size 21"W x 22"D	No Charge	

### SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10038	Ergonomic Seat Cushion Leatherette <i>Not available with seat widths less than 16" or seat depths less than 14".</i>	483.00	E2601,E2602
<input type="checkbox"/> I10461	Ergonomic Seat Cushion - Black Mesh <i>Not available with seat widths less than 16" or seat depths less than 14".</i>	493.00	E2601

## Armrest Options

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103339-99-0	4"x10" Corpus Arm Right - Dual Taper, Leatherette	No Charge	
<input type="checkbox"/> I103341-99-0	4"x12" Corpus Arm Right - Dual Taper, Leatherette	No Charge	
<input type="checkbox"/> I103343-99-0	4"x14" Corpus Arm Right - Dual Taper, Leatherette	No Charge	
<input type="checkbox"/> I103345-99-0	4"x16" Corpus Arm Right - Dual Taper, Leatherette	No Charge	
<input type="checkbox"/> I103331-99-0	4"x10" Corpus Arm Right - Dual Taper, Black Fabric	No Charge	
<input type="checkbox"/> I103333-99-0	4"x12" Corpus Arm Right - Dual Taper, Black Fabric	No Charge	
<input type="checkbox"/> I103335-99-0	4"x14" Corpus Arm Right - Dual Taper, Black Fabric	No Charge	
<input type="checkbox"/> I103337-99-0	4"x16" Corpus Arm Right - Dual Taper, Black Fabric	No Charge	

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103338-99-0	4"x10" Corpus Arm Left - Dual Taper, Leatherette	No Charge	
<input type="checkbox"/> I103340-99-0	4"x12" Corpus Arm Left - Dual Taper, Leatherette	No Charge	
<input type="checkbox"/> I103342-99-0	4"x14" Corpus Arm Left - Dual Taper, Leatherette	No Charge	
<input type="checkbox"/> I103344-99-0	4"x16" Corpus Arm Left - Dual Taper, Leatherette	No Charge	
<input type="checkbox"/> I103330-99-0	4"x10" Corpus Arm Left - Dual Taper, Black Fabric	No Charge	
<input type="checkbox"/> I103332-99-0	4"x12" Corpus Arm Left - Dual Taper, Black Fabric	No Charge	
<input type="checkbox"/> I103334-99-0	4"x14" Corpus Arm Left - Dual Taper, Black Fabric	No Charge	
<input type="checkbox"/> I103336-99-0	4"x16" Corpus Arm Left - Dual Taper, Black Fabric	No Charge	

### SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10137	Height Adjustable Armrest Assembly	289.00	E0973
<input type="checkbox"/> I10016	Height Adjustable Armrest Assembly - No Charge <i>No charge with power tilt and/or recline.</i>	No Charge	

# C300 Corpus

## Armrest Options

### SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I101003-99-0	Armrest Mount High <i>Range of adjustment in high position is 11" - 13.5" measured from the seat pan to the top of the armrest pad.</i>	No Charge	
<input type="checkbox"/> I101002-99-0	Armrest Mount Low (High Position Standard) <i>Range of adjustment in low position is 9" - 11.5" measured from the seat pan to the top of the armrest pad.</i>	No Charge	
<input type="checkbox"/> I304968-99-0	Corpus Armbar for 18" Btw the Arms - Installed	270.00	
<input type="checkbox"/> I10001	Corpus Armbar for 20" Btw the Arms - Installed	300.00	
<input type="checkbox"/> I100523-99-0	Armrest Height Ext (+2")	127.00	
<input type="checkbox"/> IARMADAP1"	Arm Trough Adaptor 1" Pair <i>For aftermarket arm troughs. Please indicate below if they should be installed or on chair sent as parts.</i>	279.00	
<input type="checkbox"/> IARMADAP7/8"	Arm Trough Adaptor 7/8" Pair <i>For aftermarket arm troughs. Please indicate below if they should be installed on chair or sent as parts.</i> Install Arm Trough Adapters: _____ Send Arm Trough Adapters as Parts: _____	279.00	
<input type="checkbox"/> I100671-99-0	Armrest Pouch-Right	193.00	
<input type="checkbox"/> I100672-99-0	Armrest Pouch-Left	193.00	
<input type="checkbox"/> I101261-99-0	Armrest Side Support Right Outer <i>For Corpus Armrest Pads only.</i>	144.00	
<input type="checkbox"/> I101262-99-0	Armrest Side Support Left Outer <i>For Corpus Armrest Pads only.</i>	144.00	
<input type="checkbox"/> I101263-99-0	Armrest Side Support Right Inner <i>Recommend 2" wider width between the armrest to ensure armrests will flip-up.</i>	144.00	
<input type="checkbox"/> I101264-99-0	Armrest Side Support Left Inner <i>Recommend 2" wider width between the armrest to ensure armrests will flip-up.</i>	144.00	
<input type="checkbox"/> I101510-99-0	Armrest Supports HD <i>Cannot be ordered with 175 degree power recline seating systems. (I10482, I10485, I10488).</i>	495.00	
<input type="checkbox"/> I10285	Mount Corpus Armrest Pads On Inside of Bracket <i>Used to achieve 14"W between the arm measurement when ordering a 14"W Backrest.</i>	No Charge	

## Legrest Options

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I102206-99-0	Manual Adjustment Legrest Elevation	No Charge	
<input type="checkbox"/> I103534-99-0	Power Articulating Center Mount Elevating Legrest	2,825.00	K0108

### SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103055-99-0	UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad. Recommended pad to be used with lower leg lengths of 12" or less.</i>	No Charge	
<input type="checkbox"/> I102982-99-0	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad.</i>	No Charge	
<input type="checkbox"/> I102978-99-0	UT Calf Support Kit (H) w Mesh Cover & Adj Hardware <i>7"W x 7"T Curved Pad.</i>	317.00	E1028
<input type="checkbox"/> I100561-99-0	Legrest Assm. 1pc Footplate (Complete)	402.00	
<input type="checkbox"/> I101276-99-0	Legrest Assm. 2pc Footplates (Complete)	No Charge	
<input type="checkbox"/> I10024	Leg Length Ext Kit (pair) 16-21" From Seat Pan <i>Standard length is 11" - 16" from seat pan to footplate. Set legs at _____" from top of seat pan to footplate.</i>	280.00	

# C300 Corpus

## Accessories

### SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10400	Bodypoint Non-Padded Hip Belt, UniTrack - Medium	No Charge	
<input type="checkbox"/> I10404	Bodypoint Padded Hip Belt UniTrack (C300) - Medium	115.00	K0108
<input type="checkbox"/> I10402	Bodypoint Non-Padded Hip Belt, UniTrack - Large	94.00	K0108
<input type="checkbox"/> I10408	Bodypoint Padded Hip Belt UniTrack - Large	158.00	K0108
<input type="checkbox"/> I100519-99-0	Positioning Belt Retractable	239.00	K0108
<input type="checkbox"/> I100600-99-0	Positioning Belt Retractable Long	259.00	K0108
<input type="checkbox"/> I104131-99-0	Headrest for Ergo Back w/ Curved Mount Brkt - Leath <i>Must also select PN I10552. This part consists of a leatherette headrest pad (10"W x 5"T) with a 10"L slotted, curved mounting bracket.</i>	268.00	E0955
<input type="checkbox"/> I103460-99-0	UT Ergo Headrest w/ Black Fabric Cover <i>Must also select PN I10552.</i>	297.00	E0955
<input type="checkbox"/> I1823249	UT Lateral Adjustment Bar for UniTrack Headrest <i>Allows for approximately 4" of lateral adjustment to the UT Ergo Headrest either right or left. NOTE: Due to the width of this item (approx 9"), it may be necessary to send other items such as back pack clips unmounted. Custom charge may apply with lower back heights.</i>	56.00	
<input type="checkbox"/> IHEADADAPTER	Universal Headrest Adapter <i>For aftermarket headrests.</i>	205.00	K0108
<input type="checkbox"/> I10552	Adj Removable - Headrest Hardware (1 Each)	216.00	E1028
<input type="checkbox"/> I103070-99-0	UT Amp. Support Pad (H) w/ Mesh Cover - (1 Each) <i>7"W x 7"T Curved Pad. Must also select PN I10553.</i>	127.00	E1020
<input type="checkbox"/> I103071-99-0	UT Amp. Support Pad (D) w/ Mesh Cover - (1 Each) <i>8"W x 8"T Flat Pad. Must also select PN I10553.</i>	127.00	E1020
<input type="checkbox"/> I10553	Adj Removable - Amp Support Hardware (1 Each) <i>Please indicate if you would like the amputation support mounted on the right or left side of the seat frame: Mount Amp Suppt on Right: ___ Mount Amp Suppt on Left: ___</i>	216.00	E1028
<input type="checkbox"/> I103077-99-0	UT Hip Suppt (B) w/ Mesh Cover - Low (Pair) <i>4.5"W x 3.5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	229.00	E0956
<input type="checkbox"/> I103079-99-0	UT Hip Suppt (A) w/ Mesh Cover - Low (Pair) <i>6"W x 5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	229.00	E0956
<input type="checkbox"/> I103094-99-0	UT Thigh/Hip Suppt (J) w/ Mesh Cover - High (Pair) <i>6.5"W x 5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	229.00	E0956
<input type="checkbox"/> I103095-99-0	UT Thigh/Hip Suppt (C) w/ Mesh Cover - High (Pair) <i>8"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	229.00	E0956
<input type="checkbox"/> I102977-99-0	UT Thigh/Hip Suppt (G) w/ Mesh Cover - High (Pair) <i>12"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	229.00	E0956
<input type="checkbox"/> I10550	Adj Removable - Thigh Support Hardware (Pair)	431.00	E1028
<input type="checkbox"/> I305266-28-90-0	Lateral Supports Std. Leatherette-ErgoBack Only <i>Must also select PN I10551.</i>	250.00	E0956
<input type="checkbox"/> I100521-99-0	Swing Away Trunk Supports Leatherette-For Ergo Back <i>Must also select PN I10551.</i>	254.00	E0956
<input type="checkbox"/> I10551	Adj Removable - Trunk Support Hardware (Pair)	431.00	E1028
<input type="checkbox"/> I100978-99-0	Abductor Kit - Leatherette <i>Must also select PN I10555.</i>	311.00	E0957
<input type="checkbox"/> I10555	Adj Removable - Abductor Hardware (1 Each)	216.00	E1028
<input type="checkbox"/> I10015	Knee Block Assembly Corpus	641.00	K0108
<input type="checkbox"/> I103651-99-0	Transfer Handles for UniTrack Accessory Rail	397.00	
<input type="checkbox"/> I100975-99-0	Push Handles for Ergo Backrest <i>Order with Ergonomic or Rectangular Backs. Custom charge may apply if ordering push handles with an articulating vent tray.</i>	345.00	
<input type="checkbox"/> I307103-00-0	Push Handles for Tubular Backrest <i>Order with Corpus Tubular Backs (PN: I10111 &amp; I10112).</i>	343.00	

# C300 Corpus

## Accessories

### SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I102059-99-0	Medical Necessities Bag <i>Includes Medical Necessity Bag Clips (Carabiner Style).</i>	229.00	
<input type="checkbox"/> I101473-99-0	Medical Necessity Bag Clips - Carabiner Style	76.00	
<input type="checkbox"/> I105144-99-0	Bodypoint Monoflex Large Chest Support Belt Kit <i>21.5"L x 4.5"W Pad Size. Item will be shipped as parts.</i>	145.00	E0960
<input type="checkbox"/> I105143-99-0	Bodypoint Monoflex Medium Chest Support Belt Kit <i>18.25"L x 4"W Pad Size. Item will be shipped as parts.</i>	145.00	E0960
<input type="checkbox"/> I105145-99-0	Bodypoint Monoflex X-Lrge Chest Support Belt Kit <i>24.5"L x 5"W Pad Size. Item will be shipped as parts.</i>	145.00	E0960
<input type="checkbox"/> I100578-99-0	Upper Extremity Support Small, Angle Adjustable <i>Must also select PN I10556.</i>	242.00	E0950
<input type="checkbox"/> I100575-99-0	Upper Extremity Support Angle Adjustable. 14" W <i>Must also select PN I10556.</i>	272.00	E0950
<input type="checkbox"/> I100576-99-0	Upper Extremity Support Angle Adjustable. 16" W <i>Must also select PN I10556.</i>	284.00	E0950
<input type="checkbox"/> I100577-99-0	Upper Extremity Support Angle Adjustable. 18" W <i>Must also select PN I10556.</i>	311.00	E0950
<input type="checkbox"/> I100571-99-0	Upper Extremity Support Joystick Cutout 16" W <i>Must also select PN I10556.</i>	328.00	E0950
<input type="checkbox"/> I100572-99-0	Upper Extremity Support Joystick Cutout 18" W <i>Must also select PN I10556.</i>	328.00	E0950
<input type="checkbox"/> I100573-99-0	Upper Extremity Support Mid Joystick Cutout 16"W <i>Must also select PN I10556.</i>	402.00	E0950
<input type="checkbox"/> I100574-99-0	Upper Extremity Support Mid Joystick Cutout 18"W <i>Must also select PN I10556.</i>	402.00	E0950
<input type="checkbox"/> I103426-99-0	Upper Extremity Support Angle Adjust - UniTrack <i>Side mounted tray that attaches to the UniTrack Accessory Rail. Dimensions are: 19.5"W x 13"D. Can only be ordered with chairs that offer the UniTrack System. Must also select PN I10556.</i>	301.00	E0950
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each)	216.00	E1028
<input type="checkbox"/> I101964-99-0	Oxygen Holder for PS or Corpus Seat	228.00	E2208
	<i>Cannot be ordered with a Corpus Tubular Backrest. (I10111 or I10112). Custom charge may be applied if ordering in combination with accessories such as medical necessities bags, attendant controls, push handles, back pack clips, etc. This item is mounted to the Ergonomic Backrest. Customization may be required if ordering on a Street.</i>		
<input type="checkbox"/> I101956-99-0	Crutch Holder for PS or Corpus Seat <i>Cannot be ordered with a Corpus Tubular Backrest. (I10111 or I10112). Accessories such as medical necessities bags and push handles may need to be sent as parts if ordering a crutch holder. This item is mounted to the ergonomic backrest. Customization may be required if ordering on a Street.</i>	228.00	E2207

## Alternative Drive Systems

### (1) CHOOSE ALTERNATIVE DRIVE INPUT:

Part Number	Description	Price	HCPCS Code
	R-net Alternative Drive Input Devices <i>R-net Alternative Drive Input Devices take the place of, or may be used in addition to, the standard R-net Joystick. They are only usable on wheelchairs that have R-net electronics. You MUST select one of each of the following items: (1) Input Device, (2) Adjustable/Removable Alt Drive Hardware, AND (3) Omni or Input Module.</i>		
<input type="checkbox"/> I10558	Gimble Joystick set up for Armrest Mount	1,412.00	E2373
<input type="checkbox"/> I10559	Gimble Joystick set up for Midline Mount	1,412.00	E2373

# C300 Corpus

## Alternative Drive Systems

### (1) CHOOSE ALTERNATIVE DRIVE INPUT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10560	Gimble Joystick set up for Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	1,412.00	E2373
<input type="checkbox"/> I10561	Gimble Joystick w/ Haims Harness - Small <i>For neck width up to 5".</i>	1,582.00	E2373
<input type="checkbox"/> I10562	Gimble Joystick w/ Haims Harness - Large <i>For neck width up to 7".</i>	1,582.00	E2373
<input type="checkbox"/> I10563	Gimble Joystick w/ Power Swing Away Chin Boom <i>Custom charge may be applied if ordering smaller backrest sizes and heights.</i>	5,265.00	E2373
<input type="checkbox"/> I10564	Mini Joystick set up for Armrest Mount	2,685.00	E2312
<input type="checkbox"/> I10565	Mini Joystick set up for Midline Mount	2,685.00	E2312
<input type="checkbox"/> I10566	Mini Joystick set up for Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	2,685.00	E2312
<input type="checkbox"/> I10567	Mini Joystick w/ Haims Harness - Small <i>For neck width up to 5".</i>	2,745.00	E2312
<input type="checkbox"/> I10568	Mini Joystick w/ Haims Harness - Large <i>For neck width up to 7".</i>	2,745.00	E2312
<input type="checkbox"/> I10569	Mini Joystick w/ Power Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	5,265.00	E2312
<input type="checkbox"/> I1821750	Tablet Control Pad <i>The Tablet Control is a unique drive input which allows a user to control a power wheelchair in a very similar manner to a mouse being controlled by a touchpad on a laptop computer. The dimensions of the Tablet Control Pad are: 12"W x 9.5"D x 1.5"H. Kit also includes: Buddy Button Switch, and Tablet Control Mounting Bracket.</i>	4,375.00	E2328
<input type="checkbox"/> I10570	Head Array - Adult 5 Switch HDR-A5	4,375.00	E2330
<input type="checkbox"/> I10571	Head Array - Adult 3 Switch HDR-A4 <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	4,175.00	E2330
<input type="checkbox"/> I10572	Head Array - Pediatric 5 Switch HDR-C5-RR	4,375.00	E2330
<input type="checkbox"/> I10573	Sip & Puff System <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	1,475.00	E2325
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	1,963.00	E2325,E2326

### SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10076	Right Side Mount For Joystick	No Charge	
<input type="checkbox"/> I10077	Left Side Mount For Joystick	No Charge	
<input type="checkbox"/> I10074	Right Side Mount For Omni	No Charge	
<input type="checkbox"/> I10075	Left Side Mount For Omni	No Charge	
<input type="checkbox"/> I10078	Mount Seat Control Panel on Right	No Charge	
<input type="checkbox"/> I10079	Mount Seat Control Panel on Left	No Charge	
<input type="checkbox"/> I10080	Mount Seat Control Panel on Backrest	No Charge	

## Alt Drive Systems Setup

### (2) CHOOSE ADJ ALT DRIVE HARDWARE (REQUIRED):

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10577	Adjustable Removable Hardware - Alt Drive	216.00	E1028

### (3) CHOOSE OMNI OR INPUT MODULE (REQUIRED):

Part Number	Description	Price	HCPCS Code
-------------	-------------	-------	------------

# C300 Corpus

## Alt Drive Systems Setup

(3) CHOOSE OMNI OR INPUT MODULE (REQUIRED):

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> ID51154IR	Omni Control Module for R-net w/ Infrared <i>Please specify which type of alternative drive input you will be using if not ordering the input device from Permobil. Input Device: _____</i>	1,496.00	
<input type="checkbox"/> I10581	Input Module (IOM) for R-net Control System <i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head Array or RIM control.</i>	1,012.00	E2351

## Alt Drive Accessories

SELECT OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> IMICROSWITCH	Microlite Switch Qty: _____	151.00	
<input type="checkbox"/> IBUDDYBUTTON	Buddy Button Qty: _____	113.00	K0108
<input type="checkbox"/> I5808	Mini Cup Switch Qty: _____	110.00	
<input type="checkbox"/> IZTS-P	Proximity Switch Qty: _____	524.00	
<input type="checkbox"/> IMINIJOY5SW	Mini Joystick 5 Switch OC	603.00	
<input type="checkbox"/> I5018	Wafer Board w/5 Switches	506.00	

Notes & Additional Instructions:

### Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, Inc. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

### Order Acknowledgement:

I, \_\_\_\_\_, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: \_\_\_\_\_

Printed Name: \_\_\_\_\_

Title: \_\_\_\_\_

Date: \_\_\_\_\_

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. MID PROFILE® Cushion is a registered trademark of the The ROHO Group. The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.