


C500 Lowrider



*** Indicates a Required Field**

Dealer Information

* Contact: _____

Dealer Code: _____

* Dealer Name: _____

Address: _____

* City: _____

* State/Zip: _____

* Phone#: _____

* Fax#: _____

PO#: _____

Email Address: _____

Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

* First Name: _____

* Last Name: _____

Diagnosis: _____

Funding Source: _____

Client Age: _____

Client Measurements

* Weight: _____

* Height: _____

A) Top of Shoulders: _____

B) Chest Depth: _____

C) Chest Width: _____

D) Seat Depth: _____

E) Top of Head: _____

F) Elbow to Hand: _____

G) Seat Pan to Elbow: _____

H) Hip Width: _____

* I) Knee to Foot: _____

Three column pricing on order/quote:

Display HCPCS Codes on order/quote:

Please send order/quote to fax#: **(800) 231-3256**

Email to: sales@permobilus.com

Permobil Inc.
6961 Eastgate Blvd.
Lebanon, TN 37090
Tel: (800) 736-0925
Fax: (800) 231-3256
www.permobil.com

Prices effective **August 9, 2010.**

C500 Lowrider

Chair Model and Colors

CHOOSE BASE:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> I104020-99-0 | C500 Lowrider R-net, 5.0mph - C500 PS2 <i>Price Includes: FWD C500 Power Base Set Up For Multiple Power Seat Functions, Shock Absorbing Suspension System, Tie-Down Hardware for Strap Systems, Charger, 17" Seat to Floor Height, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below.</i> | 12,535.00 | K0884 |
| <input type="checkbox"/> I103968-99-0 | C500 HS Lowrider R-net, 6.5mph - C500 PS2 <i>Price Includes: FWD C500 Power Base Set Up For Multiple Power Seat Functions, Shock Absorbing Suspension System, R-net ESP (Enhanced Steering Performance), High Speed Motors, 120 AMP Electronics, Tie-Down Hardware for Strap Systems, Charger, 17" Seat to Floor Height, Flat Free Tires, and Standard Positioning Belt. Must choose either Flat Free or Pneumatic Tires.</i> | 12,995.00 | K0884 |

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> I103499-99-0 | Tri-Spoke Split Rim w/ Flat Free Drive Tire 3x8" | No Charge | |
| <input type="checkbox"/> I103498-99-0 | Tri-Spoke Split Rim w/ Pneumatic Tires - 3x8" | No Charge | |

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|------------------------------------|-----------|------------|
| <input type="checkbox"/> I103816-99-0 | C500 Shroud Onyx Black | No Charge | |
| <input type="checkbox"/> I102260-99-0 | C500 Shroud Vintage Wine | No Charge | |
| <input type="checkbox"/> I102261-99-0 | C500 Shroud Crystal Blue | No Charge | |
| <input type="checkbox"/> I102262-99-0 | C500 Shroud Metallic Charcoal Gray | No Charge | |
| <input type="checkbox"/> I102263-99-0 | C500 Shroud Green Metallic | No Charge | |
| <input type="checkbox"/> I103833-99-0 | C500 Shroud Platinum Silver | No Charge | |

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> I10451 | Match Hub Cap to Shroud Color <i>If this is selected, the Tri-Spoke Hub Cap will be matched to the selected shroud color.</i> | No Charge | |
| <input type="checkbox"/> I103537-99-0 | Hub Cap Platinum Silver <i>If this is selected, the Tri-Spoke Hub Cap will be Platinum Silver.</i> | No Charge | |

Base Options

MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> I102273-99-0 | Seat Tube Fixed | No Charge | |
| <input type="checkbox"/> I103972-99-0 | Power Adjustable Seat Height 6" Travel | 3,045.00 | E2300 |

SELECT OPTIONS:

| Part Number | Description | Price | HCPCS Code |
|--|---|-----------|------------|
| <input type="checkbox"/> IM24 SLD G FT | Batteries, Group 24 Sealed LA/Gel, Installed | 775.00 | E2363 |
| <input type="checkbox"/> I102210-99-0 | Anti-Tippers for C400/C500 | 406.00 | |
| <input type="checkbox"/> I103635-99-0 | Light Kit, R-net, C400/C500 | 972.00 | |
| <input checked="" type="checkbox"/> I104082-99-0 | R-net ESP w C400/C500 Mounting Bracket <i>MUST be selected if ordering any ventilator tray (I102769-99-0, I103241-99-0, or I102768-99-0). Recommended item with switched drive inputs for improved tracking.</i> | No Charge | K0108 |
| <input type="checkbox"/> ID51113 | Programming Key for R-net - Permobil Version | 443.00 | |

Permlock

SELECT OPTIONS:

| Part Number | Description | Price | HCPCS Code |
|-------------|-------------|-------|------------|
|-------------|-------------|-------|------------|

C500 Lowrider

Permolock

SELECT OPTIONS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|----------|------------|
| <input type="checkbox"/> I102843-99-0 | Permolock C Kit - COMPLETE <i>Includes: Locking Base, Plunger, Manual Release Handle, and Vehicle Mounted Release Button. The Locking Base must be installed by a NMEDA Dealer. Please provide the NMEDA Dealer Name & Address where the locking Base will ship to</i> NAME: _____ ADDR: _____ | 2,528.00 | |
| <input type="checkbox"/> I102842-99-0 | Plunger ONLY for Permolock C <i>Price Includes: Plunger only. This item may be purchased by any Permobil Inc or NMEDA certified dealer. The standard mounting position for the button to control the plunger is at the end of the armrest pad, opposite the joystick.</i> | 1,026.00 | |

Joystick Options

CHOOSE JOYSTICK TYPE:

| Part Number | Description | Price | HCPCS Code |
|--|--|--------|------------|
| PLEASE NOTE: <i>For all remote joystick modules shown below, you must also select an expandable controller. For additional input devices, refer to Alternative Drive Systems Section.</i> | | | |
| <input type="checkbox"/> I103544-99-0 | R-net Remote Joystick w/ Mono Jacks | 539.00 | |
| <input type="checkbox"/> ID51040 | R-net Remote Joystick w/ Toggle Switches | 616.00 | |
| <input type="checkbox"/> ID51042 | R-net Remote Joystick w/ Speed Pot & Single Toggle | 743.00 | |
| <input type="checkbox"/> ID51153 | R-net Remote Color Screen Joystick w/ Mono Jacks | 850.00 | |
| <input type="checkbox"/> ID51149 | R-net Remote Color Screen JS w/ Toggle Switches | 890.00 | |
| <input type="checkbox"/> ID51190 | R-net Remote Color JS w/ Speed Pot & Single Toggle | 890.00 | |

CHOOSE JOYSTICK HANDLE TYPE:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> I10045 | Cone Shaped Joystick Knob - Standard. | No Charge | |
| <input type="checkbox"/> I100960-99-0 | Chin Cup For Joystick Knob | 96.00 | E2324 |
| <input type="checkbox"/> I100961-99-0 | "T" Handle For Joystick Knob | 141.00 | E2323 |
| <input type="checkbox"/> I100962-99-0 | "Large Ball" For Joystick Knob | 96.00 | E2323 |
| <input type="checkbox"/> I100963-99-0 | "Softball" For Joystick Knob | 129.00 | E2323 |
| <input type="checkbox"/> I100964-99-0 | "Mushroom" For Joystick Knob | 129.00 | E2323 |
| <input type="checkbox"/> I100965-99-0 | "Stick" For Joystick Knob | 129.00 | E2323 |
| <input type="checkbox"/> IPC101B | Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft | 129.00 | E2323 |
| <input type="checkbox"/> IPC102B | Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft | 129.00 | E2323 |
| <input type="checkbox"/> IPC107B | Bodypoint J/S Handle Dome Shaped, Rubber | 129.00 | E2323 |

CHOOSE JOYSTICK MOUNT:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> I10412 | Joystick Mount Right - Fixed, R-net Only | No Charge | |
| <input type="checkbox"/> I10413 | Joystick Mount Left - Fixed, R-net Only | No Charge | |
| <input type="checkbox"/> I102537-99-0 | Joystick Mount Right - Swing Away, R-net & VR2 | 297.00 | E1028 |
| <input type="checkbox"/> I102538-99-0 | Joystick Mount Left - Swing Away, R-net & VR2 | 297.00 | E1028 |
| <input type="checkbox"/> I103820-99-0 | Retractable Joystick Mnt Right - R-net, VR2 Only | 495.00 | E1028 |
| <input type="checkbox"/> I103819-99-0 | Retractable Joystick Mnt Left - R-net, VR2 Only | 495.00 | E1028 |
| <input type="checkbox"/> I101260-99-0 | Midline Control Mount - Pilot+ & R-net Only | 268.00 | E1028 |

Please indicate which armrest to mount midline below.

Mount Midline Joystick on Left Armrest: _____

Mount Midline Joystick on Right Armrest: _____

SELECT OPTIONS:

| Part Number | Description | Price | HCPCS Code |
|-------------|-------------|-------|------------|
|-------------|-------------|-------|------------|

C500 Lowrider

Joystick Options

SELECT OPTIONS:

| Part Number | Description | Price | HCPCS Code |
|--|--|----------|------------|
| <input checked="" type="checkbox"/> I10411 | Expandable Controller - R-net <i>Must be ordered with the expandable/remote joysticks (PN: I10410, I103544-99-0, ID51040, ID51108, ID51042, and I103814-99-0) OR when ordering Alternative Drive Packages.</i> | 545.00 | E2377 |
| <input checked="" type="checkbox"/> I10513 | Harness for Expandable Controller <i>Must be ordered with any Expandable Controller (I10411 or I10450)</i> | 475.00 | E2313 |
| <input type="checkbox"/> I104003-99-0 | Co-Pilot Attendant Control, Corpus R-net <i>Custom charge may apply if ordering a Corpus Tubular Backrest OR a backrest height of 21" or Lower. Not available with VR2 Controllers.</i> | 4,248.00 | E2331 |
| <input type="checkbox"/> ID50882 | Attendant Control for R-net Control System Only <i>Not available with VR2 Controllers.</i> | 1,103.00 | E2331 |
| <input type="checkbox"/> ID51111 | R-net BT Mouse Module <i>The R-net BT Mouse Module allows for mouse control through a standard R-net Joystick or other type of input device (ex: Mini Joystick w/ Omni). It is only usable as a mouse on a computer. It CANNOT be paired to work with other Bluetooth devices such as cell phones, wireless headsets, etc. This kit includes: an R-net BT Mouse Module, a 1.5m R-net Bus Cable, and a USB Style Bluetooth Adapter for the user's computer. Note: The Bluetooth Adapter may not work with all computers.</i> | 809.00 | |
| <input type="checkbox"/> I10582 | Output Module (IOM) for R-net Control System <i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). One output module is required for each separate external device to be controlled. Up to three output modules can be installed on one wheelchair.</i> | 1,012.00 | E2351 |

Seating System

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|---|-----------|------------|
| | The Seating Systems Listed Below Include: <i>Corpus Seat Frame, Solid Backrest (Standard Size) with Adjustable Recline Angle, Height Adjustable Armrest Assembly, Solid Seat Pan (17"w x 18"d), Legrest Assembly with Adjustable Elevation. Weight Limit 250lbs (310 lbs/ HD Package).</i> | | |
| <input type="checkbox"/> I10489 | Lowrider Seat w/20° Power Tilt & Manual Recline <i>Includes a manual recline actuator mounted behind the backrest.</i> | 7,195.00 | E1002 |
| <input type="checkbox"/> I10490 | Lowrider Seat w/20° Power Tilt & 120° Power Recline <i>Includes a powered recline actuator mounted behind the backrest.</i> | 11,395.00 | E1007 |

Seating System Options

CHOOSE SEAT FUNCTION CONTROL TYPE:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|--|-----------|------------|
| <input type="checkbox"/> I10418 | Single Seat Function Control Switch for R-net <i>Please choose either Push Buttons or Toggle Switches below.</i> | No Charge | |
| <input type="checkbox"/> I10419 | Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control and/or ICS Switchbox. Order with R-net Controller or remote joystick options. Please indicate below if you want push buttons or toggles on the Seat Function Control Kit.</i> | 1,495.00 | E2310 |
| <input type="checkbox"/> I10585 | Multiple Seat Function Control Kit for R-net <i>MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive control and/or ICS switch box. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the Seat Function Control Kit.</i> | 2,925.00 | E2311 |

CHOOSE PUSH BUTTONS OR TOGGLES:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> I102798-99-0 | Push Buttons on Seat Function Control Kit | No Charge | |
| <input type="checkbox"/> I102799-99-0 | Toggle Switches on Seat Function Control Kit | 249.00 | |

C500 Lowrider

Seating System Options

CHOOSE MOUNTING POSITION OF SEAT FUNCTION CONTROL KIT:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> I10626 | Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i> | No Charge | |
| <input type="checkbox"/> I103471-99-0 | Mount ICS Switchbox to End of Arm Pad - Corpus <i>ICS Switchbox will be mounted to the end of the armrest pad. If this item is selected, the ICS Switchbox will NOT swing away with the Joystick Module. Please indicate below if you would like it mounted onto the left or right arm pad.</i> Mount ICS Switchbox on Left Armpad _____ Mount ICS Switchbox on Right Armpad _____ | No Charge | |
| <input type="checkbox"/> I103470-99-0 | Mount ICS Switchbox on Backrest | 127.00 | |

SELECT OPTIONS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|----------|------------|
| <input type="checkbox"/> I10496 | Heavy Duty Package Lowrider w/ R-net <i>Includes: Anti-tippers, heavy duty armrest supports, drive inhibit for power adjustable seat height and heavy duty mount for recline. This is a required option for users between 251-310lbs.</i> | 2,590.00 | |
| <input type="checkbox"/> I104196-99-0 | ICS Alternative Switchbox for R-net Electronics <i>Not available with VR2 Electronics. Seat function interface box which accommodates up to 8 momentary switches. Provides an additional method to accessing the seat functions. You must still order the Multiple Seat Function Control Kit if ordering this part. 1/8" mini-mono switches not included.</i> | 1,002.00 | K0108 |
| <input type="checkbox"/> I104746-99-0 | Remote Stop System - R-net, Pilot+ <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Each transmitter and receiver is a serialized matched pair, meaning they must be used together. Custom charge may be required if ordering with a low backrest height or a pediatric wheelchair. If ordering Remote Stop System with a K450, push handles must also be selected.</i> | 1,002.00 | |
| <input type="checkbox"/> I10497 | Memory Seat Program, C400 or C500 Corpus w R-net <i>If this is selected, the chair will be programmed to allow the user to store and save up to 3 different memory seat positions. On the C400/C500 Corpus and Lowrider, the following seat functions are capable of remembering a seat position: Power Tilt, Power Recline, & Power Legs. The Power Adj Seat Height Function is NOT capable of recalling a saved memory seat position. Please refer to the Corpus Seat Owner's Manual for more information on saving and recalling memory seat positions.</i> | 100.00 | |

Backrest Options

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---|--|--------|------------|
| BACKREST HEIGHT NOTE: <i>Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°.</i> | | | |
| <input type="checkbox"/> I101853-99-0 | Ergo Back Fixed 14" W x 21" H - Leatherette | 852.00 | E2620 |
| <input type="checkbox"/> I101863-99-0 | Ergo Back Fixed 14"W x 26.5" High - Leatherette | 852.00 | E2620 |
| <input type="checkbox"/> I101854-99-0 | Ergo Back Sliding 16"W x 21"H - Leatherette | 852.00 | E2620 |
| <input type="checkbox"/> I101861-99-0 | Ergo Back Sliding 16"W x 26.5" H - Leatherette | 852.00 | E2620 |
| <input type="checkbox"/> I101855-99-0 | Ergo Back Sliding 18"W x 21"H - Leatherette | 852.00 | E2620 |
| <input type="checkbox"/> I101862-99-0 | Ergo Back Sliding 18"W x 26.5"H - Leatherette | 852.00 | E2620 |
| <input type="checkbox"/> I305510-00-90-0 | Ergo Back Fixed Tall 16"W x 29.5"H - Leatherette | 852.00 | E2620 |
| <input type="checkbox"/> I10455 | Ergo Back Fixed 14" W x 21" H - Black Mesh | 874.00 | E2620 |
| <input type="checkbox"/> I10456 | Ergo Back Fixed 14"W x 26.5" High - Black Mesh | 874.00 | E2620 |
| <input type="checkbox"/> I10457 | Ergo Back Sliding 16"W x 21"H - Black Mesh | 874.00 | E2620 |
| <input type="checkbox"/> I10458 | Ergo Back Sliding 16"W x 26.5" H - Black Mesh | 874.00 | E2620 |
| <input type="checkbox"/> I10459 | Ergo Back Sliding 18"W x 21"H - Black Mesh | 874.00 | E2620 |
| <input type="checkbox"/> I10460 | Ergo Back Sliding 18"W x 26.5"H - Black Mesh | 874.00 | E2620 |

C500 Lowrider

Backrest Options

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|--|-----------|------------|
| <input type="checkbox"/> I10004 | Rectangular Backrest 16"W - Plastic, Fixed <i>Choose backrest height from seat pan (19"-29"): _____ "</i> | 1,139.00 | |
| <input type="checkbox"/> I10005 | Rectangular Back 18"W - Metal Fixed <i>Choose backrest height from seat pan (19"-29"): _____ "</i> | 1,330.00 | |
| <input type="checkbox"/> I10240 | Solid Back 16"W X 21"H <i>Standard Option with Tilt and/or Recline Systems. Includes Leatherette Back Cushion ONLY.</i> | No Charge | |
| <input type="checkbox"/> I10130 | Solid Back 16"W x 26.5"H <i>Standard Option with Tilt and/or Recline Systems. Includes Leatherette Back Cushion ONLY.</i> | No Charge | |
| <input type="checkbox"/> I10241 | Solid Back 18"W x 21"H <i>Standard Option with Tilt and/or Recline Systems. Includes Leatherette Back Cushion ONLY.</i> | No Charge | |
| <input type="checkbox"/> I10131 | Solid Back 18"W x 26.5"H <i>Standard Option with Tilt and/or Recline Systems. Includes Leatherette Back Cushion ONLY.</i> | No Charge | |

SELECT OPTIONS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> I10014 | Custom Backrest Height _____" <i>Backrest height measured from seat pan and cannot be lower than 19". Price in addition to above backrest.</i> | 636.00 | |
| <input type="checkbox"/> I100972-99-0 | Split Backrest <i>Split Backrest 8" from top to allow for transfers and/or transportation. This modification cannot be made on back heights lower than 26.5". Must select either PN: I10429 or I10430 below.</i> | 544.00 | |
| <input type="checkbox"/> I10429 | Split Backrest Cushion | No Charge | |
| <input type="checkbox"/> I10430 | Do NOT Split Backrest Cushion | No Charge | |
| <input type="checkbox"/> I10284 | Flat Backrest Cushion for Ergo Back - Leatherette <i>This style of cushion is ONLY available in Leatherette. Recommended cushion when ordering lateral supports. Lateral wedges will be omitted if ordering this part.</i> | 91.00 | |
| <input type="checkbox"/> I10318 | ROHO Backrest Cushion For Ergo Back <i>Price Includes: ROHO MID PROFILE Backrest Cushion which will be matched to fit the Ergonomic Backrest that you select. Approximate height of MID PROFILE cells is 3" in their inflated state. Also includes a Mesh Cover for the ROHO Backrest Cushion, Hand Pump, and Patch Kit. This will NOT fit these part numbers: (I10004, I10005, I10111, I10112, I305510-00-90-0). Please indicate below if you would like to OMIT the Ergo Back Cushion.</i> | 636.00 | |
| <input type="checkbox"/> I10319 | ROHO Backrest Cushion For Ergo Back - Notched <i>Recommended ROHO cushion when ordering Lateral Trunk Supports. Please indicate below if you would like to OMIT the Ergo Back Cushion.</i> | 636.00 | |
| <input type="checkbox"/> I10320 | Omit Ergonomic Backrest Cushion | No Charge | |

Seat Sizing

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> I100681-99-0 | Corpus Seat 17"W x 14"D <i>An additional custom charge will be applied if ordering an 18" wide sliding backrest with this seat size (17" x 14").</i> | No Charge | |
| <input type="checkbox"/> I100684-99-0 | Corpus Seat 17"W x 16"D | No Charge | |
| <input type="checkbox"/> I100689-99-0 | Corpus Seat 17"W x 18"D | No Charge | |
| <input type="checkbox"/> I100695-99-0 | Corpus Seat 17"W x 20"D | No Charge | |
| <input type="checkbox"/> I100720-99-0 | Corpus Seat 17"W x 22"D | No Charge | |
| <input type="checkbox"/> I100690-99-0 | Corpus Seat 19"W x 18"D | No Charge | |
| <input type="checkbox"/> I100696-99-0 | Corpus Seat 19"W x 20"D | No Charge | |
| <input type="checkbox"/> I100721-99-0 | Corpus Seat 19"W x 22"D | No Charge | |
| <input type="checkbox"/> I10007 | Special Seat Size 21"W x 18"D | No Charge | |
| <input type="checkbox"/> I10008 | Special Seat Size 21"W x 20"D | No Charge | |
| <input type="checkbox"/> I10009 | Special Seat Size 21"W x 22"D | No Charge | |

C500 Lowrider

Seat Sizing

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|-------------------------------|-----------|------------|
| <input type="checkbox"/> I10010 | Special Seat Size 23"W x 18"D | No Charge | |
| <input type="checkbox"/> I10011 | Special Seat Size 23"W x 20"D | No Charge | |
| <input type="checkbox"/> I10012 | Special Seat Size 23"W x 22"D | No Charge | |

SELECT OPTIONS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|--|--------|-------------|
| <input type="checkbox"/> I10038 | Ergonomic Seat Cushion Leatherette <i>Not available with seat widths less than 16" or seat depths less than 14".</i> | 483.00 | E2601,E2602 |
| <input type="checkbox"/> I10461 | Ergonomic Seat Cushion - Black Mesh <i>Not available with seat widths less than 16" or seat depths less than 14".</i> | 493.00 | E2601 |

Armrest Options

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> I103339-99-0 | 4"x10" Corpus Arm Right - Dual Taper, Leatherette | No Charge | |
| <input type="checkbox"/> I103341-99-0 | 4"x12" Corpus Arm Right - Dual Taper, Leatherette | No Charge | |
| <input type="checkbox"/> I103343-99-0 | 4"x14" Corpus Arm Right - Dual Taper, Leatherette | No Charge | |
| <input type="checkbox"/> I103345-99-0 | 4"x16" Corpus Arm Right - Dual Taper, Leatherette | No Charge | |
| <input type="checkbox"/> I103331-99-0 | 4"x10" Corpus Arm Right - Dual Taper, Black Fabric | No Charge | |
| <input type="checkbox"/> I103333-99-0 | 4"x12" Corpus Arm Right - Dual Taper, Black Fabric | No Charge | |
| <input type="checkbox"/> I103335-99-0 | 4"x14" Corpus Arm Right - Dual Taper, Black Fabric | No Charge | |
| <input type="checkbox"/> I103337-99-0 | 4"x16" Corpus Arm Right - Dual Taper, Black Fabric | No Charge | |

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> I103338-99-0 | 4"x10" Corpus Arm Left - Dual Taper, Leatherette | No Charge | |
| <input type="checkbox"/> I103340-99-0 | 4"x12" Corpus Arm Left - Dual Taper, Leatherette | No Charge | |
| <input type="checkbox"/> I103342-99-0 | 4"x14" Corpus Arm Left - Dual Taper, Leatherette | No Charge | |
| <input type="checkbox"/> I103344-99-0 | 4"x16" Corpus Arm Left - Dual Taper, Leatherette | No Charge | |
| <input type="checkbox"/> I103330-99-0 | 4"x10" Corpus Arm Left - Dual Taper, Black Fabric | No Charge | |
| <input type="checkbox"/> I103332-99-0 | 4"x12" Corpus Arm Left - Dual Taper, Black Fabric | No Charge | |
| <input type="checkbox"/> I103334-99-0 | 4"x14" Corpus Arm Left - Dual Taper, Black Fabric | No Charge | |
| <input type="checkbox"/> I103336-99-0 | 4"x16" Corpus Arm Left - Dual Taper, Black Fabric | No Charge | |

SELECT OPTIONS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> I10137 | Height Adjustable Armrest Assembly | 289.00 | E0973 |
| <input type="checkbox"/> I10016 | Height Adjustable Armrest Assembly - No Charge <i>No charge with power tilt and/or recline.</i> | No Charge | |
| <input type="checkbox"/> I101003-99-0 | Armrest Mount High <i>Range of adjustment in high position is 11" - 13.5" measured from the seat pan to the top of the armrest pad.</i> | No Charge | |
| <input type="checkbox"/> I101002-99-0 | Armrest Mount Low (High Position Standard) <i>Range of adjustment in low position is 9" - 11.5" measured from the seat pan to the top of the armrest pad.</i> | No Charge | |
| <input type="checkbox"/> I304968-99-0 | Corpus Armbar for 18" Btw the Arms - Installed | 270.00 | |
| <input type="checkbox"/> I10001 | Corpus Armbar for 20" Btw the Arms - Installed | 300.00 | |
| <input type="checkbox"/> I100523-99-0 | Armrest Height Ext (+2") | 127.00 | |
| <input type="checkbox"/> IARMADAP1" | Arm Trough Adaptor 1" Pair <i>For aftermarket arm troughs. Please indicate below if they should be installed or on chair sent as parts.</i> | 279.00 | |

C500 Lowrider

Armrest Options

SELECT OPTIONS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|--------|------------|
| <input type="checkbox"/> IARMADAP7/8" | Arm Trough Adaptor 7/8" Pair <i>For aftermarket arm troughs. Please indicate below if they should be installed on chair or sent as parts.</i> Install Arm Trough Adapters: _____ Send Arm Trough Adapters as Parts: _____ | 279.00 | |
| <input type="checkbox"/> I100671-99-0 | Armrest Pouch-Right | 193.00 | |
| <input type="checkbox"/> I100672-99-0 | Armrest Pouch-Left | 193.00 | |
| <input type="checkbox"/> I101261-99-0 | Armrest Side Support Right Outer <i>For Corpus Armrest Pads only.</i> | 144.00 | |
| <input type="checkbox"/> I101262-99-0 | Armrest Side Support Left Outer <i>For Corpus Armrest Pads only.</i> | 144.00 | |
| <input type="checkbox"/> I101263-99-0 | Armrest Side Support Right Inner <i>Recommend 2" wider width between the armrest to ensure armrests will flip-up.</i> | 144.00 | |
| <input type="checkbox"/> I101264-99-0 | Armrest Side Support Left Inner <i>Recommend 2" wider width between the armrest to ensure armrests will flip-up.</i> | 144.00 | |
| <input type="checkbox"/> I101510-99-0 | Armrest Supports HD <i>Cannot be ordered with 175 degree power recline seating systems. (I10482, I10485, I10488).</i> | 495.00 | |

Legrest Options

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> I102206-99-0 | Manual Adjustment Legrest Elevation | No Charge | |
| <input type="checkbox"/> I103994-99-0 | Power Articulating Center Mount Elevating Legrest <i>Includes calf supports.</i> | 2,995.00 | K0108 |

SELECT OPTIONS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> I103055-99-0 | UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad. Recommended pad to be used with lower leg lengths of 12" or less.</i> | No Charge | |
| <input type="checkbox"/> I102982-99-0 | UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad.</i> | No Charge | |
| <input type="checkbox"/> I102978-99-0 | UT Calf Support Kit (H) w Mesh Cover & Adj Hardware <i>7"W x 7"T Curved Pad.</i> | 317.00 | E1028 |
| <input type="checkbox"/> I101276-99-0 | Legrest Assm. 2pc Footplates (Complete) | No Charge | |
| <input type="checkbox"/> I101277-99-0 | Legrest Assm. 1pc Footplate | 402.00 | |
| <input type="checkbox"/> I10006 | 3" Width Between Two- Piece Footplates | 504.00 | |

Not available if chair has anti-tips and lower leg length is greater than 14" (measured from seat pan to top of footplate.)

Accessories

SELECT OPTIONS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> I10400 | Bodypoint Non-Padded Hip Belt, UniTrack - Medium | No Charge | |
| <input type="checkbox"/> I10405 | Bodypoint Padded Hip Belt UniTrack (N/C) - Medium | No Charge | |
| <input type="checkbox"/> I10402 | Bodypoint Non-Padded Hip Belt, UniTrack - Large | 94.00 | K0108 |
| <input type="checkbox"/> I10408 | Bodypoint Padded Hip Belt UniTrack - Large | 158.00 | K0108 |
| <input type="checkbox"/> I100519-99-0 | Positioning Belt Retractable | 239.00 | K0108 |
| <input type="checkbox"/> I100600-99-0 | Positioning Belt Retractable Long | 259.00 | K0108 |
| <input type="checkbox"/> I104131-99-0 | Headrest for Ergo Back w/ Curved Mount Brkt - Leath <i>Must also select PN I10552. This part consists of a leatherette headrest pad (10"W x 5"T) with a 10"L slotted, curved mounting bracket.</i> | 268.00 | E0955 |

C500 Lowrider

Accessories

SELECT OPTIONS:

| Part Number | Description | Price | HCPCS Code |
|--|--|--------|------------|
| <input type="checkbox"/> I103460-99-0 | UT Ergo Headrest w/ Black Fabric Cover <i>Must also select PN I10552.</i> | 297.00 | E0955 |
| <input type="checkbox"/> I1823249 | UT Lateral Adjustment Bar for UniTrack Headrest <i>Allows for approximately 4" of lateral adjustment to the UT Ergo Headrest either right or left. NOTE: Due to the width of this item (approx 9"), it may be necessary to send other items such as back pack clips unmounted. Custom charge may apply with lower back heights.</i> | 56.00 | |
| <input type="checkbox"/> IHEADADAPTER | Universal Headrest Adapter <i>For aftermarket headrests.</i> | 205.00 | K0108 |
| <input type="checkbox"/> I10552 | Adj Removable - Headrest Hardware (1 Each) | 216.00 | E1028 |
| <input type="checkbox"/> I103070-99-0 | UT Amp. Support Pad (H) w/ Mesh Cover - (1 Each) <i>7"W x 7"T Curved Pad. Must also select PN I10553.</i> | 127.00 | E1020 |
| <input type="checkbox"/> I103071-99-0 | UT Amp. Support Pad (D) w/ Mesh Cover - (1 Each) <i>8"W x 8"T Flat Pad. Must also select PN I10553.</i> | 127.00 | E1020 |
| <input type="checkbox"/> I10553 | Adj Removable - Amp Support Hardware (1 Each) <i>Please indicate if you would like the amputation support mounted on the right or left side of the seat frame: Mount Amp Suppt on Right: ___ Mount Amp Suppt on Left: ___</i> | 216.00 | E1028 |
| <input type="checkbox"/> I103077-99-0 | UT Hip Suppt (B) w/ Mesh Cover - Low (Pair) <i>4.5"W x 3.5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i> | 229.00 | E0956 |
| <input type="checkbox"/> I103079-99-0 | UT Hip Suppt (A) w/ Mesh Cover - Low (Pair) <i>6"W x 5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i> | 229.00 | E0956 |
| <input type="checkbox"/> I103094-99-0 | UT Thigh/Hip Suppt (J) w/ Mesh Cover - High (Pair) <i>6.5"W x 5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i> | 229.00 | E0956 |
| <input type="checkbox"/> I103095-99-0 | UT Thigh/Hip Suppt (C) w/ Mesh Cover - High (Pair) <i>8"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i> | 229.00 | E0956 |
| <input type="checkbox"/> I102977-99-0 | UT Thigh/Hip Suppt (G) w/ Mesh Cover - High (Pair) <i>12"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i> | 229.00 | E0956 |
| <input type="checkbox"/> I10550 | Adj Removable - Thigh Support Hardware (Pair) | 431.00 | E1028 |
| <input type="checkbox"/> I305266-28-90-0 | Lateral Supports Std. Leatherette-ErgoBack Only <i>Must also select PN I10551.</i> | 250.00 | E0956 |
| <input type="checkbox"/> I100521-99-0 | Swing Away Trunk Supports Leatherette-For Ergo Back <i>Must also select PN I10551.</i> | 254.00 | E0956 |
| <input type="checkbox"/> I10551 | Adj Removable - Trunk Support Hardware (Pair) | 431.00 | E1028 |
| <input type="checkbox"/> I100978-99-0 | Abductor Kit - Leatherette <i>Must also select PN I10555.</i> | 311.00 | E0957 |
| <input type="checkbox"/> I10555 | Adj Removable - Abductor Hardware (1 Each) | 216.00 | E1028 |
| <input type="checkbox"/> I10015 | Knee Block Assembly Corpus | 641.00 | K0108 |
| <input type="checkbox"/> I103651-99-0 | Transfer Handles for UniTrack Accessory Rail | 397.00 | |
| <input type="checkbox"/> I100975-99-0 | Push Handles for Ergo Backrest <i>Order with Ergonomic or Rectangular Backs. Custom charge may apply if ordering push handles with an articulating vent tray.</i> | 345.00 | |
| <input type="checkbox"/> I102059-99-0 | Medical Necessities Bag <i>Includes Medical Necessity Bag Clips (Carabiner Style).</i> | 229.00 | |
| <input type="checkbox"/> I101473-99-0 | Medical Necessity Bag Clips - Carabiner Style | 76.00 | |
| <input type="checkbox"/> I101500-99-0 | Essentials Carrier for Chassis <i>It is recommended that you also select Rear Tie Down Points Extension Kit (I101515-99-0) below when selecting this item.</i> | 386.00 | |
| <input type="checkbox"/> I101515-99-0 | Rear Tie Down Points Extension Kit <i>Recommended when ordering Essentials Carrier for Chassis (I101500-99-0) in order to access rear tiedown points.</i> | 115.00 | |
| <input type="checkbox"/> I105144-99-0 | Bodypoint Monoflex Large Chest Support Belt Kit <i>21.5"L x 4.5"W Pad Size. Item will be shipped as parts.</i> | 145.00 | E0960 |

C500 Lowrider

Accessories

SELECT OPTIONS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|--------|------------|
| <input type="checkbox"/> I105143-99-0 | Bodypoint Monoflex Medium Chest Support Belt Kit <i>18.25"L x 4"W Pad Size. Item will be shipped as parts.</i> | 145.00 | E0960 |
| <input type="checkbox"/> I105145-99-0 | Bodypoint Monoflex X-Lrge Chest Support Belt Kit <i>24.5"L x 5"W Pad Size. Item will be shipped as parts.</i> | 145.00 | E0960 |
| <input type="checkbox"/> I100578-99-0 | Upper Extremity Support Small, Angle Adjustable <i>Must also select PN I10556.</i> | 242.00 | E0950 |
| <input type="checkbox"/> I100575-99-0 | Upper Extremity Support Angle Adjustable. 14" W <i>Must also select PN I10556.</i> | 272.00 | E0950 |
| <input type="checkbox"/> I100576-99-0 | Upper Extremity Support Angle Adjustable. 16" W <i>Must also select PN I10556.</i> | 284.00 | E0950 |
| <input type="checkbox"/> I100577-99-0 | Upper Extremity Support Angle Adjustable. 18" W <i>Must also select PN I10556.</i> | 311.00 | E0950 |
| <input type="checkbox"/> I100571-99-0 | Upper Extremity Support Joystick Cutout 16" W <i>Must also select PN I10556.</i> | 328.00 | E0950 |
| <input type="checkbox"/> I100572-99-0 | Upper Extremity Support Joystick Cutout 18" W <i>Must also select PN I10556.</i> | 328.00 | E0950 |
| <input type="checkbox"/> I100573-99-0 | Upper Extremity Support Mid Joystick Cutout 16"W <i>Must also select PN I10556.</i> | 402.00 | E0950 |
| <input type="checkbox"/> I100574-99-0 | Upper Extremity Support Mid Joystick Cutout 18"W <i>Must also select PN I10556.</i> | 402.00 | E0950 |
| <input type="checkbox"/> I103426-99-0 | Upper Extremity Support Angle Adjust - UniTrack <i>Side mounted tray that attaches to the UniTrack Accessory Rail. Dimensions are: 19.5"W x 13"D. Can only be ordered with chairs that offer the UniTrack System. Must also select PN I10556.</i> | 301.00 | E0950 |
| <input type="checkbox"/> I10556 | Adj Removable - Tray Hardware (1 Each) | 216.00 | E1028 |
| <input type="checkbox"/> I101964-99-0 | Oxygen Holder for PS or Corpus Seat <i>Cannot be ordered with a Corpus Tubular Backrest. (I10111 or I10112). Custom charge may be applied if ordering in combination with accessories such as medical necessities bags, attendant controls, push handles, back pack clips, etc. This item is mounted to the Ergonomic Backrest. Customization may be required if ordering on a Street.</i> | 228.00 | E2208 |
| <input type="checkbox"/> I101956-99-0 | Crutch Holder for PS or Corpus Seat <i>Cannot be ordered with a Corpus Tubular Backrest. (I10111 or I10112). Accessories such as medical necessities bags and push handles may need to be sent as parts if ordering a crutch holder. This item is mounted to the ergonomic backrest. Customization may be required if ordering on a Street.</i> | 228.00 | E2207 |

Alternative Drive Systems

(1) CHOOSE ALTERNATIVE DRIVE INPUT:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|--|----------|------------|
| | R-net Alternative Drive Input Devices <i>R-net Alternative Drive Input Devices take the place of, or may be used in addition to, the standard R-net Joystick. They are only usable on wheelchairs that have R-net electronics. You MUST select one of each of the following items: (1) Input Device, (2) Adjustable/Removable Alt Drive Hardware, AND (3) Omni or Input Module.</i> | | |
| <input type="checkbox"/> I10558 | Gimble Joystick set up for Armrest Mount | 1,412.00 | E2373 |
| <input type="checkbox"/> I10559 | Gimble Joystick set up for Midline Mount | 1,412.00 | E2373 |
| <input type="checkbox"/> I10560 | Gimble Joystick set up for Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i> | 1,412.00 | E2373 |
| <input type="checkbox"/> I10561 | Gimble Joystick w/ Haims Harness - Small <i>For neck width up to 5".</i> | 1,582.00 | E2373 |
| <input type="checkbox"/> I10562 | Gimble Joystick w/ Haims Harness - Large <i>For neck width up to 7".</i> | 1,582.00 | E2373 |

C500 Lowrider

Alternative Drive Systems

(1) CHOOSE ALTERNATIVE DRIVE INPUT:

| Part Number | Description | Price | HCPCS Code |
|-----------------------------------|---|----------|-------------|
| <input type="checkbox"/> I10563 | Gimble Joystick w/ Power Swing Away Chin Boom <i>Custom charge may be applied if ordering smaller backrest sizes and heights.</i> | 5,265.00 | E2373 |
| <input type="checkbox"/> I10564 | Mini Joystick set up for Armrest Mount | 2,685.00 | E2312 |
| <input type="checkbox"/> I10565 | Mini Joystick set up for Midline Mount | 2,685.00 | E2312 |
| <input type="checkbox"/> I10566 | Mini Joystick set up for Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i> | 2,685.00 | E2312 |
| <input type="checkbox"/> I10567 | Mini Joystick w/ Haims Harness - Small <i>For neck width up to 5".</i> | 2,745.00 | E2312 |
| <input type="checkbox"/> I10568 | Mini Joystick w/ Haims Harness - Large <i>For neck width up to 7".</i> | 2,745.00 | E2312 |
| <input type="checkbox"/> I10569 | Mini Joystick w/ Power Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i> | 5,265.00 | E2312 |
| <input type="checkbox"/> I1821750 | Tablet Control Pad <i>The Tablet Control is a unique drive input which allows a user to control a power wheelchair in a very similar manner to a mouse being controlled by a touchpad on a laptop computer. The dimensions of the Tablet Control Pad are: 12"W x 9.5"D x 1.5"H. Kit also includes: Buddy Button Switch, and Tablet Control Mounting Bracket.</i> | 4,375.00 | E2328 |
| <input type="checkbox"/> I10570 | Head Array - Adult 5 Switch HDR-A5 | 4,375.00 | E2330 |
| <input type="checkbox"/> I10571 | Head Array - Adult 3 Switch HDR-A4 <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i> | 4,175.00 | E2330 |
| <input type="checkbox"/> I10572 | Head Array - Pediatric 5 Switch HDR-C5-RR | 4,375.00 | E2330 |
| <input type="checkbox"/> I10573 | Sip & Puff System <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i> | 1,475.00 | E2325 |
| <input type="checkbox"/> I10574 | Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i> | 1,963.00 | E2325,E2326 |

SELECT OPTIONS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|--------------------------------------|-----------|------------|
| <input type="checkbox"/> I10076 | Right Side Mount For Joystick | No Charge | |
| <input type="checkbox"/> I10077 | Left Side Mount For Joystick | No Charge | |
| <input type="checkbox"/> I10074 | Right Side Mount For Omni | No Charge | |
| <input type="checkbox"/> I10075 | Left Side Mount For Omni | No Charge | |
| <input type="checkbox"/> I10078 | Mount Seat Control Panel on Right | No Charge | |
| <input type="checkbox"/> I10079 | Mount Seat Control Panel on Left | No Charge | |
| <input type="checkbox"/> I10080 | Mount Seat Control Panel on Backrest | No Charge | |

Alt Drive Systems Setup

(2) CHOOSE ADJ ALT DRIVE HARDWARE (REQUIRED):

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|---|--------|------------|
| <input type="checkbox"/> I10577 | Adjustable Removable Hardware - Alt Drive | 216.00 | E1028 |

(3) CHOOSE OMNI OR INPUT MODULE (REQUIRED):

| Part Number | Description | Price | HCPCS Code |
|------------------------------------|---|----------|------------|
| <input type="checkbox"/> ID51154IR | Omni Control Module for R-net w/ Infrared <i>Please specify which type of alternative drive input you will be using if not ordering the input device from Permobil. Input Device: _____</i> | 1,496.00 | |
| <input type="checkbox"/> I10581 | Input Module (IOM) for R-net Control System <i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head Array or RIM control.</i> | 1,012.00 | E2351 |

Alt Drive Accessories

C500 Lowrider

Alt Drive Accessories

SELECT OPTIONS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|-----------------------------|--------|------------|
| <input type="checkbox"/> IMICROSWITCH | Microlite Switch Qty: _____ | 151.00 | |
| <input type="checkbox"/> IBUDDYBUTTON | Buddy Button Qty: _____ | 113.00 | K0108 |
| <input type="checkbox"/> I5808 | Mini Cup Switch Qty: _____ | 110.00 | |
| <input type="checkbox"/> IZTS-P | Proximity Switch Qty: _____ | 524.00 | |
| <input type="checkbox"/> IMINIJOY5SW | Mini Joystick 5 Switch OC | 603.00 | |
| <input type="checkbox"/> I5018 | Wafer Board w/5 Switches | 506.00 | |

Notes & Additional Instructions:

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, Inc. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. MID PROFILE® Cushion is a registered trademark of the The ROHO Group. The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.