

C400 Corpus Junior



*** Indicates a Required Field**



Dealer Information

* **Contact:** _____

Dealer Code: _____

* **Dealer Name:** _____

Address: _____

* **City:** _____

* **State/Zip:** _____

* **Phone#:** _____

* **Fax#:** _____

PO#: _____

Email Address: _____

Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

* **First Name:** _____

* **Last Name:** _____

Diagnosis: _____

Funding Source: _____

Client Age: _____

Client Measurements

* **Weight:** _____

* **Height:** _____

A) Top of Shoulders: _____

B) Chest Depth: _____

C) Chest Width: _____

D) Seat Depth: _____


E) Top of Head: _____

F) Elbow to Hand: _____

G) Seat Pan to Elbow: _____

H) Hip Width: _____

* **I) Knee to Foot:** _____



Three column pricing on order/quote:

Display HCPCS Codes on order/quote:

Please send order/quote to fax#: **(800) 231-3256**

Email to: sales@permobilus.com

Permobil Inc.
6961 Eastgate Blvd.
Lebanon, TN 37090
Tel: (800) 736-0925
Fax: (800) 231-3256
www.permobil.com

Prices effective **March 8, 2010.**

Chair Model and Colors

Select Options:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10475	C400 Base Corpus R-net - C400 PS0 <i>Price Includes: FWD C400 Power Base Set Up For No Power Seat Functions, Shock Absorbing Suspension System, R-net Basic Controller, Tie-Down Hardware for Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below.</i>	8,995.00	K0848
<input type="checkbox"/> I10476	C400 Base Corpus R-net - C400 PS1 <i>Price Includes: FWD C400 Power Base Set Up For One Power Seat Function, Shock Absorbing Suspension System, R-net Controller, Tie Down Hardware for Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below.</i>	8,995.00	K0856
<input type="checkbox"/> I103259-99-0	C400 Base Corpus R-net - C400 PS2 <i>Price Includes: FWD C400 Power Base Set Up For Multiple Power Seat Functions, Shock Absorbing Suspension System, Tie-Down Hardware for Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below.</i>	8,995.00	K0861

Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103889-99-0	Tri Spoke Split Rim w/ FF Drive Tires - 2.5 x 8"	No Charge	
<input type="checkbox"/> I103890-99-0	Tri Spoke Split Rim w/ Pneumatic Tire - 2.5 x 8"	No Charge	

Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I102190-99-0	C400 Shroud Vintage Wine	No Charge	
<input type="checkbox"/> I103957-99-0	C400 Shroud Onyx Black	No Charge	
<input type="checkbox"/> I102194-99-0	C400 Shroud Crystal Blue	No Charge	
<input type="checkbox"/> I102191-99-0	C400 Shroud Metallic Charcoal Gray	No Charge	
<input type="checkbox"/> I102192-99-0	C400 Shroud Green Metallic	No Charge	
<input type="checkbox"/> I103958-99-0	C400 Shroud Extreme Yellow	No Charge	

Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10451	Match Hub Cap to Shroud Color <i>If this is selected, the Tri-Spoke Hub Cap will be matched to the selected shroud color.</i>	No Charge	
<input type="checkbox"/> I103537-99-0	Hub Cap Platinum Silver <i>If this is selected, the Tri-Spoke Hub Cap will be Platinum Silver.</i>	No Charge	

Base Options

Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103865-99-0	Power Adjustable Seat Height, 8" Travel (R-net) <i>To assist in transfers and improve access for function.</i>	2,595.00	E2300
<input type="checkbox"/> I102085-99-0	Seat Tube Fixed	No Charge	

Select Options:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> IM24 SLD G FT	Batteries, Group 24 Sealed LA/Gel, Installed	762.00	E2363
<input type="checkbox"/> I102210-99-0	Anti-Tippers for C400/C500	399.00	
<input type="checkbox"/> I103635-99-0	Light Kit, R-net, C400/C500	955.00	
<input checked="" type="checkbox"/> I100570-99-0	Seat Stabilizer	No Charge	

C400 Corpus Junior**MSRP Price List***Select Options:*

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103241-99-0	Articulating Vent Tray for C400/C500 <i>If selecting this part, you MUST also pick: PN I104082-99-0 - ESP Unit. Fits an LP10 Sized Vent and a Group 34 Battery. Cannot be ordered with backrest height less than 26.5", a tubular back, or a 13" seat depth. Customization may be required if ordered with any of the following items Co-Pilot, any standard attendant control, push handles, or an oxygen holder. Only available with expandable controllers (R-net Electronics.)</i>	1,295.00	E1030
<input checked="" type="checkbox"/> I104082-99-0	R-net ESP w C400/C500 Mounting Bracket <i>MUST be selected if ordering any ventilator tray (I102769-99-0, I103241-99-0, or I102768-99-0). Recommended item with switched drive inputs for improved tracking.</i>	No Charge	K0108
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	435.00	

Joystick Options*Choose One of the Following*

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard.	No Charge	
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	94.00	E2324
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	139.00	E2323
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	94.00	E2323
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	127.00	E2323
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	127.00	E2323
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	127.00	E2323
<input type="checkbox"/> IPC101B	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	127.00	E2323
<input type="checkbox"/> IPC102B	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	127.00	E2323
<input type="checkbox"/> IPC107B	Bodypoint J/S Handle Dome Shaped, Rubber	127.00	E2323

Select Options:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10412	Joystick Mount Right - Fixed, R-net Only	No Charge	
<input type="checkbox"/> I10413	Joystick Mount Left - Fixed, R-net Only	No Charge	
<input type="checkbox"/> I102537-99-0	Joystick Mount Right - Swing Away, R-net & VR2	292.00	E1028
<input type="checkbox"/> I102538-99-0	Joystick Mount Left - Swing Away, R-net & VR2	292.00	E1028
<input type="checkbox"/> I103820-99-0	Retractable Joystick Mnt Right - R-net, VR2 Only	486.00	E1028
<input type="checkbox"/> I103819-99-0	Retractable Joystick Mnt Left - R-net, VR2 Only	486.00	E1028
<input type="checkbox"/> I101260-99-0	Midline Control Mount - Pilot+ & R-net Only <i>Please indicate which armrest to mount midline below.</i>	263.00	E1028
<input type="checkbox"/> I10428	Mount Midline Joystick on Left Armrest	No Charge	
<input type="checkbox"/> I10427	Mount Midline Joystick on Right Armrest	No Charge	
<input type="checkbox"/> I10410	R-net Remote Joystick - No Charge <i>Must order expandable controller (PN: I10411) with this part. Comes standard w/ mono jacks.</i>	No Charge	
<input type="checkbox"/> I103544-99-0	R-net Remote Joystick <i>Must order expandable controller (PN: I10411) with this part. Comes standard w/ mono jacks.</i>	530.00	
<input type="checkbox"/> I103814-99-0	Compact Advanced Joystick for R-net <i>Compact sized joystick that does not contain an LCD Display. This is an expandable/remote joystick used with multiple power option bases or as an upgrade for no power option or single power option bases. Must order expandable controller (PN: I10411) with this part. Contains two mono jacks for mode & on/off. Does NOT contain a wrist support pad. If you want a wrist support pad, please select I10621, below.</i>	605.00	E2373
<input type="checkbox"/> I10621	Add Wrist Pad to Compact Joystick Advanced	195.00	
<input type="checkbox"/> ID51040	R-net Remote Joystick w/ Toggle Switches <i>Expandable/remote joystick with toggle switches. Used with multiple power option bases or as an upgrade for no power option or single power option bases. Must order expandable controller (PN: I10411) with this part.</i>	605.00	

C400 Corpus Junior

MSRP Price List

Select Options:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> ID51042	R-net Remote Joystick w/ Speed Pot & Single Toggle <i>Expandable/remote joystick with a single toggle switch for modes/profiles & a speed pot. Used with multiple power option bases or as an upgrade for no power option or single power option bases. Must order expandable controller (PN: I10411) with this part.</i>	730.00	
<input type="checkbox"/> ID51153	R-net Remote Color Screen Joystick w/ Mono Jacks <i>Must order expandable controller (PN: I10411) with this part. Comes standard w/ mono jacks.</i>	835.00	
<input type="checkbox"/> ID51149	R-net Remote Color Screen JS w/ Toggle Switches <i>Must order expandable controller (PN: I10411) with this part.</i>	875.00	
<input type="checkbox"/> ID51190	R-net Remote Color JS w/ Speed Pot & Single Toggle <i>Must order expandable controller (PN: I10411) with this part.</i>	875.00	
<input type="checkbox"/> I10411	Expandable Controller - R-net <i>Must be ordered with the expandable/remote joysticks (PN: I10410, I103544-99-0, ID51040, ID51108, ID51042, and I103814-99-0).</i>	525.00	E2377
<input type="checkbox"/> I10513	Harness for Expandable Controller <i>Must be ordered with any Expandable Controller (I10411 or I10450)</i>	445.00	E2313
<input type="checkbox"/> I104003-99-0	Co-Pilot Attendant Control, Corpus R-net <i>Custom charge may apply if ordering a Corpus Tubular Backrest OR a backrest height of 21" or Lower. Not available with VR2 Controllers.</i>	4,175.00	E2331
<input type="checkbox"/> ID50882	Attendant Control for R-net Control System Only <i>Not available with VR2 Controllers.</i>	1,084.00	E2331
<input type="checkbox"/> ID51111	R-net BT Mouse Module <i>The R-net BT Mouse Module allows for mouse control through a standard R-net Joystick or other type of input device (ex: Mini Joystick w/ Omni). It is only usable as a mouse on a computer. It CANNOT be paired to work with other Bluetooth devices such as cell phones, wireless headsets, etc. This kit includes: an R-net BT Mouse Module, a 1.5m R-net Bus Cable, and a USB Style Bluetooth Adapter for the user's computer. Note: The Bluetooth Adapter may not work with all computers.</i>	795.00	
<input type="checkbox"/> I10582	Output Module (IOM) for R-net Control System <i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). One output module is required for each separate external device to be controlled. Up to three output modules can be installed on one wheelchair.</i>	995.00	E2351

Seating System

Choose One of the Following

Part Number	Description	Price	HCPCS Code
	The Seating Systems Listed Below Include: <i>Corpus Seat Frame with 1" Depth Adjustments from 13" - 18" (& up to 22" Deep with Standard Corpus Seat Frame) Back Frame with Adjustable Recline Angle, Reclining Armrest Assembly, Solid Seat Pan, Legrest Assembly with Adjustable Legrest Elevation, Weight Limit 165 lbs for Seat Depths Less Than 16". Weight Limit of 250 lbs for Seat Depths of 16" or Greater.</i>		
<input type="checkbox"/> I10491	Corpus Junior Seat w/o Pwr Tilt or Pwr Recline <i>Includes:Solid Backrest, Height Adjustable Armrest and Seat Sizes Adjustable in 1" Increments</i>	5,689.00	K0108
<input type="checkbox"/> I10492	Corpus Junior Seat w/ 45° Pwr Tilt & Manual Recline <i>Includes:Solid Backrest, Height Adjustable Armrest and Seat Sizes Adjustable in 1" Increments</i>	6,830.00	E1002
<input type="checkbox"/> I10493	Corpus Junior Seat w 150° Power Recline <i>Includes:Solid Backrest, Height Adjustable Armrest and Seat Sizes Adjustable in 1" Increments</i>	7,356.00	E1004
<input type="checkbox"/> I10494	Corpus Jr Seat w 45° Pwr Tilt & 150° Pwr Recline <i>Includes:Solid Backrest, Height Adjustable Armrest and Seat Sizes Adjustable in 1" Increments</i>	10,695.00	E1007

Seating System Options

Select Options:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10418	Single Seat Function Control Switch for R-net <i>Please choose either Push Buttons or Toggle Switches below.</i>	No Charge	
<input type="checkbox"/> I10419	Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control and/or ICS Switchbox. Order with R-net Controller or remote joystick options. Please indicate below if you want push buttons or toggles on the Seat Function Control Kit.</i>	1,478.00	E2310
<input type="checkbox"/> I10585	Multiple Seat Function Control Kit for R-net <i>MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive control and/or ICS switch box. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the Seat Function Control Kit.</i>	2,875.00	E2311
<input type="checkbox"/> I102798-99-0	Push Buttons on Seat Function Control Kit	No Charge	
<input type="checkbox"/> I102799-99-0	Toggle Switches on Seat Function Control Kit	245.00	
<input type="checkbox"/> I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge	
<input type="checkbox"/> I103471-99-0	Mount ICS Switchbox to End of Arm Pad - Corpus <i>ICS Switchbox will be mounted to the end of the armrest pad. If this item is selected, the ICS Switchbox will NOT swing away with the Joystick Module. Please indicate below if you would like it mounted onto the left or right arm pad.</i> Mount ICS Switchbox on Left Armpad _____ Mount ICS Switchbox on Right Armpad _____	No Charge	
<input type="checkbox"/> I103470-99-0	Mount ICS Switchbox on Backrest	125.00	
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Not available with VR2 Electronics. Seat function interface box which accommodates up to 8 momentary switches. Provides an additional method to accessing the seat functions. You must still order the Multiple Seat Function Control Kit if ordering this part. 1/8" mini-mono switches not included.</i>	985.00	K0108
<input type="checkbox"/> I104746-99-0	Remote Stop System - R-net, Pilot+ <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Each transmitter and receiver is a serialized matched pair, meaning they must be used together. Custom charge may be required if ordering with a low backrest height or a pediatric wheelchair. If ordering Remote Stop System with a K450, push handles must also be selected.</i>	985.00	
<input type="checkbox"/> I10497	Memory Seat Program, C400 or C500 Corpus w R-net <i>If this is selected, the chair will be programmed to allow the user to store and save up to 3 different memory seat positions. On the C400/C500 Corpus and Lowrider, the following seat functions are capable of remembering a seat position: Power Tilt, Power Recline, & Power Legs. The Power Adj Seat Height Function is NOT capable of recalling a saved memory seat position. Please refer to the Corpus Seat Owner's Manual for more information on saving and recalling memory seat positions.</i>	98.00	

Backrest Options

Choose One of the Following

Part Number	Description	Price	HCPCS Code
	BACKREST HEIGHT NOTE: <i>Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°.</i>		
<input type="checkbox"/> I101853-99-0	Ergo Back Fixed 14" W x 21" H - Leatherette	837.00	E2620
<input type="checkbox"/> I101863-99-0	Ergo Back Fixed 14"W x 26.5" High - Leatherette	837.00	E2620
<input type="checkbox"/> I101854-99-0	Ergo Back Sliding 16"W x 21"H - Leatherette	837.00	E2620
<input type="checkbox"/> I101861-99-0	Ergo Back Sliding 16"W x 26.5" H - Leatherette	837.00	E2620
<input type="checkbox"/> I10455	Ergo Back Fixed 14" W x 21" H - Black Mesh	859.00	E2620
<input type="checkbox"/> I10456	Ergo Back Fixed 14"W x 26.5" High - Black Mesh	859.00	E2620

C400 Corpus Junior**MSRP Price List**

Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10457	Ergo Back Sliding 16"W x 21"H - Black Mesh	859.00	E2620
<input type="checkbox"/> I10458	Ergo Back Sliding 16"W x 26.5" H - Black Mesh	859.00	E2620
<input type="checkbox"/> I10004	Rectangular Backrest 16"W - Plastic, Fixed <i>Choose backrest height from seat pan (19"-29"): _____ "</i>	1,119.00	
<input type="checkbox"/> I10240	Solid Back 16"W X 21"H <i>Standard Option with Tilt and/or Recline Systems. Includes Leatherette Back Cushion ONLY.</i>	No Charge	
<input type="checkbox"/> I10130	Solid Back 16"W x 26.5"H <i>Standard Option with Tilt and/or Recline Systems. Includes Leatherette Back Cushion ONLY.</i>	No Charge	
<input type="checkbox"/> I10111	Tubular Backrest Frame Corpus, 16" Wide <i>Backrest canes to mount aftermarket backrest systems. For tilt only systems. Measurement between the armrests is 17.5". Seat pan to top of canes is 22".</i>	357.00	

Select Options:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10014	Custom Backrest Height _____" <i>Backrest height measured from seat pan and cannot be lower than 19". Price in addition to above backrest.</i>	625.00	
<input type="checkbox"/> I100972-99-0	Split Backrest <i>Split Backrest 8" from top to allow for transfers and/or transportation. This modification cannot be made on back heights lower than 26.5". Must select either PN: I10429 or I10430 below.</i>	535.00	
<input type="checkbox"/> I10429	Split Backrest Cushion	No Charge	
<input type="checkbox"/> I10430	Do NOT Split Backrest Cushion	No Charge	
<input type="checkbox"/> I10284	Flat Backrest Cushion for Ergo Back - Leatherette <i>This style of cushion is ONLY available in Leatherette. Recommended cushion when ordering lateral supports. Lateral wedges will be omitted if ordering this part.</i>	89.00	
<input type="checkbox"/> I10318	ROHO Backrest Cushion For Ergo Back <i>Price Includes: ROHO MID PROFILE Backrest Cushion which will be matched to fit the Ergonomic Backrest that you select. Approximate height of MID PROFILE cells is 3" in their inflated state. Also includes a Mesh Cover for the ROHO Backrest Cushion, Hand Pump, and Patch Kit. This will NOT fit these part numbers: (I10004, I10005, I10111, I10112, I305510-00-90-0). Please indicate below if you would like to OMIT the Ergo Back Cushion.</i>	625.00	
<input type="checkbox"/> I10319	ROHO Backrest Cushion For Ergo Back - Notched <i>Recommended ROHO cushion when ordering Lateral Trunk Supports. Please indicate below if you would like to OMIT the Ergo Back Cushion.</i>	625.00	
<input type="checkbox"/> I10320	Omit Ergonomic Backrest Cushion	No Charge	
<input type="checkbox"/> I10278	Notch Backrest Cushion for Lateral Supports <i>Backrest Cushion will be notched 2" in on either side to allow for closer adjustment of trunk supports.</i>	185.00	

Seat Sizing

Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10264	Corpus Junior Seat Width - Standard <i>Corpus seat frame is 17" wide at the outside and 15" wide between the seat rails. Actual hip width can be accommodated with size appropriate cushion and/or hip guides.</i>	No Charge	

Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10265	Corpus Junior Seat 13" Deep	No Charge	
<input type="checkbox"/> I10266	Corpus Junior Seat 14" Deep	No Charge	
<input type="checkbox"/> I10267	Corpus Junior Seat 15" Deep	No Charge	
<input type="checkbox"/> I10268	Corpus Junior Seat 16" Deep	No Charge	
<input type="checkbox"/> I10269	Corpus Junior Seat 17" Deep	No Charge	

C400 Corpus Junior**MSRP Price List**

Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10270	Corpus Junior Seat 18" Deep	No Charge	

Select Options:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10038	Ergonomic Seat Cushion Leatherette <i>Not available with seat widths less than 16" or seat depths less than 14".</i>	475.00	E2601,E2602
<input type="checkbox"/> I10461	Ergonomic Seat Cushion - Black Mesh <i>Not available with seat widths less than 16" or seat depths less than 14".</i>	485.00	E2601
<input type="checkbox"/> I10276	Miniflex Seat Cushion Leatherette <i>Approximately 2" thick. Not available with seat widths or depths greater than 14".</i>	475.00	

Armrest Options

Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103339-99-0	4"x10" Corpus Arm Right - Dual Taper, Leatherette	No Charge	
<input type="checkbox"/> I103341-99-0	4"x12" Corpus Arm Right - Dual Taper, Leatherette	No Charge	
<input type="checkbox"/> I103343-99-0	4"x14" Corpus Arm Right - Dual Taper, Leatherette	No Charge	
<input type="checkbox"/> I103331-99-0	4"x10" Corpus Arm Right - Dual Taper, Black Fabric	No Charge	
<input type="checkbox"/> I103333-99-0	4"x12" Corpus Arm Right - Dual Taper, Black Fabric	No Charge	
<input type="checkbox"/> I103335-99-0	4"x14" Corpus Arm Right - Dual Taper, Black Fabric	No Charge	

Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103338-99-0	4"x10" Corpus Arm Left - Dual Taper, Leatherette	No Charge	
<input type="checkbox"/> I103340-99-0	4"x12" Corpus Arm Left - Dual Taper, Leatherette	No Charge	
<input type="checkbox"/> I103342-99-0	4"x14" Corpus Arm Left - Dual Taper, Leatherette	No Charge	
<input type="checkbox"/> I103330-99-0	4"x10" Corpus Arm Left - Dual Taper, Black Fabric	No Charge	
<input type="checkbox"/> I103332-99-0	4"x12" Corpus Arm Left - Dual Taper, Black Fabric	No Charge	
<input type="checkbox"/> I103334-99-0	4"x14" Corpus Arm Left - Dual Taper, Black Fabric	No Charge	

Select Options:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10137	Height Adjustable Armrest Assembly	284.00	E0973
<input type="checkbox"/> I10016	Height Adjustable Armrest Assembly - No Charge <i>No charge with power tilt and/or recline.</i>	No Charge	
<input type="checkbox"/> I101003-99-0	Armrest Mount High <i>Range of adjustment in high position is 11" - 13.5" measured from the seat pan to the top of the armrest pad.</i>	No Charge	
<input type="checkbox"/> I101002-99-0	Armrest Mount Low (High Position Standard) <i>Range of adjustment in low position is 9" - 11.5" measured from the seat pan to the top of the armrest pad.</i>	No Charge	
<input type="checkbox"/> IARMADAP1"	Arm Trough Adaptor 1" Pair <i>For aftermarket arm troughs. Please indicate below if they should be installed on chair sent as parts.</i>	274.00	
<input type="checkbox"/> IARMADAP7/8"	Arm Trough Adaptor 7/8" Pair <i>For aftermarket arm troughs. Please indicate below if they should be installed on chair or sent as parts.</i>	274.00	
Install Arm Trough Adapters: _____ Send Arm Trough Adapters as Parts: _____			
<input type="checkbox"/> I100671-99-0	Armrest Pouch-Right	190.00	
<input type="checkbox"/> I100672-99-0	Armrest Pouch-Left	190.00	
<input type="checkbox"/> I101261-99-0	Armrest Side Support Right Outer <i>For Corpus Armrest Pads only.</i>	142.00	

C400 Corpus Junior**MSRP Price List***Select Options:*

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I101262-99-0	Armrest Side Support Left Outer <i>For Corpus Armrest Pads only.</i>	142.00	
<input type="checkbox"/> I101263-99-0	Armrest Side Support Right Inner <i>Recommend 2" wider width between the armrest to ensure armrests will flip-up.</i>	142.00	
<input type="checkbox"/> I101264-99-0	Armrest Side Support Left Inner <i>Recommend 2" wider width between the armrest to ensure armrests will flip-up.</i>	142.00	
<input type="checkbox"/> I10285	Mount Corpus Armrest Pads On Inside of Bracket <i>Used to achieve 14"W between the arm measurement when ordering a 14"W Backrest.</i>	No Charge	

Legrest Options*Choose One of the Following*

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103867-99-0	Power Articulating Center Mount Elevating Legrest	2,749.00	K0108
<input type="checkbox"/> I102212-99-0	Manual Adjustment Legrest Elevation <i>Includes Calf Supports.</i>	No Charge	

Select Options:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103055-99-0	UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad. Recommended pad to be used with lower leg lengths of 12" or less.</i>	No Charge	
<input type="checkbox"/> I102982-99-0	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad.</i>	No Charge	
<input type="checkbox"/> I102978-99-0	UT Calf Support Kit (H) w Mesh Cover & Adj Hardware <i>7"W x 7"T Curved Pad.</i>	312.00	E1028
<input type="checkbox"/> I100561-99-0	Legrest Assm. 1pc Footplate (Complete)	395.00	
<input type="checkbox"/> I102051-99-0	Legrest Assm. 2pc Narrow Footplates (Complete)	No Charge	
<input type="checkbox"/> I10006	3" Width Between Two- Piece Footplates <i>Not available if chair has anti-tips and lower leg length is greater than 14" (measured from seat pan to top of footplate.)</i>	495.00	
<input type="checkbox"/> I10271	Corpus Junior Lower Leg Length <i>Modification to allow for lower leg lengths of 8" - 10". Standard length is 11"-15". Measurement is taken from seat pan to the top of the footplates. Please specify lower leg length measurement in inches: _____</i>	No Charge	
<input type="checkbox"/> IFT202SR	BP Ankle Huggers, Med Side Release, Adj (Pair) <i>This item will be sent as a part with the wheelchair. It will not be installed.</i>	225.00	K0108

Accessories*Select Options:*

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10400	Bodypoint Non-Padded Hip Belt, UniTrack - Medium	No Charge	
<input type="checkbox"/> I10405	Bodypoint Padded Hip Belt UniTrack (N/C) - Medium	No Charge	
<input type="checkbox"/> I10623	Bodypoint Padded Hip Belt UniTrack - Extra Small	113.00	
<input type="checkbox"/> I10402	Bodypoint Non-Padded Hip Belt, UniTrack - Large	92.00	K0108
<input type="checkbox"/> I10408	Bodypoint Padded Hip Belt UniTrack - Large	155.00	K0108
<input type="checkbox"/> I100519-99-0	Positioning Belt Retractable	235.00	K0108
<input type="checkbox"/> I100749-99-0	Miniflex Headrest Removable - Leatherette <i>Must also select PN I10552.</i>	263.00	E0955
<input type="checkbox"/> I104131-99-0	Headrest for Ergo Back w/ Curved Mount Brkt - Leath <i>Must also select PN I10552. This part consists of a leatherette headrest pad (10"W x 5"T) with a 10"L slotted, curved mounting bracket.</i>	263.00	E0955

C400 Corpus Junior

MSRP Price List

Select Options:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103460-99-0	UT Ergo Headrest w/ Black Fabric Cover <i>Must also select PN I10552.</i>	292.00	E0955
<input type="checkbox"/> I1823249	UT Lateral Adjustment Bar for UniTrack Headrest <i>Allows for approximately 4" of lateral adjustment to the UT Ergo Headrest either right or left. NOTE: Due to the width of this item (approx 9"), it may be necessary to send other items such as back pack clips unmounted. Custom charge may apply with lower back heights.</i>	55.00	
<input type="checkbox"/> IHEADADAPTER	Universal Headrest Adapter <i>For aftermarket headrests.</i>	197.00	K0108
<input type="checkbox"/> I10552	Adj Removable - Headrest Hardware (1 Each)	212.00	E1028
<input type="checkbox"/> I103070-99-0	UT Amp. Support Pad (H) w/ Mesh Cover - (1 Each) <i>7"W x 7"T Curved Pad. Must also select PN I10553.</i>	125.00	E1020
<input type="checkbox"/> I103071-99-0	UT Amp. Support Pad (D) w/ Mesh Cover - (1 Each) <i>8"W x 8"T Flat Pad. Must also select PN I10553.</i>	125.00	E1020
<input type="checkbox"/> I10553	Adj Removable - Amp Support Hardware (1 Each) <i>Please indicate if you would like the amputation support mounted on the right or left side of the seat frame: Mount Amp Suppt on Right: _____ Mount Amp Suppt on Left: _____</i>	212.00	E1028
<input type="checkbox"/> I103077-99-0	UT Hip Suppt (B) w/ Mesh Cover - Low (Pair) <i>4.5"W x 3.5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	225.00	E0956
<input type="checkbox"/> I103079-99-0	UT Hip Suppt (A) w/ Mesh Cover - Low (Pair) <i>6"W x 5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	225.00	E0956
<input type="checkbox"/> I103094-99-0	UT Thigh/Hip Suppt (J) w/ Mesh Cover - High (Pair) <i>6.5"W x 5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	225.00	E0956
<input type="checkbox"/> I103095-99-0	UT Thigh/Hip Suppt (C) w/ Mesh Cover - High (Pair) <i>8"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	225.00	E0956
<input type="checkbox"/> I102977-99-0	UT Thigh/Hip Suppt (G) w/ Mesh Cover - High (Pair) <i>12"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	225.00	E0956
<input type="checkbox"/> I10550	Adj Removable - Thigh Support Hardware (Pair)	424.00	E1028
<input type="checkbox"/> I10272	Fixed Trunk Supports Jr Corpus - Leatherette <i>Must also select PN I10551.</i>	246.00	E0956
<input type="checkbox"/> I10274	Swing Away Trunk Supports Jr Corpus - Leatherette <i>Must also select PN I10551.</i>	265.00	E0956
<input type="checkbox"/> I10551	Adj Removable - Trunk Support Hardware (Pair)	424.00	E1028
<input type="checkbox"/> I100978-99-0	Abductor Kit - Leatherette <i>Must also select PN I10555.</i>	306.00	E0957
<input type="checkbox"/> I10555	Adj Removable - Abductor Hardware (1 Each)	212.00	E1028
<input type="checkbox"/> I10015	Knee Block Assembly Corpus	630.00	K0108
<input type="checkbox"/> I103651-99-0	Transfer Handles for UniTrack Accessory Rail	390.00	
<input type="checkbox"/> I100975-99-0	Push Handles for Ergo Backrest <i>Order with Ergonomic or Rectangular Backs. Custom charge may apply if ordering push handles with an articulating vent tray.</i>	339.00	
<input type="checkbox"/> I307103-00-0	Push Handles for Tubular Backrest <i>Order with Corpus Tubular Backs (PN: I10111 & I10112).</i>	337.00	
<input type="checkbox"/> I102059-99-0	Medical Necessities Bag <i>Includes Medical Necessity Bag Clips (Carabiner Style).</i>	225.00	
<input type="checkbox"/> I101473-99-0	Medical Necessity Bag Clips - Carabiner Style	75.00	
<input type="checkbox"/> I101500-99-0	Essentials Carrier for Chassis	379.00	
<input type="checkbox"/> ISH212	BP Shoulder Harness Medium, Front Pull <i>This item will be sent as a part with the wheelchair. It will not be installed.</i>	196.00	E0960
<input type="checkbox"/> I100578-99-0	Upper Extremity Support Small, Angle Adjustable <i>Must also select PN I10556.</i>	238.00	E0950
<input type="checkbox"/> I100575-99-0	Upper Extremity Support Angle Adjustable. 14" W <i>Must also select PN I10556.</i>	267.00	E0950

C400 Corpus Junior**MSRP Price List***Select Options:*

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I100576-99-0	Upper Extremity Support Angle Adjustable. 16" W <i>Must also select PN I10556.</i>	279.00	E0950
<input type="checkbox"/> I100571-99-0	Upper Extremity Support Joystick Cutout 16" W <i>Must also select PN I10556.</i>	322.00	E0950
<input type="checkbox"/> I100573-99-0	Upper Extremity Support Mid Joystick Cutout 16"W <i>Must also select PN I10556.</i>	395.00	E0950
<input type="checkbox"/> I103426-99-0	Upper Extremity Support Angle Adjust - UniTrack <i>Side mounted tray that attaches to the UniTrack Accessory Rail. Dimensions are: 19.5"W x 13"D. Can only be ordered with chairs that offer the UniTrack System. Must also select PN I10556.</i>	296.00	E0950
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each)	212.00	E1028
<input type="checkbox"/> I101964-99-0	Oxygen Holder for PS or Corpus Seat <i>Cannot be ordered with a Corpus Tubular Backrest. (I10111 or I10112). Custom charge may be applied if ordering in combination with accessories such as medical necessities bags, attendant controls, push handles, back pack clips, etc. This item is mounted to the Ergonomic Backrest.</i>	224.00	
<input type="checkbox"/> I101956-99-0	Crutch Holder for PS or Corpus Seat <i>Cannot be ordered with a Corpus Tubular Backrest. (I10111 or I10112). Accessories such as medical necessities bags and push handles may need to be sent as parts if ordering a crutch holder. This item is mounted to the ergonomic backrest.</i>	224.00	

Alternative Drive Systems*Choose One of the Following*

Part Number	Description	Price	HCPCS Code
	R-net Alternative Drive Input Devices <i>R-net Alternative Drive Input Devices take the place of, or may be used in addition to, the standard R-net Joystick. They are only usable on wheelchairs that have R-net electronics. You MUST select one of each of the following items: (1) Input Device, (2) Expandable Controller, (3) Adjustable/Removable Alt Drive Hardware, AND (4) Omni or Input Module.</i>		
	(1) INPUT DEVICES FOR ALTERNATIVE DRIVE CONTROL		
<input type="checkbox"/> I10558	Gimble Joystick set up for Armrest Mount	1,388.00	E2373
<input type="checkbox"/> I10559	Gimble Joystick set up for Midline Mount	1,388.00	E2373
<input type="checkbox"/> I10560	Gimble Joystick set up for Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	1,388.00	E2373
<input type="checkbox"/> I10561	Gimble Joystick w/ Haims Harness - Small <i>For neck width up to 5".</i>	1,555.00	E2373
<input type="checkbox"/> I10562	Gimble Joystick w/ Haims Harness - Large <i>For neck width up to 7".</i>	1,555.00	E2373
<input type="checkbox"/> I10563	Gimble Joystick w/ Power Swing Away Chin Boom <i>Custom charge may be applied if ordering smaller backrest sizes and heights.</i>	5,174.00	E2399
<input type="checkbox"/> I10564	Mini Joystick set up for Armrest Mount	2,035.00	E2312
<input type="checkbox"/> I10565	Mini Joystick set up for Midline Mount	2,035.00	E2312
<input type="checkbox"/> I10566	Mini Joystick set up for Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	2,035.00	E2312
<input type="checkbox"/> I10567	Mini Joystick w/ Haims Harness - Small <i>For neck width up to 5".</i>	2,135.00	E2312
<input type="checkbox"/> I10568	Mini Joystick w/ Haims Harness - Large <i>For neck width up to 7".</i>	2,135.00	E2312
<input type="checkbox"/> I10569	Mini Joystick w/ Power Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	5,174.00	E2399

C400 Corpus Junior**MSRP Price List**

Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I1821750	Tablet Control Pad <i>The Tablet Control is a unique drive input which allows a user to control a power wheelchair in a very similar manner to a mouse being controlled by a touchpad on a laptop computer. The dimensions of the Tablet Control Pad are: 12"W x 9.5"D x 1.5"H. Kit also includes: Buddy Button Switch, and Tablet Control Mounting Bracket.</i>	3,614.00	E2328
<input type="checkbox"/> I10570	Head Array - Adult 5 Switch HDR-A5	3,614.00	E2330
<input type="checkbox"/> I10571	Head Array - Adult 3 Switch HDR-A4 <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	3,614.00	E2330
<input type="checkbox"/> I10572	Head Array - Pediatric 5 Switch HDR-C5-RR	3,614.00	E2330
<input type="checkbox"/> I10573	Sip & Puff System <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	1,409.00	E2325
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	1,929.00	E2325,E2326

Select Options:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10076	Right Side Mount For Joystick	No Charge	
<input type="checkbox"/> I10077	Left Side Mount For Joystick	No Charge	
<input type="checkbox"/> I10074	Right Side Mount For Omni	No Charge	
<input type="checkbox"/> I10075	Left Side Mount For Omni	No Charge	
<input type="checkbox"/> I10078	Mount Seat Control Panel on Right	No Charge	
<input type="checkbox"/> I10079	Mount Seat Control Panel on Left	No Charge	
<input type="checkbox"/> I10080	Mount Seat Control Panel on Backrest	No Charge	

Alt Drive Systems Setup

Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10576	(2) EXPANDABLE CONTROLLER (REQUIRED) Expandable Controller - Alternative Drive <i>Must also select Harness for Expandable Controller (I10513), located in "Joystick Options".</i>	525.00	E2377

Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10577	(3) ADJUSTABLE REMOVABLE HARDWARE (REQUIRED) Adjustable Removable Hardware - Alt Drive	212.00	E1028

Choose One of the Following

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> ID51154IR	(4) OMNI PLUS CONTROL MODULE (REQUIRED) Omni Control Module for R-net w/ Infrared <i>Please specify which type of alternative drive input you will be using if not ordering the input device from Permobil. Input Device: _____</i>	1,470.00	
<input type="checkbox"/> I10581	Input Module (IOM) for R-net Control System <i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head Array or RIM control.</i>	995.00	E2351

Alt Drive Accessories

Select Options:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> IMICROSWITCH	Microlite Switch	148.00	

C400 Corpus Junior**MSRP Price List***Select Options:*

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> IBUDDYBUTTON	Buddy Button	111.00	K0108
<input type="checkbox"/> I5808	Mini Cup Switch	108.00	
<input type="checkbox"/> IZTS-P	Proximity Switch	515.00	
<input type="checkbox"/> IMINIJOY5SW	Mini Joystick 5 Switch OC	593.00	
<input type="checkbox"/> I5018	Wafer Board w/5 Switches	497.00	

Notes & Additional Instructions:

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, Inc. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. MID PROFILE® Cushion is a registered trademark of the The ROHO Group. *The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.*