


C400 VS Junior



** Indicates a Required Field*

Dealer Information

* Contact: _____

Dealer Code: _____

* Dealer Name: _____

Address: _____

* City: _____

* State/Zip: _____

* Phone#: _____

* Fax#: _____

PO#: _____

Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

* First Name: _____

* Last Name: _____

Diagnosis: _____

Funding Source: _____

Client Age: _____

Client Measurements

* **Weight:** _____

* **Height:** _____

A) Top of Shoulders: _____

B) Chest Depth: _____

C) Chest Width: _____

D) Seat Depth: _____

E) Top of Head: _____

F) Elbow to Hand: _____

G) Seat Pan to Elbow: _____

H) Hip Width: _____

*I) Knee to Foot: _____

Three column pricing on order/quote: _____

Display HCPCS Codes on order/quote: _____

Please send order/quote to fax#: **(800) 231-3256**

Prices effective July 1, 2009.

Permobil Inc.
6961 Eastgate Blvd.
Lebanon, TN 37090
Tel: (800) 736-0925
Fax: (800) 231-3256
www.permobil.com

Chair Model and Colors

Select Options:

Part Number	Description	Price	HCPCS Code
I104057-99-0	C400 Base VS Jr R-net - C400 PS2 <i>Price Includes: FWD C400 Power Base Set Up For Multiple Power Seat Functions, Shock Absorbing Suspension System, Tie-Down Hardware for Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below.</i>	9,995.00 <input checked="" type="checkbox"/>	K0861

Required Items

Part Number	Description	Price	HCPCS Code
I103889-99-0	Tri Spoke Split Rim w/ FF Drive Tires - 2.5 x 8"	No Charge <input type="checkbox"/>	
I103890-99-0	Tri Spoke Split Rim w/ Pneumatic Tire - 2.5 x 8"	No Charge <input type="checkbox"/>	

Required Items

Part Number	Description	Price	HCPCS Code
I102190-99-0	C400 Shroud Vintage Wine	No Charge <input type="checkbox"/>	
I103957-99-0	C400 Shroud Onyx Black	No Charge <input type="checkbox"/>	
I102194-99-0	C400 Shroud Crystal Blue	No Charge <input type="checkbox"/>	
I102191-99-0	C400 Shroud Metallic Charcoal Gray	No Charge <input type="checkbox"/>	
I102192-99-0	C400 Shroud Green Metallic	No Charge <input type="checkbox"/>	
I103958-99-0	C400 Shroud Extreme Yellow	No Charge <input type="checkbox"/>	

Required Items

Part Number	Description	Price	HCPCS Code
I10451	Match Hub Cap to Shroud Color <i>If this is selected, the Tri-Spoke Hub Cap will be matched to the selected shroud color.</i>	No Charge <input type="checkbox"/>	
I103537-99-0	Hub Cap Platinum Silver <i>If this is selected, the Tri-Spoke Hub Cap will be Platinum Silver.</i>	No Charge <input type="checkbox"/>	

Base Options

Required Items

Part Number	Description	Price	HCPCS Code
I103547-99-0	Power Adjustable Seat Height 10" Travel - C400 <i>This item must be selected. The lowest seat to floor height may be 18.5", 19.5" or 20.5"; however, this is a setting that must be adjusted prior to chair shipping from Permobil. Please indicate below what seat to floor height you would like.</i>	2,495.00 <input checked="" type="checkbox"/>	E2300

Required Items

Part Number	Description	Price	HCPCS Code
I10504	Set Seat Height at 18.5"	No Charge <input type="checkbox"/>	
I10505	Set Seat Height at 19.5"	No Charge <input type="checkbox"/>	
I10506	Set Seat Height at 20.5"	No Charge <input type="checkbox"/>	

Select Options:

Part Number	Description	Price	HCPCS Code
IM24 SLD G FT	Batteries, Group 24 Sealed LA/Gel, Installed	762.00 <input type="checkbox"/>	E2363
I103635-99-0	Light Kit, R-net, C400/C500	955.00 <input type="checkbox"/>	
I309297-28-0	Anti-Tippers for C400 VS Jr & C500 VS <i>Anti-Tippers are required for any user weight above 220lbs. They are included in the High Speed Package for the C500 VS Base (I104060-99-0).</i>	399.00 <input type="checkbox"/>	
I104082-99-0	R-net ESP w C400/C500 Mounting Bracket <i>MUST be selected if ordering any ventilator tray (I102769-99-0, I103241-99-0, or I102768-99-0). Recommended item with switched drive inputs for improved tracking.</i>	1,250.00 <input type="checkbox"/>	K0108
ID51113	Programming Key for R-net - Permobil Version	435.00 <input type="checkbox"/>	

Permolock

Select Options:

Part Number	Description	Price	HCPCS Code
I102843-99-0	Permolock C Kit - COMPLETE <i>Includes: Locking Base, Plunger, and Manual Release Handle. The Locking Base must be installed by a NMEDA Dealer. Please provide the NMEDA Dealer Name & Address where the locking Base will ship to.</i> NAME: _____ ADDR: _____	2,485.00 <input type="checkbox"/>	
I102842-99-0	Plunger ONLY for Permolock C <i>Price Includes: Plunger only. This item may be purchased by any Permobil Inc or NMEDA certified dealer. The standard mounting position for the button to control the plunger is at the end of the armrest pad, opposite the joystick.</i>	1,008.00 <input type="checkbox"/>	
I103242-99-0	3.6m Plunger Sw Ext Cable for VS with Permolock C <i>This Cable MUST be selected when a Permolock C is ordered with a VS seating system.</i>	No Charge <input type="checkbox"/>	

Joystick Options

Required Items

Part Number	Description	Price	HCPCS Code
I10045	Cone Shaped Joystick Knob - Standard.	No Charge <input type="checkbox"/>	
I100960-99-0	Chin Cup For Joystick Knob	94.00 <input type="checkbox"/>	E2324
I100961-99-0	"T" Handle For Joystick Knob	139.00 <input type="checkbox"/>	E2323
I100962-99-0	"Large Ball" For Joystick Knob	94.00 <input type="checkbox"/>	E2323
I100963-99-0	"Softball" For Joystick Knob	127.00 <input type="checkbox"/>	E2323
I100964-99-0	"Mushroom" For Joystick Knob	127.00 <input type="checkbox"/>	E2323
I100965-99-0	"Stick" For Joystick Knob	127.00 <input type="checkbox"/>	E2323
IPC101B	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	127.00 <input type="checkbox"/>	E2323
IPC102B	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	127.00 <input type="checkbox"/>	E2323
IPC107B	Bodypoint J/S Handle Dome Shaped, Rubber	127.00 <input type="checkbox"/>	E2323

Select Options:

Part Number	Description	Price	HCPCS Code
I10412	Joystick Mount Right - Fixed, R-net Only	No Charge <input type="checkbox"/>	
I10413	Joystick Mount Left - Fixed, R-net Only	No Charge <input type="checkbox"/>	
I102537-99-0	Joystick Mount Right - Swing Away, R-net & VR2	292.00 <input type="checkbox"/>	E1028
I102538-99-0	Joystick Mount Left - Swing Away, R-net & VR2	292.00 <input type="checkbox"/>	E1028
I103820-99-0	Retractable Joystick Mnt Right - R-net, VR2 Only	486.00 <input type="checkbox"/>	E1028
I103819-99-0	Retractable Joystick Mnt Left - R-net, VR2 Only	486.00 <input type="checkbox"/>	E1028
I101260-99-0	Midline Control Mount - Pilot+ & R-net Only <i>Please indicate which armrest to mount midline below.</i>	263.00 <input type="checkbox"/>	E1028
I10428	Mount Midline Joystick on Left Armrest	No Charge <input type="checkbox"/>	
I10427	Mount Midline Joystick on Right Armrest	No Charge <input type="checkbox"/>	
I103544-99-0	R-net Remote Joystick <i>Must order expandable controller (PN: I10411) with this part. Comes standard w/ mono jacks.</i>	530.00 <input type="checkbox"/>	
ID51040	R-net Remote Joystick w/ Toggle Switches <i>Expandable/remote joystick with toggle switches. Used with multiple power option bases or as an upgrade for no power option or single power option bases. Must order expandable controller (PN: I10411) with this part.</i>	605.00 <input type="checkbox"/>	
ID51042	R-net Remote Joystick w/ Speed Pot & Single Toggle <i>Expandable/remote joystick with a single toggle switch for modes/profiles & a speed pot. Used with multiple power option bases or as an upgrade for no power option or single power option bases. Must order expandable controller (PN: I10411) with this part.</i>	730.00 <input type="checkbox"/>	
ID51153	R-net Remote Color Screen Joystick w/ Mono Jacks <i>Must order expandable controller (PN: I10411) with this part. Comes standard w/ mono jacks.</i>	835.00 <input type="checkbox"/>	

C400 VS Junior**MSRP Price List***Select Options:*

Part Number	Description	Price	HCPCS Code
ID51190	R-net Remote Color JS w/ Speed Pot & Single Toggle <i>Must order expandable controller (PN: I10411) with this part.</i>	875.00 <input type="checkbox"/>	
ID51149	R-net Remote Color Screen JS w/ Toggle Switches <i>Must order expandable controller (PN: I10411) with this part.</i>	875.00 <input type="checkbox"/>	
I10411	Expandable Controller - R-net <i>Must be ordered with the expandable/remote joysticks (PN: I10410, I103544-99-0, ID51040, ID51108, ID51042, and I103814-99-0).</i>	525.00 <input checked="" type="checkbox"/>	E2377
I10513	Harness for Expandable Controller <i>Must be ordered with any Expandable Controller (I10411 or I10450)</i>	445.00 <input checked="" type="checkbox"/>	E2313
ID50882	Attendant Control for R-net Control System Only <i>Not available with VR2 Controllers.</i>	1,084.00 <input type="checkbox"/>	E2331
I104002-99-0	Co-Pilot Attendant Control, VS R-net <i>Custom charge may apply if ordering a backrest height of 21" or Lower.</i>	4,175.00 <input type="checkbox"/>	E2331
ID51111	R-net BT Mouse Module <i>The R-net BT Mouse Module allows for mouse control through a standard R-net Joystick or other type of input device (ex: Mini Joystick w/ Omni). It is only usable as a mouse on a computer. It CANNOT be paired to work with other Bluetooth devices such as cell phones, wireless headsets, etc. This kit includes: an R-net BT Mouse Module, a 1.5m R-net Bus Cable, and a USB Style Bluetooth Adapter for the user's computer. Note: The Bluetooth Adapter may not work with all computers.</i>	795.00 <input type="checkbox"/>	
I10582	Output Module (IOM) for R-net Control System <i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). One output module is required for each separate external device to be controlled. Up to three output modules can be installed on one wheelchair.</i>	995.00 <input type="checkbox"/>	E2351

Seating System*Required Items*

Part Number	Description	Price	HCPCS Code
	VS Jr Seat Note: <i>You must select both part numbers listed below. The combined price includes: Seat Frame with Adjustable Depth, Fixed Chest Bar, Solid Seat Pan, Solid Backrest, and Height Adjustable Armrests. The weight capacity of the VS Jr Seat is 165 lbs.</i>		
I10508	Power Tilt & Power Recline - VS Junior Seat <i>This item must be selected. Includes 25° Power Tilt & 175° Power Recline.</i>	9,855.00 <input checked="" type="checkbox"/>	E1007
I104133-99-0	Power Standing Function - VS Junior Seat <i>This item must be selected.</i>	4,295.00 <input checked="" type="checkbox"/>	E2301

Seating System Options*Required Items*

Part Number	Description	Price	HCPCS Code
I10585	Multiple Seat Function Control Kit for R-net <i>MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive control and/or ICS switch box. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the Seat Function Control Kit.</i>	2,875.00 <input checked="" type="checkbox"/>	E2311

Select Options:

Part Number	Description	Price	HCPCS Code
I102798-99-0	Push Buttons on Seat Function Control Kit	No Charge <input type="checkbox"/>	
I102799-99-0	Toggle Switches on Seat Function Control Kit	245.00 <input type="checkbox"/>	
I103469-99-0	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge <input type="checkbox"/>	

Select Options:

Part Number	Description	Price	HCPCS Code
I103471-99-0	Mount ICS Switchbox to End of Arm Pad - Corpus <i>ICS Switchbox will be mounted to the end of the armrest pad. If this item is selected, the ICS Switchbox will NOT swing away with the Joystick Module. Please indicate below if you would like it mounted onto the left or right arm pad.</i> Mount ICS Switchbox on Left Armpad _____ Mount ICS Switchbox on Right Armpad _____	No Charge	<input type="checkbox"/>
I103470-99-0	Mount ICS Switchbox on Backrest	125.00	<input type="checkbox"/>
I104746-99-0	Remote Stop System - R-net, Pilot+ <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Each transmitter and receiver is a serialized matched pair, meaning they must be used together. Custom charge may be required if ordering with a low backrest height or a pediatric wheelchair. If ordering Remote Stop System with a K450, push handles must also be selected.</i>	985.00	<input type="checkbox"/>

Backrest Options

Required Items

Part Number	Description	Price	HCPCS Code
	BACKREST HEIGHT NOTE: <i>Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°.</i>		
I100794-99-0	Ergo Back 14"W x 26.5"T - Leatherette, For Stander	837.00	<input type="checkbox"/> E2620
I100795-99-0	Ergo Back 16"W x 26.5"T - Leatherette, For Stander	837.00	<input type="checkbox"/> E2620
I100796-99-0	Ergo Back 18"W x 26.5"T - Leatherette, For Stander	837.00	<input type="checkbox"/> E2620
I10464	Ergo Back 14"W x 26.5"T - Black Mesh, For Stander	847.00	<input type="checkbox"/> E2620
I10465	Ergo Back 16"W x 26.5"T - Black Mesh, For Stander	847.00	<input type="checkbox"/> E2620
I10466	Ergo Back 18"W x 26.5"T - Black Mesh, For Stander	847.00	<input type="checkbox"/> E2620
I100434-99-0	Omit Backrest Kit 14" W Between Armrests	No Charge	<input type="checkbox"/>
I100435-99-0	Omit Backrest Kit 16" W Between Armrests	No Charge	<input type="checkbox"/>
I100436-99-0	Omit Backrest Kit 18" W Between Armrests	No Charge	<input type="checkbox"/>
I10242	Solid Back for Standing Seat 16"W x 26.5"H	No Charge	<input type="checkbox"/>
I10243	Solid Back for Standing Seat - 18"W x 26.5"H <i>Standard Option with Tilt and/or Recline Systems. Includes Leatherette back cushion ONLY.</i>	No Charge	<input type="checkbox"/>

Select Options:

Part Number	Description	Price	HCPCS Code
I10047	21" Back Height (Ergo Only)	No Charge	<input type="checkbox"/>
I10014	Custom Backrest Height _____" <i>Backrest height measured from seat pan and cannot be lower than 19". Price in addition to above backrest.</i>	625.00	<input type="checkbox"/>
I100972-99-0	Split Backrest <i>Split Backrest 8" from top to allow for transfers and/or transportation. This modification cannot be made on back heights lower than 26.5". Must select either PN: I10429 or I10430 below.</i>	535.00	<input type="checkbox"/>
I10429	Split Backrest Cushion	No Charge	<input type="checkbox"/>
I10430	Do NOT Split Backrest Cushion	No Charge	<input type="checkbox"/>
I10284	Flat Backrest Cushion for Ergo Back - Leatherette <i>This style of cushion is ONLY available in Leatherette. Recommended cushion when ordering lateral supports. Lateral wedges will be omitted if ordering this part.</i>	89.00	<input type="checkbox"/>
I101101-99-0	Backrest Reinforcement	260.00	<input type="checkbox"/>
I1824006	Sliding Backrest Kit for VS & VS Jr Seats <i>Provides 4" of mechanical shear reduction. Custom charge may be applied if ordering this item with any backrest shorter than 26.5".</i>	1,270.00	<input type="checkbox"/>

Seat Sizing

Required Items

Part Number	Description	Price	HCPCS Code
I103368-99-0	VS Jr Seat 17"W x 14"D	No Charge	<input type="checkbox"/>
I103369-99-0	VS Jr Seat 17"W x 15"D	No Charge	<input type="checkbox"/>
I103370-99-0	VS Jr Seat 17"W x 16"D	No Charge	<input type="checkbox"/>
I103371-99-0	VS Jr Seat 17"W x 17"D	No Charge	<input type="checkbox"/>
I103372-99-0	VS Jr Seat 17"W x 18"D	No Charge	<input type="checkbox"/>

Select Options:

Part Number	Description	Price	HCPCS Code
I10038	Ergonomic Seat Cushion Leatherette	475.00	<input type="checkbox"/> E2601,E2602
	<i>Not available with seat widths less than 16" or seat depths less than 14".</i>		
I10461	Ergonomic Seat Cushion - Black Mesh	485.00	<input type="checkbox"/> E2601
	<i>Not available with seat widths less than 16" or seat depths less than 14".</i>		

Armrest Options

Required Items

Part Number	Description	Price	HCPCS Code
I103339-99-0	4"x10" Corpus Arm Right - Dual Taper, Leatherette	No Charge	<input type="checkbox"/>
I103341-99-0	4"x12" Corpus Arm Right - Dual Taper, Leatherette	No Charge	<input type="checkbox"/>
I103343-99-0	4"x14" Corpus Arm Right - Dual Taper, Leatherette	No Charge	<input type="checkbox"/>
I103345-99-0	4"x16" Corpus Arm Right - Dual Taper, Leatherette	No Charge	<input type="checkbox"/>
I103331-99-0	4"x10" Corpus Arm Right - Dual Taper, Black Fabric	No Charge	<input type="checkbox"/>
I103333-99-0	4"x12" Corpus Arm Right - Dual Taper, Black Fabric	No Charge	<input type="checkbox"/>
I103335-99-0	4"x14" Corpus Arm Right - Dual Taper, Black Fabric	No Charge	<input type="checkbox"/>
I103337-99-0	4"x16" Corpus Arm Right - Dual Taper, Black Fabric	No Charge	<input type="checkbox"/>

Required Items

Part Number	Description	Price	HCPCS Code
I103338-99-0	4"x10" Corpus Arm Left - Dual Taper, Leatherette	No Charge	<input type="checkbox"/>
I103340-99-0	4"x12" Corpus Arm Left - Dual Taper, Leatherette	No Charge	<input type="checkbox"/>
I103342-99-0	4"x14" Corpus Arm Left - Dual Taper, Leatherette	No Charge	<input type="checkbox"/>
I103344-99-0	4"x16" Corpus Arm Left - Dual Taper, Leatherette	No Charge	<input type="checkbox"/>
I103330-99-0	4"x10" Corpus Arm Left - Dual Taper, Black Fabric	No Charge	<input type="checkbox"/>
I103332-99-0	4"x12" Corpus Arm Left - Dual Taper, Black Fabric	No Charge	<input type="checkbox"/>
I103334-99-0	4"x14" Corpus Arm Left - Dual Taper, Black Fabric	No Charge	<input type="checkbox"/>
I103336-99-0	4"x16" Corpus Arm Left - Dual Taper, Black Fabric	No Charge	<input type="checkbox"/>

Required Items

Part Number	Description	Price	HCPCS Code
I101003-99-0	Armrest Mount High	No Charge	<input type="checkbox"/>
	<i>Range of adjustment in high position is 11" - 13.5" measured from the seat pan to the top of the armrest pad.</i>		
I101002-99-0	Armrest Mount Low (High Position Standard)	No Charge	<input type="checkbox"/>
	<i>Range of adjustment in low position is 9" - 11.5" measured from the seat pan to the top of the armrest pad.</i>		

Select Options:

Part Number	Description	Price	HCPCS Code
I10016	Height Adjustable Armrest Assembly - No Charge	No Charge	<input type="checkbox"/>
	<i>No charge with power tilt and/or recline.</i>		
I304969-99-0	18" Between Armrests w/ 16" Backrest	265.00	<input type="checkbox"/>
I100523-99-0	Armrest Height Ext (+2")	125.00	<input type="checkbox"/>

C400 VS Junior**MSRP Price List***Select Options:*

Part Number	Description	Price	HCPCS Code
IARMADAP1"	Arm Trough Adaptor 1" Pair <i>For aftermarket arm troughs. Please indicate below if they should be installed or on chair sent as parts.</i>	274.00 <input type="checkbox"/>	
IARMADAP7/8"	Arm Trough Adaptor 7/8" Pair <i>For aftermarket arm troughs. Please indicate below if they should be installed on chair or sent as parts.</i>	274.00 <input type="checkbox"/>	
Install Arm Trough Adapters: _____			
Send Arm Trough Adapters as Parts: _____			
I100671-99-0	Armrest Pouch-Right	190.00 <input type="checkbox"/>	
I100672-99-0	Armrest Pouch-Left	190.00 <input type="checkbox"/>	
I101261-99-0	Armrest Side Support Right Outer <i>For Corpus Armrest Pads only.</i>	142.00 <input type="checkbox"/>	
I101262-99-0	Armrest Side Support Left Outer <i>For Corpus Armrest Pads only.</i>	142.00 <input type="checkbox"/>	
I101263-99-0	Armrest Side Support Right Inner <i>Recommend 2" wider width between the armrest to ensure armrests will flip-up.</i>	142.00 <input type="checkbox"/>	
I101264-99-0	Armrest Side Support Left Inner <i>Recommend 2" wider width between the armrest to ensure armrests will flip-up.</i>	142.00 <input type="checkbox"/>	

Legrest Options*Required Items*

Part Number	Description	Price	HCPCS Code
I10502	Power Articulating Center Mount Elevating Legrests <i>This item must be selected.</i>	2,749.00 <input checked="" type="checkbox"/>	K0108

Required Items

Part Number	Description	Price	HCPCS Code
I103351-99-0	Stand and Drive Legrest Assembly - Jr <i>This item must be selected.</i>	685.00 <input checked="" type="checkbox"/>	

Required Items

Part Number	Description	Price	HCPCS Code
I103353-99-0	Kneeblock Upholstery - Leatherette	No Charge <input type="checkbox"/>	
I103352-99-0	Kneeblock Upholstery - Black Mesh	147.00 <input type="checkbox"/>	

Select Options:

Part Number	Description	Price	HCPCS Code
I103055-99-0	UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad. Recommended pad to be used with lower leg lengths of 12" or less.</i>	No Charge <input type="checkbox"/>	
I103394-99-0	UT Calf Support Kit (J) w Mesh Cover & Adj Hdw - VS	No Charge <input type="checkbox"/>	

Accessories*Select Options:*

Part Number	Description	Price	HCPCS Code
I10400	Bodypoint Non-Padded Hip Belt, UniTrack - Medium	No Charge <input type="checkbox"/>	
I10405	Bodypoint Padded Hip Belt UniTrack (N/C) - Medium	No Charge <input type="checkbox"/>	
I10402	Bodypoint Non-Padded Hip Belt, UniTrack - Large	92.00 <input type="checkbox"/>	K0108
I10408	Bodypoint Padded Hip Belt UniTrack - Large	155.00 <input type="checkbox"/>	K0108
I104131-99-0	Headrest for Ergo Back w/ Curved Mount Brkt - Leath <i>Must also select PN I10552. This part consists of a leatherette headrest pad (10"W x 5"T) with a 10"L slotted, curved mounting bracket.</i>	263.00 <input type="checkbox"/>	E0955

C400 VS Junior

MSRP Price List

Select Options:

Part Number	Description	Price	HPCPS Code
I103460-99-0	UT Ergo Headrest w/ Black Fabric Cover <i>Must also select PN I10552.</i>	292.00 <input type="checkbox"/>	E0955
I1823249	UT Lateral Adjustment Bar for UniTrack Headrest <i>Allows for approximately 4" of lateral adjustment to the UT Ergo Headrest either right or left. NOTE: Due to the width of this item (approx 9"), it may be necessary to send other items such as back pack clips unmounted. Custom charge may apply with lower back heights.</i>	55.00 <input type="checkbox"/>	
IHEADADAPTER	Universal Headrest Adapter <i>For aftermarket headrests.</i>	197.00 <input type="checkbox"/>	K0108
I10552	Adj Removable - Headrest Hardware (1 Each)	212.00 <input type="checkbox"/>	E1028
I103077-99-0	UT Hip Suppt (B) w/ Mesh Cover - Low (Pair) <i>4.5"W x 3.5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	225.00 <input type="checkbox"/>	E0956
I103079-99-0	UT Hip Suppt (A) w/ Mesh Cover - Low (Pair) <i>6"W x 5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	225.00 <input type="checkbox"/>	E0956
I103094-99-0	UT Thigh/Hip Suppt (J) w/ Mesh Cover - High (Pair) <i>6.5"W x 5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	225.00 <input type="checkbox"/>	E0956
I103095-99-0	UT Thigh/Hip Suppt (C) w/ Mesh Cover - High (Pair) <i>8"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	225.00 <input type="checkbox"/>	E0956
I102977-99-0	UT Thigh/Hip Suppt (G) w/ Mesh Cover - High (Pair) <i>12"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	225.00 <input type="checkbox"/>	E0956
I10550	Adj Removable - Thigh Support Hardware (Pair)	424.00 <input type="checkbox"/>	E1028
I10272	Fixed Trunk Supports Jr Corpus - Leatherette <i>Must also select PN I10551.</i>	246.00 <input type="checkbox"/>	E0956
I305266-28-90-0	Lateral Supports Std. Leatherette-ErgoBack Only <i>Must also select PN I10551.</i>	246.00 <input type="checkbox"/>	E0956
I10274	Swing Away Trunk Supports Jr Corpus - Leatherette <i>Must also select PN I10551.</i>	265.00 <input type="checkbox"/>	E0956
I100521-99-0	Swing Away Trunk Supports Leatherette-For Ergo Back <i>Must also select PN I10551.</i>	250.00 <input type="checkbox"/>	E0956
I10551	Adj Removable - Trunk Support Hardware (Pair)	424.00 <input type="checkbox"/>	E1028
I101099-99-0	Chest Bar 14" Adjustable Leatherette	245.00 <input type="checkbox"/>	
I307649-99-60-0	Chest Bar 14" Adjustable Black Mesh	265.00 <input type="checkbox"/>	
I101095-99-0	Chest Bar 16" Adjustable Leatherette	245.00 <input type="checkbox"/>	
I306688-99-60-0	Chest Bar 16" Adjustable, Black Mesh	265.00 <input type="checkbox"/>	
I101097-99-0	Chest Bar 18" Adjustable Leatherette	245.00 <input type="checkbox"/>	
I307650-99-60-0	Chest Bar 18" Adjustable Black Mesh	265.00 <input type="checkbox"/>	
I100978-99-0	Abductor Kit - Leatherette <i>Must also select PN I10555.</i>	306.00 <input type="checkbox"/>	E0957
I10555	Adj Removable - Abductor Hardware (1 Each)	212.00 <input type="checkbox"/>	E1028
I103651-99-0	Transfer Handles for UniTrack Accessory Rail	390.00 <input type="checkbox"/>	
I101100-99-0	Push Handles for Standing Seat	339.00 <input type="checkbox"/>	
I101473-99-0	Back Pack Clips Carabiner Style	75.00 <input type="checkbox"/>	
I102059-99-0	Medical Necessities Bag <i>Includes Back Pack Clips (Carabiner Style).</i>	265.00 <input type="checkbox"/>	
I102697001331	Medical Necessities Bag - Case Logic	265.00 <input type="checkbox"/>	
I101500-99-0	Essentials Carrier for Chassis	379.00 <input type="checkbox"/>	
I100578-99-0	Upper Extremity Support Small, Angle Adjustable <i>Must also select PN I10556.</i>	238.00 <input type="checkbox"/>	E0950
I100575-99-0	Upper Extremity Support Angle Adjustable. 14" W <i>Must also select PN I10556.</i>	267.00 <input type="checkbox"/>	E0950
I100576-99-0	Upper Extremity Support Angle Adjustable. 16" W <i>Must also select PN I10556.</i>	279.00 <input type="checkbox"/>	E0950
I100577-99-0	Upper Extremity Support Angle Adjustable. 18" W <i>Must also select PN I10556.</i>	306.00 <input type="checkbox"/>	E0950

C400 VS Junior**MSRP Price List**

Select Options:

Part Number	Description	Price	HCPCS Code
I100571-99-0	Upper Extremity Support Joystick Cutout 16" W <i>Must also select PN I10556.</i>	322.00 <input type="checkbox"/>	E0950
I100572-99-0	Upper Extremity Support Joystick Cutout 18" W <i>Must also select PN I10556.</i>	322.00 <input type="checkbox"/>	E0950
I100574-99-0	Upper Extremity Support Mid Joystick Cutout 18"W <i>Must also select PN I10556.</i>	395.00 <input type="checkbox"/>	E0950
I10556	Adj Removable - Tray Hardware (1 Each)	212.00 <input type="checkbox"/>	E1028

Alternative Drive Systems

Required Items

Part Number	Description	Price	HCPCS Code
	R-net Alternative Drive Input Devices <i>R-net Alternative Drive Input Devices take the place of, or may be used in addition to, the standard R-net Joystick. They are only usable on wheelchairs that have R-net electronics. You MUST select one of each of the following items: (1) Input Device, (2) Expandable Controller, (3) Adjustable/Removable Alt Drive Hardware, AND (4) Omni or Input Module.</i>		
	(1) INPUT DEVICES FOR ALTERNATIVE DRIVE CONTROL		
I10558	Gimble Joystick set up for Armrest Mount	1,388.00 <input type="checkbox"/>	E2373
I10559	Gimble Joystick set up for Midline Mount	1,388.00 <input type="checkbox"/>	E2373
I10560	Gimble Joystick set up for Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	1,388.00 <input type="checkbox"/>	E2373
I10561	Gimble Joystick w/ Haims Harness - Small <i>For neck width up to 5".</i>	1,555.00 <input type="checkbox"/>	E2373
I10562	Gimble Joystick w/ Haims Harness - Large <i>For neck width up to 7".</i>	1,555.00 <input type="checkbox"/>	E2373
I10563	Gimble Joystick w/ Power Swing Away Chin Boom <i>Custom charge may be applied if ordering smaller backrest sizes and heights.</i>	5,174.00 <input type="checkbox"/>	E2399
I10564	Mini Joystick set up for Armrest Mount	2,035.00 <input type="checkbox"/>	E2312
I10565	Mini Joystick set up for Midline Mount	2,035.00 <input type="checkbox"/>	E2312
I10566	Mini Joystick set up for Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	2,035.00 <input type="checkbox"/>	E2312
I10567	Mini Joystick w/ Haims Harness - Small <i>For neck width up to 5".</i>	2,135.00 <input type="checkbox"/>	E2312
I10568	Mini Joystick w/ Haims Harness - Large <i>For neck width up to 7".</i>	2,135.00 <input type="checkbox"/>	E2312
I10569	Mini Joystick w/ Power Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	5,174.00 <input type="checkbox"/>	E2399
I1821750	Tablet Control Pad <i>The Tablet Control is a unique drive input which allows a user to control a power wheelchair in a very similar manner to a mouse being controlled by a touchpad on a laptop computer. The dimensions of the Tablet Control Pad are: 12"W x 9.5"D x 1.5"H. Kit also includes: Buddy Button Switch, and Tablet Control Mounting Bracket.</i>	3,614.00 <input type="checkbox"/>	E2328
I10570	Head Array - Adult 5 Switch HDR-A5	3,614.00 <input type="checkbox"/>	E2330
I10571	Head Array - Adult 3 Switch HDR-A4 <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	3,614.00 <input type="checkbox"/>	E2330
I10572	Head Array - Pediatric 5 Switch HDR-C5-RR	3,614.00 <input type="checkbox"/>	E2330
I10573	Sip & Puff System <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	1,409.00 <input type="checkbox"/>	E2325
I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	1,929.00 <input type="checkbox"/>	E2325,E2326

Select Options:

Part Number	Description	Price	HCPCS Code
-------------	-------------	-------	------------

C400 VS Junior**MSRP Price List***Select Options:*

Part Number	Description	Price	HCPCS Code
I10076	Right Side Mount For Joystick	No Charge	<input type="checkbox"/>
I10077	Left Side Mount For Joystick	No Charge	<input type="checkbox"/>
I10074	Right Side Mount For Omni	No Charge	<input type="checkbox"/>
I10075	Left Side Mount For Omni	No Charge	<input type="checkbox"/>
I10078	Mount Seat Control Panel on Right	No Charge	<input type="checkbox"/>
I10079	Mount Seat Control Panel on Left	No Charge	<input type="checkbox"/>
I10080	Mount Seat Control Panel on Backrest	No Charge	<input type="checkbox"/>

Alt Drive Systems Setup*Required Items*

Part Number	Description	Price	HCPCS Code
I10576	(2) EXPANDABLE CONTROLLER (REQUIRED) Expandable Controller - Alternative Drive	525.00	<input type="checkbox"/> E2377

Required Items

Part Number	Description	Price	HCPCS Code
I10577	(3) ADJUSTABLE REMOVABLE HARDWARE (REQUIRED) Adjustable Removable Hardware - Alt Drive	212.00	<input type="checkbox"/> E1028

Required Items

Part Number	Description	Price	HCPCS Code
ID51154IR	(4) OMNI PLUS CONTROL MODULE (REQUIRED) Omni Control Module for R-net w/ Infrared	1,470.00	<input type="checkbox"/>
I10581	<i>Please specify which type of alternative drive input you will be using if not ordering the input device from Permobil. Input Device: _____</i> Input Module (IOM) for R-net Control System	995.00	<input type="checkbox"/> E2351
	<i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head Array or RIM control.</i>		

Alt Drive Accessories*Select Options:*

Part Number	Description	Price	HCPCS Code
IMICROSWITCH	Microlite Switch	148.00	<input type="checkbox"/>
IBUDDYBUTTON	Buddy Button	111.00	<input type="checkbox"/> K0108
I5808	Mini Cup Switch	108.00	<input type="checkbox"/>
IZTS-P	Proximity Switch	515.00	<input type="checkbox"/>
IMINIJOY5SW	Mini Joystick 5 Switch OC	593.00	<input type="checkbox"/>
I5018	Wafer Board w/5 Switches	497.00	<input type="checkbox"/>

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, Inc. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. MID PROFILE® Cushion is a registered trademark of the The ROHO Group.

The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.