

# Chairman HD3



**\* Indicates a Required Field**

## Dealer Information

\* **Contact:** \_\_\_\_\_  
Dealer Code: \_\_\_\_\_  
\* **Dealer Name:** \_\_\_\_\_  
Address: \_\_\_\_\_  
\* **City:** \_\_\_\_\_  
\* **State/Zip:** \_\_\_\_\_  
\* **Phone#:** \_\_\_\_\_  
\* **Fax#:** \_\_\_\_\_  
PO#: \_\_\_\_\_  
Email Address: \_\_\_\_\_


## Client Information

*Permobil recommends that the client is evaluated by a certified rehab specialist.*

\* **First Name:** \_\_\_\_\_  
\* **Last Name:** \_\_\_\_\_  
Diagnosis: \_\_\_\_\_  
Funding Source: \_\_\_\_\_  
Client Age: \_\_\_\_\_

## Client Measurements

\* **Weight:** \_\_\_\_\_  
\* **Height:** \_\_\_\_\_  
A) Top of Shoulders: \_\_\_\_\_  
B) Chest Depth: \_\_\_\_\_  
C) Chest Width: \_\_\_\_\_  
D) Seat Depth: \_\_\_\_\_  
E) Top of Head: \_\_\_\_\_  
F) Elbow to Hand: \_\_\_\_\_  
G) Seat Pan to Elbow: \_\_\_\_\_  
H) Hip Width: \_\_\_\_\_  
\* **I) Knee to Foot:** \_\_\_\_\_



Three column pricing on order/quote:

Display HCPCS Codes on order/quote:

Please send order/quote to fax#: **(800) 231-3256**

Email to: [sales@permobilus.com](mailto:sales@permobilus.com)

**Permobil Inc.**  
6961 Eastgate Blvd.  
Lebanon, TN 37090  
Tel: (800) 736-0925  
Fax: (800) 231-3256  
[www.permobil.com](http://www.permobil.com)

Prices effective **March 8, 2010.**

## Chair Model and Colors

Select Options:

| Part Number                           | Description   | Price     | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> I101233-99-0 | Chairman HD3 Base W/ Lights<br><i>Price includes: FWD Power Base Set Up for Multiple Power Options, Tie Down Hardware for Strap Systems. Lights, Seat Stabilizers, HD Seat Post, Anti-tippers, Dual Rear Casters, Long Positioning Belt, and Charger.</i> | 11,595.00 |            |

Choose One of the Following

| Part Number                           | Description                           | Price     | HCPCS Code |
|---------------------------------------|---------------------------------------|-----------|------------|
| <input type="checkbox"/> I310158-FF   | Flat Free Drive Tires for HD3 or Robo | No Charge |            |
| <input type="checkbox"/> I310158-99-0 | Pneumatic Drive Tires for HD3 or Robo | No Charge |            |

## Base Options

Choose One of the Following

| Part Number                           | Description  | Price     | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> I100500-99-0 | Power Adjustable Seat Height, HD3 8" Travel<br><i>To assist in transfers. Drive inhibited when raised.</i> | 2,895.00  | E2300      |
| <input type="checkbox"/> I100702-99-0 | Seat Tube Fixed HD3  | No Charge |            |

Select Options:

| Part Number                            | Description   | Price  | HCPCS Code |
|--|---|--------|------------|
| <input type="checkbox"/> IM24 SLD G FT | Batteries, Group 24 Sealed LA/Gel, Installed        | 762.00 | E2363      |
| <input type="checkbox"/> I602307-99-0  | P+ PP1A Programmer - VSI, VR2 or Pilot+ Controllers | 595.00 |            |

## Joystick Options

Choose One of the Following

| Part Number                           | Description                                    | Price     | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> I10045       | Cone Shaped Joystick Knob - Standard.          | No Charge |            |
| <input type="checkbox"/> I100960-99-0 | Chin Cup For Joystick Knob                     | 94.00     | E2324      |
| <input type="checkbox"/> I100961-99-0 | "T" Handle For Joystick Knob                   | 139.00    | E2323      |
| <input type="checkbox"/> I100962-99-0 | "Large Ball" For Joystick Knob                 | 94.00     | E2323      |
| <input type="checkbox"/> I100963-99-0 | "Softball" For Joystick Knob                   | 127.00    | E2323      |
| <input type="checkbox"/> I100964-99-0 | "Mushroom" For Joystick Knob                   | 127.00    | E2323      |
| <input type="checkbox"/> I100965-99-0 | "Stick" For Joystick Knob                      | 127.00    | E2323      |
| <input type="checkbox"/> IPC101B      | Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft | 127.00    | E2323      |
| <input type="checkbox"/> IPC102B      | Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft | 127.00    | E2323      |
| <input type="checkbox"/> IPC107B      | Bodypoint J/S Handle Dome Shaped, Rubber       | 127.00    | E2323      |

Select Options:

| Part Number                           | Description   | Price     | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> I10127       | Fixed Joystick Mount - Right  | No Charge |            |
| <input type="checkbox"/> I10136       | Fixed Joystick Mount - Left   | No Charge |            |
| <input type="checkbox"/> I101243-99-0 | Joystick Mount Right - Swing Away   | 292.00    | E1028      |
| <input type="checkbox"/> I101242-99-0 | Joystick Mount Left - Swing Away  | 292.00    | E1028      |
| <input type="checkbox"/> I100819-99-0 | Retractable Joystick Mnt Right  | 486.00    | E1028      |
| <input type="checkbox"/> I100824-99-0 | Retractable Joystick Mnt Left   | 486.00    | E1028      |
| <input type="checkbox"/> I101260-99-0 | Midline Control Mount - Pilot+ & R-net Only<br><i>Please indicate which armrest to mount midline below.</i> | 263.00    | E1028      |
| <input type="checkbox"/> I10428       | Mount Midline Joystick on Left Armrest  | No Charge |            |
| <input type="checkbox"/> I10427       | Mount Midline Joystick on Right Armrest   | No Charge |            |

## Chairman HD3

## MSRP Price List

Select Options:

| Part Number                           | Description   | Price     | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> ID49670/3    | Pilot+ series Joystick w/ 7 Buttons - No Charge<br><i>Comes standard when ordering zero or one power option.</i>  | No Charge |            |
| <input type="checkbox"/> I10207       | Pilot+ Series Joystick w/ 7 Buttons<br><i>Expandable/remote joystick used with multiple power option bases or as an upgrade for no power option or single power option bases. MUST order expandable controller (PN: I10450) with this part.</i>                           | 530.00    |            |
| <input type="checkbox"/> ID49880      | Pilot+ Joystick M/J<br><i>Expandable/remote joystick used with multiple power option bases or as an upgrade for no power option or single power option bases. MUST order expandable controller (PN: I10450) with this part. 1/8" mini-mono jacks for power and mode.</i>  | 530.00    |            |
| <input type="checkbox"/> ID50002      | Pilot+Series Joystick T/S<br><i>Expandable/remote joystick used with multiple power option bases or as an upgrade for no power option or single power option bases. MUST order expandable controller (PN: I10450) with this part. Toggle switches for power and mode.</i> | 530.00    |            |
| <input type="checkbox"/> I10450       | Expandable Controller - Pilot Plus<br><i>MUST be ordered with the expandable/remote joysticks (PN: I10207, ID49880, or ID50002).</i>  | 525.00    | E2377      |
| <input type="checkbox"/> ID49682      | P+DAC Attendant Control - Pilot+ Controller Only<br><i>Not available with VR2 controllers. Will be sent as part if ordering sling back upholstery.</i>  | 1,136.00  | E2331      |
| <input type="checkbox"/> I10513       | Harness for Expandable Controller<br><i>Must be ordered with any Expandable Controller (I10411 or I10450)</i>   | 445.00    | E2313      |
| <input type="checkbox"/> I103181-99-0 | Co-Pilot Attendant Control - Corpus<br><i>Custom charge may apply if ordering a Corpus Tubular Backrest OR a backrest height of 21" or Lower. Not available with VR2 Controllers.</i>   | 2,626.00  | E2331      |

## Seating System

Choose One of the Following

| Part Number                           | Description  | Price     | HCPCS Code |
|---------------------------------------|--|-----------|------------|
|                                       | The Seating System Listed Below Includes:<br><i>HD3 Corpus Seat Frame, Back Frame with Power Recline, Reclining Armrest Assembly, Solid Seat Pan (19"w x 19"d), Legrest Assembly with Adjustable Elevation. Weight Limit 440lbs.</i> |           |            |
| <input type="checkbox"/> I100501-99-0 | HD3 Seat w/ NO Tilt & 120° Power Recline   | 8,834.00  | K0108      |
| <input type="checkbox"/> I100473-99-0 | HD3 Seat w/ 30° Power Tilt & 120° Power Recline<br><i>Drive inhibited when tilted past 7 degrees.</i>  | 11,395.00 | E1007      |

## Seating System Options

Select Options:

| Part Number                           | Description   | Price     | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> I10128       | Single Seat Function Control Switch   | No Charge |            |
| <input type="checkbox"/> I102057-99-0 | Multiple Seat Function Control Kit<br><i>MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive controls and/or membrane switch box. Order with Pilot+ or remote joystick options. Not available with VSI Controllers.</i> | 2,875.00  | E2311      |
| <input type="checkbox"/> I10129       | Single Seat Function Control Kit<br><i>For operation of a single seat function through the wheelchair electronics via the drive controls and/or membrane switch box. Order with Pilot+ or remote joystick options. Not available with VSI Controllers.</i>  | 1,478.00  | E2310      |
| <input type="checkbox"/> I100524-99-0 | Toggle Switch Kit For Seat Function Control Kit<br><i>Add on part for either the multiple seat function control kit or single seat function control kit in order.</i>   | 128.00    | K0108      |

## Chairman HD3

## MSRP Price List

### Select Options:

| Part Number                           | Description  | Price     | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> ISW-JOY      | Seat Control Joystick<br><i>Multi directional toggle to operate up to 4 power seat functions. Provides an additional method to accessing the seat functions. You must still order the Multiple Seat Function Control Kit if selecting this part. Please specify where to mount membrane switchbox: Right Arm Pad: _____ Left Armrest Pad: _____</i>  | 1,451.00  | K0108      |
| <input type="checkbox"/> I10253       | Mount Joystick on Left   | No Charge |            |
| <input type="checkbox"/> I10254       | Mount Joystick on Right  | No Charge |            |
| <input type="checkbox"/> ISW-ALT      | Altswitch Alternating Seat Control<br><i>Seat function interface box which accommodates 4 single momentary switches. 1/8" mini-mono switches not included. Provides an additional method to accessing the seat functions. You must still order the Multiple Seat Function Control Kit if selecting this part.</i>  | 985.00    | K0108      |
| <input type="checkbox"/> I104746-99-0 | Remote Stop System - R-net, Pilot+<br><i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Each transmitter and receiver is a serialized matched pair, meaning they must be used together. Custom charge may be required if ordering with a low backrest height or a pediatric wheelchair. If ordering Remote Stop System with a K450, push handles must also be selected.</i> | 985.00    |            |

## Backrest Options

### Choose One of the Following

| Part Number | Description  | Price | HCPCS Code |
|-------------|--|-------|------------|
|             | BACKREST HEIGHT NOTE:<br><i>Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°.</i> |       |            |

### Select Options:

| Part Number                     | Description  | Price     | HCPCS Code |
|---------------------------------|--|-----------|------------|
| <input type="checkbox"/> I10023 | HD 3 Backrest Fixed, Rectangular - Leatherette<br><i>HD3 Backrest has a metal pan and is 28" High from the seat pan. The width is 1" narrower than the seat width.</i> | No Charge |            |
| <input type="checkbox"/> I10014 | Custom Backrest Height _____"<br><i>Backrest height measured from seat pan and cannot be lower than 19". Price in addition to above backrest.</i>                      | 625.00    |            |

## Seat Sizing

### Choose One of the Following

| Part Number                           | Description         | Price     | HCPCS Code |
|---------------------------------------|---------------------|-----------|------------|
| <input type="checkbox"/> I101019-99-0 | HD Seat 19"W x 15"D | No Charge |            |
| <input type="checkbox"/> I101020-99-0 | HD Seat 19"W x 17"D | No Charge |            |
| <input type="checkbox"/> I101021-99-0 | HD Seat 19"W x 19"D | No Charge |            |
| <input type="checkbox"/> I101022-99-0 | HD Seat 19"W x 21"D | No Charge |            |
| <input type="checkbox"/> I101023-99-0 | HD Seat 19"W x 23"D | No Charge |            |
| <input type="checkbox"/> I101024-99-0 | HD Seat 21"W x 15"D | No Charge |            |
| <input type="checkbox"/> I101025-99-0 | HD Seat 21"W x 17"D | No Charge |            |
| <input type="checkbox"/> I101026-99-0 | HD Seat 21"W x 19"D | No Charge |            |
| <input type="checkbox"/> I101027-99-0 | HD Seat 21"W x 21"D | No Charge |            |
| <input type="checkbox"/> I101028-99-0 | HD Seat 21"W x 23"D | No Charge |            |
| <input type="checkbox"/> I101029-99-0 | HD Seat 23"W x 15"D | No Charge |            |
| <input type="checkbox"/> I101030-99-0 | HD Seat 23"W x 17"D | No Charge |            |
| <input type="checkbox"/> I101031-99-0 | HD Seat 23"W x 19"D | No Charge |            |
| <input type="checkbox"/> I101032-99-0 | HD Seat 23"W x 21"D | No Charge |            |
| <input type="checkbox"/> I101033-99-0 | HD Seat 23"W x 23"D | No Charge |            |

**Chairman HD3****MSRP Price List***Select Options:*

| Part Number                     | Description  | Price  | HCPCS Code  |
|---------------------------------|--|--------|-------------|
| <input type="checkbox"/> I10038 | Ergonomic Seat Cushion Leatherette<br><i>Not available with seat widths less than 16" or seat depths less than 14".</i>  | 475.00 | E2601,E2602 |
| <input type="checkbox"/> I10461 | Ergonomic Seat Cushion - Black Mesh<br><i>Not available with seat widths less than 16" or seat depths less than 14".</i> | 485.00 | E2601       |

**Armrest Options***Choose One of the Following*

| Part Number                     | Description  | Price     | HCPCS Code |
|---------------------------------|--|-----------|------------|
| <input type="checkbox"/> I10137 | Height Adjustable Armrest Assembly   | 284.00    | E0973      |
| <input type="checkbox"/> I10016 | Height Adjustable Armrest Assembly - No Charge<br><i>No charge with power tilt and/or recline.</i> | No Charge |            |

*Select Options:*

| Part Number                           | Description  | Price     | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> I103344-99-0 | 4"x16" Corpus Arm Left - Dual Taper, Leatherette   | No Charge |            |
| <input type="checkbox"/> I103345-99-0 | 4"x16" Corpus Arm Right - Dual Taper, Leatherette  | No Charge |            |
| <input type="checkbox"/> I103336-99-0 | 4"x16" Corpus Arm Left - Dual Taper, Black Fabric  | No Charge |            |
| <input type="checkbox"/> I103337-99-0 | 4"x16" Corpus Arm Right - Dual Taper, Black Fabric   | No Charge |            |
| <input type="checkbox"/> I100523-99-0 | Armrest Height Ext (+2")   | 125.00    |            |
| <input type="checkbox"/> IARMADAP1"   | Arm Trough Adaptor 1" Pair<br><i>For aftermarket arm troughs. Please indicate below if they should be installed on chair sent as parts.</i>  | 274.00    |            |
| <input type="checkbox"/> IARMADAP7/8" | Arm Trough Adaptor 7/8" Pair<br><i>For aftermarket arm troughs. Please indicate below if they should be installed on chair or sent as parts.</i><br>Install Arm Trough Adapters: _____<br>Send Arm Trough Adapters as Parts: _____ | 274.00    |            |
| <input type="checkbox"/> I100671-99-0 | Armrest Pouch-Right  | 190.00    |            |
| <input type="checkbox"/> I100672-99-0 | Armrest Pouch-Left   | 190.00    |            |
| <input type="checkbox"/> I101261-99-0 | Armrest Side Support Right Outer<br><i>For Corpus Armrest Pads only.</i>   | 142.00    |            |
| <input type="checkbox"/> I101262-99-0 | Armrest Side Support Left Outer<br><i>For Corpus Armrest Pads only.</i>  | 142.00    |            |
| <input type="checkbox"/> I101263-99-0 | Armrest Side Support Right Inner<br><i>Recommend 2" wider width between the armrest to ensure armrests will flip-up.</i>   | 142.00    |            |
| <input type="checkbox"/> I101264-99-0 | Armrest Side Support Left Inner<br><i>Recommend 2" wider width between the armrest to ensure armrests will flip-up.</i>  | 142.00    |            |

**Legrest Options***Choose One of the Following*

| Part Number                           | Description   | Price     | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> I100503-99-0 | Power Articulating Center Mount Elevating Legrests<br><i>Elevation limited to 65 degrees. Includes Calf Supports.</i> | 2,749.00  | K0108      |
| <input type="checkbox"/> I100502-99-0 | Manual Adjustment Legrest Elevation<br><i>Includes Calf Supports.</i>   | No Charge |            |

*Select Options:*

| Part Number                           | Description   | Price     | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> I100515-99-0 | Calf Supports Leatherette<br><i>Not available with lower leg lengths of 9" or less</i>                  | No Charge |            |
| <input type="checkbox"/> I10462       | UT Calf Support (H) w/ Mesh Cover & Adj Hardware<br><i>7"W x 7"T Curved Pad. No Charge on HD3 ONLY.</i> | No Charge |            |

**Chairman HD3**

**MSRP Price List**

Select Options:

| Part Number                           | Description              | Price     | HCPCS Code |
|---------------------------------------|--------------------------|-----------|------------|
| <input type="checkbox"/> I100475-99-0 | Footplate assembly--1 pc | 395.00    |            |
| <input type="checkbox"/> I100474-99-0 | Footplate assembly--2 pc | No Charge |            |

**Accessories**

Select Options:

| Part Number                              | Description  | Price     | HCPCS Code |
|--|--|-----------|------------|
| <input type="checkbox"/> I100600-99-0    | Positioning Belt Retractable Long  | 255.00    | K0108      |
| <input type="checkbox"/> I100507-99-0    | Positioning Belt Long<br><i>Included in Basic Equipment Package.</i>   | No Charge |            |
| <input type="checkbox"/> I104131-99-0    | Headrest for Ergo Back w/ Curved Mount Brkt - Leath<br><i>Must also select PN I10552. This part consists of a leatherette headrest pad (10"W x 5"T) with a 10"L slotted, curved mounting bracket.</i>  | 263.00    | E0955      |
| <input type="checkbox"/> I103460-99-0    | UT Ergo Headrest w/ Black Fabric Cover<br><i>Must also select PN I10552.</i>   | 292.00    | E0955      |
| <input type="checkbox"/> I1823249        | UT Lateral Adjustment Bar for UniTrack Headrest<br><i>Allows for approximately 4" of lateral adjustment to the UT Ergo Headrest either right or left. NOTE: Due to the width of this item (approx 9"), it may be necessary to send other items such as back pack clips unmounted. Custom charge may apply with lower back heights.</i> | 55.00     |            |
| <input type="checkbox"/> IHEADADAPTER    | Universal Headrest Adapter<br><i>For aftermarket headrests.</i>  | 197.00    | K0108      |
| <input type="checkbox"/> I10552          | Adj Removable - Headrest Hardware (1 Each)   | 212.00    | E1028      |
| <input type="checkbox"/> I103070-99-0    | UT Amp. Support Pad (H) w/ Mesh Cover - (1 Each)<br><i>7"W x 7"T Curved Pad. Must also select PN I10553.</i>   | 125.00    | E1020      |
| <input type="checkbox"/> I103071-99-0    | UT Amp. Support Pad (D) w/ Mesh Cover - (1 Each)<br><i>8"W x 8"T Flat Pad. Must also select PN I10553.</i>   | 125.00    | E1020      |
| <input type="checkbox"/> I10553          | Adj Removable - Amp Support Hardware (1 Each)<br><i>Please indicate if you would like the amputation support mounted on the right or left side of the seat frame: Mount Amp Suppt on Right: _____ Mount Amp Suppt on Left: _____</i>   | 212.00    | E1028      |
| <input type="checkbox"/> I103077-99-0    | UT Hip Suppt (B) w/ Mesh Cover - Low (Pair)<br><i>4.5"W x 3.5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>   | 225.00    | E0956      |
| <input type="checkbox"/> I103079-99-0    | UT Hip Suppt (A) w/ Mesh Cover - Low (Pair)<br><i>6"W x 5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>   | 225.00    | E0956      |
| <input type="checkbox"/> I103094-99-0    | UT Thigh/Hip Suppt (J) w/ Mesh Cover - High (Pair)<br><i>6.5"W x 5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>  | 225.00    | E0956      |
| <input type="checkbox"/> I103095-99-0    | UT Thigh/Hip Suppt (C) w/ Mesh Cover - High (Pair)<br><i>8"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>  | 225.00    | E0956      |
| <input type="checkbox"/> I102977-99-0    | UT Thigh/Hip Suppt (G) w/ Mesh Cover - High (Pair)<br><i>12"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>   | 225.00    | E0956      |
| <input type="checkbox"/> I10550          | Adj Removable - Thigh Support Hardware (Pair)  | 424.00    | E1028      |
| <input type="checkbox"/> I305266-28-90-0 | Lateral Supports Std. Leatherette-ErgoBack Only<br><i>Must also select PN I10551.</i>  | 246.00    | E0956      |
| <input type="checkbox"/> I100521-99-0    | Swing Away Trunk Supports Leatherette-For Ergo Back<br><i>Must also select PN I10551.</i>  | 250.00    | E0956      |
| <input type="checkbox"/> I10551          | Adj Removable - Trunk Support Hardware (Pair)  | 424.00    | E1028      |
| <input type="checkbox"/> I103651-99-0    | Transfer Handles for UniTrack Accessory Rail   | 390.00    |            |
| <input type="checkbox"/> I100975-99-0    | Push Handles for Ergo Backrest<br><i>Order with Ergonomic or Rectangular Backs. Custom charge may apply if ordering push handles with an articulating vent tray.</i>   | 339.00    |            |
| <input type="checkbox"/> I102059-99-0    | Medical Necessities Bag<br><i>Includes Medical Necessity Bag Clips (Carabiner Style).</i>  | 225.00    |            |
| <input type="checkbox"/> I101473-99-0    | Medical Necessity Bag Clips - Carabiner Style  | 75.00     |            |

**Chairman HD3**

**MSRP Price List**

Select Options:

| Part Number                           | Description  | Price  | HCPCS Code |
|---------------------------------------|--|--------|------------|
| <input type="checkbox"/> I101964-99-0 | Oxygen Holder for PS or Corpus Seat<br><i>Cannot be ordered with a Corpus Tubular Backrest. (I10111 or I10112). Custom charge may be applied if ordering in combination with accessories such as medical necessities bags, attendant controls, push handles, back pack clips, etc. This item is mounted to the Ergonomic Backrest.</i> | 224.00 |            |

**Alternative Drive Systems**

Choose One of the Following

| Part Number                       | Description   | Price    | HCPCS Code  |
|-----------------------------------|---|----------|-------------|
|                                   | <b>Pilot Plus Alternative Drive Input Devices</b><br><i>Pilot + Alt Drive Input Devices take the place of, or may be used in addition to, the standard Pilot + Joystick. They are only usable on wheelchairs that have Pilot + electronics. You MUST select one of each of the following items: (1) Input Device, (2) Expandable Controller, (3) Adjustable/Removable Alt Drive Hardware, AND (4) Omni Plus Control Module.</i> |          |             |
|                                   | <b>(1) INPUT DEVICES FOR ALTERNATIVE DRIVE CONTROL</b>  |          |             |
| <input type="checkbox"/> I10558   | Gimble Joystick set up for Armrest Mount  | 1,388.00 | E2373       |
| <input type="checkbox"/> I10559   | Gimble Joystick set up for Midline Mount  | 1,388.00 | E2373       |
| <input type="checkbox"/> I10560   | Gimble Joystick set up for Swing Away Chin Boom<br><i>Custom charge may be applied if ordering lower backrest heights.</i>  | 1,388.00 | E2373       |
| <input type="checkbox"/> I10561   | Gimble Joystick w/ Haims Harness - Small<br><i>For neck width up to 5".</i>   | 1,555.00 | E2373       |
| <input type="checkbox"/> I10562   | Gimble Joystick w/ Haims Harness - Large<br><i>For neck width up to 7".</i>   | 1,555.00 | E2373       |
| <input type="checkbox"/> I10563   | Gimble Joystick w/ Power Swing Away Chin Boom<br><i>Custom charge may be applied if ordering smaller backrest sizes and heights.</i>  | 5,174.00 | E2399       |
| <input type="checkbox"/> I10564   | Mini Joystick set up for Armrest Mount  | 2,035.00 | E2312       |
| <input type="checkbox"/> I10565   | Mini Joystick set up for Midline Mount  | 2,035.00 | E2312       |
| <input type="checkbox"/> I10566   | Mini Joystick set up for Swing Away Chin Boom<br><i>Custom charge may be applied if ordering lower backrest heights.</i>  | 2,035.00 | E2312       |
| <input type="checkbox"/> I10567   | Mini Joystick w/ Haims Harness - Small<br><i>For neck width up to 5".</i>   | 2,135.00 | E2312       |
| <input type="checkbox"/> I10568   | Mini Joystick w/ Haims Harness - Large<br><i>For neck width up to 7".</i>   | 2,135.00 | E2312       |
| <input type="checkbox"/> I10569   | Mini Joystick w/ Power Swing Away Chin Boom<br><i>Custom charge may be applied if ordering lower backrest heights.</i>  | 5,174.00 | E2399       |
| <input type="checkbox"/> I1821750 | Tablet Control Pad<br><i>The Tablet Control is a unique drive input which allows a user to control a power wheelchair in a very similar manner to a mouse being controlled by a touchpad on a laptop computer. The dimensions of the Tablet Control Pad are: 12"W x 9.5"D x 1.5"H. Kit also includes: Buddy Button Switch, and Tablet Control Mounting Bracket.</i>   | 3,614.00 | E2328       |
| <input type="checkbox"/> I10570   | Head Array - Adult 5 Switch HDR-A5  | 3,614.00 | E2330       |
| <input type="checkbox"/> I10571   | Head Array - Adult 3 Switch HDR-A4<br><i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>   | 3,614.00 | E2330       |
| <input type="checkbox"/> I10572   | Head Array - Pediatric 5 Switch HDR-C5-RR   | 3,614.00 | E2330       |
| <input type="checkbox"/> I10573   | Sip & Puff System<br><i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>  | 1,409.00 | E2325       |
| <input type="checkbox"/> I10574   | Sip & Puff System with Fiber Optic Opti-Stop<br><i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>   | 1,929.00 | E2325,E2326 |

Select Options:

| Part Number                     | Description                   | Price     | HCPCS Code |
|---------------------------------|-------------------------------|-----------|------------|
| <input type="checkbox"/> I10076 | Right Side Mount For Joystick | No Charge |            |

**Chairman HD3****MSRP Price List***Select Options:*

| Part Number                     | Description                          | Price     | HCPCS Code |
|---------------------------------|--------------------------------------|-----------|------------|
| <input type="checkbox"/> I10077 | Left Side Mount For Joystick         | No Charge |            |
| <input type="checkbox"/> I10074 | Right Side Mount For Omni            | No Charge |            |
| <input type="checkbox"/> I10075 | Left Side Mount For Omni             | No Charge |            |
| <input type="checkbox"/> I10078 | Mount Seat Control Panel on Right    | No Charge |            |
| <input type="checkbox"/> I10079 | Mount Seat Control Panel on Left     | No Charge |            |
| <input type="checkbox"/> I10080 | Mount Seat Control Panel on Backrest | No Charge |            |

**Alt Drive Systems Setup***Choose One of the Following*

| Part Number                     | Description   | Price  | HCPCS Code |
|---------------------------------|---|--------|------------|
| <input type="checkbox"/> I10576 | (2) EXPANDABLE CONTROLLER (REQUIRED)<br>Expandable Controller - Alternative Drive<br><i>Must also select Harness for Expandable Controller (I10513), located in "Joystick Options".</i> | 525.00 | E2377      |

*Choose One of the Following*

| Part Number                     | Description   | Price  | HCPCS Code |
|---------------------------------|---|--------|------------|
| <input type="checkbox"/> I10577 | (3) ADJUSTABLE REMOVABLE HARDWARE (REQUIRED)<br>Adjustable Removable Hardware - Alt Drive | 212.00 | E1028      |

*Choose One of the Following*

| Part Number                     | Description  | Price    | HCPCS Code |
|---------------------------------|--|----------|------------|
| <input type="checkbox"/> I10073 | (4) OMNI PLUS CONTROL MODULE (REQUIRED)<br>Omni Plus Control Module for Pilot Plus<br><i>Please specify input method if not ordering from Permobil. Input Device:</i><br>_____ | 1,470.00 | E2399      |

**Alt Drive Accessories***Select Options:*

| Part Number                           | Description  | Price    | HCPCS Code  |
|---------------------------------------|--|----------|-------------|
| <input type="checkbox"/> IMICROSWITCH | Microlite Switch   | 148.00   |             |
| <input type="checkbox"/> IBUDDYBUTTON | Buddy Button   | 111.00   | K0108       |
| <input type="checkbox"/> I5808        | Mini Cup Switch  | 108.00   |             |
| <input type="checkbox"/> IZTS-P       | Proximity Switch   | 515.00   |             |
| <input type="checkbox"/> IMINIJOY5SW  | Mini Joystick 5 Switch OC  | 593.00   |             |
| <input type="checkbox"/> I5018        | Wafer Board w/5 Switches   | 497.00   |             |
| <input type="checkbox"/> ID00005      | Dual Control Switch Box  | 859.00   |             |
| <input type="checkbox"/> ID49871      | Auxiliary Control Module   | 877.00   | E2351       |
| <input type="checkbox"/> I10082       | Chin Control Unit Stand Alone<br><i>Includes: P &amp; G Chin Control Unit and swing away mount off the backrest. (Do not need Omni Plus Control Module).</i> | 1,215.00 | E2373,E1028 |

Notes & Additional Instructions:

**Chair Order Policy:**

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, Inc. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

**Order Acknowledgement:**

I, \_\_\_\_\_, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: \_\_\_\_\_

Printed Name: \_\_\_\_\_

Title: \_\_\_\_\_

Date: \_\_\_\_\_

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. MID PROFILE® Cushion is a registered trademark of the The ROHO Group. *The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.*