


C300 Corpus Junior



Dealer Information

Contact: _____
 Dealer Code: _____
 Dealer Name: _____
 Address: _____
 City: _____
 Province/Postal Code: _____
 Phone#: _____
 Fax#: _____
 PO#: _____

Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

First Name: _____
 Last Name: _____
 Diagnosis: _____
 Age: _____

Client Measurements

Weight: _____
Height: _____
 A) Top of Shoulders: _____
 B) Chest Depth: _____
 C) Chest Width: _____
 D) Seat Depth: _____
 E) Top of Head: _____
 F) Elbow to Hand: _____
 G) Seat Pan to Elbow: _____
 H) Hip Width: _____
 I) Knee to Foot: _____

Display ADP Codes on order/quote: _____

Please send order/quote to fax#: **(800) 231-3256**

Canadian prices effective **October 28, 2008** and payable in Canadian dollars. All specifications and prices are subject to change without notice. **Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty.**

Permobil Ltd.
 Toronto Canada
 C/O Permobil Inc
 6961 Eastgate Blvd.
 Lebanon, TN 37090
 Tel: (800) 736-0925
 Fax: (800) 231-3256
 www.permobil.com

Base

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
I103260-99-0	C300 Base Corpus R-net - C300 PS2 <i>Price Includes: FWD C300 Power Base Set Up For Multiple Power Seat Functions, Shock Absorbing Suspension System, Tie-Down Hardware for Strap Systems, Anti-Tippers, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below.</i>	6,895.00 <input type="checkbox"/>	WPA31P010

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I103499-99-0	Tri-Spoke Split Rim w/ Flat Free Drive Tire 3x8"	No Charge <input type="checkbox"/>	
I103498-99-0	Tri-Spoke Split Rim w/ Pneumatic Tires - 3x8"	No Charge <input type="checkbox"/>	

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I103576-99-0	C300/K300 Shroud - Onyx Black	No Charge <input type="checkbox"/>	
I101381-99-0	C300/K300 Shroud - Vintage Wine	No Charge <input type="checkbox"/>	
I101383-99-0	C300/K300 Shroud - Crystal Blue	No Charge <input type="checkbox"/>	
I101384-99-0	C300/K300 Shroud - Metallic Charcoal Gray	No Charge <input type="checkbox"/>	
I103577-99-0	C300/K300 Shroud - Platinum Silver	No Charge <input type="checkbox"/>	
I101385-99-0	C300/K300 Shroud - Green Metallic	No Charge <input type="checkbox"/>	
I103575-99-0	C300/K300 Shroud - Titan Blue	No Charge <input type="checkbox"/>	
I103578-99-0	C300/K300 Shroud - Pop Star Pink	No Charge <input type="checkbox"/>	
I103574-99-0	C300/K300 Shroud - Chili Pepper Red	No Charge <input type="checkbox"/>	

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I10451	Match Hub Cap to Shroud Color <i>If this is selected, the Tri-Spoke Hub Cap will be matched to the selected shroud color.</i>	No Charge <input type="checkbox"/>	
I103537-99-0	Hub Cap Platinum Silver <i>If this is selected, the Tri-Spoke Hub Cap will be Platinum Silver.</i>	No Charge <input type="checkbox"/>	

Base Options

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I103779-99-0	Power Adjustable Seat Height - 8" Travel <i>To assist in transfers and improve access for function.</i>	2,495.00 <input type="checkbox"/>	
I101994-99-0	Seat Tube Fixed	No Charge <input type="checkbox"/>	

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
IM34 SLDG	Batteries - Group 34, 12V 60 Amp Hours, Installed	762.00 <input type="checkbox"/>	
I102768-99-0	Fixed Vent Shelf for C300 <i>If selecting this part, you MUST also pick: PN I0036-7002A - ESP Unit. Fits an LP10 Sized Vent & a U1 Battery. Must be ordered with an expandable controller (ex. R-net Electronics).</i>	395.00 <input type="checkbox"/>	WEPY
I102769-99-0	Articulating Vent Tray for C300 <i>If selecting this part, you MUST also pick: PN I0036-7002A - ESP Unit. Fits an LP10 Sized Vent and a Group 34 Battery. Cannot be ordered with any backrest height less than 26.5", a tubular back, or a 13" seat depth. Customization may be required if ordered with any of the following items Co-Pilot, any standard attendant control, push handles, or an oxygen holder. Only available with expandable controllers (R-net Electronics.)</i>	1,295.00 <input type="checkbox"/>	WEPY
I0036-7002A	Enhanced Steering Perform Unit (ESP), R-net <i>MUST be selected if ordering any ventilator tray (I102769-99-0, I103241-99-0, or I102768-99-0).</i>	1,250.00 <input type="checkbox"/>	WEPQ
I103564-99-0	Light Kit for C300 w/ R-net Electronics <i>Shrouds with lights cannot be ordered with VR2 Electronics.</i>	955.00 <input type="checkbox"/>	

Base Options

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
ID51113	Programming Key for R-net - Permobil Version	435.00 <input type="checkbox"/>	

Joystick Options

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I10045	Cone Shaped Joystick Knob - Standard.	No Charge <input type="checkbox"/>	
I100960-99-0	Chin Cup For Joystick Knob	94.00 <input type="checkbox"/>	
I100961-99-0	"T" Handle For Joystick Knob	139.00 <input type="checkbox"/>	
I100962-99-0	"Large Ball" For Joystick Knob	94.00 <input type="checkbox"/>	
I100963-99-0	"Softball" For Joystick Knob	127.00 <input type="checkbox"/>	
I100964-99-0	"Mushroom" For Joystick Knob	127.00 <input type="checkbox"/>	
I100965-99-0	"Stick" For Joystick Knob	127.00 <input type="checkbox"/>	

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
I102537-99-0	Joystick Mount Right - Swing Away, R-net & VR2	325.00 <input type="checkbox"/>	WEPM
I102538-99-0	Joystick Mount Left - Swing Away, R-net & VR2	325.00 <input type="checkbox"/>	WEPM
I103820-99-0	Retractable Joystick Mnt Right - R-net, VR2 Only	375.00 <input type="checkbox"/>	WEPM
I103819-99-0	Retractable Joystick Mnt Left - R-net, VR2 Only	375.00 <input type="checkbox"/>	WEPM
I101260-99-0	Midline Control Mount - Pilot+ & R-net Only <i>Please indicate which armrest to mount midline below.</i>	325.00 <input type="checkbox"/>	WE3D
I10428	Mount Midline Joystick on Left Armrest	No Charge <input type="checkbox"/>	
I10427	Mount Midline Joystick on Right Armrest	No Charge <input type="checkbox"/>	
I103612-99-0	VR2 Basic Joystick	No Charge <input type="checkbox"/>	
I103814-99-0	<i>Cannot be selected with any Corpus Seat that has any power seat functions. If selecting any power seat function, you must upgrade to a R-net Joystick (ex: I103544-99-0).</i> Compact Advanced Joystick for R-net	1,650.00 <input type="checkbox"/>	
I103544-99-0	<i>Compact sized joystick that does not contain an LCD Display. This is an expandable/remote joystick used with multiple power option bases or as an upgrade for no power option or single power option bases. Must order expandable controller (PN: I10411) with this part. Contains two mono jacks for mode & on/off.</i> R-net Remote Joystick w/ Mono Jacks	1,495.00 <input type="checkbox"/>	
ID51040	<i>Expandable/remote joystick used with multiple power option bases or as an upgrade for no power option or single power option bases. Must order expandable controller (PN: I10411) with this part. Comes standard w/ mono jacks.</i> R-net Remote Joystick w/ Toggle Switches	1,595.00 <input type="checkbox"/>	WEX1
ID51042	<i>Expandable/remote joystick with toggle switches. Used with multiple power option bases or as an upgrade for no power option or single power option bases. Must order expandable controller (PN: I10411) with this part.</i> R-net Remote Joystick w/ Speed Pot & Single Toggle	1,595.00 <input type="checkbox"/>	WEX1
I10411	<i>Expandable/remote joystick with a single toggle switch for modes/profiles & a speed pot. Used with multiple power option bases or as an upgrade for no power option or single power option bases. Must order expandable controller (PN: I10411) with this part.</i> Expandable Controller - R-net	No Charge <input type="checkbox"/>	
I104003-99-0	<i>Must be ordered with the expandable/remote joysticks (PN: I103544-99-0, ID51040, ID51108, ID51042, and I103814-99-0).</i> Co-Pilot Attendant Control, Corpus R-net	4,175.00 <input type="checkbox"/>	
ID50882	<i>Custom charge may apply if ordering a Corpus Tubular Backrest OR a backrest height of 21" or Lower. Not available with VR2 Controllers.</i> Attendant Control for R-net Control System Only	1,084.00 <input type="checkbox"/>	
	<i>Not available with VR2 controllers.</i>		

Joystick Options

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
I10582	Output Module (IOM) for R-net Control System <i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). One output module is required for each separate external device to be controlled. Up to three output modules can be installed on one wheelchair.</i>	1,450.00 <input type="checkbox"/>	

Seating System

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
	The Seating Systems Listed Below Include: <i>Corpus Seat Frame with 1" Depth Adjustments from 13" - 18" (& up to 22" Deep with Standard Corpus Seat Frame) Back Frame with Adjustable Recline Angle, Reclining Armrest Assembly, Solid Seat Pan, Legrest Assembly with Adjustable Legrest Elevation, Weight Limit 165 lbs for Seat Depths Less Than 16". Weight Limit of 250 lbs for Seat Depths of 16" or Greater.</i>		
I10256	Corpus Junior Seat w/o Power Tilt or Recline	5,689.00 <input type="checkbox"/>	WEPN
I10257	Corpus Junior Seat w/ 45° CG Power Tilt <i>Includes: Solid Backrest, Height Adjustable Armrest and Seat Sizes Adjustable in 1" Increments</i>	6,830.00 <input type="checkbox"/>	SEC1P0005
I10258	Corpus Junior Seat w/ 150° Power Recline <i>Includes: Solid Backrest, Height Adjustable Armrest and Seat Sizes Adjustable in 1" Increments</i>	7,655.00 <input type="checkbox"/>	SEC1P0010
I10259	Corpus Jr Seat w 45° Pwr Tilt and 150° Pwr Recline <i>Includes: Solid Backrest, Height Adjustable Armrest and Seat Sizes Adjustable in 1" Increments</i>	10,695.00 <input type="checkbox"/>	SEC1P0015

Seating System Options

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
I10418	Single Seat Function Control Switch for R-net <i>Please choose either Push Buttons or Toggle Switches below.</i>	No Charge <input type="checkbox"/>	
I10419	Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control and/or ICS Switchbox. Order with R-net Controller or remote joystick options. Please indicate below if you want push buttons or toggles on the Seat Function Control Kit.</i>	1,378.00 <input type="checkbox"/>	
I10585	Multiple Seat Function Control Kit for R-net <i>MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive control and/or ICS switch box. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the Seat Function Control Kit.</i>	2,295.00 <input type="checkbox"/>	SEC1P1005
I102798-99-0	Push Buttons on Seat Function Control Kit	No Charge <input type="checkbox"/>	
I102799-99-0	Toggle Switches on Seat Function Control Kit	245.00 <input type="checkbox"/>	
I103469-99-0	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge <input type="checkbox"/>	
I103471-99-0	Mount ICS Switchbox to End of Arm Pad - Corpus <i>ICS Switchbox will be mounted to the end of the armrest pad. If this item is selected, the ICS Switchbox will NOT swing away with the Joystick Module. Please indicate below if you would like it mounted onto the left or right arm pad.</i>	No Charge <input type="checkbox"/>	
I10453	Mount ICS Switchbox on Left Armpad _____	No Charge <input type="checkbox"/>	
I10452	Mount ICS Switchbox on Right Armpad _____	No Charge <input type="checkbox"/>	
I103470-99-0	Mount ICS Switchbox on Backrest	125.00 <input type="checkbox"/>	
I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Not available with VR2 Electronics. Seat function interface box which accommodates up to 8 momentary switches. Provides an additional method to accessing the seat functions. You must still order the Multiple Seat Function Control Kit if ordering this part. 1/8" mini-mono switches not included.</i>	1,084.00 <input type="checkbox"/>	

Backrest Options

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
	BACKREST HEIGHT NOTE: <i>Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°.</i>		
I101853-99-0	Ergo Back Fixed 14" W x 21" H - Leatherette	837.00 <input type="checkbox"/>	SEB1P0005
I101863-99-0	Ergo Back Fixed 14"W x 26.5" High - Leatherette	837.00 <input type="checkbox"/>	SEB1P0005
I101854-99-0	Ergo Back Sliding 16"W x 21"H - Leatherette	837.00 <input type="checkbox"/>	SEB1P0010
I101861-99-0	Ergo Back Sliding 16"W x 26.5" H - Leatherette	837.00 <input type="checkbox"/>	SEB1P0010
I10455	Ergo Back Fixed 14" W x 21" H - Black Mesh	859.00 <input type="checkbox"/>	SEB1P0005
I10456	Ergo Back Fixed 14"W x 26.5" High - Black Mesh	859.00 <input type="checkbox"/>	SEB1P0005
I10457	Ergo Back Sliding 16"W x 21"H - Black Mesh	859.00 <input type="checkbox"/>	SEB1P0010
I10458	Ergo Back Sliding 16"W x 26.5" H - Black Mesh	859.00 <input type="checkbox"/>	SEB1P0010
I10004	Rectangular Backrest 16"W - Plastic, Fixed <i>Choose backrest height from seat pan: 19"-29" ____"</i>	1,119.00 <input type="checkbox"/>	SEB1P0015
I10240	Solid Back 16"W X 21"H <i>Standard Option with Tilt and/or Recline Systems. Includes Leatherette Back Cushion ONLY.</i>	No Charge <input type="checkbox"/>	
I10130	Solid Back 16"W x 26.5"H <i>Standard Option with Tilt and/or Recline Systems. Includes Leatherette Back Cushion ONLY.</i>	No Charge <input type="checkbox"/>	

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
I10014	Custom Backrest Height _____" <i>Backrest height measured from seat pan and cannot be lower than 19". Price in addition to above backrest.</i>	625.00 <input type="checkbox"/>	
I100972-99-0	Split Backrest <i>Split Backrest 8" from top to allow for transfers and/or transportation. This modification cannot be made on back heights lower than 26.5". Must select either PN: I10429 or I10430 below.</i>	535.00 <input type="checkbox"/>	
I10429	Split Backrest Cushion	No Charge <input type="checkbox"/>	
I10430	Do NOT Split Backrest Cushion	No Charge <input type="checkbox"/>	
I10284	Flat Backrest Cushion for Ergo Back - Leatherette <i>This style of cushion is ONLY available in Leatherette. Recommended cushion when ordering lateral supports. Lateral wedges will be omitted if ordering this part.</i>	89.00 <input type="checkbox"/>	
I10318	ROHO Backrest Cushion For Ergo Back <i>Price Includes: ROHO MID PROFILE Backrest Cushion which will be matched to fit the Ergonomic Backrest that you select. Approximate height of MID PROFILE cells is 3" in their inflated state. Also includes a Mesh Cover for the ROHO Backrest Cushion, Hand Pump, and Patch Kit. This will NOT fit these part numbers: (I10004, I10005, I10111, I10112, I305510-00-90-0). Please indicate below if you would like to OMIT the Ergo Back Cushion.</i>	745.00 <input type="checkbox"/>	
I10319	ROHO Backrest Cushion For Ergo Back - Notched <i>Recommended ROHO cushion when ordering Lateral Trunk Supports. Please indicate below if you would like to OMIT the Ergo Back Cushion.</i>	745.00 <input type="checkbox"/>	
I10320	Omit Ergonomic Backrest Cushion	No Charge <input type="checkbox"/>	
I10278	Notch Backrest Cushion for Lateral Supports <i>Backrest Cushion will be notched 2" in on either side to allow for closer adjustment of trunk supports.</i>	185.00 <input type="checkbox"/>	

Seat Sizing

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I10264	Corpus Junior Seat Width - Standard <i>Corpus seat frame is 17" wide at the outside and 15" wide between the seat rails. Actual hip width can be accommodated with size appropriate cushion and/or hip guides.</i>	No Charge <input type="checkbox"/>	

Seat Sizing

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I10265	Corpus Junior Seat 13" Deep	No Charge <input type="checkbox"/>	
I10266	Corpus Junior Seat 14" Deep	No Charge <input type="checkbox"/>	
I10267	Corpus Junior Seat 15" Deep	No Charge <input type="checkbox"/>	
I10268	Corpus Junior Seat 16" Deep	No Charge <input type="checkbox"/>	
I10269	Corpus Junior Seat 17" Deep	No Charge <input type="checkbox"/>	
I10270	Corpus Junior Seat 18" Deep	No Charge <input type="checkbox"/>	

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
I10038	Ergonomic Seat Cushion Leatherette	475.00 <input type="checkbox"/>	SES1P0005
	<i>Not available with seat widths less than 16" or seat depths less than 14".</i>		
I10461	Ergonomic Seat Cushion - Black Mesh	485.00 <input type="checkbox"/>	SES1P0005
	<i>Not available with seat widths less than 16" or seat depths less than 14".</i>		
I10276	Miniflex Seat Cushion Leatherette	475.00 <input type="checkbox"/>	
	<i>Approximately 2" thick. Not available with seat widths or depths greater than 14".</i>		
I10277	Miniflex Seat Cushion Fabric	485.00 <input type="checkbox"/>	
	<i>Approximately 2" thick. Not available with seat widths or depths greater than 14".</i>		

Armrest Options

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I103339-99-0	4"x10" Corpus Arm Right - Dual Taper, Leatherette	No Charge <input type="checkbox"/>	SEAND0010
I103341-99-0	4"x12" Corpus Arm Right - Dual Taper, Leatherette	No Charge <input type="checkbox"/>	SEAND0010
I103343-99-0	4"x14" Corpus Arm Right - Dual Taper, Leatherette	No Charge <input type="checkbox"/>	SEAND0010
I103331-99-0	4"x10" Corpus Arm Right - Dual Taper, Black Fabric	No Charge <input type="checkbox"/>	SEAND0010
I103333-99-0	4"x12" Corpus Arm Right - Dual Taper, Black Fabric	No Charge <input type="checkbox"/>	SEAND0010
I103335-99-0	4"x14" Corpus Arm Right - Dual Taper, Black Fabric	No Charge <input type="checkbox"/>	SEAND0010

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I103338-99-0	4"x10" Corpus Arm Left - Dual Taper, Leatherette	No Charge <input type="checkbox"/>	SEAND0010
I103340-99-0	4"x12" Corpus Arm Left - Dual Taper, Leatherette	No Charge <input type="checkbox"/>	SEAND0010
I103342-99-0	4"x14" Corpus Arm Left - Dual Taper, Leatherette	No Charge <input type="checkbox"/>	SEAND0010
I103330-99-0	4"x10" Corpus Arm Left - Dual Taper, Black Fabric	No Charge <input type="checkbox"/>	SEAND0010
I103332-99-0	4"x12" Corpus Arm Left - Dual Taper, Black Fabric	No Charge <input type="checkbox"/>	SEAND0010
I103334-99-0	4"x14" Corpus Arm Left - Dual Taper, Black Fabric	No Charge <input type="checkbox"/>	SEAND0010

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
I10137	Height Adjustable Armrest Assembly	284.00 <input type="checkbox"/>	
I10016	Height Adjustable Armrest Assembly - No Charge	No Charge <input type="checkbox"/>	
	<i>No charge with power tilt and/or recline.</i>		
I101003-99-0	Armrest Mount High	No Charge <input type="checkbox"/>	
	<i>Range of adjustment in high position is 10" - 12.5" measured from the seat pan to the top of the armrest pad.</i>		
I101002-99-0	Armrest Mount Low (High Position Standard)	No Charge <input type="checkbox"/>	
	<i>Range of adjustment in low position is 8" - 10.5" measured from the seat pan to the top of the armrest pad.</i>		
IARMADAP1"	Arm Trough Adaptor 1" Pair	305.00 <input type="checkbox"/>	SEAND0025
	<i>For aftermarket arm troughs.</i>		

Armrest Options

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
IARMADAP7/8"	Arm Trough Adaptor 7/8" Pair <i>For aftermarket arm troughs.</i>	305.00 <input type="checkbox"/>	SEAND0025
I100671-99-0	Armrest Pouch-Right	190.00 <input type="checkbox"/>	
I100672-99-0	Armrest Pouch-Left	190.00 <input type="checkbox"/>	
I101261-99-0	Armrest Side Support Right Outer <i>For Corpus Armrest Pads only.</i>	142.00 <input type="checkbox"/>	
I101262-99-0	Armrest Side Support Left Outer <i>For Corpus Armrest Pads only.</i>	142.00 <input type="checkbox"/>	
I101263-99-0	Armrest Side Support Right Inner <i>Recommend 2" wider width between the armrest to ensure armrests will flip-up.</i>	142.00 <input type="checkbox"/>	
I101264-99-0	Armrest Side Support Left Inner <i>Recommend 2" wider width between the armrest to ensure armrests will flip-up.</i>	142.00 <input type="checkbox"/>	
I10285	Mount Corpus Armrest Pads On Inside of Bracket <i>Used to achieve 14"W between the arm measurement when ordering a 14"W Backrest.</i>	No Charge <input type="checkbox"/>	

Legrest Options

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I102206-99-0	Manual Adjustment Legrest Elevation	No Charge <input type="checkbox"/>	
I103534-99-0	Power Articulating Center Mount Elevating Legrest	2,825.00 <input type="checkbox"/>	SEF1P0005
Select Options:			
Part Number	Description	MSRP - CAD	ADP Code
I103055-99-0	UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad. Recommended pad to be used with lower leg lengths of 12" or less.</i>	65.00 <input type="checkbox"/>	SEFND0010
I102982-99-0	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad.</i>	65.00 <input type="checkbox"/>	SEFND0010
I102978-99-0	UT Calf Support Kit (H) w Mesh Cover & Adj Hardware <i>7"W x 7"T Curved Pad.</i>	312.00 <input type="checkbox"/>	SEFND0010
I100561-99-0	Legrest Assm. 1pc Footplate (Complete)	395.00 <input type="checkbox"/>	SEFND0015
I102051-99-0	Legrest Assm. 2pc Narrow Footplates (Complete)	No Charge <input type="checkbox"/>	
I10271	Corpus Junior Lower Leg Length <i>Modification to allow for lower leg lengths of 8" - 11". Standard length is 12-17". Measurement is taken from seat pan to the top of the footplates. Please specify lower leg length measurement in inches: _____</i>	No Charge <input type="checkbox"/>	

Accessories

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
I10400	Bodypoint Non-Padded Hip Belt, UniTrack - Medium	No Charge <input type="checkbox"/>	
I10404	Bodypoint Padded Hip Belt UniTrack (C300) - Medium	119.00 <input type="checkbox"/>	
I10402	Bodypoint Non-Padded Hip Belt, UniTrack - Large	92.00 <input type="checkbox"/>	
I10408	Bodypoint Padded Hip Belt UniTrack - Large	155.00 <input type="checkbox"/>	
I100519-99-0	Positioning Belt Retractable	235.00 <input type="checkbox"/>	
I100748-99-0	Miniflex Headrest Removable - Fabric <i>Must also select PN I10552.</i>	410.00 <input type="checkbox"/>	
I100749-99-0	Miniflex Headrest Removable - Leatherette <i>Must also select PN I10552.</i>	410.00 <input type="checkbox"/>	
I100508-99-0	Headrest-Removable Fabric-For Ergo Back <i>Must also select PN I10552.</i>	410.00 <input type="checkbox"/>	

Accessories

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
I100509-99-0	Headrest-Removable Leatherette-For Ergo Back <i>Must also select PN I10552.</i>	410.00 <input type="checkbox"/>	SEH1P0005
I103460-99-0	UT Ergo Headrest w/ Black Fabric Cover <i>Must also select PN I10552.</i>	504.00 <input type="checkbox"/>	
I1823249	UT Lateral Adjustment Bar for UniTrack Headrest <i>Allows for approximately 4" of lateral adjustment to the UT Ergo Headrest either right or left. NOTE: Due to the width of this item (approx 9"), it may be necessary to send other items such as back pack clips unmounted. Custom charge may apply with lower back heights.</i>	55.00 <input type="checkbox"/>	
IHEADADAPTER	Universal Headrest Adapter <i>For aftermarket headrests.</i>	197.00 <input type="checkbox"/>	
I10552	Adj Removable - Headrest Hardware (1 Each)	110.00 <input type="checkbox"/>	SEHND0020
I103070-99-0	UT Amp. Support Pad (H) w/ Mesh Cover - (1 Each) <i>7"W x 7"T Curved Pad. Must also select PN I10553.</i>	310.00 <input type="checkbox"/>	SEFND0005
I103071-99-0	UT Amp. Support Pad (D) w/ Mesh Cover - (1 Each) <i>8"W x 8"T Flat Pad. Must also select PN I10553.</i>	310.00 <input type="checkbox"/>	SEFND0005
I10553	Adj Removable - Amp Support Hardware (1 Each) <i>Please indicate if you would like the amputation support mounted on the right or left side of the seat frame: Mount Amp Suppt on Right: _____ Mount Amp Suppt on Left: _____</i>	110.00 <input type="checkbox"/>	
I103077-99-0	UT Hip Suppt (B) w/ Mesh Cover - Low (Pair) <i>4.5"W x 3.5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	410.00 <input type="checkbox"/>	SELND0015
I103079-99-0	UT Hip Suppt (A) w/ Mesh Cover - Low (Pair) <i>6"W x 5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	410.00 <input type="checkbox"/>	SELND0015
I103094-99-0	UT Thigh/Hip Suppt (J) w/ Mesh Cover - High (Pair) <i>6.5"W x 5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	410.00 <input type="checkbox"/>	SELND0015
I103095-99-0	UT Thigh/Hip Suppt (C) w/ Mesh Cover - High (Pair) <i>8"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	410.00 <input type="checkbox"/>	SELND0015
I102977-99-0	UT Thigh/Hip Suppt (G) w/ Mesh Cover - High (Pair) <i>12"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	410.00 <input type="checkbox"/>	SELND0015
I10550	Adj Removable - Thigh Support Hardware (Pair)	110.00 <input type="checkbox"/>	
I10273	Fixed Trunk Supports Jr Corpus - Fabric <i>Must also select PN I10551.</i>	425.00 <input type="checkbox"/>	
I10272	Fixed Trunk Supports Jr Corpus - Leatherette <i>Must also select PN I10551.</i>	410.00 <input type="checkbox"/>	
I10275	Swing Away Trunk Supports Jr Corpus - Fabric <i>Must also select PN I10551.</i>	425.00 <input type="checkbox"/>	
I10274	Swing Away Trunk Supports Jr Corpus - Leatherette <i>Must also select PN I10551.</i>	410.00 <input type="checkbox"/>	
I10551	Adj Removable - Trunk Support Hardware (Pair)	110.00 <input type="checkbox"/>	
I100978-99-0	Abductor Kit - Leatherette <i>Must also select PN I10555.</i>	410.00 <input type="checkbox"/>	SESND1043
I10555	Adj Removable - Abductor Hardware (1 Each)	110.00 <input type="checkbox"/>	
I10015	Knee Block Assembly Corpus	630.00 <input type="checkbox"/>	
I103651-99-0	Transfer Handles for UniTrack Accessory Rail	390.00 <input type="checkbox"/>	
I100975-99-0	Push Handles for Ergo Backrest <i>Order with Ergonomic or Rectangular Backs. Custom charge may apply if ordering push handles with an articulating vent tray.</i>	339.00 <input type="checkbox"/>	
I307103-00-0	Push Handles for Tubular Backrest <i>Order with Corpus Tubular Backs (PN: I10111 & I10112).</i>	337.00 <input type="checkbox"/>	
I101473-99-0	Back Pack Clips Carabiner Style	75.00 <input type="checkbox"/>	
I102059-99-0	Medical Necessities Bag <i>Includes Back Pack Clips (Carabiner Style).</i>	265.00 <input type="checkbox"/>	
I102697001331	Medical Necessities Bag - Case Logic	265.00 <input type="checkbox"/>	

Accessories

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
I100578-99-0	Upper Extremity Support Small, Angle Adjustable <i>Must also select PN I10556.</i>	356.00 <input type="checkbox"/>	SETND0010
I100575-99-0	Upper Extremity Support Angle Adjustable. 14" W <i>Must also select PN I10556.</i>	356.00 <input type="checkbox"/>	SETND0010
I100576-99-0	Upper Extremity Support Angle Adjustable. 16" W <i>Must also select PN I10556.</i>	356.00 <input type="checkbox"/>	SETND0010
I100571-99-0	Upper Extremity Support Joystick Cutout 16" W <i>Must also select PN I10556.</i>	356.00 <input type="checkbox"/>	SETND0010
I100573-99-0	Upper Extremity Support Mid Joystick Cutout 16"W <i>Must also select PN I10556.</i>	356.00 <input type="checkbox"/>	SETND0010
I103426-99-0	Upper Extremity Support Angle Adjust - UniTrack <i>Side mounted tray that attaches to the UniTrack Accessory Rail. Dimensions are: 19.5"W x 13"D. Can only be ordered with chairs that offer the UniTrack System. Must also select PN I10556.</i>	356.00 <input type="checkbox"/>	SETND0010
I10556	Adj Removable - Tray Hardware (1 Each)	110.00 <input type="checkbox"/>	
I101964-99-0	Oxygen Holder for PS, Corpus or Standing Seat <i>Cannot be ordered with a Corpus Tubular Backrest. (I10111 or I10112). Custom charge may be applied if ordering in combination with accessories such as medical necessities bags, attendant controls, push handles, back pack clips, etc.</i>	224.00 <input type="checkbox"/>	WEPX
I101956-99-0	Crutch Holder for PS, Corpus or Standing Seat <i>Cannot be ordered with a Corpus Tubular Backrest. (I10111 or I10112). Accessories such as medical necessities bags and push handles may need to be sent as parts if ordering a crutch holder.</i>	224.00 <input type="checkbox"/>	

Alternative Drive Systems

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
	R-net Alternative Drive Input Devices <i>R-net Alternative Drive Input Devices take the place of, or may be used in addition to, the standard R-net Joystick. They are only usable on wheelchairs that have R-net electronics. You MUST select one of each of the following items: (1) Input Device, (2) Expandable Controller, (3) Adjustable/Removable Alt Drive Hardware, AND (4) Omni or Input Module.</i>		
	(1) INPUT DEVICES FOR ALTERNATIVE DRIVE CONTROL		
I10558	Gimble Joystick set up for Armrest Mount	1,388.00 <input type="checkbox"/>	WEX2
I10559	Gimble Joystick set up for Midline Mount	1,388.00 <input type="checkbox"/>	WEX2
I10560	Gimble Joystick set up for Swing Away Chin Boom	1,388.00 <input type="checkbox"/>	WEX2
I10561	Gimble Joystick w/ Haims Harness - Small <i>Custom charge may be applied if ordering lower backrest heights.</i> <i>For neck width up to 5".</i>	1,555.00 <input type="checkbox"/>	WEX2
I10562	Gimble Joystick w/ Haims Harness - Large <i>For neck width up to 7".</i>	1,555.00 <input type="checkbox"/>	WEX2
I10563	Gimble Joystick w/ Power Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	5,174.00 <input type="checkbox"/>	WEX2
I10564	Mini Joystick set up for Armrest Mount	2,035.00 <input type="checkbox"/>	WEX2
I10565	Mini Joystick set up for Midline Mount	2,035.00 <input type="checkbox"/>	WEX2
I10566	Mini Joystick set up for Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	2,035.00 <input type="checkbox"/>	WEX2
I10567	Mini Joystick w/ Haims Harness - Small <i>For neck width up to 5".</i>	2,135.00 <input type="checkbox"/>	WEX2
I10568	Mini Joystick w/ Haims Harness - Large <i>For neck width up to 7".</i>	2,135.00 <input type="checkbox"/>	WEX2
I10569	Mini Joystick w/ Power Swing Away Chin Boom <i>Custom charge may be applied if ordering lower backrest heights.</i>	5,174.00 <input type="checkbox"/>	

Alternative Drive Systems

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I1821750	Tablet Control Pad <i>The Tablet Control is a unique drive input which allows a user to control a power wheelchair in a very similar manner to a mouse being controlled by a touchpad on a laptop computer. The dimensions of the Tablet Control Pad are: 12"W x 9.5"D x 1.5"H. Kit also includes: Buddy Button Switch, and Tablet Control Mounting Bracket.</i>	4,085.00 <input type="checkbox"/>	WEX3
I10570	Head Array - Adult 5 Switch HDR-A5	4,075.00 <input type="checkbox"/>	WEX4
I10571	Head Array - Adult 3 Switch HDR-A4 <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	4,075.00 <input type="checkbox"/>	WEX4
I10572	Head Array - Pediatric 5 Switch HDR-C5-RR	4,075.00 <input type="checkbox"/>	WEX4
I10573	Sip & Puff System <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	1,409.00 <input type="checkbox"/>	WEX5
I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	1,929.00 <input type="checkbox"/>	WEX5

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
I10076	Right Side Mount For Joystick	No Charge <input type="checkbox"/>	
I10077	Left Side Mount For Joystick	No Charge <input type="checkbox"/>	
I10074	Right Side Mount For Omni	No Charge <input type="checkbox"/>	
I10075	Left Side Mount For Omni	No Charge <input type="checkbox"/>	
I10078	Mount Seat Control Panel on Right	No Charge <input type="checkbox"/>	
I10079	Mount Seat Control Panel on Left	No Charge <input type="checkbox"/>	
I10080	Mount Seat Control Panel on Backrest	No Charge <input type="checkbox"/>	

Alt Drive Systems Setup

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I10576	(2) EXPANDABLE CONTROLLER (REQUIRED) Expandable Controller - Alternative Drive	No Charge <input type="checkbox"/>	

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
I10577	(3) ADJUSTABLE REMOVABLE HARDWARE (REQUIRED) Adjustable Removable Hardware - Alt Drive	110.00 <input type="checkbox"/>	

Must select one of the following:

Part Number	Description	MSRP - CAD	ADP Code
ID50881	(4) OMNI OR INPUT MODULE (REQUIRED) <i>Cannot select Input Module with Sip-n-Puff Input Devices (PN I10573 or I10574). Omni Control Module for R-net</i>	2,055.00 <input type="checkbox"/>	
I10581	<i>Please specify which type of alternative drive input you will be using if not ordering the input device from Permobil. Input Device: _____</i> Input Module (IOM) for R-net Control System <i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head Array or RIM control.</i>	1,450.00 <input type="checkbox"/>	

Alt Drive Accessories

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
-------------	-------------	------------	----------

Alt Drive Accessories

Select Options:

Part Number	Description	MSRP - CAD	ADP Code
IMICROSWITCH	Microlite Switch	148.00 <input type="checkbox"/>	
IBUDDYBUTTON	Buddy Button	111.00 <input type="checkbox"/>	
I5808	Mini Cup Switch	108.00 <input type="checkbox"/>	
IZTS-P	Proximity Switch	515.00 <input type="checkbox"/>	
IMINIJOY5SW	Mini Joystick 5 Switch OC	875.00 <input type="checkbox"/>	
I5018	Wafer Board w/5 Switches	497.00 <input type="checkbox"/>	

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code
- The client's measurements
- The client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil Ltd is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____