



permobil



## EMPOWERMENT RULES.

With its new clean design and R-Net electronics, there is no protesting the fact the C500 VS is the best seat in town. Perfect for both indoor and outdoor use, the C500 VS has enough power and torque to navigate through rough terrain. Yet the independent suspension of the C500 VS provides a smooth ride. Users can maneuver in tight spaces with a 26" turning radius, perfect for making your way through a crowd. Plus, the seating system with three programmable memory positions allows you to stand upright from any position as well as to drive while standing—an ideal feature when leading the latest rally for your rights.

**C500**  
VS STANDER



### POWER FUNCTIONS



R-Net electronics with ICS




Power sit to stand from any position. Drive at reduced speed in standing position.

### SHROUD COLORS



Arctic White



Crystal Blue



Vintage Wine



Metallic Green



Platinum Silver



Black Onyx



Power adjustable recline operates independently from power adjustable legrests



Power adjustable tilt to 25°



Power adjustable legrests operate independently from power adjustable recline

### SPECIFICATIONS

Maximum Speed	5 MPH	Gradient Capacity	12°
Maximum Speed w/ ESP	6.5 MPH	Power Seat Elevator	10"
Driving Range	19 - 25 Miles	Seat Widths	17" - 21" (adjustable)
Max User Weight	220 lbs (180° recline)	Seat Depths	14" - 22" (adjustable in 1" increments)
Max User Weight	221 - 265 lbs (140° recline)	Backrest Heights	21" & 26"
Ground Clearance	3"	Seat To Floor Height	18.5" - 21"
Battery Capacity	Group 24	Electronics (standard)	R-Net
Width	25.75" - 28.5"	3 memory seat positions available on all actuators	
Length	38"	Alternative Drive Controls	Yes
Recommended Battery	Group 24		



THE POWER OF MOBILITY

Permobil Inc.  
 6961 Eastgate Boulevard  
 Lebanon, Tennessee 37090  
 800.736.0925  
[www.permobil.com](http://www.permobil.com)

C500VS050508

