



## THE NEED FOR SPEED

Whether racing to a meeting or crossing campus to get to class, Permobil understands the demands of daily life and is proud to introduce the C300 High-Speed (HS) Package. For proven performance and reliability, the C300 HS utilizes the same motor and gearbox drive train as the C350. With a 120A power module and a top-speed of 6.5 MPH, the C300 HS offers added speed without sacrificing power or torque. Plus, the Enhanced Steering Performance (ESP) system, which operates in combination with the R-Net drive electronics, provides directional stability and ensures safe travel at high speeds...perfect for keeping pace in an ever-changing world.

# C300





**POWER FUNCTIONS**



**C300 HIGH-SPEED PACKAGE INCLUDES:**

- 120 Amp R-Net Power Module
- ESP (Enhanced Steering Performance)
- Powerful C350 6.5 MPH Motor Package

Designed to move at the speed of life



Powerful 120 Amp R-Net Electronics



Top speed of 6.5 MPH



High-speed package designation



ESP standard for added safety & stability

**SHROUD COLORS**

- Crystal Blue
- Vintage Wine
- Metallic Green
- Platinum Silver
- Black Onyx
- Chili Pepper Red
- Popstar Pink
- Titan Blue

**SPECIFICATIONS**

Maximum Speed	6.5 MPH	Gradient Capacity	10° Up 6° Dn
Driving Range	15-19 Miles	Power Seat Elevator	8"
Max User Weight	265 lbs	Seat Widths	17" - 23" (Adjustable)
Ground Clearance	3"	Seat Depths	16" - 22" (Adjustable)
Battery Capacity	60Ah	Backrest Heights	21.5" - 29.5"
Base Width	24.25"	Seat To Floor Height	17.5" - 19.5"
Base Length	43"	Electronics (standard)	R-Net
Recommended Battery	Group 34	Alternative Drive Controls	Yes w/R-Net



THE POWER OF MOBILITY

Permobil Inc.  
 6961 Eastgate Boulevard  
 Lebanon, Tennessee 37090  
 800.736.0925  
 www.permobil.com