



**C350**



THE POWER OF MOBILITY

# C350



C350 Corpus™

C350 PS

## A class of its own.

Building on the success of the C300 series, Permobil introduces the rear-wheel drive C350 power base. This cool new base accommodates a variety of power seating systems from the Corpus to the PS. Defying classification, the C350 is somewhere between a low-rider with a seat-to-floor height of 18.25" and a hot rod with two different speed packages of up to 6.5 MPH. Quick access to the battery allows users to easily check under the hood, and classic paint colors truly showcase this ride. Available in Black Onyx, Crystal Blue, Extreme Yellow, Platinum Silver and Vintage Wine, the C350 comes with standard matching tri-spoke wheel covers to complete the look.



THE POWER OF MOBILITY

[www.permobil.com](http://www.permobil.com)



**C350 Corpus™**

The Corpus seating system is now available on the all new rear-wheel drive C350 base. The Corpus seating system allows you to recline and even elevate the seat a full eight inches. Customization of the Corpus seating system is almost limitless with depth, width, height, and angle of multiple seat components being easily adjustable. Ergonomically pleasing with power tilt and recline, the Corpus is the world's most popular powered seating system.



**SHROUD COLORS**

	
Crystal Blue	Extreme Yellow
	
Black Onyx	Platinum Silver
	
Vintage Wine	

**POWER FUNCTIONS**



R-Net electronics with ICS Intelligent control system. Two programmable tilt and power legrest positions.



4-wheel adjustable independent suspension



Tri-spoke split rims with color matching hubcaps



Easy access external circuit breaker switch

Power tilt and recline



Power seat elevator



Power legrest





### C350 PS

Adjustability and versatility define the PS. The telescoping frame and expanding seat pan allow for full range of growth. The PS is constructed of a durable tubular material, designed to support after-market seating options. The adjustability is unmatched – 5 degrees at the footplate, armrests that flip up and out of the way, and footplate hangers that swing away or inside.



### POWER FUNCTIONS



VR2 electronics



Power seat elevator

\*shown with optional after-market seat cushion



Power tilt

### C350 CORPUS™ SPECIFICATIONS

Battery Capacity	Group 34
Base Width	24.25"
Base Length w/o Front Riggings	36.5"
Drive Wheel Suspension	Yes
Battery Charger Capacity	8A
Weight of Wheelchair w/Batteries	350 lbs.
Seat Elevator	Yes/8"

Drive Motor	4 Pole
Seat Widths	17" - 23" (Adjustable)
Seat Depths	16" - 22" (Adjustable)
Backrest Heights	14" - 29"
Seat To Floor Height w/ Tilt	18.25"
Power Tilt Angle	45°
Power Recline Angle	150°

Electronics (Standard)	R-Net
Electronics (Upgraded)	R-Net
Alternative Drive Controls	Yes
Effect of Driving	
When Tilted/Reclined	1/2 Speed @ 135° back angle
Color Options	Crystal Blue
	Extreme Yellow      Black Onyx
	Platinum Silver      Vintage Wine



THE POWER OF MOBILITY

Permobil Inc.

6961 Eastgate Boulevard • Lebanon, Tennessee 37090 • 800.736.0925

www.permobil.com