



THE POWER OF MOBILITY™

An almost imperceptible Co-pilot.



The most important thing when assistants operate a wheelchair is that it should be easy to control for those who are not sitting in the chair. Permobil's new Co-pilot for assistant steered wheelchairs works completely different from earlier generations of attendant controls.

The Co-pilot's great advantage is that you hardly need to use muscle power, and that steering the wheelchair is very intuitive; as easy as steering a child's stroller.

The control system responds instantly when you regulate speed or direction. The Co-pilot also has the unique feature allowing the assistant the ability to control seat functions and speed.

Lightly Grip the control bar and the wheelchair will move in the direction that you want to go. A light touch to the right, and the chair turns without effort to the right. If you walk faster, the chair increases its speed. Proportional control makes driving easier and safer for the assistant.



Ergonomically Designed Handle folds neatly out of the way when not in use.



Ergonomic positioning of the handle makes operation easier and provides better comfort for the assistant.

For more information and pricing contact your local Permobil Sales Representative or call 800-736-0925. Or visit our website at www.permobil.com