



READY. SET. ROLL!

With Permobil's new K300 PS Junior, cool kids of all sizes can get in on the action. The K300 comes with a standard adult power base, and the expandable seat pan allows for a custom fit no matter how fast a child grows. Features such as EZ adjust armrests, the power seat elevator and power tilt functions, R-Net electronics and the Uni-track system provide for maximum comfort and allow for an infinite number of adjustments and add-on accessories. Plus, the hot new colors—Chili Pepper Red, Black Onyx, Platinum Silver, Pop Star Pink and Titan Blue—and the tri-spoke wheel covers make any child the envy of the block.

K300

permobil

THE POWER OF MOBILITY



POWER FUNCTIONS



Power seat elevator and power tilt

SHROUD COLORS



R-Net electronics with ICS



Adjustable seat pan fits kids of all sizes



Push handles



Uni-track system allows for extensive customization



One-step removable armrest



EZ adjust armrests easily adjust height angle and rotation



Easy touch swing-away footrests with adjustable high strength multi-axis footplates. Available with 70° or 90° front hanger.



Optional stroller handle with manual tilt feature easily adjusts handle height without tools



Stroller handle folds down for ease of transport

SPECIFICATIONS

Maximum Speed	5 MPH	Seat Elevator	Yes
Max User Weight	165 lbs	Seat Widths	11" - 16" (Adjustable) in 1" increments
Ground Clearance	3"	Seat Depths	10" - 18" (Adjustable)
Battery Capacity	Group 34	Backrest Heights	16" - 22"
Minimum Width	24.25"	Backrest Angle Adjustments	85° - 120°
Base Length	34.5"	Seat To Floor Height w/ Tilt	16.5"
Drive Wheel Suspension	Yes/8"	Tilt Angle	45° Power/Manual
Battery Charger Capacity	8A	Electronics (standard)	R-Net
Weight of Wheelchair w/Batteries	265 lbs	Alternative Drive Controls	Yes w/R-Net
Uphill Angle	10°	Effect of Driving When Tilted	1/2 Speed @ 20°
Downhill Angle	6°	Color Options	Chili Pepper Red
Side Slope	10°		Crystal Blue
			Black Onyx
		Platinum Silver	
		Popstar Pink	
		Titan Blue	



THE POWER OF MOBILITY

Permobil Inc.
6961 Eastgate Boulevard
Lebanon, Tennessee 37090
800.736.0925
www.permobil.com