

permobil



FREEDOM OF EXPRESSION

As unique as you are, the Street is perfect for those living in an urban environment. The Street has a top speed of 7.5 mph making it easy to zip through the city. Designed to flexibly handle narrow and awkward spaces, the 4-wheel adjustable independent suspension system provides a smooth, comfortable ride. The suspension is also visible, giving the Street a stylish and powerful look. Plus, the center of gravity tilt system with recline keeps you balanced over your center of gravity for a new level of confidence. With standard rear-wheel drive, the Street separates you from the rest of the crowd.

Street
CORPUS

* Shown with custom painted shroud



POWER FUNCTIONS



Power adjustable tilt and recline

SHROUD COLORS



Blue/Green Pearl



Black/Gold



Copper Metallic



R-Net electronics with ICS



Power adjustable tilt



Power seat height adjustment



Power adjustable recline

SPECIFICATIONS

Maximum Speed	7.5 MPH	Gradient Capacity	10°
Driving Range	15-22 Miles	Power Seat Elevator	10"
Max User Weight	250 lbs	Seat Widths	17" - 23"
Ground Clearance	2.4"	Seat Depths	14" - 22"
Battery Capacity	Group 24	Backrest Heights	16" - 26.5"
Width	24.75"	Seat To Floor Height	20"
Length	32.5"	Electronics (standard)	R-Net
Recommended Battery	Group 24	Alternative Drive Controls	Yes

* Shown with custom painted shroud



THE POWER OF MOBILITY

Permobil Inc.
6961 Eastgate Boulevard
Lebanon, Tennessee 37090
800.736.0925
www.permobil.com