

UniTrack

S Y S T E M



THE POWER OF MOBILITY



UniTrack

S Y S T E M



Simple by design.

Permobil's new UniTrack system combines a few modular parts that work together to easily provide you with the custom fit you demand. Simply choose the correct pad and adjustment link and position them exactly where needed. The adjustment and customization possibilities are endless. Kind of makes you wonder why no one else has thought of this.

That's because no one else is Permobil.



THE POWER OF MOBILITY

www.permobil.com

* Only available on the
Corpus™ seating system.



UNITRACK SYSTEM

The UniTrack system offers infinite adjustments.



ADJUSTMENT LINKS

The adjustment links allow you to rotate the pads and lock them into place. The adjustment increments are so small that it provides almost infinite adjustments.



QUICK RELEASE FASTENER

It only takes a half turn of the quick release lever to make adjustments to the position and angle of the support pads. Positioned to prevent snags and hang ups on objects.



D
8"W x 8"T



H
7"W x 7"T



G
12"W x 3.5"T



C
8"W x 3.5"T



J
5"W x 6.5"T



A
6"W x 5"T



B
3.5"W x 4.5"T

SUPPORT PADS

The support pads are available in many shapes and sizes with breathable covers to offer additional comfort. Each pad has only one part number in order to avoid confusion. The pads are interchangeable and can be used as thigh support pads or calf support pads, depending on your needs.

UNITRACK KITS



I102977-99-0
UT Thigh/Hip Support Kit (G)
w/ Mesh Cover - High



I102978-99-0
UT Calf Support Kit (H)
w/ Mesh Cover & Adj Hardware



I102982-99-0
UT Calf Support Kit (J)
w/ Mesh Cover & Adj Hardware



I103055-99-0
UT Calf Support Kit (B)
w/ Mesh Cover & Adj Hardware



I103070-99-0
UT Amp. Support Kit (H)
w/ Mesh Cover - 1 Each



I103071-99-0
UT Amp. Support Kit (D)
w/ Mesh Cover - 1 Each



I103077-99-0
UT Thigh/Hip Support Kit (B)
w/ Mesh Cover - Low



I103079-99-0
UT Thigh/Hip Support Kit (A)
w/ Mesh Cover - Low



I103094-99-0
UT Thigh/Hip Support Kit (J)
w/ Mesh Cover - High



I103095-99-0
UT Thigh/Hip Support Kit (C)
w/ Mesh Cover - High



THE POWER OF MOBILITY

Permobil Inc.
6961 Eastgate Boulevard
Lebanon, Tennessee 37090
800.736.0925

www.permobil.com