

permobil



## REBEL WITH A CAUSE

With two different speed packages of up to 6.5 MPH, the C350 Corpus is perfect for speedy getaways. Quick access to the battery allows users to easily check under the hood, and the C350 Corpus offers a seat-to-floor height of 18.25." Available in Black Onyx, Crystal Blue, Extreme Yellow, Platinum Silver and Vintage Wine, the C350 Corpus comes with standard matching tri-spoke wheel covers to help complete the look. Plus, the Corpus seating system features built-in, adjustable lumbar and trunk supports, and the modular design of the Corpus allows for growth and added seat functions. The C350 Corpus is ideal for life in the fast lane.

**C350**  
CORPUS



**POWER FUNCTIONS**



Power adjustable tilt and recline



R-Net electronics with ICS



4-wheel adjustable independent suspension



Power seat height adjustment




Power adjustable legrests

**SHROUD COLORS**



Crystal Blue




Vintage Wine



Platinum Silver



Xtreme Yellow



Black Onyx

**SPECIFICATIONS**

Battery Capacity	Group 34	Seat To Floor Height w/ Tilt	18.25"
Base Width	25"	Ground Clearance	3.1"
Base Length w/o Front Riggings	36.5"	Power Tilt Angle	45°
Drive Wheel Suspension	Yes	Power Recline Angle	150°
Battery Charger Capacity	8A	Electronics (Standard)	R-Net
Weight of Wheelchair w/Batteries	313 lbs.	Alternative Drive Controls	Yes
Max User Weight	300 lbs.Seat	Effect of Driving	1/2 Speed @ 20° tilt angle
Power Seat Elevator	8"	When Tilted/Reclined	1/2 Speed @ 120° back angle
Seat Widths	17" - 23" (Adjustable)	Color Options	Crystal Blue
Seat Depths	16" - 22" (Adjustable)		Xtreme Yellow      Black Onyx
Backrest Heights	18" - 29"		Platinum Silver      Vintage Wine



THE POWER OF MOBILITY

**Permobil Inc.**  
**6961 Eastgate Boulevard**  
**Lebanon, Tennessee 37090**  
**800.736.0925**  
**www.permobil.com**