



SWEET AND SASSY

With sleek good looks, the C300 PS offers the perfect blend of stability and agility—all in one economical package. The adjustable suspension and increased ground clearance of up to three inches ensure a smooth ride in the quiet, but powerful C300 PS. The seating system features a sliding, telescoping frame and expandable seat pan that allows for full range of inherent growth. The PS seating system also provides unmatched adjustability with 5° at the footplate, armrests that both flip up and move out of the way, and footplate hangers that swing away or inside. Constructed of universal one-inch tubing, the C300 PS is designed to support aftermarket options and endless accessories.

C300

PS



POWER FUNCTIONS



Power adjustable tilt



VR2 electronics standard with PS seating system



Expandable seat pan grows with you



Power seat elevator



Simple flip-up and removeable armrest

SHROUD COLORS

-  Crystal Blue
-  Vintage Wine
-  Metallic Green
-  Platinum Silver
-  Black Onyx
-  Chili Pepper Red
-  Popstar Pink
-  Titan Blue

SPECIFICATIONS

Maximum Speed	5 MPH	Power Seat Elevator	8"
Driving Range	15-19 Miles	Seat Widths	16" - 22" (Adjustable)
Max User Weight	300 lbs	Seat Depths	16" - 22" (Adjustable)
Ground Clearance	3.1"	Backrest Heights	22"
Battery Capacity	Group 34	Seat To Floor Height	16.75" - 19"
Width	24.25" - 28.5"	Electronics (standard)	VR2 Basic R-Net electronics optional
Length	43"	Alternative Drive Controls	Yes w/R-Net
Recommended Battery	Group 34	Effect of Driving	1/2 Speed @ 20° tilt angle
Gradient Capacity	10° Up 6° Dn		



THE POWER OF MOBILITY

Permobil Inc.
6961 Eastgate Boulevard
Lebanon, Tennessee 37090
800.736.0925
www.permobil.com