

TECHNICAL SPECIFICATIONS • C400 CORPUS JUNIOR

MEASUREMENT AND PERFORMANCE

Length	.34 in
Minimum transport length	31 in
Width	24 in
Height (floor to top of standard back)	41.5 in
Seat to floor height	19.25 in
Seat width	17 in
Seat depth	13-18-22 in
Backrest height	21/25.5 in
Armrest height	5-12.5 in
Armrest width	3-4 in
Backrest angle – adjustable	90°-150°
Leg rest angle – adjustable	90°-170°
Seat elevator	8 in
Max. user weight	250 lbs
Weight incl. batteries	350 lbs
Range	19-25 miles
Max. speed (forward)	4 mph
Turning space 180°	45 in
Ground clearance	2.8 in
Vertical obstacle climb	2.5 in
Uphill angle	10°
Downhill angle	10°
Front tire size	2.5 x 8
Rear tire size	200 x 50
Motor type	4 pole
Battery capacity	2x72AH
Battery needed (x2)	GR 24
Charging time	approx. 8 hours
Electronics available	R-Net, VR2

- Compact front wheel drive base for active users
- Small footprint combines maneuverability with advanced obstacle climbing capabilities
- Rugged, advanced suspension system
- Switched ESP (Enhanced Steering Performance) option
- Our highly regarded Corpus seat in pediatric sizing
- Seat depth grows 13-18" in 1" increments, option to grow to 22"
- Standard and optional accessories sized for smaller users
- Modular design allows for easily added seat functions
- Ergonomic back with adjustable lumbar and trunk supports
- Adjustable suspension on compact base
- Increased ground clearance close to 3 inches
- Powerful and quiet four pole motors
- Side entry battery access

STANDARD OPTIONS

Color options: Vintage Wine, Crystal Blue, Metallic Green, Black Onyx, Extreme Yellow
 Choice of seat frame width/depth: 17x13/14/15/16/17/18 in
 Choice of back width/height: 14x21/25.5 in, 16x21/25,5/29.5 in
 Choice of armrest size width/length: 4x10/12/14 in
 Choice of lower leg length: 8-12/12-17 in

ACCESSORY OPTIONS

Small lateral and hip guides
 Miniflex or Ergo cushions
 Optional joystick knobs
 Bodypoint Positioning belt
 Removable leatherette or fabric headrest in two sizes
 Back Pack Hooks
 Upper extremity support - swing away mounting
 Calf support
 Foot straps
 Vent Tray – Fixed or Articulating
 Bodypoint Ankle Huggers
 Bodypoint Chest Harness



Permobil, Inc.
 6961 Eastgate Blvd.
 Lebanon, TN 37090
 Tel: 800-736-0925 Fax: 800-231-3256
www.permobil.com