

permobil



MAKING THE GRADE

The C400 VS Jr. redefines the way standing is incorporated into a child's daily routine and allows the user to drive while in a standing position. Standing can be achieved from either a sitting or a horizontal position, enabling the child to move about throughout the day and remain comfortable. Designed specifically for children, the C400 VS Jr. features a front-wheel drive base and a rugged suspension system. The small footprint of the C400 VS Jr. combines maneuverability with advanced obstacle climbing capabilities. Plus, the modular design of the Corpus seat ensures optimal fit with room for growth. With the C400 VS Jr., every child will be at the head of the class.

C400
VS JR. STANDER



POWER FUNCTIONS



R-Net electronics with ICS

Power sit to stand from any position. Drive at reduced speed in standing position.

SHROUD COLORS



Crystal Blue




Vintage Wine



Metallic Green




Xtreme Yellow



Black Onyx



Power adjustable recline operates independently from power adjustable legrests



Power adjustable tilt



Power adjustable legrests operate independently from power adjustable recline

SPECIFICATIONS

Maximum Speed	4 MPH	Power Seat Elevator	10"
Driving Range	19 - 25 Miles	Seat Widths	17"
Max User Weight	165 lbs	Seat Depths	14" - 18" (adjustable in 1" increments)
Ground Clearance	2.8"	Backrest Heights	21" - 26.5"
Battery Capacity	Group 24	Seat To Floor Height	18.5" - 20.5"
Width	24"	Electronics (standard)	R-Net
Length	40.5"	3 memory seat positions available on all actuators	
Recommended Battery	Group 24	Alternative Drive Controls	Yes
Gradient Capacity	10°		



THE POWER OF MOBILITY

Permobil Inc.
6961 Eastgate Boulevard
Lebanon, Tennessee 37090
800.736.0925
www.permobil.com