

permobil



EMPOWERMENT RULES.

With its new clean design and R-Net electronics, there is no protesting the fact the C500 VS is the best seat in town. Perfect for both indoor and outdoor use, the C500 VS has enough power and torque to navigate through rough terrain. Yet the independent suspension of the C500 VS provides a smooth ride. Users can maneuver in tight spaces with a 26" turning radius, perfect for making your way through a crowd. Plus, the seating system allows you to stand upright from any position as well as to drive while standing—an ideal feature when leading the latest rally for your rights.

C500
VS STANDER

POWER FUNCTIONS



R-Net electronics with ICS for advanced programmability.




Power sit to stand from any position. Allows driving at reduced speed

SHROUD COLORS



Arctic White



Crystal Blue




Vintage Wine



Metallic Green



Platinum Silver



Black Onyx



Power adjustable recline



Power adjustable tilt



Power adjustable legrests operate independently from power recline

SPECIFICATIONS

Maximum Speed	5.0 MPH	Drive Wheel Suspension	Yes
High Speed	6.5 MPH	Power Seat Elevator	10"
Driving Range	19 - 25 Miles	Seat Widths	17" - 19" (Adjustable)
Max User Weight	265 lbs*	Seat Depths	16" - 22" (Adjustable)
Ground Clearance	2.75"	Backrest Heights	21" - 26.5"
Battery Capacity	Group 24	Seat To Floor Height	18.5" - 20.5"
Width	25.75" - 28.5"	Electronics (standard)	R-Net
Length	38"	Alternative Drive Controls	Yes w/R-Net
Gradient Capacity	12°		

* Back angle limited to 140° for users over 220 lbs. Tilt, table to stand not allowed for users over 220 lbs.



THE POWER OF MOBILITY

Permobil Inc.
6961 Eastgate Boulevard
Lebanon, Tennessee 37090
800.736.0925
www.permobil.com