

POWER FUNCTIONS



Power adjustable tilt



R-Net electronics with ICS



Power adjustable recline



Power seat height adjustment




Power adjustable legrests


SHROUD COLORS



Arctic White




Crystal Blue



Vintage Wine



Metallic Green



Platinum Silver



Black Onyx

SPECIFICATIONS

Maximum Speed	6.5 MPH	Seat Widths	17" - 23" (Adjustable)
Driving Range	19 - 25 Miles	Seat Depths	14" - 22" (Adjustable)
Max User Weight	250 lbs	Backrest Heights	16" - 22"
Ground Clearance	2.75"	Seat To Floor Height	19.5" - 23.5"
Battery Capacity	Group 24	Electronics (standard)	R-Net
Width	25.75" - 28.5"	Alternative Drive Controls	Yes w/R-Net
Length	38"	Color Options	Arctic White Black Onyx
Recommended Battery	Group 24		Crystal Blue Metallic Green
Gradient Capacity	12°		Platinum Silver Vintage Wine
Power Seat Elevator	8"		



THE POWER OF MOBILITY

Permobil Inc.
6961 Eastgate Boulevard
Lebanon, Tennessee 37090
800.736.0925
www.permobil.com