



POWER FUNCTIONS




R-Net electronics with ICS



Power seat elevator and power tilt

SHROUD COLORS



Crystal Blue



Platinum Silver



Black Onyx



Chili Pepper Red



Popstar Pink



Titan Blue



Adjustable seat pan fits kids of all sizes



Push handles



Uni-track system allows for extensive customization



One-step removable armrest




EZ adjust armrests easily adjust height angle and rotation



Easy touch swing-away footrests with adjustable high strength multi-axis footplates. Available with 70° or 90° front hanger.



Optional stroller handle with manual tilt feature easily adjusts handle height without tools



Stroller handle folds down for ease of transport

SPECIFICATIONS

Maximum Speed	5 MPH	Seat Elevator	Yes
Max User Weight	165 lbs	Seat Widths	11" - 16" (Adjustable) in 1" increments
Ground Clearance	3"	Seat Depths	10" - 18" (Adjustable)
Battery Capacity	Group 34	Backrest Heights	16" - 22"
Minimum Width	24.25"	Backrest Angle Adjustments	85° - 120°
Base Length	34.5"	Seat To Floor Height w/ Tilt	16.5"
Drive Wheel Suspension	Yes/8"	Tilt Angle	45° Power/Manual
Battery Charger Capacity	8A	Electronics (standard)	R-Net
Weight of Wheelchair w/Batteries	265 lbs	Alternative Drive Controls	Yes w/R-Net
Uphill Angle	10°	Effect of Driving When Tilted	1/2 Speed @ 20°
Downhill Angle	6°	Color Options	Chili Pepper Red Crystal Blue
Side Slope	10°		Black Onyx Platinum Silver
			Popstar Pink Titan Blue



THE POWER OF MOBILITY

Permobil Inc.
6961 Eastgate Boulevard
Lebanon, Tennessee 37090
800.736.0925
www.permobil.com