


permobil

POWER WHEELCHAIRS



READY SET
ROLL!

K300
PS Jr.



POWERFUL SUSPENSION

Provides a smoother ride, even over rough terrain, to help you handle the adventures ahead.



HIGHLY CONFIGURABLE

Allows you to maintain proper positioning as you grow.



MULTI-MOTION ARMREST

Perfectly positioned armrest with unique rotation and flip back movements that allow for easier transfers.

COLOR OPTIONS



Volcano Red



Pop Star Pink



Cobalt Blue



Purple Rain



Phantom Grey



Onyx Black

SPECIFICATIONS

Maximum User Weight	165 lbs
Maximum Speed	5 mph, 6.5 mph
Range	16 mi*
Base Width	24"
Base Length	40"
Minimum Turning Radius	25"
Ground Clearance	3"
Recommended Battery Type	Group 34
Weight Including Batteries	272 lbs
Drive Electronics	90A R-Net
Maximum Slope	10° Up, 6° Down†
Seat-to-Floor Height	16.5" (elevator) 16.5", 17.25", 18" (seat tube)

Optional Seat Elevator	8"
Tilt Options	0° – 45° (power) 0° – 45° (manual)
Manual Recline	85° – 120°
Legrest Elevation	70° (swing-away) 90° – 150° (center mount elevating)
Backrest Height	16.5" – 22.5"
Seat Depth	10" – 18"
Seat Width	11" – 16"
Armrest Pad Lengths	10", 13"

*Actual driving range will vary based on driving and battery conditions.
†Refer to Owner Manual for safe driving practices.