

**permobil**  
POWER WHEELCHAIRS

ADVENTURE  
AWAITS...

**K450**  
MX



### EXPLORE MORE

Uniquely engineered seating system that brings the world to your fingertips.

### KID-TOUGH TIRES

Large flat tires improve traction and reduce maintenance.

### ALL-WHEEL SUSPENSION

Suspension system designed for increased balance, improved ride and safety.

High or low, explore more.



### SPECIFICATIONS

Maximum User Weight	125 lbs
Maximum Speed	5 mph
Range	17 mi*
Base Width	25"
Base Length	35"
Minimum Turning Radius	28"
Ground Clearance	2.5"
Recommended Battery Type	Group 34
Weight Including Batteries	337 lbs
Drive Electronics	90A R-Net
Maximum Slope	10° Up, 10° Down†
Seat-to-Floor Height	3" – 26" (elevator)
Seat Elevator	8"
Power Tilt	0° – 45°
Manual Recline	85° – 105°
Manual Legrest Elevation	90° – 125°
Backrest Height	14" – 22.5"
Seat Depth	10" – 18"
Seat Width	10", 12", 14", 16"
Armrest Pad Lengths	8", 12.5"

\*Actual driving range will vary based on driving and battery conditions.  
†Refer to Owner Manual for safe driving practices.

### COLOR OPTIONS



Chili  
Pepper  
Red



Pop Star  
Pink



Extreme  
Yellow



Titan  
Blue



Platinum  
Silver



Onyx  
Black