

permobil

TILITE MANUAL WHEELCHAIRS



TWIST

PLAY. GIGGLE. GROW.

Designed from the ground up for kiddos, the TWIST encourages exploration, independence, and, yes, playfulness. There are no parts or growth kits to buy. Instead, the TWIST grows with you.

TWIST



GROWTH BUILT IN.

Other chairs say they grow, but new parts or growth kits are required. With the TWIST, growth is built in. Sized from 8" x 8" up to 15" x 16", the TWIST offers 2" of growth in seat width and 3" of growth in seat depth.

CUSTOM CONFIGURATION.

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimize your riding experience. This includes choices of wheels, handrims, tires, wheel locks, color, and more. TiLite offers more options than any other manufacturer!

TECHNICAL SPECIFICATIONS

Frame Features	Aluminum Mono-Tube Rigid Growable Full Adjustability	Camber	0°, 2°, 4°, 6°, 8°, or 12°
Total Chair Width*	19"	Center of Gravity Adjustment	6"
Seat Width	8" – 15" (2" built-in Seat Width Adj.)	Transport Options	Transit Tie-Down
Seat Depth	8" – 16" (3" built-in Seat Depth Adj.)	Transport Weight*	12 lb.
Front Seat Height	13" – 20" (up to 2.25" of built-in Front Seat Height Adj.)	Weight Limit	165 lb.
Rear Seat Height	12" – 19.5" (up to 3.5" of built-in Rear Seat Height Adj.)	HCPCS Code	E1235, E1237
Front Frame Angle	70°, 80°, or 90°	Warranty	Five Years on Frame One Year on TiLite Components
Seat Back Type	Folding Seat Back Adj. Height 8.5" – 21" Adj. Angle 80° – 101°		

**Weights and widths based on 12" x 12" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weights. For reference, a pair of standard 22" rear wheels weighs 7 lb.*


COLOR OPTIONS

HIGH GLOSS PAINT



TiLite Red Burgundy Neon Pink Super Wet Black

MATTE PAINT



Deep Blue Ice Titanium Grey Black Lava

CANDY PAINT



Ruby Red Azure Blue Electric Plum

PEARLESCENT PAINT



Midnight Blue Pearl White Pearl

METALLIC PAINT



Salsa Red Metallic Tangerine Metallic Acid Green Metallic Ocean Blue Metallic

ANODIZED COLORS



Red Pepper Burnt Orange Sapphire Ultraviolet Fuchsia Black