

# F3 Corpus



*\* Indicates a Required Field*

## Dealer Information

\* Contact: \_\_\_\_\_

Dealer Code: \_\_\_\_\_

\* Dealer Name: \_\_\_\_\_

Address: \_\_\_\_\_

\* City: \_\_\_\_\_

\* State: \_\_\_\_\_

\* Phone#: \_\_\_\_\_

\* Fax#: \_\_\_\_\_

PO#: \_\_\_\_\_

Email Address: \_\_\_\_\_

## Client Information

*Permobil recommends that the client is evaluated by a certified rehab specialist.*

\* First Name: \_\_\_\_\_

\* Last Name: \_\_\_\_\_

Diagnosis: \_\_\_\_\_

Funding Source: \_\_\_\_\_

Client Age: \_\_\_\_\_

## Client Measurements

\* Weight: \_\_\_\_\_

\* Height: \_\_\_\_\_

A) Top of Shoulders: \_\_\_\_\_

B) Chest Depth: \_\_\_\_\_

C) Chest Width: \_\_\_\_\_

D) Seat Depth: \_\_\_\_\_


E) Top of Head: \_\_\_\_\_

F) Elbow to Hand: \_\_\_\_\_

G) Seat Pan to Elbow: \_\_\_\_\_

H) Hip Width: \_\_\_\_\_

\*I) Knee to Foot: \_\_\_\_\_





**Chair Model and Colors**

**CHOOSE BASE:**

Part Number	Description	
<input checked="" type="checkbox"/> I108398-99-0	F3 Base Corpus <i>Price Includes: FWD F3 Power Base, ComfortRide Fully Independent Adjustable Suspension, ESP (Enhanced Steering Performance), Tie-Down Hardware for Strap Systems, Charger, Flat Free or Pneumatic Tires, a Standard Positioning Belt and Permobil Connect. Max Speed = 10kmh, Weight Capacity = 136kg, STFH = 17.5".</i>	

**SUPPORT WHEEL OPTIONS:**

Part Number	Description	Price AUD
<input checked="" type="checkbox"/> I107330-99-0	Support Wheels for F3	No Charge
<input type="checkbox"/> I11313	Omit Front Support Wheels <i>User Weight must not exceed 220lbs. Max speed will be limited to 5mph on the F3 Corpus and 6mph on the F5 Corpus. Additional restrictions may apply based on chair configuration.</i>	No Charge

**CHOOSE DRIVE TIRES:**

Part Number	Description	Price AUD
<input type="checkbox"/> I108414-99-0	Flat Free Tire 3"W x 14"T Black - Standard Tread <i>Designed for a smooth ride on indoor and outdoor surfaces. Includes flat free standard tread caster wheels.</i>	No Charge
<input type="checkbox"/> I108411-99-0	Pneumatic Tire 3"W x 14"T Black - Standard Tread <i>Includes flat free standard tread caster wheels.</i>	No Charge
<input type="checkbox"/> I11333	Flat Free Tire 3"W x 14"T Black - Aggressive Tread <i>Includes flat free aggressive tread caster wheels and mud guards.</i>	

**ADDITIONAL WHEEL OPTIONS:**

Part Number	Description	Price AUD
<input type="checkbox"/> I108412-99-0	Pneumatic/Snow Tire 3"W x 14"T Black - Aggr Tread <i>Snow tires are designed for use in cold weather conditions. The nylon compound in these tires remains flexible at low temperatures, but wears more easily. Therefore, the snow tires (sold as wheel assemblies) are not recommended for year-long use and must be ordered along with standard drive wheels. Does not include caster wheels.</i>	
<input type="checkbox"/> I11334	Mud Guard Kit - F3	

**CHOOSE SHROUD COLOR:**

Part Number	Description	Price AUD
<input type="checkbox"/> I108940-99-0	F3 Color Kit - Aqua	No Charge
<input type="checkbox"/> I108941-99-0	F3 Color Kit - Sky	No Charge
<input type="checkbox"/> I108942-99-0	F3 Color Kit - Violet	No Charge
<input type="checkbox"/> I108943-99-0	F3 Color Kit - Flamingo	No Charge
<input type="checkbox"/> I108944-99-0	F3 Color Kit - Chili	No Charge
<input type="checkbox"/> I108945-99-0	F3 Color Kit - Canary	No Charge
<input type="checkbox"/> I108946-99-0	F3 Color Kit - Spring	No Charge
<input type="checkbox"/> I108947-99-0	F3 Color Kit - Arctic	No Charge
<input type="checkbox"/> I108948-99-0	F3 Color Kit - Sunset	No Charge
<input type="checkbox"/> I108949-99-0	F3 Color Kit - Midnight	No Charge
<input type="checkbox"/> I108950-99-0	F3 Color Kit - Gold	No Charge
<input type="checkbox"/> I108951-99-0	F3 Color Kit - Bronze	No Charge
<input type="checkbox"/> I108952-99-0	F3 Color Kit - Platinum	No Charge
<input type="checkbox"/> I108953-99-0	F3 Color Kit - Carbon Fiber	
<input type="checkbox"/> I108954-99-0	F3 Color Kit - Mossy Oak Break-Up Country	
<input type="checkbox"/> I108955-99-0	F3 Color Kit - CDG Poly	
<input type="checkbox"/> I108957-99-0	F3 Color Kit - Fireball	
<input type="checkbox"/> I108958-99-0	F3 Color Kit - Water Lily	



**Chair Model and Colors**

**CHOOSE HUB CAP KIT:**

Part Number	Description
<input type="checkbox"/> I11343	Hub Cap Color: <i>Please write in Hub Cap Color below. Hub Cap Color will be matched to the Shroud Color if no selection is made.</i> _____
	Hydrographic Hub Cap Color: <i>Please write in Hub Cap Color below. Hub Cap Color will be matched to the Shroud Color if no selection is made.</i> _____

**OPTIONAL ITEMS:**

Part Number	Description	Price AUD
<input checked="" type="checkbox"/> I10901	Australian Freight & Handling Fee	

**Base Options**

**MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:**

Part Number	Description	Price AUD
<input type="checkbox"/> I11110	Active Height Power Adj. Seat Height - 12" Travel <i>Must also select Power Tilt when selecting this option.</i>	
<input type="checkbox"/> I11117	Seat Tube Fixed	No Charge

**SET SEAT HEIGHT:**

Part Number	Description	Price AUD
<input type="checkbox"/> I107435-99-0	Set Seat Height to 17.5"	No Charge
<input type="checkbox"/> I107433-99-0	Set Seat Height at 18.5" <i>If ordering Power Adjustable Seat Height, 1" of elevator travel is used to set seat at 18.5".</i>	No Charge
<input type="checkbox"/> I107434-99-0	Set Seat Height at 19.5" <i>If ordering Power Adjustable Seat Height, 2" of elevator travel is used to set seat at 19.5".</i>	No Charge

**MUST CHOOSE BATTERIES:**

Part Number	Description
<input checked="" type="checkbox"/> I108999-99-0	<b>BATTERY OPTIONS:</b> Permobil Batteries, Grp 24 AGM GEL, Pair <i>For maximum range. Includes Permobil VoltPro charger optimized for Permobil Batteries.</i>

**OPTIONAL ITEMS:**

Part Number	Description	Price AUD
<input type="checkbox"/> I107459-99-0	Essentials Box with Reflectors <i>Please Note: The total backrest angle will be limited to 150° to avoid contact between backrest components and the Essentials Box. New design available Fall 2020!</i>	
<input type="checkbox"/> I108417-99-0	Essentials Box with Lights <i>Please Note: The total backrest angle will be limited to 150° to avoid contact between backrest components and the Essentials Box. Must also select Light Kit. New design available Fall 2020!</i>	
<input type="checkbox"/> I108418-99-0	LED Light Kit for F3	
<input type="checkbox"/> I10715	Compact Laptop Ventilator Tray - Corpus Seat <i>Accommodates Pulmonetics LTV Series and similar sized Laptop Ventilators. The tray is 13"T, adjustable in width from 8" - 13", and adjustable in depth from 2" - 6.5". Overall backrest angle limited to 150° when mounting this ventilator tray. Custom Option. Additional drive restrictions may apply.</i>	
<input type="checkbox"/> I10787	Battery Tray for Laptop Ventilator Tray <i>Measures 11 3/8" W x 6 1/2" H x 1 5/8" D. Accommodates Pulmonetic Systems Lithium-Ion Battery Pack, CareFusion SprintPack Battery System, or Richardson PowerTech Vent Power Center. Custom Option.</i>	
<input type="checkbox"/> I10788	O2 Holder for Laptop Ventilator Tray <i>Mounts to Left or Right Side of Laptop Vent Tray. O2 Holder is 7" Tall. Width infinitely adjustable from 3.5" - 7". Depth infinitely adjustable from 3.75" - 5". Custom Option.</i>	
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	



**Joystick Options**

**PLEASE NOTE:**

*For additional input devices, refer to Alternative Drive Systems section.*

I107779-99-0

Color Joystick Module, (CJSM)

*This R-net color joystick features: A Large High-Resolution 3.5" VGA color display, control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.*

I107778-99-0

Permobil Joystick Module

*Designed exclusively for Permobil, this R-net color joystick features: a 2.8" color display, control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.*

I107829-99-0

3.5mm Y-Splitter Cable for PJSM/CJSM

*For use with the PJSM or CJSM Mode/Profile jack. Allows 2 mechanical switches to control 4 separate functions. Requires PC Programming to assign the short hold and long hold commands to each switch. Switches sold separately.*

I105246-99-0

R-net LED Remote Joystick - No Display

No Charge

*This Joystick does NOT have a Display Screen, does NOT contain mono jack ports, and does NOT include Bluetooth or Assignable Buttons. PC Programming required and profiles limited to five.*

I105700-99-0

R-net Remote Color JS w/ Lights, Heavy Duty Gimbal

*Heavy-Duty modification enables gimbal to hold up to extreme operating forces. Reinforced housing is attached to a standard Color LCD R-net Joystick. 3" long joystick knob has 3/4" diameter. Not compatible with alternative joystick handles.*

**CHOOSE JOYSTICK HANDLE TYPE:**

Part Number	Description	Price AUD
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard	No Charge
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	
<input type="checkbox"/> I105442-99-0	1.75" Black Foam Ball For Joystick Knob	
<input type="checkbox"/> I105846-99-0	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	
<input type="checkbox"/> I105847-99-0	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	
<input type="checkbox"/> I105848-99-0	Bodypoint J/S Handle Dome Shaped, Rubber	

**CHOOSE JOYSTICK MOUNT:**

Part Number	Description	Price AUD
<input type="checkbox"/> I10664	Joystick Mount Right - Fixed	No Charge
<input type="checkbox"/> I10665	Joystick Mount Left - Fixed	No Charge
<input type="checkbox"/> I102537-99-0	Joystick Mount Right - Swing Away	
<input type="checkbox"/> I102538-99-0	Joystick Mount Left - Swing Away	
<input type="checkbox"/> I108037-99-0	Slimline Retractable Joystick Mount - Right	
<input type="checkbox"/> I108038-99-0	Slimline Retractable Joystick Mount - Left	
<input type="checkbox"/> I11232	Adjustable Height Panel Bracket <i>Provides 1.77" of vertical adjustment below the armrest and 25° of angle adjustment. Compatible with all Permobil joystick mounts except for the M1 Fixed Joystick Mount.</i>	
<input type="checkbox"/> I105210-99-0	Mount Standard Joystick on Corpus Back Hinge <i>Standard R-net Joystick is mounted to backrest hinge. Can be chosen instead of PN I105408-99-0. Mount on Left: _____ Mount on Right: _____</i>	
<input type="checkbox"/> I105225-99-0	Co-Pilot Attendant Control R-net - Corpus	
<input type="checkbox"/> I105408-99-0	Attendant Control R-net	



## Joystick Options

### OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I107350-99-0	<p>R-net Bluetooth iDevice Module</p> <p><i>Allows input device commands (such as joystick nudges or switch clicks) to be used to operate an iPhone, iPad or iPod. By assigning each input to an iOS function, users can make phone calls, send texts and emails and access applications. Kit includes a Bluetooth iDevice Module, an R-net Bus Cable, and two mono jack cables for direct switch inputs (switches not included). Must be ordered with a Color JSM or Omni for mode display.</i></p>	
<input type="checkbox"/> I107367-99-0	<p>R-net Bluetooth Mouse Module</p> <p><i>Enables mouse control through any R-net input device. It can operate the mouse on a computer and on select Android devices (with R-net for Android app). Kit includes a Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click (switches not included) and a USB BT Adapter for user's computer. Note: BT Adapter may not work with all computers. Must be ordered with a Color JSM or Omni for mode display.</i></p>	
<input type="checkbox"/> I107368-99-0	<p>Add a second Bluetooth Mouse Module</p> <p><i>Allows R-net to simultaneously pair to a second compatible computer or Android device. Each module is accessible through a separate R-net mode. Kit includes a unique Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click and a USB BT Adapter.</i></p>	
<input type="checkbox"/> I10582	<p>Output Module (IOM) for R-net Control System</p> <p><i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). Each module can provide up to 4 switch closures. May require specialty cable or wiring. Up to three output modules can be installed on one wheelchair.</i></p>	

## Seating System

### CHOOSE POSTERIOR TILT/RECLINE PACKAGE:

Part Number	Description
	<p>PLEASE NOTE:</p> <p><i>A half speed limitation is imposed when the Total Backrest Angle (TBA) &gt; 150°. TBA = Tilt Angle + Backrest Recline Angle. Example 1: 45° Tilt + 100° Recline = 145° TBA (Full Speed). Example 2: 35° Tilt + 120° Recline = 155° TBA (Half Speed).</i></p>
<input type="checkbox"/> I11130	<p>Corpus Seat w/ No Tilt &amp; Manual Recline</p> <p><i>Manual Recline angle limited to 120°.</i></p>
<input type="checkbox"/> I11131	<p>Corpus Seat w/ No Tilt &amp; 180° Power Recline</p>
<input type="checkbox"/> I11132	<p>Corpus Seat w/ 50° Tilt &amp; Manual Recline</p> <p><i>Manual Recline angle limited to 120°.</i></p>
<input type="checkbox"/> I11133	<p>Corpus Seat w/ 50° Tilt &amp; 180° Power Recline</p>

### ACTIVE REACH PACKAGES:

Part Number	Description	Price AUD
	<p>PLEASE NOTE:</p> <p><i>ActiveReach should be recommended by a licensed clinical professional and/or ATP. While operating ActiveReach, positioning equipment such as lap belts and chest supports should be considered to meet the user's specific safety needs. Note: A low speed limitation is imposed from 5° of ActiveReach and Drive inhibit will occur at 10° or greater.</i></p>	
<input type="checkbox"/> I11141	<p>5° Active Reach</p> <p><i>When ordered with a Power Adjustable Seat Height, tilt range will be -5° to +50°. When ordered without a Power Adjustable Seat Height, tilt range will be -5° to +45°. If this option is selected with a fixed seat, the minimum seat to floor height at 0° of tilt will be approximately 19.5". Power Legrest recommended for full functionality. Minimum legrest angle will be ~ 97° when manual legrest is selected.</i></p>	No Charge
<input type="checkbox"/> I107090-99-0	<p>10° Active Reach</p> <p><i>Please Note: This option is only available on configurations ordered with Power Adjustable Seat Height and Power Tilt. A chest support belt (see Accessories section) is strongly recommended when selecting this option. Power Legrest recommended for full functionality. Minimum legrest angle will be ~ 103° when manual legrest is selected.</i></p>	No Charge





**Seating System**

**ACTIVE REACH PACKAGES:**

Part Number	Description	Price AUD
<input type="checkbox"/> I107091-99-0	20° Active Reach <i>Please Note: This option is only available on configurations ordered with Power Adjustable Seat Height, Power Tilt, Power Recline and Power Legrest Elevation. Must select Corpus VS Power Elevating Legrest, Corpus VS Footplates and Corpus VS Knee Support from VS Legrest Options section. A chest support belt (see Accessories section) is strongly recommended when selecting this option. Not compatible with retractable positioning belts.</i>	
<input type="checkbox"/> I111100	30° Active Reach w/Transfer <i>Please Note: This option is only available on configurations ordered with Power Adjustable Seat Height, Power Tilt, Power Recline and Power Legrest Elevation. Must select Corpus VS Power Elevating Legrest, Corpus VS Footplates and Corpus VS Knee Support from VS Legrest Options section. A chest support belt (see Accessories section) is strongly recommended when selecting this option. Not compatible with retractable positioning belts.</i>	

**Seating System Options**

<input checked="" type="checkbox"/> I10585	Seat Function Control Kit for R-net <i>Please indicate below if you would like push buttons or toggles on the ICS Switchbox.</i>	
--	---	--

**CHOOSE PUSH BUTTONS OR TOGGLES:**

Part Number	Description	Price AUD
<input type="checkbox"/> I102798-99-0	Push Buttons on ICS Switchbox	No Charge
<input type="checkbox"/> I102799-99-0	Toggle Switches on ICS Switchbox	
<input type="checkbox"/> I105409-99-0	Do Not Install ICS Switchbox, Send as Part	No Charge

**CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:**

Part Number	Description	Price AUD
<input type="checkbox"/> I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge
<input type="checkbox"/> I104651-99-0	Mount ICS Switchbox onto Armrest <i>Mount on Left: _____ Mount on Right: _____</i>	No Charge
<input type="checkbox"/> I106256-99-0	Mount ICS Switchbox on Backrest <i>Mount on Left: _____ Mount on Right: _____</i>	

**OPTIONAL ITEMS:**

Part Number	Description	Price AUD
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Seat function interface box that accommodates up to 8 momentary switches. Provides an additional method for accessing the seat functions. 1/8" mini-mono switches not included.</i>	
<input type="checkbox"/> I108526-99-0	Remote Stop System - AUS/NZ/EU <i>Consists of a 433MHz Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair.</i>	
<input type="checkbox"/> I10750	Memory Seat Program, Corpus w/ R-net	
<input type="checkbox"/> I106209-99-0	Legrest Movement Linked to Recline Electronically <i>Power Legrest will raise when Power Recline moves backward and will lower when back moves forward. If Power Legrest is operated, Power Recline will NOT automatically reposition. Feature works through ICS Switchbox and Driver Control.</i>	
<input type="checkbox"/> I10898	Enable Independent Repositioning Mode <i>IRM allows for sequential movement of power seat functions into a posterior position (35° tilt, 115° recline, 135° legrest) and returns to upright position (5° tilt, 100° recline, 110° legrest). This can be accessed using a standard drive control or a single switch. On the Permobil Joystick Module, IRM is accessible through the mode/profile jack (long switch hold). Please see the Alt Drive Accessories section for switch options.</i>	No Charge

**Backrest Options**



**Backrest Options**

**CHOOSE FIXED OR SLIDING BACKREST KIT:**

Part Number	Description
<input type="checkbox"/> I105067-99-0	Sliding Backrest Kit - 5" Mechanical Shear
<input type="checkbox"/> I105066-99-0	Fixed Backrest Kit - No Mechanical Shear

**Price AUD**  
No Charge

**CHOOSE BACKREST TYPE:**

Part Number	Description
<b>BACKREST HEIGHT NOTE:</b>	
<i>Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Corpus Ergo Backs include lateral wedges and foam lumbar pads. Solid Backs do not include foam lumbar pads.</i>	
<input type="checkbox"/> I106545-99-0	Corpus Ergo Back - 14"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>
<input type="checkbox"/> I106546-99-0	Corpus Ergo Back - 16"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>
<input type="checkbox"/> I106547-99-0	Corpus Ergo Back - 18"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>
<input type="checkbox"/> I106548-99-0	Corpus Ergo Back - 20"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>

**PLEASE NOTE:**

*Taller backrests listed below include: lateral wedges, foam lumbar pads\*, UT rail mounted at the top for lateral adjustment of the headrest, and split backrest feature with height adjustment in 1" increments. Please indicate below where you want the height set at so that a properly sized backrest cushion can be provided. \* Foam lumbar pads not included with Solid Backs.*

<input type="checkbox"/> I106532-99-0	Corpus Ergo Back - 14"W x 23" - 28"T (C3GERGBCA)
<input type="checkbox"/> I106533-99-0	Corpus Ergo Back - 16"W x 23" - 28"T (C3GERGBCA)
<input type="checkbox"/> I106534-99-0	Corpus Ergo Back - 18"W x 23" - 28"T (C3GERGBCA)
<input type="checkbox"/> I106535-99-0	Corpus Ergo Back - 20"W x 23" - 28"T (C3GERGBCA)

**MUST ALSO SELECT BACK HEIGHT:**

- Set Corpus Backrest at 23" \_\_\_\_\_
- Set Corpus Backrest at 24" \_\_\_\_\_
- Set Corpus Backrest at 25" \_\_\_\_\_
- Set Corpus Backrest at 26" \_\_\_\_\_
- Set Corpus Backrest at 27" \_\_\_\_\_
- Set Corpus Backrest at 28" \_\_\_\_\_

**CHOOSE BACKREST CUSHION:**

Part Number	Description
<input type="checkbox"/> I11335	Letherette Cushion - Corpus Ergo Back <i>Includes Dual Density Backrest Foam.</i>
<input type="checkbox"/> I11336	Stretch-Air Cushion - Corpus Ergo Back <i>Includes Dual Density Backrest Foam with Stretch-Air breathable fabric cover.</i>

**Price AUD**  
No Charge  
No Charge

**CHOOSE ONE OF THE FOLLOWING:**

Part Number	Description
<input type="checkbox"/> I109029-99-0	BodiLink Premium 3" x 5" Lateral Supports - Pair <i>Must also select I109031-99-0. Includes foam wedge for additional support. Only compatible with Corpus Ergo backrest shells.</i>
<input type="checkbox"/> I109030-99-0	BodiLink Premium 4" x 6" Lateral Supports - Pair <i>Must also select I109031-99-0. Includes foam wedge for additional support. Only compatible with Corpus Ergo backrest shells.</i>
<input type="checkbox"/> I109031-99-0	Swing-Away Adj BodiLink Lateral Hdwr Corpus Pair

**Price AUD**

**OPTIONAL ITEMS:**

Part Number	Description	Price AUD
-------------	-------------	-----------



**Backrest Options**

**OPTIONAL ITEMS:**

Part Number	Description	Price AUD
<input type="checkbox"/> I105782-99-0	<p><b>Tubular Backrest Frame - Corpus (1" Canes)</b></p> <p><i>Provides 1" diameter, vertical canes for mounting aftermarket backrests. Frame can be adjusted to support backrest widths from 14" to 20". Canes are attached with 2" long UniTrack linkage disks that can rotate to provide space for backrest hardware. The tops of the canes can be adjusted up to a height of 27" above the seat pan. Depending on the aftermarket backrest, several inches of seat depth may be lost when mounting this system.</i></p>	
<input type="checkbox"/> I106072-99-0	<p><b>Direct Backrest Frame - Corpus</b></p> <p><i>Used to mount ROHO Agility, Jay 3, Matrx Elite, Varilite Icon and Ride Java backrests ranging in width from 14" to 20". Slotted brackets and 1" UniTrack linkage disks provide a close, adjustable mount. Backrest pitch can be tuned by rotating the linkage disks. Some seat depth may be lost due to the position and orientation of the backrest.</i></p>	
<input type="checkbox"/> I10725	<p><b>Split Backrest Cushion - Corpus</b></p> <p><i>Custom Option.</i></p>	
<input type="checkbox"/> I10786	<p><b>Tailored Backrest Cushion - Corpus</b></p> <p><i>Sides of cushion are trimmed to provide additional clearance for lateral trunk supports. Must also select a backrest cushion (I10659 or I10660) when ordering this custom option.</i></p>	

**Seat Sizing**

**CHOOSE SEAT WIDTH:**

Part Number	Description	Price AUD
<input type="checkbox"/> I107023-99-0	Set Seat Width at 17": _____	No Charge
<input type="checkbox"/> I107024-99-0	Set Seat Width at 19": _____	No Charge
<input type="checkbox"/> I107025-99-0	Set Seat Width at 21": _____	No Charge
<input type="checkbox"/> I107026-99-0	Set Seat Width at 23": _____	No Charge

**CHOOSE SEAT DEPTH:**

Part Number	Description	Price AUD
<b>PLEASE NOTE:</b>		
<i>There are three seat depth kits, each providing three seat depth settings: (1) 14" - 16"; (2) 17" - 19" and (3) 20" - 22". Please select the required seat depth and the appropriate kit will be provided.</i>		
<input type="checkbox"/> I107013-99-0	Seat Depth 14"	No Charge
<input type="checkbox"/> I107014-99-0	Seat Depth 15"	No Charge
<input type="checkbox"/> I107015-99-0	Seat Depth 16"	No Charge
<input type="checkbox"/> I107016-99-0	Seat Depth 17"	No Charge
<input type="checkbox"/> I107017-99-0	Seat Depth 18"	No Charge
<input type="checkbox"/> I107018-99-0	Seat Depth 19"	No Charge
<input type="checkbox"/> I107019-99-0	Seat Depth 20"	No Charge
<input type="checkbox"/> I107020-99-0	Seat Depth 21"	No Charge
<input type="checkbox"/> I107021-99-0	Seat Depth 22"	No Charge

**CHOOSE SEAT CUSHION:**

Part Number	Description	Price AUD
<b>NOTE:</b>		
<i>See supplemental order form for ROHO seat cushion options. A Corpus II Ergonomic Seat Cushion will be provided for select seat sizes.</i>		
<input type="checkbox"/> I11337	Corpus Ergo Cushion - Leatherette Cover	
<input type="checkbox"/> I11338	Corpus Ergo Cushion - Stretch-Air Cover	
<input type="checkbox"/> I10965	<p><b>Corpus Ergo Air Cushion (C3GRHOSCA)</b></p> <p><i>Features ROHO Dry Floatation technology. The removable ROHO cell pad provides adjustable skin protection while the contoured foam and thigh support wedges help maintain proper positioning and stability. Includes removable stretch fabric cover. Available in widths of 17", 19" and 21" and depths of 17", 19" and 21", as well as 17" x 23". Not Available in 21" x 17"</i></p>	





## Seat Sizing

### CHOOSE SEAT CUSHION:

Part Number	Description	Price AUD
<input type="checkbox"/> I11113	Incontinence Cover for Ergo Cushion <i>Made from a stretch Dartex fabric and suitable for all Ergo cushions including the Corpus Ergo Air. The cover is cut to the contours of the Ergo cushion and is designed to be used as an outer layer of protection. Cover size will be matched to cushion size. Sent as a part.</i>	
<input type="checkbox"/> I105118-99-0	Omit Seat Cushion	No Charge

## Armrest Options

### CHOOSE ARMBAR KIT AND MOUNTING POSITION:

Part Number	Description	Price AUD
PLEASE NOTE:		
<i>Narrow Armbar Kit can be adjusted from 15"W to 19"W between the armrest pads without any additional parts. Wide Armbar Kit can be adjusted from 19"W to 23"W between the armrest pads without any additional parts. Please select both Armbar Kit and Width Setting below.</i>		
<input type="checkbox"/> I106543-99-0	Narrow Armbar Kit, Inf Adjustable 15"W - 19"W	No Charge
<input type="checkbox"/> I106544-99-0	Wide Armrest Kit, Inf Adjustable 19"W - 23"W <i>Required when selecting a 20" wide backrest.</i>	No Charge
Set the Width Between the Armrests to: _____ <i>Must be wider than backrest width to ensure proper fit. Please round to nearest 1".</i>		

### CHOOSE ARMREST HEIGHT MOUNTING POSITION:

Part Number	Description	Price AUD
<input type="checkbox"/> I105069-99-0	Height Adjustable Armrests - NC <i>No charge when ordering power tilt and/or recline.</i>	No Charge
<input type="checkbox"/> I105081-99-0	Height Adjustable Armrests SET ARMREST HEIGHT AT: _____ <i>Corpus Armrests are infinitely adjustable from 8" - 15", measured from the seat pan to the top of the armrest pad. No additional parts are needed to get this full range of adjustment. In one inch increments (ex: 8", 9", etc.), please indicate where you would like the armrests set.</i>	

### CHOOSE ARMREST SWIVEL FIXED OR ENABLED:

Part Number	Description	Price AUD
<input type="checkbox"/> I105071-99-0	Left Armrest Bracket - Fixed <i>Medial swivel cannot be enabled when selecting this option. Must be selected when ordering arm trough adapters.</i>	No Charge
<input type="checkbox"/> I105070-99-0	Right Armrest Bracket - Fixed <i>Medial swivel cannot be enabled when selecting this option. Must be selected when ordering arm trough adapters.</i>	No Charge
<input type="checkbox"/> I105073-99-0	Left Armrest Bracket - Medial Swivel Enabled	No Charge
<input type="checkbox"/> I105072-99-0	Right Armrest Bracket - Medial Swivel Enabled	No Charge

### CHOOSE ARMREST PAD RIGHT:

Part Number	Description	Price AUD
<input type="checkbox"/> I105100-99-0	4"x10" Corpus Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I105104-99-0	4"x13" Corpus Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I105105-99-0	4"x16" Corpus Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I105999-99-0	4"x18" Corpus Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I107212-99-0	4"x10" Corpus Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/> I106047-99-0	4"x13" Corpus Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/> I107214-99-0	4"x16" Corpus Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/> I105995-99-0	4"x18" Corpus Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/> I107591-99-0	4"x10" Corpus Arm w/ Permobil Gel Pad - Right	
<input type="checkbox"/> I107593-99-0	4"x13" Corpus Arm w/ Permobil Gel Pad - Right	
<input type="checkbox"/> I107595-99-0	4"x16" Corpus Arm w/ Permobil Gel Pad - Right	



**Armrest Options**

**CHOOSE ARMREST PAD RIGHT:**

Part Number	Description	Price AUD
<input type="checkbox"/> I107597-99-0	4"x18" Corpus HD Arm w/ Permobil Gel Pad - Right	
<input type="checkbox"/> I108861-99-0	Moldable Forearm Support - Right <i>11.3"L x 5"W. Not compatible with armrest pouch, trays (except Small Tray) or the VS Chest Support.</i>	
<input type="checkbox"/> I108863-99-0	Contoured Hand Pad Medium - Right <i>6.4"L x 5.5" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.</i>	
<input type="checkbox"/> I108865-99-0	Contoured Hand Pad Large - Right <i>7.7"L x 6.1" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.</i>	

**CHOOSE ARMREST PAD LEFT:**

Part Number	Description	Price AUD
<input type="checkbox"/> I105098-99-0	4"x10" Corpus Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I105102-99-0	4"x13" Corpus Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I105107-99-0	4"x16" Corpus Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I105998-99-0	4"x18" Corpus Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I107213-99-0	4"x10" Corpus Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/> I106046-99-0	4"x13" Corpus Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/> I107215-99-0	4"x16" Corpus Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/> I105994-99-0	4"x18" Corpus Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/> I107590-99-0	4"x10" Corpus Arm w/ Permobil Gel Pad - Left	
<input type="checkbox"/> I107592-99-0	4"x13" Corpus Arm w/ Permobil Gel Pad - Left	
<input type="checkbox"/> I107594-99-0	4"x16" Corpus Arm w/ Permobil Gel Pad - Left	
<input type="checkbox"/> I107596-99-0	4"x18" Corpus HD Arm w/ Permobil Gel Pad - Left	
<input type="checkbox"/> I108860-99-0	Moldable Forearm Support - Left <i>11.3"L x 5"W. Not compatible with armrest pouch, trays (except Small Tray) or the VS Chest Support.</i>	
<input type="checkbox"/> I108862-99-0	Contoured Hand Pad Medium - Left <i>6.4"L x 5.5" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.</i>	
<input type="checkbox"/> I108864-99-0	Contoured Hand Pad Large - Left <i>7.7"L x 6.1" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.</i>	

**OPTIONAL ITEMS:**

Part Number	Description	Price AUD
<input type="checkbox"/> I105463-99-0	Corpus Arm Trough Adapter 1" - Pair <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked (I10645 &amp; I10646) when selecting this option.</i>	
<input type="checkbox"/> I105462-99-0	Corpus Arm Trough Adapter 7/8" - Pair <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked (I10645 &amp; I10646) when selecting this option.</i> Install Arm Trough Adapters: _____ Send Arm Trough Adapters as Parts: _____	
<input type="checkbox"/> I107339-99-0	Armrest Pouch - Right	
<input type="checkbox"/> I107338-99-0	Armrest Pouch - Left	
<input type="checkbox"/> I107840-99-0	UT Armrest Side Support - Right Outer <i>Includes 2"T x 6.5"W tapered pad.</i>	
<input type="checkbox"/> I107863-99-0	UT Armrest Side Support - Left Outer <i>Includes 2"T x 6.5"W tapered pad.</i>	



**Armrest Options**

**OPTIONAL ITEMS:**

Part Number	Description	Price AUD
<input type="checkbox"/> I107864-99-0	UT Armrest Side Support - Right Inner <i>Includes 2"T x 6.5"W tapered pad. May require an increase to the width between the armrests to ensure armrests will flip-up.</i>	
<input type="checkbox"/> I107865-99-0	UT Armrest Side Support - Left Inner <i>Includes 2"T x 6.5"W tapered pad. May require an increase to the width between the armrests to ensure armrests will flip-up.</i>	
<input type="checkbox"/> I107839-99-0	UT Elbow Support - Right <i>Must also select PN I10881.</i>	
<input type="checkbox"/> I107838-99-0	UT Elbow Support - Left <i>Must also select PN I10882.</i>	
<input type="checkbox"/> I10881	Adjustable Elbow Support Hardware - Right	
<input type="checkbox"/> I10882	Adjustable Elbow Support Hardware - Left	

**Legrest Options**

**CORPUS LEGREST OPTIONS:**

Part Number	Description	Price AUD
<input type="checkbox"/> I107763-99-0	CHOOSE LEGREST TYPE: Corpus Manual Elevating Legrest	No Charge
<input type="checkbox"/> I107757-99-0	Corpus Power Elevating Legrest - 170° <i>Includes mechanical articulation.</i>	
<input type="checkbox"/> I105437-99-0	CHOOSE FOOTPLATE STYLE: 2pc Wide Corpus Footplates - 9"D x 7.5"W	No Charge
<input type="checkbox"/> I105436-99-0	2pc Narrow Corpus Footplates - 9"D x 6"W	
<input type="checkbox"/> I106760-99-0	1pc Narrow Corpus Footpl - 9.75"D x 12.5"W <i>Footplate is 12.5"W in front, 8"W in back.</i>	
<input type="checkbox"/> I105252-99-0	CALF SUPPORT OPTIONS: UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad.</i>	No Charge
<input type="checkbox"/> I105258-99-0	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad.</i>	No Charge
<input type="checkbox"/> I105250-99-0	UT Calf Support Kit (A) w Mesh Cover & Adj Hardware <i>5.5"W x 6"T Pad.</i>	No Charge
<input type="checkbox"/> I105256-99-0	UT Calf Support Kit (H) w Mesh Cover & Adj Hardware <i>7"W x 7"T Curved Pad.</i>	
<input type="checkbox"/> I105254-99-0	UT Calf Support Kit (D) w Mesh Cover & Adj Hardware <i>8"W x 8"T Pad.</i>	

**CORPUS VS LEGREST OPTIONS (FOR SELECT ACTIVE REACH PACKAGES):**

Part Number	Description	Price AUD
<input type="checkbox"/> I107758-99-0	CORPUS VS LEGREST AND KNEE SUPPORT: <i>NOTE: The Corpus VS Legrest is only available with 20°, 30° and 45° of Active Reach.</i> Corpus VS Power Elevating Legrest - 180° <i>Includes power articulation.</i>	
<input type="checkbox"/> I107038-99-0	Knee Support - Corpus VS <i>Select when ordering Corpus VS Legrest. Must also select I11129.</i>	
<input type="checkbox"/> I11129	Adj Removable Knee Support Hardware (1 Each)	
<input type="checkbox"/> I106969-99-0	CHOOSE FOOTPLATE STYLE: 1pc Corpus VS Footplate - 9"D x 16.5"W	No Charge



**Legrest Options**

**CORPUS VS LEGREST OPTIONS (FOR SELECT ACTIVE REACH PACKAGES):**

Part Number	Description	Price AUD
<input type="checkbox"/> I106971-99-0	2pc Wide Corpus VS Footplates - 9"D x 7.5"W	
<input type="checkbox"/> I106967-99-0	2pc Narrow Corpus VS Footplates - 9"D x 6"W	
<b>CALF SUPPORT OPTIONS:</b>		
<input type="checkbox"/> I107211-99-0	VS Calf Support Kit (J) w Mesh Cover & Adj Hardware 5"W x 6.5"T Pad. Compatible with footplate heights of 10" - 18".	No Charge
<input type="checkbox"/> I106190-99-0	VS Calf Support Kit (A) w Mesh Cover & Adj Hardware 5.5"W x 6"T Pad. Compatible with footplate heights of 10" - 18".	No Charge
<input type="checkbox"/> I106188-99-0	VS Calf Support Kit (H) w Mesh Cover & Adj Hardware 7"W x 7"T Curved Pad. Compatible with footplate heights of 11" - 18".	
<input type="checkbox"/> I106186-99-0	VS Calf Support Kit (D) w Mesh Cover & Adj Hardware 8"W x 8"T Pad. Compatible with footplate heights of 12" - 18".	

**FOOTPLATE HEIGHT POSITION:**

Part Number	Description
	SET FOOTPLATE HEIGHT AT: _____ <i>This is measured from the SEAT PAN to the top of the FOOTPLATE, when the legrest assembly is at 90°. DO NOT INCLUDE CUSHION THICKNESS IN YOUR MEASUREMENT. In one inch increments, please indicate desired footplate height. NOTE: Corpus footplate height range is 6"-17" and Corpus VS height range is 9"-17". Height may need to be reduced for shipping. Calf Pads are not available at heights less than 9".</i>

**OPTIONAL ITEMS:**

Part Number	Description	Price AUD
<input type="checkbox"/> I106530-99-0	Stealth Gel Pads - Narrow Footplates (Pair)	
<input type="checkbox"/> I106225-99-0	Stealth Gel Pads - Wide Footplates (Pair)	

**Accessories**

**MUST CHOOSE POSITIONING BELT:**

Part Number	Description	Price AUD
	NOTE: <i>See supplemental order form for more Bodypoint options.</i>	
<input checked="" type="checkbox"/> I104518-99-0	Bodypoint Non-Padded Hip Belt - Medium	No Charge

**CHOOSE ONE OF THE FOLLOWING:**

Part Number	Description	Price AUD
<input type="checkbox"/> I107470-99-0	Permobil 7"W Headrest w/ Link Hardware <i>Must also select I107480-99-0 or I107481-99-0.</i>	
<input type="checkbox"/> I107469-99-0	Permobil 10"W Headrest w/ Link Hardware <i>Must also select I107480-99-0 or I107481-99-0.</i>	
<input type="checkbox"/> I107468-99-0	Permobil 14"W Headrest w/ Link Hardware <i>Must also select I107480-99-0 or I107481-99-0.</i>	
<input type="checkbox"/> I107480-99-0	Adj/Removable Headrest Link Hardware - Standard	
<input type="checkbox"/> I108721-99-0	Adj/Removable BodiLink Headrest Hardware Max8 <i>Must be selected when ordering any of the Swing Away Chin Boom assemblies.</i>	
<input type="checkbox"/> I111305	Swing Away Egg Switch Assembly <i>Mount on Upper Left: _____ Mount on Upper Right: _____</i>	
<input type="checkbox"/> I111306	Swing Away Microlight Assembly <i>Mount on Upper Left: _____ Mount on Upper Right: _____</i>	
<input type="checkbox"/> I111307	Swing Away Facial Pad Assembly 4"L x 2"T <i>Mount on Upper Left: _____ Mount on Upper Right: _____</i>	



**Accessories**

**CHOOSE ONE OF THE FOLLOWING:**

Part Number	Description	Price AUD
<input type="checkbox"/> I11308	Swing Away Facial Pad Assembly 6"L x 2"T <i>Mount on Upper Left: _____ Mount on Upper Right: _____</i>	
<input type="checkbox"/> I11276	Provide Headrest Receiver for AGILITY Backrest <i>Select for AGILITY backs and third party backs ordered off chair.</i>	No Charge
<input type="checkbox"/> I106331-99-0	Permobil 10"W Headrest, Leatherette <i>Not available with ROHO AGILITY Back Systems. Must also select PN I10552.</i>	
<input type="checkbox"/> I10552	Adj/Removable Headrest Hardware	
<input type="checkbox"/> I107499-99-0	Universal Headrest Adapter <i>Please Note: Recline is limited to 160° or less (depending on backrest height) when ordering a Universal Headrest Adapter on a 23" or 24" tall Corpus backrest.</i>	

**OPTIONAL ITEMS:**

Part Number	Description	Price AUD
<input type="checkbox"/> I103070-99-0	UT Amp. Support Pad (H) w/ Mesh Cover - (1 Each) <i>7"W x 7"T Curved Pad. Must also select PN I10553.</i>	
<input type="checkbox"/> I103071-99-0	UT Amp. Support Pad (D) w/ Mesh Cover - (1 Each) <i>8"W x 8"T Flat Pad. Must also select PN I10553.</i>	
<input type="checkbox"/> I10553	Adj Removable - Amp Support Hardware (1 Each) <i>Please indicate if you would like the amputation support mounted on the right or left side of the seat frame: Mount Amp Suppt on Left: _____ Mount Amp Suppt on Right: _____</i>	
<input type="checkbox"/> I109040-99-0	BodiLink Premium LPTS Pad (B), Stretch-Air - Pair <i>4.5"W x 3.5"T Pad w/additional foam insert and Stretch-Air cover. Must also select BodiLink LPTS Mounting Hardware.</i>	
<input type="checkbox"/> I109043-99-0	BodiLink Premium LPTS Pad (A), Stretch-Air - Pair <i>6"W x 5.5"T Pad w/additional foam insert and Stretch-Air cover. Must also select BodiLink LPTS Mounting Hardware.</i>	
<input type="checkbox"/> I109041-99-0	BodiLink Premium LPTS Pad (C), Stretch-Air - Pair <i>8"W x 3.5"T Pad w/additional foam insert and Stretch-Air cover. Must also select BodiLink LPTS Mounting Hardware.</i>	
<input type="checkbox"/> I109042-99-0	BodiLink Premium LPTS Pad (G), Stretch-Air - Pair <i>12"W x 3.5"T Pad w/additional foam insert and Stretch-Air cover. Must also select BodiLink LPTS Mounting Hardware.</i>	
<input type="checkbox"/> I109044-99-0	Adj/Removable BodiLink LTPS Hdwr Med Corpus Pair <i>Medium TT Power Mount hardware can be extended upward to clear cushions of 3-5.5" thick.</i>	
<input type="checkbox"/> I109194-99-0	Adj/Removable Lateral Thigh/Hip Hardware - Pair	
<input type="checkbox"/> I106762-99-0	Abductor Kit Corpus - Leatherette <i>Must also select PN I10555.</i>	
<input type="checkbox"/> I10555	Adj Removable - Abductor Hardware (1 Each)	
<input type="checkbox"/> I105260-99-0	Transfer Handles, Corpus Seat	
<input type="checkbox"/> I105292-99-0	Push Handles, Corpus Seat	
<input type="checkbox"/> I106210-99-0	RAM Self-Leveling Cup Holder	
<input type="checkbox"/> I106806-99-0	RAM X-Grip Phone Holder	
<input type="checkbox"/> I107100-99-0	RAM X-Grip 10" Tablet Holder w/Swing Away Mount <i>Compatible with most 10" screen tablets with or without a protective case. Includes adjustable swing away mount (I103476-99-0). Note: Customization may be required if ordered along with an Omni.</i>	
<input type="checkbox"/> I107222-99-0	USB Charger <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets. The charger features an isolated circuit design for improved safety; and by pulling power from the ICS, prevents battery drainage when the chair is powered off. Includes mounting hardware and ICS connector. Keep in mind charger may need to be relocated with some configurations. Mount on left armrest: _____ Mount on right armrest: _____</i>	





## Accessories

### OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I105805-99-0	Medical Necessities Bag Hooks <i>Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with 23" or 24" tall Ergo backrest. Custom required if ordered with a ROHO AGILITY backrest. May ship as a part at shorter backrest heights.</i>	
<input type="checkbox"/> I102059-99-0	Medical Necessities Bag <i>Includes hardware to attach bag to backrest. Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with a 23" or 24" tall Ergo backrest. Custom required if ordered with a ROHO AGILITY backrest.</i>	
<input type="checkbox"/> I105262-99-0	Joystick Cutout Tray - 16" W <i>Must also select PN I10556.</i>	
<input type="checkbox"/> I105263-99-0	Joystick Cutout Tray - 18" W <i>Must also select PN I10556.</i>	
<input type="checkbox"/> I10843	Small Angle Adjustable Tray <i>Must also select PN I10556.</i>	
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each) <i>Mount on Left: _____ Mount on Right: _____</i>	
<input type="checkbox"/> I108230-99-0	Permobil Adjustable Tray - 13"W	
<input type="checkbox"/> I108228-99-0	Permobil Adjustable Tray - 15"W	
<input type="checkbox"/> I108227-99-0	Permobil Adjustable Tray - 17"W	
<input type="checkbox"/> I108229-99-0	Permobil Adjustable Tray - 19"W	
<input type="checkbox"/> I108232-99-0	Permobil Tray Hardware Right Side Adjustable	
<input type="checkbox"/> I108233-99-0	Permobil Tray Hardware Left Side Adjustable	
<input type="checkbox"/> I105224-99-0	Oxygen Holder, Corpus	
<input type="checkbox"/> I105223-99-0	Crutch Holder, Corpus	

## Alternative Drive Systems

### ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10751	Mini-Joystick for Omni - (9-pin connector) <i>Operating Force = 13 g Neutral to Full Deflection = 5mm</i>	
<input type="checkbox"/> I103893-99-0	Compact Joystick for Omni - (9-pin connector) <i>Operating Force = 220 g Neutral to Full Deflection = 28mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	
<input type="checkbox"/> I103894-99-0	Compact Joystick LITE for Omni - (9-pin connector) <i>Operating Force = 50 g Neutral to Full Deflection = 17mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	

### ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description
<input type="checkbox"/> I107946-99-0	Swing Away Chin Boom Assm for Compact Lite - Left
<input type="checkbox"/> I107948-99-0	Swing Away Chin Boom Assm for Compact Lite - Right
<input type="checkbox"/> I107857-99-0	Swing Away Chin Boom Assm for Mini Joystick - Left
<input type="checkbox"/> I107947-99-0	Swing Away Chin Boom Assm for Mini Joystick - Right
<input type="checkbox"/> I104690-99-0	Alt. JS Swing Away Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>



**Alternative Drive Systems**

**ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:**

<b>Part Number</b>	<b>Description</b>	<b>Price AUD</b>
<input type="checkbox"/> I107135-99-0	Alt JS/Omni Combo Swing Away Mount Kit - Armrest <i>Used to mount Omni and Alternative Proportional Joystick to the same armrest. Includes Swing Away Omni/Tablet Mount, with a longer mounting rod, two PERMOfix modules, one 3" rod and one 4" rod.</i> Mount in front of Left Arm: _____ Mount in front of Right Arm: _____	
<input type="checkbox"/> I104143-99-0	Alt. JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	
<input type="checkbox"/> I10752	Haims Harness for Chin Control Mount, Large <i>For neck width up to 19cm.</i>	
<input type="checkbox"/> I10753	Haims Harness for Chin Control Mount, Small <i>For neck width up to 14cm.</i>	

**ALT PROPORTIONAL JOYSTICK ACCESSORIES:**

<b>Part Number</b>	<b>Description</b>	<b>Price AUD</b>
<input type="checkbox"/> I104904-99-0	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid arms. Switches set-up for Mode &amp; On/Off function. Switch operation force = 130 g</i> Mount at front of Joystick: _____ Mount at rear of Joystick: _____ Mount at left of Joystick: _____ Mount at right of Joystick: _____	
<input type="checkbox"/> I104641-99-0	Dual Switches for Compact Joystick - Gooseneck Mnt <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches set-up for Mode &amp; On/Off function. Switch operation force = 130 g</i>	

**SIP & PUFF SYSTEMS:**

<b>Part Number</b>	<b>Description</b>	<b>Price AUD</b>
<input type="checkbox"/> I10573	Sip & Puff System for Omni <i>Requires Omni Control Module. Includes breath tube kit.</i>	
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Includes breath tube kit.</i>	
<input type="checkbox"/> I106807-99-0	Sip & Puff Kit for Total Control Head Array <i>For use with Total Control Head Array. Requires Omni Control Module. Includes breath tube kit.</i>	
<input type="checkbox"/> I10762	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573, I10574 or I106807-99-0.</i>	

**PROPORTIONAL FOOT CONTROL:**

<b>Part Number</b>	<b>Description</b>	<b>Price AUD</b>
<input type="checkbox"/> I10767	Proportional Foot Control <i>Not available with Active Reach Packages.</i> Mount on Left: _____ Mount on Right: _____	
<input type="checkbox"/> I10766	Adjustable Mounting Hrdw - Foot Control <i>Required when selecting PN I10767.</i>	

**Total Control Head Array**

**CONFIGURE HEAD ARRAY:**

<b>Part Number</b>	<b>Description</b>	<b>Price AUD</b>
<input type="checkbox"/> I106805-99-0	Permobil Head Array - 3-Switch System <i>Features an occipital pad sensor and two stainless steel arms configured with proximity sensors. The arms are mounted to the upper attachment points on the head array. Longer arms are also included to extend the sensor adjustment range. ESP recommended. Choose Omni Control Module separately.</i>	



**Total Control Head Array**

**CONFIGURE HEAD ARRAY:**

Part Number	Description	Price AUD
<input type="checkbox"/> I11072	Program Head Array as a 5-Switch System <i>Select additional arms/sensors below and indicate the intended use of each sensor in the Notes &amp; Comments section. Choose Omni Control Module separately.</i>	No Charge
CHOOSE OCCIPITAL PAD SIZE (REQUIRED):		
<input type="checkbox"/> I105920-99-0	Small Occipital Pad for Head Array - 4.25"W x 3"H	No Charge
<input type="checkbox"/> I105921-99-0	Large Occipital Pad for Head Array - 6"W x 4"H	No Charge
CHOOSE OCCIPITAL PAD MOUNT (REQUIRED):		
<input type="checkbox"/> I10847	Fixed Mount for Occipital Pad - Head Array	No Charge
<input type="checkbox"/> I105936-99-0	Adj. Swivel Mount for Head Array Occipital Pad <i>Enables Occipital Pad to be angled up/down and left/right independent of the lateral arms.</i>	No Charge
CHOOSE HEAD ARRAY MOUNT HARDWARE (REQUIRED):		
<input type="checkbox"/> I105927-99-0	Stealth TWB480 Backrest Mounting Hdwr with Links <i>Please select the length of the vertical stem below: Standard length (10"): _____ Short length (6"): _____</i>	
<input type="checkbox"/> I105941-99-0	Stealth TWB480 Backrest Mtng Hdwr with Flipback	
<input type="checkbox"/> I105925-99-0	Stealth HMO475 Backr Mtng Hdwr w 17mm Ball	
SELECT ADDITIONAL LEFT ARM PACKAGE (OPTIONAL):		
<input type="checkbox"/> I10976	Add Left Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	
<input type="checkbox"/> I10977	Add Left Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	
<input type="checkbox"/> I10978	Add Left Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	
SELECT ADDITIONAL RIGHT ARM PACKAGE (OPTIONAL):		
<input type="checkbox"/> I10979	Add Right Lateral Arm w Proximity Sensor <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	
<input type="checkbox"/> I10980	Add Right Lateral Arm w Egg Switch <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	
<input type="checkbox"/> I10981	Add Right Lateral Arm w Pad (No Sensor or Switch) <i>Choose Lateral Arm Adjustment Range: Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>	

**Alt Drive Systems Setup**

**CHOOSE OMNI (REQUIRED):**

Part Number	Description	Price AUD
<input type="checkbox"/> I108500-99-0	Omni2 Control Module for R-net - Corpus <i>Two-section design features 3.5" VGA high resolution display, USB charger port, built in IR control and integrated Bluetooth technology for control of up to 4 Bluetooth devices. Must also select mounting hardware below. Please specify which type of alternative drive input device you will be using. Input Device: _____</i>	Price AUD
<input type="checkbox"/> I103476-99-0	Adjustable Swing-Away Hardware - Omni/Tablet <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	Price AUD

**CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:**

Part Number	Description	Price AUD
-------------	-------------	-----------



**Alt Drive Systems Setup**

**CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:**

Part Number	Description	Price AUD
<input type="checkbox"/> I10080	Mount Seat Control Switchbox on Backrest	No Charge

**Alt Drive Accessories**

**PERMOFIX MOUNTING SOLUTIONS:**

Part Number	Description	Price AUD
<input type="checkbox"/> I106655-99-0	PERMOfix UniTrack Fixed Mount Kit <i>Includes two base joints, two 4" rods and one 8" rod.</i>	
<input type="checkbox"/> I106656-99-0	PERMOfix UniTrack Removable Mount Kit <i>Includes quick release mount and thigh support hardware, two base joints and three different length rods.</i>	
<input type="checkbox"/> I106657-99-0	PERMOfix UniTrack Swing Away Mount Kit <i>Includes PERMOfix Swing Away module, two base joints, one 2" rod and one 4" rod.</i>	
<input type="checkbox"/> I105907-99-0	Proximity Sensor + Swivel Assm - Not Mounted <i>Only available with the Total Control Head Array. Swivel assembly allows sensor to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array. Since Proximity Sensors require power, they CANNOT be used as On/Off switches.</i>	
<input type="checkbox"/> I105919-99-0	Egg Switch + Swivel Mount Assm - Not Mounted <i>Swivel assembly allows switch to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	
<input type="checkbox"/> I105912-99-0	Solid Lateral Pad+Swivel Joint - Not Mounted <i>Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	
<input type="checkbox"/> I106675-99-0	Oval Plate + Swivel Assm - No Switch - Not Mounted <i>Blank plate provides mounting option for Micro Lite, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	

**OPTIONAL ITEMS:**

Part Number	Description	Price AUD
<input type="checkbox"/> I101769-99-0	Microlight Switch Qty: _____	
<input type="checkbox"/> I101767-99-0	Buddy Button Qty: _____	
<input type="checkbox"/> I106293-99-0	Egg Switch Black Qty: _____	
<input type="checkbox"/> I101771-99-0	Mini Cup Switch Qty: _____	
<input type="checkbox"/> I107500-99-0	Proximity Switch Qty: _____	
<input type="checkbox"/> I107495-99-0	Mini Joystick 5 Switch OC for Omni - (9-pin)	
<input type="checkbox"/> I10772	Dome Switch, Backrest Mnt w Gooseneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____</i> <i>Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights. Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 150mm rod, one 200mm rod and one 225mm gooseneck switch with 1.2m cable and mounting clamp.</i>	
<input type="checkbox"/> I10846	Dome Switch, Backrest Swg Awy Mnt w Gsneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____</i> <i>w 90° Swing-Away Function. Mounts at top of back. Custom charge applies if ordered with short back. Kit Includes: One PERMOfix joint with Flat Plate, one PERMOfix joint with 45° swivel adjustment, two 6" rods, one 90° Swing-Away joint, one 9"</i>	
<input type="checkbox"/> I105935-99-0	3.5mm Ext. Cable for Sensor/Switch, 0.3m Long	
<input type="checkbox"/> I105931-99-0	3.5mm Ext. Cable for Sensor/Switch, 1m Long	
<input type="checkbox"/> I105932-99-0	3.5mm Ext. Cable for Sensor/Switch, 2m Long	
<input type="checkbox"/> I105933-99-0	3.5mm Ext. Cable for Sensor/Switch, 3m Long	
<input type="checkbox"/> I105934-99-0	3.5mm Y-splitter Cable for Sensor/Switch <i>Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.</i>	



Additional Order Instructions:

**Chair Order Policy:**

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

**Order Acknowledgement:**

I, \_\_\_\_\_, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: \_\_\_\_\_

Printed Name: \_\_\_\_\_

Title: \_\_\_\_\_

Date: \_\_\_\_\_

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty.