

# K450 MX



*\* Indicates a Required Field*

## Dealer Information

\* Contact: \_\_\_\_\_

Dealer Code: \_\_\_\_\_

\* Dealer Name: \_\_\_\_\_

Address: \_\_\_\_\_

\* City: \_\_\_\_\_

\* State: \_\_\_\_\_

\* Phone#: \_\_\_\_\_

\* Fax#: \_\_\_\_\_

PO#: \_\_\_\_\_

Email Address: \_\_\_\_\_

## Client Information

*Permobil recommends that the client is evaluated by a certified rehab specialist.*

\* First Name: \_\_\_\_\_

\* Last Name: \_\_\_\_\_

Diagnosis: \_\_\_\_\_

Funding Source: \_\_\_\_\_

Client Age: \_\_\_\_\_

## Client Measurements

\* Weight: \_\_\_\_\_

\* Height: \_\_\_\_\_

A) Top of Shoulders: \_\_\_\_\_

B) Chest Depth: \_\_\_\_\_

C) Chest Width: \_\_\_\_\_

D) Seat Depth: \_\_\_\_\_


E) Top of Head: \_\_\_\_\_

F) Elbow to Hand: \_\_\_\_\_

G) Seat Pan to Elbow: \_\_\_\_\_

H) Hip Width: \_\_\_\_\_

\*I) Knee to Foot: \_\_\_\_\_





## Chair Model and Colors

|                                     |  | Price AUD |
|-------------------------------------|--|-----------|
| <input checked="" type="checkbox"/> | I104524-99-0 K450 Base w/ R-net 8kmh<br><i>Price Includes: K450 Power Base, Tie-Down Hardware for Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. Max. Speed = 8 km/h, Weight Capacity = 60 kg, STFH = 75 mm - 635 mm.</i> |           |

### CHOOSE DRIVE TIRES:

|                          | Part Number  | Description                  | Price AUD |
|--------------------------|--------------|------------------------------|-----------|
| <input type="checkbox"/> | I106204-99-0 | Flat Free Casters - K450     | No Charge |
| <input type="checkbox"/> | I106205-99-0 | Flat Free Drive Tires - K450 | No Charge |
| <input type="checkbox"/> | I106206-99-0 | Pneumatic Casters - K450     | No Charge |
| <input type="checkbox"/> | I106207-99-0 | Pneumatic Drive Tires - K450 | No Charge |

### CHOOSE SHROUD COLOR:

|                          | Part Number  | Description                       | Price AUD |
|--------------------------|--------------|-----------------------------------|-----------|
| <input type="checkbox"/> | I104535-99-0 | Shroud for K450 - Onyx Black      | No Charge |
| <input type="checkbox"/> | I104536-99-0 | Shroud for K450 - Platinum Silver | No Charge |
| <input type="checkbox"/> | I104537-99-0 | Shroud for K450 - Pop Star Pink   | No Charge |
| <input type="checkbox"/> | I104539-99-0 | Shroud for K450 - Extreme Yellow  | No Charge |
| <input type="checkbox"/> | I106367-99-0 | Shroud for K450 - Volcano Red     | No Charge |
| <input type="checkbox"/> | I106368-99-0 | Shroud for K450 - Cobalt Blue     | No Charge |
| <input type="checkbox"/> | I107234-99-0 | Shroud for K450 - Purple Rain     | No Charge |

### CHOOSE ONE OF THE FOLLOWING:

|                                     | Part Number | Description                       | Price AUD |
|-------------------------------------|-------------|-----------------------------------|-----------|
| <input checked="" type="checkbox"/> | I10901      | Australian Freight & Handling Fee |           |

## Base Options

### MUST CHOOSE SEAT ELEVATOR:

|                                     | Part Number | Description  | Price AUD |
|-------------------------------------|-------------|--|-----------|
| <input checked="" type="checkbox"/> | I10510      | Power Adjustable Seat Height - 8" Travel<br><i>This item must be selected.</i> |           |

### MUST CHOOSE BATTERIES:

|                                     | Part Number  | Description  | Price AUD |
|-------------------------------------|--------------|--|-----------|
| <input checked="" type="checkbox"/> | I107032-99-0 | Batteries, Grp 34 Gel, Pair<br><i>Includes Permobil Volt charger for non-Permobil gel batteries.</i> |           |

### OPTIONAL ITEMS:

|                          | Part Number  | Description   | Price AUD |
|--------------------------|--------------|---|-----------|
| <input type="checkbox"/> | I104550-99-0 | Support Wheels for K450   |           |
| <input type="checkbox"/> | I104540-99-0 | Light Kit - K450 R-net  |           |
| <input type="checkbox"/> | I104553-99-0 | WC 19 Transit Kit - K450<br><i>This kit includes additional anchor points for tie down straps, seat stabilizing bars and hardware, and a WC-19 Approved integrated lap and shoulder belt. This kit has been approved according to ANSI/RESNA WC19. Cannot be ordered with a 10"W Seat Size (I104317-99-0). An ICS Switchbox must be installed on the chair in order to access WC19 mode.</i>  |           |
| <input type="checkbox"/> | I104554-99-0 | Fixed Ventilator Tray - K450<br><i>Tray is 9.75"D x 13"W and increases the overall length of a K450 by approximately 11". The ventilator tray is designed to accommodate an LP10 style ventilator and one U1 Battery. A custom charge will need to be applied if ordering vent tray with any of the following items: WC-19 Transit Kit (I104553-99-0), Attendant Control (ID50882), Push Handles (I104576-99-0) or Bag Hooks (I105805-99-0). Please Note: Must also select Support Wheels (I104550-99-0).</i> |           |
| <input type="checkbox"/> | ID51113      | Programming Key for R-net - Permobil Version  |           |



## Joystick Options

**PLEASE NOTE:**

*For additional input devices, refer to Alternative Drive Systems section.*

I108319-99-0

Permobil Joystick Module

*Designed exclusively for Permobil, this R-net color joystick features: control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.*

I105246-99-0

R-net LED Remote Joystick - No Display

No Charge

*This Joystick does NOT have a Display Screen, does NOT contain mono jack ports, and does NOT include Bluetooth or Assignable Buttons. PC Programming required and profiles limited to five.*

**CHOOSE JOYSTICK HANDLE TYPE:**

| Part Number                           | Description                                    | Price AUD |
|---------------------------------------|--|-----------|
| <input type="checkbox"/> I10045       | Cone Shaped Joystick Knob - Standard           | No Charge |
| <input type="checkbox"/> I100960-99-0 | Chin Cup For Joystick Knob                     |           |
| <input type="checkbox"/> I100961-99-0 | "T" Handle For Joystick Knob                   |           |
| <input type="checkbox"/> I100962-99-0 | "Large Ball" For Joystick Knob                 |           |
| <input type="checkbox"/> I100963-99-0 | "Softball" For Joystick Knob                   |           |
| <input type="checkbox"/> I100964-99-0 | "Mushroom" For Joystick Knob                   |           |
| <input type="checkbox"/> I100965-99-0 | "Stick" For Joystick Knob                      |           |
| <input type="checkbox"/> I105442-99-0 | 1.75" Black Foam Ball For Joystick Knob        |           |
| <input type="checkbox"/> I105846-99-0 | Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft |           |
| <input type="checkbox"/> I105847-99-0 | Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft |           |
| <input type="checkbox"/> I105848-99-0 | Bodypoint J/S Handle Dome Shaped, Rubber       |           |

**CHOOSE JOYSTICK MOUNT:**

| Part Number                           | Description   | Price AUD |
|---------------------------------------|---|-----------|
| <input type="checkbox"/> I10412       | Joystick Mount Right - Fixed, R-net Only  | No Charge |
| <input type="checkbox"/> I10413       | Joystick Mount Left - Fixed, R-net Only   | No Charge |
| <input type="checkbox"/> I102537-99-0 | Joystick Mount Right - Swing Away   |           |
| <input type="checkbox"/> I102538-99-0 | Joystick Mount Left - Swing Away  |           |
| <input type="checkbox"/> ID50882-A    | Attendant Control for R-net Control System Only<br><i>On the K450, the mounting location for the Attendant Control is off of the Push Handle Bracket. If ordering an attendant control, you must also order I104576-99-0 (Push Handles - MX Seat).</i>  |           |
| <input type="checkbox"/> I107350-99-0 | R-net Bluetooth iDevice Module<br><i>Allows input device commands (such as joystick nudges or switch clicks) to be used to operate an iPhone, iPad or iPod. By assigning each input to an iOS function, users can make phone calls, send texts and emails and access applications. Kit includes a Bluetooth iDevice Module, an R-net Bus Cable, and two mono jack cables for direct switch inputs (switches not included). Must be ordered with a Color JSM or Omni for mode display.</i> |           |
| <input type="checkbox"/> I107367-99-0 | R-net Bluetooth Mouse Module<br><i>Enables mouse control through any R-net input device. It can operate the mouse on a computer and on select Android devices (with R-net for Android app). Kit includes a Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click (switches not included) and a USB BT Adapter for user's computer. Note: BT Adapter may not work with all computers. Must be ordered with a Color JSM or Omni for mode display.</i>              |           |
| <input type="checkbox"/> I107368-99-0 | Add a second Bluetooth Mouse Module<br><i>Allows R-net to simultaneously pair to a second compatible computer or Android device. Each module is accessible through a separate R-net mode. Kit includes a unique Bluetooth Mouse Module, an R-net Bus Cable, two mono jack cables for L/R click and a USB BT Adapter.</i>  |           |

## Seating System

**CHOOSE ONE OF THE FOLLOWING:**

| Part Number | Description | Price AUD |
|-------------|-------------|-----------|
|-------------|-------------|-----------|



## Seating System

### CHOOSE ONE OF THE FOLLOWING:

| Part Number                                      | Description  | Price AUD |
|--|--|-----------|
| <input checked="" type="checkbox"/> I104324-99-0 | MX Power Seating System<br><i>Price Includes: MX Seating System, Adjustable Seat Frame, Power Seat to Floor System, Adjustable Backrest Angle, and One Piece Footplate. Weight Limit 57kg.</i> |           |

### CHOOSE ONE OF THE FOLLOWING:

| Part Number                           | Description                                 | Price AUD |
|---------------------------------------|---|-----------|
| <input type="checkbox"/> I104529-99-0 | No Tilt<br><i>Weight Limit 57kg.</i>        | No Charge |
| <input type="checkbox"/> I104551-99-0 | 45° Power Tilt<br><i>Weight Limit 57kg.</i> |           |

## Seating System Options

### CHOOSE SEAT FUNCTION CONTROL TYPE:

| Part Number                                | Description   | Price AUD |
|--|---|-----------|
| <input checked="" type="checkbox"/> I10585 | Seat Function Control Kit for R-net<br><i>Please indicate below if you would like push buttons or toggles on the ICS Switchbox.</i> |           |

### CHOOSE PUSH BUTTONS OR TOGGLES:

| Part Number                           | Description                               | Price AUD |
|---------------------------------------|---|-----------|
| <input type="checkbox"/> I104135-99-0 | Push Buttons on Seat Function Control Kit | No Charge |
| <input type="checkbox"/> I104136-99-0 | Toggles on Seat Function Control Kit      |           |

### CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

| Part Number                           | Description  | Price AUD |
|---------------------------------------|--|-----------|
| <input type="checkbox"/> I10626       | Mount ICS Switchbox onto R-net Joystick Bracket<br><i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>                            | No Charge |
| <input type="checkbox"/> I104651-99-0 | Mount ICS Switchbox onto Armrest<br><i>Mount on Left: _____ Mount on Right: _____</i>  | No Charge |
| <input type="checkbox"/> I11050       | Mount ICS Switchbox on Backrest<br><i>On the K450, the backrest mounting location for the ICS Switchbox is on the Push Handle bracket; therefore, you must also order I104576-99-0 (Push Handles - MX).<br/>Mount on Left: _____ Mount on Right: _____</i> |           |

### OPTIONAL ITEMS:

| Part Number                           | Description  | Price AUD |
|---------------------------------------|--|-----------|
| <input type="checkbox"/> I105409-99-0 | Do Not Install ICS Switchbox, Send as Part   | No Charge |
| <input type="checkbox"/> I104196-99-0 | ICS Alternative Switchbox for R-net Electronics<br><i>Seat function interface box that accommodates up to 8 momentary switches. Provides an additional method for accessing the seat functions. 1/8" mini-mono switches not included.</i>  |           |
| <input type="checkbox"/> I108526-99-0 | Remote Stop System - AUS/NZ/EU<br><i>Consists of a 433MHz Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair.</i>  |           |
| <input type="checkbox"/> I10515       | Memory Seat Program for K450 MX<br><i>If this is selected, the chair will be programmed to allow the user to store and save up to 3 different memory seat positions through the ICS Switchbox. On the K450 MX, ALL seat functions are capable of recalling a saved seat position. Please refer to the K450 Owner's Manual for more information on saving and recalling a memory seat position.</i> | No Charge |

## Backrest Options

### CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price AUD |
|-------------|-------------|-----------|
|-------------|-------------|-----------|



## Backrest Options

### CHOOSE ONE OF THE FOLLOWING:

| Part Number  | Description   |           |
|--|---|-----------|
| <b>MX STYLE BACKRESTS:</b>   |   |           |
| <i>Use this section to select MX Style Backrests. Back heights measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Please be sure to pick both a backrest shell size and a backrest cushion style (Leatherette or Black Mesh Fabric). Note: A 16"W MX Backrest is NOT available.</i> |   |           |
| <input type="checkbox"/> I104321-99-0  | MX Backrest - 10"W x 14"H<br><i>Backrest can be adjusted to either 14" or 16" high. If ordering a headrest, backrest will be mounted in the 16" high position unless client measurements require a 14" backrest height, in which case, custom modifications will be done.</i> | No Charge |
| <input type="checkbox"/> I104322-99-0  | MX Backrest - 12"W x 16"H<br><i>Backrest can be adjusted to either 16" or 18" high.</i>   | No Charge |
| <input type="checkbox"/> I104323-99-0  | MX Backrest - 14"W x 20"H<br><i>Backrest can be adjusted to either 20" or 22" high.</i>   | No Charge |

### CHOOSE ONE OF THE FOLLOWING:

| Part Number                               | Description                                       |           |
|---|---|-----------|
| <b>Please Select Cushion Style Below:</b> |   |           |
| <input type="checkbox"/> I104309-99-0     | MX/MF Backrest Cushion, 10"W x 14"T - Leatherette | No Charge |
| <input type="checkbox"/> I104311-99-0     | MX/MF Backrest Cushion, 12"W x 16"T - Leatherette | No Charge |
| <input type="checkbox"/> I104313-99-0     | MX/MF Backrest Cushion, 14"W x 20"T - Leatherette | No Charge |
| <input type="checkbox"/> I104308-99-0     | MX/MF Backrest Cushion, 10"W x 14"T - Black Mesh  | No Charge |
| <input type="checkbox"/> I104310-99-0     | MX/MF Backrest Cushion, 12"W x 16"T - Black Mesh  | No Charge |
| <input type="checkbox"/> I104312-99-0     | MX/MF Backrest Cushion, 14"W x 20"T - Black Mesh  | No Charge |

### CHOOSE ONE OF THE FOLLOWING:

| Part Number   | Description   |           |
|---|---|-----------|
| <b>PS JR TUBULAR STYLE BACKRESTS</b>  |   |           |
| <i>Use this section to select PS Jr Tubular Style Backrest Canes. Note: A 10"W Tubular Backrest is NOT available.</i> |   |           |
| <input type="checkbox"/> I103717-99-0   | PS Jr Back Canes, Height & Angle Adjustable<br><i>Must also select PS Jr Style Armrests if ordering this backrest. Choose PS Jr Style Backrest Type (Metal Back or Omit Back) below. Backrest width should match seat width. Back height is adjustable from 16.5"T - 22.5"T, measured from the seat pan to the top of the back.</i> | No Charge |

### CHOOSE ONE OF THE FOLLOWING:

| Part Number                           | Description                 | Price AUD |
|---------------------------------------|-----------------------------|-----------|
| <input type="checkbox"/> I103719-99-0 | PS Jr Solid Backrest - 12"W | No Charge |
| <input type="checkbox"/> I103721-99-0 | PS Jr Solid Backrest - 14"W | No Charge |
| <input type="checkbox"/> I103723-99-0 | PS Jr Solid Backrest - 16"W | No Charge |
| <input type="checkbox"/> I10449       | PS Jr - Omit Backrest       | No Charge |

## Seat Sizing

### CHOOSE ONE OF THE FOLLOWING:

| Part Number                           | Description   | Price AUD |
|---------------------------------------|---|-----------|
| <input type="checkbox"/> I104317-99-0 | MX Seat Size 10"W x 10"D<br><i>Seat depth can be adjusted to 10", 11", or 12". If ordering this seat size, the lowest allowed seat position when using the seat to floor function will be approximately 3.5" from the floor to the top of the seat pan.</i>   | No Charge |
| <input type="checkbox"/> I104318-99-0 | MX Seat Size 12"W x 12"D<br><i>Seat depth can be adjusted to 12", 13", or 14". If ordering this seat size with MX Arms, the lowest allowed seat position when using the seat to floor function will be approximately 3.5" from the floor to the top of the seat pan. If ordering Flip Back Arms, the lowest allowed seat position will be 8".</i> | No Charge |





## Seat Sizing

### CHOOSE ONE OF THE FOLLOWING:

| Part Number                           | Description   | Price AUD |
|---------------------------------------|---|-----------|
| <input type="checkbox"/> I104319-99-0 | MX Seat Size 14"W x 14"D<br><i>Seat depth can be adjusted to 14", 15", or 16". If ordering this seat size with MX Arms, the lowest allowed seat position when using the seat to floor function will be approximately 3.5" from the floor to the top of the seat pan. If ordering Flip Back Arms, the lowest allowed seat position will be 8".</i> | No Charge |
| <input type="checkbox"/> I104320-99-0 | MX Seat Size 16"W x 16"D<br><i>Seat depth can be adjusted to 16", 17", or 18". If ordering this seat size with MX Arms, the lowest allowed seat position when using the seat to floor function will be approximately 10" from the floor to the top of the seat pan. If ordering Flip Back Arms, the lowest allowed seat position will be 16".</i> | No Charge |

### OPTIONAL ITEMS:

| Part Number                              | Description   |
|--|---|
| <input type="checkbox"/> I306384-99-90-0 | Seat Cushions Below Are Approximately 2" Thick:<br>Leatherette Seat Cush MX - 10"W x 10"D (250x250) |
| <input type="checkbox"/> I306385-99-90-0 | Leatherette Seat Cush MX - 12"W x 12"D (300x300)  |
| <input type="checkbox"/> I306386-99-90-0 | Leatherette Seat Cush MX - 14"W x 14"D (350x350)  |
| <input type="checkbox"/> I104703-99-0    | Leatherette Seat Cush MX - 16"W x 16"D (400x400)  |
| <input type="checkbox"/> I104305-99-0    | Black Mesh Seat Cushion MX/MF - 10"W x 10"D   |
| <input type="checkbox"/> I104306-99-0    | Black Mesh Seat Cushion MX/MF - 12"W x 12"D   |
| <input type="checkbox"/> I104307-99-0    | Black Mesh Seat Cushion MX/MF - 14"W x 14"D   |
| <input type="checkbox"/> I104702-99-0    | Black Mesh Seat Cush MX - 16"W x 16"D (400x400)   |

## Armrest Options

### CHOOSE ARMREST STYLE:

| Part Number                     | Description   | Price AUD |
|---------------------------------|---|-----------|
| <input type="checkbox"/> I10512 | Height Adjustable Armrest Assembly MX Seat<br><i>If you are not selecting Tilt, you MUST order this item.</i>   |           |
| <input type="checkbox"/> I10514 | Height Adjustable Armrest Assembly MX Seat - NC<br><i>No Charge if ordering Power Tilt (I104551-99-0)</i>   | No Charge |
| <input type="checkbox"/> I10824 | Flip Back, Height Adjustable Armrests<br><i>Mounted off the backrest canes. Includes: 45° of angle adjustability, rotational adjustment, 3" offset adjustment both medially and laterally, and height adjustment of 6" - 14" from the seat pan. In one inch increments, please indicate where you would like the armrests set.</i><br>SET ARMRESTS AT (6" - 14"): _____ |           |
| <input type="checkbox"/> I10825 | Flip Back, Height Adjustable Armrests - NC<br><i>No charge when ordering power tilt. In one inch increments, please indicate where you would like the armrests set.</i><br>ARMRESTS AT (6" - 14"): _____  | No Charge |

### MX ARMREST PADS:

| Part Number                           | Description                                  | Price AUD |
|---------------------------------------|--|-----------|
| <input type="checkbox"/> I104303-99-0 | MX Armrest Pad Short Leatherette 8"L x 3.5"W | No Charge |
| <input type="checkbox"/> I104304-99-0 | MX Armrest Pad Long Leatherette 12.5"Lx3.5"W | No Charge |

### FLIP BACK ARMREST PAD RIGHT:

| Part Number                           | Description                                    | Price AUD |
|---------------------------------------|--|-----------|
| <input type="checkbox"/> I105100-99-0 | 4"x10" Corpus Armrest Pad Right - Leatherette  | No Charge |
| <input type="checkbox"/> I105104-99-0 | 4"x13" Corpus Armrest Pad Right - Leatherette  | No Charge |
| <input type="checkbox"/> I107212-99-0 | 4"x10" Corpus Armrest Pad Right - Black Fabric | No Charge |
| <input type="checkbox"/> I106047-99-0 | 4"x13" Corpus Armrest Pad Right - Black Fabric | No Charge |
| <input type="checkbox"/> I104951-99-0 | 3.5"x8" Unitrack Armrest Pad Right (C)         | No Charge |
| <input type="checkbox"/> I104950-99-0 | 3.5"x12.5" Unitrack Armrest Pad Right (G)      | No Charge |



## Armrest Options

### FLIP BACK ARMREST PAD RIGHT:

| Part Number                           | Description                                   | Price AUD |
|---------------------------------------|---|-----------|
| <input type="checkbox"/> I107591-99-0 | 4"x10" Corpus Arm w/ Permobil Gel Pad - Right |           |
| <input type="checkbox"/> I107593-99-0 | 4"x13" Corpus Arm w/ Permobil Gel Pad - Right |           |

### FLIP BACK ARMREST PAD LEFT:

| Part Number                           | Description                                   | Price AUD |
|---------------------------------------|---|-----------|
| <input type="checkbox"/> I105098-99-0 | 4"x10" Corpus Armrest Pad Left - Leatherette  | No Charge |
| <input type="checkbox"/> I105102-99-0 | 4"x13" Corpus Armrest Pad Left - Leatherette  | No Charge |
| <input type="checkbox"/> I107213-99-0 | 4"x10" Corpus Armrest Pad Left - Black Fabric | No Charge |
| <input type="checkbox"/> I106046-99-0 | 4"x13" Corpus Armrest Pad Left - Black Fabric | No Charge |
| <input type="checkbox"/> I104952-99-0 | 3.5"x8" Unitrack Armrest Pad Left (C)         | No Charge |
| <input type="checkbox"/> I104949-99-0 | 3.5"x12.5" Unitrack Armrest Pad Left (G)      | No Charge |
| <input type="checkbox"/> I107590-99-0 | 4"x10" Corpus Arm w/ Permobil Gel Pad - Left  |           |
| <input type="checkbox"/> I107592-99-0 | 4"x13" Corpus Arm w/ Permobil Gel Pad - Left  |           |

## Legrest Options

### CHOOSE ONE OF THE FOLLOWING:

| Part Number                           | Description                            | Price AUD |
|---------------------------------------|--|-----------|
| <input type="checkbox"/> I104314-99-0 | MX Footplate - Small (11.5" W x 6.5"D) | No Charge |
| <input type="checkbox"/> I104315-99-0 | MX Footplate - Large (14" W x 8"D)     | No Charge |

### OPTIONAL ITEMS:

| Part Number                           | Description                                    | Price AUD |
|---------------------------------------|--|-----------|
| <input type="checkbox"/> I103427-99-0 | Footstrap Kit                                  |           |
| <input type="checkbox"/> I103104-99-0 | Footstrap Kit - Padded                         |           |
| <input type="checkbox"/> I104338-99-0 | Small Anti-Skid Overlay for Footplate - K450   |           |
| <input type="checkbox"/> I104337-99-0 | Large Anti-Skid Overlay for Footplate - K450   |           |
| <input type="checkbox"/> I107503-99-0 | BP Ankle Huggers, Med Side Release, Adj (Pair) |           |

*This item will be sent as a part with the wheelchair. It will not be installed.*

## Accessories

### MUST CHOOSE POSITIONING BELT:

| Part Number                                      | Description                             | Price AUD |
|--|---|-----------|
| <input checked="" type="checkbox"/> I104518-99-0 | Bodypoint Non-Padded Hip Belt - Medium  | No Charge |
| <input type="checkbox"/> I105849-99-0            | Bodypoint Padded Hip Belt - Extra Small |           |

### CHOOSE ONE OF THE FOLLOWING:

| Part Number                           | Description   | Price AUD |
|---------------------------------------|---|-----------|
| <input type="checkbox"/> I100749-99-0 | MF Headrest for Koala/PS Jr - Leatherette<br><i>Must also select PN I10552.</i>   |           |
| <input type="checkbox"/> I104339-99-0 | MF Headrest for MX Back - Leatherette<br><i>Must also select PN I10552. Also includes a UT Bracket for lateral adjustment of the headrest. Please indicate below which size you need.</i><br>UT Brack for Lat Head Adj (10"W MX Back): ____<br><i>Allows for Lateral Adjustment of the headrest.</i><br>UT Brack for Lat Head Adj(12-14"W MX Back): ____<br><i>Allows for Lateral Adjustment of the headrest.</i> |           |
| <input type="checkbox"/> I10552       | Adj/Removable Headrest Hardware   |           |



## Accessories

### CHOOSE ONE OF THE FOLLOWING:

| Part Number                           | Description  | Price AUD |
|---------------------------------------|--|-----------|
| <input type="checkbox"/> I107499-99-0 | Universal Headrest Adapter<br><i>Please Note: Recline is limited to 160° or less (depending on backrest height) when ordering a Universal Headrest Adapter on a 23" or 24" tall Corpus backrest.</i> |           |

### OPTIONAL ITEMS:

| Part Number                           | Description   | Price AUD |
|---------------------------------------|---|-----------|
| <input type="checkbox"/> I106104-99-0 | UT Hip Suppt (B) w/ Mesh Cover - High (Pair)<br>4.5"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.  |           |
| <input type="checkbox"/> I106103-99-0 | UT Hip Suppt (A) w/ Mesh Cover - High (Pair)<br>6"W x 5.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.  |           |
| <input type="checkbox"/> I104575-99-0 | UT Trunk Suppt (B) for MX Back - 3.5"T x 4.5"D<br>Must also select the Swing Away Trunk Support Hardware (I10551).  |           |
| <input type="checkbox"/> I104574-99-0 | UT Trunk Suppt (F) for MX Back - 3.5"T x 8.5"D<br>Must also select the Swing Away Trunk Support Hardware (I10551).  |           |
| <input type="checkbox"/> I10550       | Adj Removable - Thigh Support Hardware (Pair)   |           |
| <input type="checkbox"/> I10551       | Adj Removable - Trunk Support Hardware (Pair)   |           |
| <input type="checkbox"/> I106210-99-0 | RAM Self-Leveling Cup Holder  |           |
| <input type="checkbox"/> I106806-99-0 | RAM X-Grip Phone Holder   |           |
| <input type="checkbox"/> I107100-99-0 | RAM X-Grip 10" Tablet Holder w/Swing Away Mount<br><i>Compatible with most 10" screen tablets with or without a protective case. Includes adjustable swing away mount (I103476-99-0). Note: Customization may be required if ordered along with an Omni.</i>  |           |
| <input type="checkbox"/> I107191-99-0 | USB Charger - MX<br><i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets. The charger features an isolated circuit design for improved safety; and by pulling power from the ICS, prevents battery drainage when the chair is powered off. Includes mounting hardware and ICS connector.</i><br>Mount on left armrest: _____ Mount on right armrest: _____ |           |
| <input type="checkbox"/> I104576-99-0 | Push Handles - K450<br><i>Mounted off the chassis. Height and angle adjustable. Must be selected if ordering bag clips or an attendant control or if mounting the ICS Switchbox on the backrest.</i>  |           |
| <input type="checkbox"/> I105805-99-0 | Medical Necessities Bag Hooks<br><i>Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with 23" or 24" tall Ergo backrest. Custom required if ordered with a ROHO AGILITY backrest. May ship as a part at shorter backrest heights.</i>   |           |
| <input type="checkbox"/> I100598-99-0 | Chest Harness<br><i>Chest harness kit includes a standard positioning belt.</i>   |           |
| <input type="checkbox"/> ISH212       | BP Shoulder Harness Medium, Front Pull<br><i>This item will be sent as a part with the wheelchair. It will not be installed.</i>  |           |
| <input type="checkbox"/> I100578-99-0 | Small Angle Adjustable Tray<br><i>Must also select PN I10556.</i>   |           |
| <input type="checkbox"/> I10556       | Adj Removable - Tray Hardware (1 Each)  |           |
| <input type="checkbox"/> I108231-99-0 | Permobil Adjustable Tray - 11"W   |           |
| <input type="checkbox"/> I108230-99-0 | Permobil Adjustable Tray - 13"W   |           |
| <input type="checkbox"/> I108232-99-0 | Permobil Tray Hardware Right Side Adjustable  |           |
| <input type="checkbox"/> I108233-99-0 | Permobil Tray Hardware Left Side Adjustable   |           |

## Alternative Drive Systems

### ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

| Part Number                     | Description  | Price AUD |
|---------------------------------|--|-----------|
| <input type="checkbox"/> I10751 | Mini-Joystick for Omni - (9-pin connector)<br><i>Operating Force = 13 g Neutral to Full Deflection = 5mm</i> |           |





## Alternative Drive Systems

### ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

| Part Number                           | Description  | Price AUD |
|---------------------------------------|--|-----------|
| <input type="checkbox"/> I103893-99-0 | Compact Joystick for Omni - (9-pin connector)<br><i>Operating Force = 220 g Neutral to Full Deflection = 28mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>     |           |
| <input type="checkbox"/> I103894-99-0 | Compact Joystick LITE for Omni - (9-pin connector)<br><i>Operating Force = 50 g Neutral to Full Deflection = 17mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i> |           |

### ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

| Part Number                           | Description   |
|---------------------------------------|---|
|                                       | <b>NOTE:</b><br><i>Please see the PERMOfix Mounting Solutions, within the Alt Drive Accessories section, for additional mounting options. The Swing Away Chin Boom Assemblies are Compatible with Permobil Headrests I107470-99-0, I107469-99-0 or I107468-99-0 and must be selected with I108721-99-0.</i>                                 |
| <input type="checkbox"/> I104690-99-0 | Alt. JS Swing Away Mount Kit - Armrest<br><i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>   |
| <input type="checkbox"/> I107135-99-0 | Alt JS/Omni Combo Swing Away Mount Kit - Armrest<br><i>Used to mount Omni and Alternative Proportional Joystick to the same armrest. Includes Swing Away Omni/Tablet Mount, with a longer mounting rod, two PERMOfix modules, one 3" rod and one 4" rod.</i><br><i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i> |
| <input type="checkbox"/> I104143-99-0 | Alt. JS Adjustable Mount Kit - Chin Boom<br><i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>  |
| <input type="checkbox"/> I10752       | Haims Harness for Chin Control Mount, Large<br><i>For neck width up to 19cm.</i>  |
| <input type="checkbox"/> I10753       | Haims Harness for Chin Control Mount, Small<br><i>For neck width up to 14cm.</i>  |

### ALT PROPORTIONAL JOYSTICK ACCESSORIES:

| Part Number                           | Description   | Price AUD |
|---------------------------------------|---|-----------|
| <input type="checkbox"/> I104904-99-0 | Dual Switches for Compact Joystick - Fixed Mount<br><i>Includes two pushbutton switches mounted on rigid arms. Switches set-up for Mode &amp; On/Off function. Switch operation force = 130 g</i><br><i>Mount at front of Joystick: _____ Mount at rear of Joystick: _____</i><br><i>Mount at left of Joystick: _____ Mount at right of Joystick: _____</i> |           |
| <input type="checkbox"/> I104641-99-0 | Dual Switches for Compact Joystick - Gooseneck Mnt<br><i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches set-up for Mode &amp; On/Off function. Switch operation force = 130 g</i>   |           |

### SIP & PUFF SYSTEMS:

| Part Number                           | Description  | Price AUD |
|---------------------------------------|--|-----------|
| <input type="checkbox"/> I10573       | Sip & Puff System for Omni<br><i>Requires Omni Control Module. Includes breath tube kit.</i>   |           |
| <input type="checkbox"/> I10574       | Sip & Puff System with Fiber Optic Opti-Stop<br><i>Requires Omni Control Module. Includes breath tube kit.</i>                                       |           |
| <input type="checkbox"/> I106807-99-0 | Sip & Puff Kit for Total Control Head Array<br><i>For use with Total Control Head Array. Requires Omni Control Module. Includes breath tube kit.</i> |           |
| <input type="checkbox"/> I10762       | Adjustable Mounting Hrdw - Sip & Puff<br><i>Required when selecting PN I10573, I10574 or I106807-99-0.</i>   |           |

## Total Control Head Array



## Total Control Head Array

### CONFIGURE HEAD ARRAY:

| Part Number                                     | Description  | Price AUD |
|---|--|-----------|
| <input type="checkbox"/> I106805-99-0           | Permobil Head Array - 3-Switch System<br><i>Features an occipital pad sensor and two stainless steel arms configured with proximity sensors. The arms are mounted to the upper attachment points on the head array. Longer arms are also included to extend the sensor adjustment range. ESP recommended. Choose Omni Control Module separately.</i> |           |
| <input type="checkbox"/> I11072                 | Program Head Array as a 5-Switch System<br><i>Select additional arms/sensors below and indicate the intended use of each sensor in the Notes &amp; Comments section. Choose Omni Control Module separately.</i>  | No Charge |
| CHOOSE OCCIPITAL PAD SIZE (REQUIRED):           |  |           |
| <input type="checkbox"/> I105920-99-0           | Small Occipital Pad for Head Array - 4.25"W x 3"H  | No Charge |
| <input type="checkbox"/> I105921-99-0           | Large Occipital Pad for Head Array - 6"W x 4"H   | No Charge |
| CHOOSE OCCIPITAL PAD MOUNT (REQUIRED):          |  |           |
| <input type="checkbox"/> I10847                 | Fixed Mount for Occipital Pad - Head Array   | No Charge |
| <input type="checkbox"/> I105936-99-0           | Adj. Swivel Mount for Head Array Occipital Pad<br><i>Enables Occipital Pad to be angled up/down and left/right independent of the lateral arms.</i>  | No Charge |
| CHOOSE HEAD ARRAY MOUNT HARDWARE (REQUIRED):    |  |           |
| <input type="checkbox"/> I105927-99-0           | Stealth TWB480 Backrest Mounting Hdwr with Links<br><i>Please select the length of the vertical stem below:<br/>Standard length (10"): _____ Short length (6"): _____</i>  |           |
| <input type="checkbox"/> I105941-99-0           | Stealth TWB480 Backrest Mtng Hdwr with Flipback  |           |
| <input type="checkbox"/> I105925-99-0           | Stealth HMO475 Backr Mtng Hdwr w 17mm Ball   |           |
| SELECT ADDITIONAL LEFT ARM PACKAGE (OPTIONAL):  |  |           |
| <input type="checkbox"/> I10976                 | Add Left Lateral Arm w Proximity Sensor<br><i>Choose Lateral Arm Adjustment Range:<br/>Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>  |           |
| <input type="checkbox"/> I10977                 | Add Left Lateral Arm w Egg Switch<br><i>Choose Lateral Arm Adjustment Range:<br/>Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>  |           |
| <input type="checkbox"/> I10978                 | Add Left Lateral Arm w Pad (No Sensor or Switch)<br><i>Choose Lateral Arm Adjustment Range:<br/>Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>   |           |
| SELECT ADDITIONAL RIGHT ARM PACKAGE (OPTIONAL): |  |           |
| <input type="checkbox"/> I10979                 | Add Right Lateral Arm w Proximity Sensor<br><i>Choose Lateral Arm Adjustment Range:<br/>Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>   |           |
| <input type="checkbox"/> I10980                 | Add Right Lateral Arm w Egg Switch<br><i>Choose Lateral Arm Adjustment Range:<br/>Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>   |           |
| <input type="checkbox"/> I10981                 | Add Right Lateral Arm w Pad (No Sensor or Switch)<br><i>Choose Lateral Arm Adjustment Range:<br/>Short (5.5" - 10"): _____ Long (5.5" - 14"): _____</i>  |           |

## Alt Drive Systems Setup

### CHOOSE OMNI (REQUIRED):

| Part Number                           | Description  | Price AUD |
|---------------------------------------|--|-----------|
| <input type="checkbox"/> I108501-99-0 | Omni2 Control Module for R-net - K450<br><i>Two-section design features 3.5" VGA high resolution display, USB charger port, built in IR control and integrated Bluetooth technology for control of up to 4 Bluetooth devices. Must also select mounting hardware below. Please specify which type of alternative drive input device you will be using. Input Device: _____</i> |           |



## Alt Drive Systems Setup

### CHOOSE OMNI (REQUIRED):

| Part Number                           | Description  | Price AUD |
|---------------------------------------|--|-----------|
| <input type="checkbox"/> I103476-99-0 | Adjustable Swing-Away Hardware - Omni/Tablet<br>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____ |           |

### CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

| Part Number                     | Description                              | Price AUD |
|---------------------------------|--|-----------|
| <input type="checkbox"/> I10080 | Mount Seat Control Switchbox on Backrest | No Charge |

## Alt Drive Accessories

### PERMOFIX MOUNTING SOLUTIONS:

| Part Number                           | Description   | Price AUD |
|---------------------------------------|---|-----------|
| <input type="checkbox"/> I106655-99-0 | PERMOfix UniTrack Fixed Mount Kit<br><i>Includes two base joints, two 4" rods and one 8" rod.</i>   |           |
| <input type="checkbox"/> I106656-99-0 | PERMOfix UniTrack Removable Mount Kit<br><i>Includes quick release mount and thigh support hardware, two base joints and three different length rods.</i>   |           |
| <input type="checkbox"/> I106657-99-0 | PERMOfix UniTrack Swing Away Mount Kit<br><i>Includes PERMOfix Swing Away module, two base joints, one 2" rod and one 4" rod.</i>   |           |
| <input type="checkbox"/> I105907-99-0 | Proximity Sensor + Swivel Assm - Not Mounted<br><i>Only available with the Total Control Head Array. Swivel assembly allows sensor to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array. Since Proximity Sensors require power, they CANNOT be used as On/Off switches.</i> |           |
| <input type="checkbox"/> I105919-99-0 | Egg Switch + Swivel Mount Assm - Not Mounted<br><i>Swivel assembly allows switch to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>  |           |
| <input type="checkbox"/> I105912-99-0 | Solid Lateral Pad+Swivel Joint - Not Mounted<br><i>Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>   |           |
| <input type="checkbox"/> I106675-99-0 | Oval Plate + Swivel Assm - No Switch - Not Mounted<br><i>Blank plate provides mounting option for Micro Lite, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>                                    |           |

### OPTIONAL ITEMS:

| Part Number                           | Description  | Price AUD |
|---------------------------------------|--|-----------|
| <input type="checkbox"/> I101769-99-0 | Microlight Switch Qty: _____   |           |
| <input type="checkbox"/> I101767-99-0 | Buddy Button Qty: _____  |           |
| <input type="checkbox"/> I106293-99-0 | Egg Switch Black Qty: _____  |           |
| <input type="checkbox"/> I101771-99-0 | Mini Cup Switch Qty: _____   |           |
| <input type="checkbox"/> I107500-99-0 | Proximity Switch Qty: _____  |           |
| <input type="checkbox"/> I107495-99-0 | Mini Joystick 5 Switch OC for Omni - (9-pin)   |           |
| <input type="checkbox"/> I10772       | Dome Switch, Backrest Mnt w Gooseneck Qty: _____<br><i>Mount on Left: _____ Mount on Right: _____</i><br><i>Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights. Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 150mm rod, one 200mm rod and one 225mm gooseneck switch with 1.2m cable and mounting clamp.</i> |           |
| <input type="checkbox"/> I10846       | Dome Switch, Backrest Swg Awy Mnt w Gsneck Qty: _____<br><i>Mount on Left: _____ Mount on Right: _____</i><br><i>w 90° Swing-Away Function. Mounts at top of back. Custom charge applies if ordered with short back. Kit Includes: One PERMOfix joint with Flat Plate, one PERMOfix joint with 45° swivel adjustment, two 6" rods, one 90° Swing-Away joint, one 9"</i>  |           |
| <input type="checkbox"/> I105935-99-0 | 3.5mm Ext. Cable for Sensor/Switch, 0.3m Long  |           |
| <input type="checkbox"/> I105931-99-0 | 3.5mm Ext. Cable for Sensor/Switch, 1m Long  |           |
| <input type="checkbox"/> I105932-99-0 | 3.5mm Ext. Cable for Sensor/Switch, 2m Long  |           |



**Alt Drive Accessories**

**OPTIONAL ITEMS:**

| Part Number                           | Description  | Price AUD |
|---------------------------------------|--|-----------|
| <input type="checkbox"/> I105933-99-0 | 3.5mm Ext. Cable for Sensor/Switch, 3m Long  |           |
| <input type="checkbox"/> I105934-99-0 | 3.5mm Y-splitter Cable for Sensor/Switch<br><i>Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.</i> |           |

Additional Order Instructions:

**Chair Order Policy:**

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

**Order Acknowledgement:**

I, \_\_\_\_\_, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: \_\_\_\_\_

Printed Name: \_\_\_\_\_

Title: \_\_\_\_\_

Date: \_\_\_\_\_

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty.