



# UPDATE

The label on the charger has been updated. Please refer to the new label image in place of the label image in the user manual. All associated instructions still apply.

Old label



New label



Please read and follow all instructions and warnings in all manuals supplied with your power wheelchair and its accessories. Incorrect use may both injure the user and damage the wheelchair. Please note the following additional warning:

Make sure the charger is properly set to charge the batteries installed in your power wheelchair. If the wrong charger setting is used, the batteries can be exposed to overcharging which increases the risk of overheating. If the batteries overheat, they can leak hydrogen gas and acid.