

**Account information**

Date	Account Number
Client: First Name	Last Name
Phone Number	Fax Number
Email Address	Primary Contact
P.O. Number (Enter "Quote only" if just requesting a quote)	

**Ship To**

First Name	Last Name
Address	
City	State/Zip Code

**Overall desired cushion width:**

Back: inches

Seat: inches

**As narrow as possible?** Yes No

*The cushion will be cut as narrow as possible (ANAP) (which may exceed overall cushion width)*

**OK to wedge to desired width?** Yes No

\* "OK to wedge". The mounting surface will be cut to the desired overall width.

**Desired wheelchair back-to-seat angle:**

Specific degrees  
(90°- 135° for full system or 90° for a single cushion)

**Desired cushion height:**

(Back height from seat pan to top of back)


Back: inches

**Check one of the following:**

Provide quote only\*

Begin fabrication

\*Quote pricing valid for six months from date quoted.



# Tru-Shape Itemized

Order form & retail price list

Effective August 16, 2021

**3 Pieces to a successful order**

1. Completed order form
2. 3D scan or digital file
3. Electronic purchase order or P.O. number

*\*All 3 pieces must be submitted together  
Please review remake policy on page 7*

**Permobil**

Email: orders.shape@permobil.com

*\*Please note: Adobe Acrobat Reader 6.0 or newer required to view and complete form. It is recommended that you print a copy for your records as older and free versions of Adobe may prevent you from saving your progress before submitting.*

**Internal Lateral Dimensions:**

inches

**Internal Adductor Dimensions:**

inches

**Chair Dimensions:**

*(Reference purposes only; measurements are not used when designing the system)*

Width: inches

Height: inches

©2021 Permobil • 08/12/21  
All specifications and pricing are subject to change without notice.

Back & Seat cushion—Itemized						
			HCPCS	Part #	MSRP	
<b>Back/Seat cushion—COMPLETE</b> (choose foam and all desired options, below)						
Foam type (choose one)	Standard <i>(price includes cover and spray coating)</i>	Back	E2617	6364910	\$1,581	
		Seat	E2609	6364915	\$1,581	
	Firm <i>(price includes cover and spray coating)</i>	Back	E2617	6364911	\$1,674	
		Seat	E2609	6364916	\$1,674	
	Hybrid <i>(price includes cover and spray coating)</i> Inches of firm foam from front of seat or from bottom of back	Back	E2617	6364911	\$1,674	
		Seat	E2609	6364916	\$1,674	
<b>Back/Seat cushion—Unsprayed</b>						
Foam type (choose one)	Standard <i>(price does not include spray coating or cover)</i>	Back	K0669	6364900	\$1,179	
		Seat	K0669	6364905	\$1,179	
	Firm <i>(price does not include spray coating or cover)</i>	Back	K0669	6364901	\$1,194	
		Seat	K0669	6364906	\$1,194	
	Hybrid <i>(price does not include spray coating or cover)</i> Inches of firm foam from front of seat or from bottom of back	Back	K0669	6364901	\$1,194	
		Seat	K0669	6364906	\$1,194	
<b>Fabrication options for back &amp; seat cushions</b>						
***CHECK ALL DESIRED OPTIONS***				Part #	MSRP	
<b>1. Moisture resistant spray sealant</b>						
Spray Sealant		Back	Seat	6364350/ 6364430	\$237*	
<b>2. Removable fabric cover options*</b>						
1st Removable cover (choose fabric type) <i>(Cover pattern included with order)</i>			Back	6364335	\$305*	
Neoprene                      Spacer fabric			Seat	6364335	\$305*	
Breathe 4™ (Lightweight, breathable, four-way stretch)						
COMFORT-TEK® (For fluid protection and an easily cleaned surface)						
STRETCH-AIR® (For client comfort and heat dissipation)						
2nd Removable cover (choose fabric type)			Back	6364340	\$239	
Neoprene                      Spacer fabric                      Breathe 4™			Seat	6364410	\$239	
COMFORT-TEK®                      STRETCH-AIR®						
Replacement pattern from existing cover				6364305	\$99	
<b>3. Relief options</b> (Relief of Seating Surface—Must relieve BOTH seat and back cushion the same amount in order for contours to match.)						
Surface relief	1/8"	1/4"	1/2"	3/4"	6364422	\$112
<b>4. Soft spot options</b>						
Soft spot overlay (1" Visco foam; recommend a minimum of a 1/8" surface relief)				6364325	\$350/ea.	
Soft spot (1" Visco foam) <i>(must be noted on digital file)</i>		Back:	Seat:	Total quantity:	6364311	\$249/ea.
Standard IT area		As digitized	No soft spot			

\* Included with complete cushions

5. Lateral/adductor reinforcement					
		Back	Seat	Part #	MSRP
	Left Thinned	Back	Seat		N/C
	Left Reinforced	Back	Seat		N/C
	Left Double Reinforcement	Back	Seat	6364426/ 6364435	\$73
	Right Thinned	Back	Seat		N/C
	Right Reinforced	Back	Seat		N/C
	Right Double Reinforcement	Back	Seat	6364427/ 6364515	\$73
6. Swing-away laterals/adductors					
	Left (price for 1 lateral; includes custom cover)	Back	Seat	6364355	\$859
	Right (price for 1 lateral; includes custom cover)	Back	Seat	6364355	\$859
	Swing-away lateral hardware (Includes hardware for both sides, must be chosen if using our hardware.)	Back	Seat	471C49= PBBKIT	\$427
7. Belt guides/notches					
	Pass-through belt guides - must be marked on digital file		Back	6364416	\$192
			Seat	6364416	\$192
	Standard notches	2" Long, 2.5" Deep, 0.5" in front of parting angle		6364415	\$110
	Specific notches Long: _____ Location: _____ in front of parting angle Deep: _____		Back	6364415	\$110
			Seat	6364415	\$110
8. Ventilation					
Back	Ventilation holes (breathable Spacer fabric or Breathe 4 covers recommended) Evenly spaced      Where marked on digital file			6364330	\$355
Seat	Ventilation holes (breathable Spacer fabric or Breathe 4 covers recommended) Evenly spaced      Where marked on digital file			6364330	\$355
9. Back profile					
	Round top corners (3" radius)				N/C
	Taper profile				N/C
10. Clearance for back canes (*For Full Systems ONLY, if adding clearance for back canes, ALSO fill out rear seat notches in section 16)					
	Width between canes:	(inches)			
	Depth of canes:	(inches)			
11. Mounting options					
	Plywood—1/2" thickness	Back	Seat	6364320/ 6364385	\$132
	ABS plastic—1/4" thickness	Back	Seat	6364345/ 6364420	\$149
	ABS plastic—3/8" thickness	Back	Seat	6364360/ 6364390	\$149
	Velcro—Two ea. 2"-wide strips	Back	Seat	6364346/ 6364421	\$32
	T-nut package (flat fee for all requested t-nuts including headrest)	Back	Seat	6364912	\$75

12. T-nut configuration for Back				
T-nut size	1/4 - 20	10-32	5mm	
Pattern: (Choose One)				
<p>A. Specific</p>		<p>B. Three columns (Default)</p>		<p>C. Four corners</p>
Number of columns			Standard placement of T-nuts: 1" from the top, 1" from the sides, and 1" from between columns and rows	
Number of rows			Standard placement of T-nuts: 1" from the top, 1" from the sides, and 1" from between columns and rows	
A	= from top (minimum of 1")			
B	= from sides (minimum of 1")			
C	= from bottom (minimum of 1")			
D	= between columns (minimum of 1")			
E	= between rows (minimum of 1")			
13. T-nut configuration for Head Supports				
T-nut size	1/4 - 20	10-32	5mm	
Pattern:				
<p>Other</p>		<p>5 x 5 1" Spacing (Default)</p>		<p>Whitmeyer</p>
F	= number of rows			
G	= number of columns			
H	= down from top			
I	= between columns			
J	= between rows			
Location down from the top		(inches) (default is 3")		
Position:	Center	Offset to client's left	(inches)	Offset to client's right (inches)
Universal Adapter Plate (including mounting hardware) <i>cannot be paired with Whitmeyer pattern</i> (allows you to mount most mfg. head supports to cushions)			Part #	MSRP
			471C49=HR-HRAP	\$57
14. T-nut configuration for Seat				
T-nut size	1/4 - 20	10-32		
Pattern: (Choose One)				
<p>A. Specific</p>		<p>B. Three columns (Default)</p>		<p>C. Four corners</p>
Number of columns		Standard placement of T-nuts: 1" from the rear, 1" from the sides, and 1" from between columns and rows		Standard placement of T-nuts: 1" from the rear, 1" from the sides, and 1" from between columns and rows
Number of rows				
A	= from rear (minimum of 1")			
B	= from sides (minimum of 1")			
C	= from front (minimum of 1")			
D	= between columns (minimum of 1")			
E	= between rows (minimum of 1")			

15. Undercutting front of seat cushion					
Standard					N/C
Specific	A. Behind trimline B. Below leg trough			6364413	\$108
16. Clearance for seat rails					
Width between rails:	(inches)				
Depth of rails:	(inches)				
17. Corner notches					
Notch front seat corners	Width	Length	Depth*	Front 	Rear 
Notch rear seat corners	Width	Length			
*If not full depth					
18. Leg length discrepancy modification					
Carve front edge of seat as scanned					N/C
Specific					N/C
If specific cutback, choose one of the contouring options below:  = material to be cut back					
Cut back right only		Cut back left only		Cut back both	
Standard A		inches			N/C
Cut back entire side, to center of medial knee support					
Standard B		inches			N/C
Cut back entire side, but leave medial knee support					
Standard C		inches			N/C
Cut back entire side, including medial knee support					
Standard D		R inches	L inches		N/C
Cut back leg trough only, leaving both adductors and medial knee support full length					

**Back & Seat cushion—Itemized** continued from the previous page

19. Mounting hardware			
	Number required to mount cushion	Part #	MSRP
<b>Knob release, adjustable mounting kits (WC-20 Crash Test Approved)</b>			
Seat Mount Kit for 1" Tubing	Includes: (4) Seat Clamps & (4) Seat Forks	471C49=OBSK88	\$558
Back Mount Kit for 1" Tubing	Includes: (4) Seat Clamps, (2) Left Forks, (2) Right Forks	471C49=OBBK88	\$571
Seat & Back Mount Kit for 1" Tubing	Includes: Back-to-Seat Hinge, (2) Seat Clamps, (2) Seat Forks, (2) Back Clamps, (1) Back Left Fork, (1) Back Right Fork	471C49=OBSBK88	\$707
<b>Quick release, non-adjustable mounting kits</b>			
Flush Mount Kit - 7/8"	1	471C49=NT-K-FH78	\$245
Flush Mount Kit - 1"	1	471C49=NT-K-FH88	\$245
1" Drop Hook with Stops Kit - 1"	1	471C49=NT-K-DH188	\$264
2" Drop Hook with Stops Kit - 1"	1	471C49=NT-K-DH288	\$276
<b>Adjustable J and L bracket kits</b>			
Adjustable Hook Hardware Kit - 5 Hole L Bracket/4 Hole J Bracket for 1" Tubing	Complete kit for seat and back on 1" tubing	470C49=KAH5L4J88	\$210
Adjustable Hook Hardware Kit - 5 Hole L Bracket/2 Hole J Bracket for 1" Tubing	Complete kit for seat and back on 1" tubing	470C49=KAH5L2J88	\$210
<b>Miscellaneous</b>		<b>Part #</b>	<b>MSRP</b>
<b>Quantity:</b>			
Black patch paint—Latex free*		TS-PATCHPAINT	\$19
Sewing pattern only		6364336	\$167
Replacement pattern from existing cover		6364305	\$99

\* For touch-up of minor cuts, holes, or abrasions in vinyl coating that may occur during normal use. Shipped as repair material in individual box. Add \$20 hazardous material handling charge.

**Specific instructions not otherwise covered in order form.** Note: We encourage you to provide additional instructions and sketches to ensure a successful fabrication. Examples: soft spot location, special t-nut patterns, soft spot overlays, leg length discrepancies, double reinforcements, swing away laterals, etc.

**3 Pieces to a successful order (All 3 pieces must be submitted together)**

1. Completed order form
2. 3D scan or digital file
3. Electronic purchase order or P.O. number

**Please review remake policy on next page**

## TRU-SHAPE remake and warranty policy

### Remake policy

A remake is defined as a product that did not meet the customer's expectations and/or did not meet the user's needs due to an error in the manufacturing process or error in shape translation.

- Within the first 180 days of the original invoice date: If, for one of the reasons listed above, the seating system outcome is not satisfactory, Permobil will work to modify and adjust the product to remedy the deficiency. Errors made by the customer during the scanning/ordering process are not considered a remake. If these modifications and adjustments are unsuccessful, the customer will be contacted by Permobil for advisement if a remake is appropriate and, if so, will initiate the Remake Request Process.
- A Remake Request Form, completed by the ordering ATP, is required to begin the remake process.
- Any approved remake request for an improperly fit product that is 6-12 months from the original invoice date will result in a charge equal to 50% of the original invoice amount. This remake policy does not extend beyond 12 months from the original invoice date.
- Return authorization must be obtained prior to returning any product. Permobil will not process any item without a valid return authorization number.
- Permobil will determine if the product in question needs to be returned for review. If Permobil determines the product does not need to be returned, the customer must fill out and return a Certificate of Destruction to continue with the remake process.
- The customer is responsible for all return shipping costs associated with any remake.

### Warranty policy

Listed below are the applicable warranty periods from the original invoice date:

- Tru-Shape seat and back: 2 years
- Covers: 6 months
- Hardware: 1 year

### What is covered?

Defects related to materials and workmanship.

### What is not covered?

Any seating system that has been misused, damaged by an accident or damage from an "act of Nature" (e.g., flood, tornado, earthquake, fire, etc.). Any substitutions of the supplied accessories or modification to the seating system may void the warranty.