

K450 MX



*** Indicates a Required Field**

Dealer Information

Order ID: _____

*ATP: _____

Dealer Code: _____

*Dealer Name: _____

Address: _____

*City: _____

* State/Zip: _____

* Phone#: _____

PO#: _____

Email Address: _____

Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

* First Name: _____

* Last Name: _____

Diagnosis: _____

Funding Source: _____

Client Age: _____

Client Measurements

* Weight: _____

*Height: _____

A) Top of Shoulders: _____

B) Chest Depth: _____

C) Chest Width: _____

D) Seat Depth: _____


E) Top of Head: _____

F) Elbow to Hand: _____

G) Seat Pan to Elbow: _____

H) Hip Width: _____

*I) Knee to Foot: _____



Three column pricing on order/quote:

Display HCPCS Codes on order/quote:

Please send order/quote to fax#: **(800) 231-3256**

Email to: sales@permobil.com

Permobil Inc.
300 Duke Dr.
Lebanon, TN 37090
Tel: (800) 736-0925
Fax: (800) 231-3256
www.permobil.com

Prices effective **January 1, 2022**

K450 MX

Chair Model and Colors

CHOOSE BASE:

Part Number	Description	Price	HCPCS Code
<input checked="" type="checkbox"/> I104524-99-0	K450 Base w/ R-net 5.2mph - K450 PS2 <i>Price Includes: K450 Power Base Set Up For Multiple Power Seat Functions, Tie-Down Hardware for Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. 5.0mph, Weight Capacity = 125lbs, STFH = 3"-25".</i>	16,091.00	K0891

CHOOSE DRIVE TIRES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I106204-99-0	Flat Free Casters - K450	No Charge	
<input type="checkbox"/> I106205-99-0	Flat Free Drive Tires - K450	No Charge	
<input type="checkbox"/> I106206-99-0	Pneumatic Casters - K450	No Charge	
<input type="checkbox"/> I106207-99-0	Pneumatic Drive Tires - K450	No Charge	

CHOOSE SHROUD COLOR:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104535-99-0	Shroud for K450 - Onyx Black	No Charge	
<input type="checkbox"/> I104536-99-0	Shroud for K450 - Platinum Silver	No Charge	
<input type="checkbox"/> I104537-99-0	Shroud for K450 - Pop Star Pink	No Charge	
<input type="checkbox"/> I104539-99-0	Shroud for K450 - Extreme Yellow	No Charge	
<input type="checkbox"/> I106367-99-0	Shroud for K450 - Volcano Red	No Charge	
<input type="checkbox"/> I106368-99-0	Shroud for K450 - Cobalt Blue	No Charge	
<input type="checkbox"/> I107234-99-0	Shroud for K450 - Purple Rain	No Charge	

Base Options

MUST CHOOSE SEAT ELEVATOR:

Part Number	Description	Price	HCPCS Code
<input checked="" type="checkbox"/> I10510	Power Adjustable Seat Height - 8" Travel <i>This item must be selected.</i>	3,369.00	E2300

MUST CHOOSE BATTERIES:

Part Number	Description	Price	HCPCS Code
<input checked="" type="checkbox"/> I107032-99-0	Batteries, Grp 34 Gel, Pair <i>Includes Permobil VoltPro charger for non-Permobil gel batteries.</i>	998.00	E2359

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104550-99-0	Support Wheels for K450	500.00	
<input type="checkbox"/> I104540-99-0	Light Kit - K450 R-net	1,280.00	
<input type="checkbox"/> I104553-99-0	WC 19 Transit Kit - K450 <i>This kit includes additional anchor points for tie down straps, seat stabilizing bars and hardware, and a WC-19 Approved integrated lap and shoulder belt. This kit has been approved according to ANSI/RESNA WC19. An ICS Switchbox must be installed on the chair in order to access WC19 mode.</i>	1,997.00	
<input type="checkbox"/> I104554-99-0	Fixed Ventilator Tray - K450 <i>Tray is 9.75"D x 13"W and increases the overall length of a K450 by approximately 11". The ventilator tray is designed to accommodate an LP10 style ventilator and one U1 Battery. A custom charge will be applied if ordering vent tray with any of the following items: WC-19 Transit Kit (I104553-99-0), Attendant Control (ID50882), Push Handles (I104576-99-0) or Bag Hooks (I105805-99-0). Please Note: Must also select Support Wheels (I104550-99-0).</i>	661.00	E1029
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	547.00	

Joystick Options

CHOOSE JOYSTICK TYPE:

Part Number	Description	Price	HCPCS Code
	PLEASE NOTE: <i>For additional input devices, please refer to the Alternative Drive Systems Section.</i>		

K450 MX

Joystick Options

CHOOSE JOYSTICK TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I108319-99-0	Permobil Joystick Module <i>Designed exclusively for Permobil, this R-net color joystick features: control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.</i>	1,125.00	
<input type="checkbox"/> I105246-99-0	R-net LED Remote Joystick - No Display <i>This Joystick does NOT have a Display Screen, does NOT contain mono jack ports, and does NOT include Bluetooth or Assignable Buttons. PC Programming required and profiles limited to five.</i>	No Charge	

CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard	No Charge	
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	132.00	E2324
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	186.00	E2323
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	132.00	E2323
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	174.00	E2323
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	174.00	E2323
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	174.00	E2323
<input type="checkbox"/> I105442-99-0	1.75" Black Foam Ball For Joystick Knob	132.00	E2323
<input type="checkbox"/> I105846-99-0	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	174.00	E2323
<input type="checkbox"/> I105847-99-0	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	174.00	E2323
<input type="checkbox"/> I105848-99-0	Bodypoint J/S Handle Dome Shaped, Rubber	174.00	E2323

CHOOSE JOYSTICK MOUNT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10412	Joystick Mount Right - Fixed, R-net Only	No Charge	
<input type="checkbox"/> I10413	Joystick Mount Left - Fixed, R-net Only	No Charge	
<input type="checkbox"/> I102537-99-0	Joystick Mount Right - Swing Away	389.00	E1028
<input type="checkbox"/> I102538-99-0	Joystick Mount Left - Swing Away	389.00	E1028

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10411	Expandable Controller - R-net <i>This item is required when ordering power tilt, power recline, and power elevating legrests together OR when ordering an Alternative Drive Control.</i>	727.00	E2377
<input type="checkbox"/> I10513	Harness for Expandable Controller <i>Must be ordered with any Expandable Controller.</i>	645.00	E2313
<input type="checkbox"/> ID50882-A	Attendant Control for R-net Control System Only <i>On the K450, the mounting location for the Attendant Control is off of the Push Handle Bracket. If ordering an attendant control, you must also order I104576-99-0 (Push Handles - MX Seat).</i>	1,427.00	E2331

Seating System

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input checked="" type="checkbox"/> I104324-99-0	MX Power Seating System <i>Price Includes: MX Seating System, Adjustable Seat Frame, Power Seat to Floor System, Adjustable Backrest Angle, and One Piece Footplate. Weight Limit 125lbs.</i>	8,869.00	K0108

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104529-99-0	No Tilt <i>Weight Limit 125lbs.</i>	No Charge	
<input type="checkbox"/> I104551-99-0	45° Power Tilt <i>Weight Limit 125lbs.</i>	6,991.00	E1002

Seating System Options

CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	Price	HCPCS Code
-------------	-------------	-------	------------

K450 MX

Seating System Options

CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	Price	HCPCS Code
<input checked="" type="checkbox"/> I10585	Multiple Seat Function Control Kit for R-net <i>For operation of two or more seat functions through the wheelchair electronics via the drive control. Order with R-net electronics or remote joystick options.</i>	3,803.00	E2311

CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104135-99-0	Push Buttons on ICS Switchbox (3m)	No Charge	
<input type="checkbox"/> I104136-99-0	Toggle Switches on ICS Switchbox (3m)	340.00	

CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge	
<input type="checkbox"/> I104651-99-0	Mount ICS Switchbox onto Armrest <i>Mount on Left: _____ Mount on Right: _____</i>	No Charge	
<input type="checkbox"/> I11050	Mount ICS Switchbox on Backrest <i>On the K450, the backrest mounting location for the ICS Switchbox is on the Push Handle bracket; therefore, you must also order I104576-99-0 (Push Handles - MX). Mount on Left: _____ Mount on Right: _____</i>	163.00	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105409-99-0	Do Not Install ICS Switchbox, Send as Part	No Charge	
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Seat function interface box that accommodates up to 8 momentary switches. Provides an additional method for accessing the seat functions. 1/8" mini-mono switches not included.</i>	1,329.00	
<input type="checkbox"/> I106743-99-0	Remote Stop System - All Electronics <i>Consists of a 315MHz Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Custom charges may be required when ordering with a PS/PS Jr seat, a low backrest height or a pediatric wheelchair. If ordering on a K450, push handles must also be selected.</i>	1,315.00	E2321
<input type="checkbox"/> I10515	Memory Seat Program for K450 MX <i>If this is selected, the chair will be programmed to allow the user to store and save up to 3 different memory seat positions through the ICS Switchbox. On the K450 MX, ALL seat functions are capable of recalling a saved seat position. Please refer to the K450 Owner's Manual for more information on saving and recalling a memory seat position.</i>	No Charge	

Backrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
MX STYLE BACKRESTS: <i>Use this section to select MX Style Backrests. Back heights measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Please be sure to pick both a backrest shell size and a backrest cushion style (Leatherette or Black Mesh Fabric). Note: A 16"W MX Backrest is NOT available.</i>			
<input type="checkbox"/> I104321-99-0	MX Backrest - 10"W x 14"H <i>Backrest can be adjusted to either 14" or 16" high. If ordering a headrest, backrest will be mounted in the 16" high position unless client measurements require a 14" high backrest, in which case, custom modifications will be made.</i>	No Charge	
<input type="checkbox"/> I104322-99-0	MX Backrest - 12"W x 16"H <i>Backrest can be adjusted to either 16" or 18" high.</i>	No Charge	
<input type="checkbox"/> I104323-99-0	MX Backrest - 14"W x 20"H <i>Backrest can be adjusted to either 20" or 22" high.</i>	No Charge	

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104309-99-0	MX/MF Backrest Cushion, 10"W x 14"T - Leatherette	No Charge	
<input type="checkbox"/> I104311-99-0	MX/MF Backrest Cushion, 12"W x 16"T - Leatherette	No Charge	

K450 MX

Backrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104313-99-0	MX/MF Backrest Cushion, 14"W x 20"T - Leatherette	No Charge	
<input type="checkbox"/> I104308-99-0	MX/MF Backrest Cushion, 10"W x 14"T - Black Mesh	No Charge	
<input type="checkbox"/> I104310-99-0	MX/MF Backrest Cushion, 12"W x 16"T - Black Mesh	No Charge	
<input type="checkbox"/> I104312-99-0	MX/MF Backrest Cushion, 14"W x 20"T - Black Mesh	No Charge	

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
PS JR TUBULAR STYLE BACKRESTS			
<i>Use this section to select PS Jr Tubular Style Backrest Canes. Note: A 10"W Tubular Backrest is NOT available.</i>			
<input type="checkbox"/> I103717-99-0	PS Jr Back Canes, Height & Angle Adjustable	No Charge	
<i>Must also select PS Jr Style Armrests if ordering this backrest. Choose PS Jr Style Backrest Type (Metal Back or Omit Back) below. Backrest width should match seat width. Back height is adjustable from 16.5"T - 22.5"T, measured from the seat pan to the top of the back.</i>			

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103719-99-0	PS Jr Solid Backrest - 12"W	No Charge	
<input type="checkbox"/> I103721-99-0	PS Jr Solid Backrest - 14"W	No Charge	
<input type="checkbox"/> I103723-99-0	PS Jr Solid Backrest - 16"W	No Charge	
<input type="checkbox"/> I10449	PS Jr - Omit Backrest	No Charge	

Seat Sizing

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104317-99-0	MX Seat Size 10"W x 10"D <i>Seat depth can be adjusted to 10", 11", or 12". If ordering this seat size, the lowest allowed seat position when using the seat to floor function will be approximately 3.5" from the floor to the top of the seat pan.</i>	No Charge	
<input type="checkbox"/> I104318-99-0	MX Seat Size 12"W x 12"D <i>Seat depth can be adjusted to 12", 13", or 14". If ordering this seat size with MX Arms, the lowest allowed seat position when using the seat to floor function will be approximately 3.5" from the floor to the top of the seat pan. If ordering Flip Back Arms, the lowest allowed seat position will be 8".</i>	No Charge	
<input type="checkbox"/> I104319-99-0	MX Seat Size 14"W x 14"D <i>Seat depth can be adjusted to 14", 15", or 16". If ordering this seat size with MX Arms, the lowest allowed seat position when using the seat to floor function will be approximately 3.5" from the floor to the top of the seat pan. If ordering Flip Back Arms, the lowest allowed seat position will be 8".</i>	No Charge	
<input type="checkbox"/> I104320-99-0	MX Seat Size 16"W x 16"D <i>Seat depth can be adjusted to 16", 17", or 18". If ordering this seat size with MX Arms, the lowest allowed seat position when using the seat to floor function will be approximately 10" from the floor to the top of the seat pan. If ordering Flip Back Arms, the lowest allowed seat position will be 16".</i>	No Charge	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
Seat Cushions Below Are Approximately 2" Thick:			
<input type="checkbox"/> I306384-99-90-0	Leatherette Seat Cush MX - 10"W x 10"D (250x250)	621.00	
<input type="checkbox"/> I306385-99-90-0	Leatherette Seat Cush MX - 12"W x 12"D (300x300)	621.00	
<input type="checkbox"/> I306386-99-90-0	Leatherette Seat Cush MX - 14"W x 14"D (350x350)	621.00	
<input type="checkbox"/> I104703-99-0	Leatherette Seat Cush MX - 16"W x 16"D (400x400)	621.00	
<input type="checkbox"/> I104305-99-0	Black Mesh Seat Cushion MX/MF - 10"W x 10"D	625.00	
<input type="checkbox"/> I104306-99-0	Black Mesh Seat Cushion MX/MF - 12"W x 12"D	625.00	
<input type="checkbox"/> I104307-99-0	Black Mesh Seat Cushion MX/MF - 14"W x 14"D	625.00	
<input type="checkbox"/> I104702-99-0	Black Mesh Seat Cush MX - 16"W x 16"D (400x400)	621.00	

Armrest Options

K450 MX

Armrest Options

CHOOSE ARMREST STYLE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10512	Height Adjustable Armrest Assembly MX Seat <i>If you are not selecting Tilt, you MUST order this item.</i>	360.00	E0973
<input type="checkbox"/> I10514	Height Adjustable Armrest Assembly MX Seat - NC <i>No Charge if ordering Power Tilt (I104551-99-0)</i>	No Charge	
<input type="checkbox"/> I10824	Flip Back, Height Adjustable Armrests <i>Mounted off the backrest canes. Includes: 45° of angle adjustability, rotational adjustment, 3" offset adjustment both medially and laterally, and height adjustment of 6" - 14" from the seat pan. In one inch increments, please indicate where you would like the armrests set.</i> <i>SET ARMRESTS AT (6" - 14"): _____</i>	405.00	E0973
<input type="checkbox"/> I10825	Flip Back, Height Adjustable Armrests - NC <i>No charge when ordering power tilt. In one inch increments, please indicate where you would like the armrests set.</i> <i>SET ARMRESTS AT (6" - 14"): _____</i>	No Charge	

MX ARMREST PADS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104303-99-0	MX Armrest Pad Short Leatherette 8"L x 3.5"W	No Charge	
<input type="checkbox"/> I104304-99-0	MX Armrest Pad Long 12.5"Lx3.5"W	No Charge	

FLIP BACK ARMREST PAD RIGHT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105100-99-0	4"x10" Corpus Armrest Pad Right - Leatherette	No Charge	
<input type="checkbox"/> I105104-99-0	4"x13" Corpus Armrest Pad Right - Leatherette	No Charge	
<input type="checkbox"/> I107212-99-0	4"x10" Corpus Armrest Pad Right - Black Fabric	No Charge	
<input type="checkbox"/> I106047-99-0	4"x13" Corpus Armrest Pad Right - Black Fabric	No Charge	
<input type="checkbox"/> I104951-99-0	3.5"x8" Unitrack Armrest Pad Right (C)	No Charge	
<input type="checkbox"/> I104950-99-0	3.5"x12.5" Unitrack Armrest Pad Right (G)	No Charge	
<input type="checkbox"/> I107591-99-0	4"x10" Corpus Arm w/ Permobil Gel Pad - Right	182.00	
<input type="checkbox"/> I107593-99-0	4"x13" Corpus Arm w/ Permobil Gel Pad - Right	187.00	

FLIP BACK ARMREST PAD LEFT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105098-99-0	4"x10" Corpus Armrest Pad Left - Leatherette	No Charge	
<input type="checkbox"/> I105102-99-0	4"x13" Corpus Armrest Pad Left - Leatherette	No Charge	
<input type="checkbox"/> I107213-99-0	4"x10" Corpus Armrest Pad Left - Black Fabric	No Charge	
<input type="checkbox"/> I106046-99-0	4"x13" Corpus Armrest Pad Left - Black Fabric	No Charge	
<input type="checkbox"/> I104952-99-0	3.5"x8" Unitrack Armrest Pad Left (C)	No Charge	
<input type="checkbox"/> I104949-99-0	3.5"x12.5" Unitrack Armrest Pad Left (G)	No Charge	
<input type="checkbox"/> I107590-99-0	4"x10" Corpus Arm w/ Permobil Gel Pad - Left	182.00	
<input type="checkbox"/> I107592-99-0	4"x13" Corpus Arm w/ Permobil Gel Pad - Left	187.00	

Legrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104314-99-0	MX Footplate - Small (11.5" W x 6.5"D)	No Charge	
<input type="checkbox"/> I104315-99-0	MX Footplate - Large (14" W x 8"D)	No Charge	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103427-99-0	Footstrap Kit	208.00	
<input type="checkbox"/> I103104-99-0	Footstrap Kit - Padded	230.00	
<input type="checkbox"/> I104338-99-0	Small Anti-Skid Overlay for Footplate - K450	170.00	
<input type="checkbox"/> I104337-99-0	Large Anti-Skid Overlay for Footplate - K450	170.00	
<input type="checkbox"/> I107503-99-0	BP Ankle Huggers, Med Side Release, Adj (Pair) <i>This item will be sent as a part with the wheelchair. It will not be installed.</i>	300.00	

Accessories

K450 MX

Accessories

<input type="checkbox"/>	I104518-99-0	Bodypoint Non-Padded Hip Belt - Medium	No Charge
<input type="checkbox"/>	I105849-99-0	Bodypoint Padded Hip Belt - Extra Small	159.00
<input type="checkbox"/>	I104520-99-0	Bodypoint Padded Hip Belt - Medium	159.00
<input type="checkbox"/>	I100519-99-0	Positioning Belt Retractable	316.00
<input type="checkbox"/>	I108267-99-0	Bodypoint Padded Hip Belt "Airplane" Latch - Med	213.00

HEADREST OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/>	I100749-99-0 MF Headrest for Koala/PS Jr - Leatherette <i>Must also select PN I10552.</i>	356.00	E0955
<input type="checkbox"/>	I104339-99-0 MF Headrest for MX Back - Leatherette <i>Must also select PN I10552. Also includes a UniTrack Bracket for Lateral Adjustment of the headrest. Please select appropriate UT Bracket below:</i> UT Bracket for Lat Head Adj (10"W MX Back): _____ UT Brack for Lat Head Adj (12-14"W MX Back): _____	356.00	E0955
<input type="checkbox"/>	I10552 Adj/Removable Headrest Hardware	291.00	E1028
<input type="checkbox"/>	I107499-99-0 Universal Headrest Adapter <i>Please Note: Recline is limited to 160° or less (depending on backrest height) when ordering a Universal Headrest Adapter on a 23" or 24" tall Corpus backrest.</i>	287.00	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/>	I106104-99-0 UT Hip Suppt (B) w/ Mesh Cover - High (Pair) <i>4.5"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	300.00	E0953
<input type="checkbox"/>	I106103-99-0 UT Hip Suppt (A) w/ Mesh Cover - High (Pair) <i>6"W x 5.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	300.00	E0953
<input type="checkbox"/>	I10550 Adj Removable - Thigh Support Hardware (Pair)	573.00	E1028
<input type="checkbox"/>	I106210-99-0 RAM Self-Leveling Cup Holder	88.00	
<input type="checkbox"/>	I106806-99-0 RAM X-Grip Phone Holder	88.00	
<input type="checkbox"/>	I107100-99-0 RAM X-Grip 10" Tablet Holder w/Swing Away Mount <i>Compatible with most 10" screen tablets with or without a protective case. Includes adjustable swing away mount (I103476-99-0). Customization may be required if ordered with an Omni.</i>	444.00	
<input type="checkbox"/>	I107191-99-0 USB Charger - MX <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets. The charger features an isolated circuit design for improved safety; and by pulling power from the ICS, prevents battery drainage when the chair is powered off. Includes mounting hardware and ICS connector.</i> <i>Mount on left armrest: _____ Mount on right armrest: _____</i>	220.00	
<input type="checkbox"/>	I104576-99-0 Push Handles - K450 <i>Mounted off the chassis. Height and angle adjustable. Must be selected if ordering bag clips or an attendant control, or if mounting the ICS Switchbox on the backrest.</i>	466.00	
<input type="checkbox"/>	I105805-99-0 Medical Necessities Bag Hooks <i>Please Note: Recline is limited to 165° or less (depending on backrest height) when ordered with 23" or 24" tall Ergo backrest. Custom required if ordered with a ROHO AGILITY backrest. May ship as a part at shorter backrest heights.</i>	105.00	
<input type="checkbox"/>	I108276-99-0 Bodypoint PivotFit Shoulder Harness - Medium	275.00	E0960
<input type="checkbox"/>	I100578-99-0 Small Angle Adjustable Tray <i>Must also select PN I10556.</i>	322.00	E0950
<input type="checkbox"/>	I10556 Adj Removable - Tray Hardware (1 Each)	291.00	
<input type="checkbox"/>	I108231-99-0 Permobil Adjustable Tray - 11"W	270.00	E0950
<input type="checkbox"/>	I108230-99-0 Permobil Adjustable Tray - 13"W	270.00	E0950
<input type="checkbox"/>	I108232-99-0 Permobil Tray Hardware Right Side Adjustable	198.00	
<input type="checkbox"/>	I108233-99-0 Permobil Tray Hardware Left Side Adjustable	198.00	

Alternative Drive Systems

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/>	I10751 Mini-Joystick for Omni - (9-pin connector) <i>Operating Force = 0.5 oz (13 g) Neutral to Full Deflection = 0.2" (5mm)</i>	3,447.00	E2312

K450 MX

Alternative Drive Systems

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103893-99-0	Compact Joystick for Omni - (9-pin connector) <i>Operating Force = 7.8 oz (220 g) Neutral to Full Deflection = 1.1" (28mm) Joystick includes two 1/8" mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,686.00	E2373
<input type="checkbox"/> I103894-99-0	Compact Joystick LITE for Omni - (9-pin connector) <i>Operating Force = 1.8 oz (50 g) Neutral to Full Deflection = 0.7" (17mm) Joystick includes two 1/8" mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,907.00	E2373

ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description	Price	HCPCS Code
	NOTE: <i>Please see the PERMOfix Mounting Solutions, within the Alt Drive Accessories section, for additional mounting options. The Swing Away Chin Boom Assemblies are Compatible with Permobil Headrests I107470-99-0, I107469-99-0 or I107468-99-0 and must be selected with I108721-99-0.</i>		
<input type="checkbox"/> I104690-99-0	Alt JS Swing Away Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	295.00	E1028
<input type="checkbox"/> I107135-99-0	Alt JS / Omni Combo Swing Away Mount Kit - Armrest <i>Used to mount Omni and Alternative Proportional Joystick to the same armrest. Includes Swing Away Omni/Tablet Mount, with a longer mounting rod, two PERMOfix modules, one 3" rod and one 4" rod. Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	558.00	E1028
<input type="checkbox"/> I104143-99-0	Alt JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	920.00	E1028
<input type="checkbox"/> I10752	Haims Harness for Chin Control Mount, Large <i>For neck width up to 7½".</i>	654.00	E1028
<input type="checkbox"/> I10753	Haims Harness for Chin Control Mount, Small <i>For neck width up to 5½".</i>	654.00	E1028

ALT PROPORTIONAL JOYSTICK ACCESSORIES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104904-99-0	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid arms. Switches set-up for Mode & On/Off function. Switch operation force = 4.5oz (130 g) Mount at front of Joystick: _____ Mount at rear of Joystick: _____ Mount at left of Joystick: _____ Mount at right of Joystick: _____</i>	274.00	
<input type="checkbox"/> I104641-99-0	Dual Switches for Compact Joystick -Gooseneck Mount <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches set-up for Mode & On/Off function. Switch operation force = 4.5oz (130 g)</i>	493.00	

SIP & PUFF SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10573	Sip & Puff System for Omni <i>Requires Omni Control Module. Includes breath tube kit.</i>	1,844.00	E2325,E2326
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Includes breath tube kit.</i>	2,454.00	E2325,E2326
<input type="checkbox"/> I106807-99-0	Sip & Puff Kit for Total Control Head Array <i>For use with Total Control Head Array. Requires Omni Control Module. Includes breath tube kit.</i>	951.00	E2325,E2326
<input type="checkbox"/> I10762	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573, I10574 or I106807-99-0.</i>	295.00	

Head Array Options

TOTAL CONTROL HEAD ARRAY:

Part Number	Description	Price	HCPCS Code
-------------	-------------	-------	------------

K450 MX

Head Array Options

TOTAL CONTROL HEAD ARRAY:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I106805-99-0	Permobil Head Array - 3-Switch System <i>Features an occipital pad sensor and two lateral arms configured with proximity sensors. The arms are mounted to the upper attachment points on the head array. Must select hardware below to complete assembly. ESP recommended. Choose Omni Control Module separately.</i>	4,589.00	E2330
<input type="checkbox"/> I11575	Swing Away Hardware for Side Sensor - Left	250.00	E1028
<input type="checkbox"/> I11576	Swing Away Hardware for Side Sensor - Right	250.00	E1028
<input type="checkbox"/> I11072	Program Head Array as a 5-Switch System <i>Select additional arms/sensors below and indicate the intended use of each sensor in the Notes & Comments section. Choose Omni Control Module separately.</i>	No Charge	
CHOOSE HEAD ARRAY MOUNT HARDWARE (REQUIRED):			
<input type="checkbox"/> I105927-99-0	Stealth TWB480 Backrest Mounting Hdwr with Links	305.00	E1028
<input type="checkbox"/> I105941-99-0	Stealth TWB480 Backrest Mnt Hdwr with Flipback	492.00	E1028
SELECT ADDITIONAL LEFT ARM PACKAGE (OPTIONAL):			
<input type="checkbox"/> I10976	Add Left Lateral Arm w Proximity Sensor	655.00	E1028
<input type="checkbox"/> I10977	Add Left Lateral Arm w Egg Switch	536.00	E1028
<input type="checkbox"/> I10978	Add Left Lateral Arm w Pad (No Sensor or Switch)	437.00	E1028
SELECT ADDITIONAL RIGHT ARM PACKAGE (OPTIONAL):			
<input type="checkbox"/> I10979	Add Right Lateral Arm w Proximity Sensor	655.00	E1028
<input type="checkbox"/> I10980	Add Right Lateral Arm w Egg Switch	536.00	E1028
<input type="checkbox"/> I10981	Add Right Lateral Arm w Pad (No Sensor or Switch)	437.00	E1028

Alt Drive Systems Setup

CHOOSE OMNI (REQUIRED):

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I108501-99-0	Omni2 Control Module for R-net - K450 <i>Two-section design features 3.5" VGA high resolution display, USB charger port, built in IR control and integrated Bluetooth technology for control of up to 4 Bluetooth devices. Must also select mounting hardware below. Please specify which type of alternative drive input device you will be using. Input Device: _____</i>	1,962.00	
<input type="checkbox"/> I103476-99-0	Adjustable Swing-Away Hardware - Omni/Tablet <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	295.00	

CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10080	Mount Seat Control Switchbox on Backrest	No Charge	

Alt Drive Accessories

PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I106655-99-0	PERMOfix UniTrack Fixed Mount Kit <i>Includes two base joints, two 4" rods and one 8" rod.</i>	297.00	
<input type="checkbox"/> I106656-99-0	PERMOfix UniTrack Removable Mount Kit <i>Includes quick release mount and thigh support hardware, two base joints and three different length rods.</i>	522.00	
<input type="checkbox"/> I106657-99-0	PERMOfix UniTrack Swing Away Mount Kit <i>Includes PERMOfix Swing Away module, two base joints, one 2" rod and one 4" rod.</i>	748.00	
<input type="checkbox"/> I105907-99-0	Proximity Sensor + Swivel Assembly <i>Only available with the Total Control Head Array. Swivel assembly allows sensor to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array. Since Proximity Sensors require power, they CANNOT be used as On/Off switches.</i>	378.00	E1028
<input type="checkbox"/> I105919-99-0	Egg Switch + Swivel Mount Assembly <i>Swivel assembly allows switch to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	255.00	
<input type="checkbox"/> I105912-99-0	Solid Oval Pad + Swivel Assembly <i>Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	154.00	

K450 MX

Alt Drive Accessories

PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I106675-99-0	Oval Plate + Swivel Assembly - No Switch <i>Blank plate provides mounting option for Micro Lite, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	135.00	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I101769-99-0	Microlight Switch Qty: _____	206.00	
<input type="checkbox"/> I101767-99-0	Buddy Button Qty: _____	174.00	
<input type="checkbox"/> I106293-99-0	Egg Switch Black Qty: _____	174.00	
<input type="checkbox"/> I101771-99-0	Mini Cup Switch Qty: _____	174.00	
<input type="checkbox"/> I107500-99-0	Proximity Switch Qty: _____	688.00	
<input type="checkbox"/> I107495-99-0	Mini Joystick 5 Switch OC for Omni - (9-pin)	747.00	
<input type="checkbox"/> I10772	Dome Switch, Backrest Mnt w Gooseneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____ Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 6" rod, one 8" rod and one 9" gooseneck switch with 50" cable and mounting clamp.</i>	545.00	
<input type="checkbox"/> I10846	Dome Switch, Backrest Swg Awy Mnt w Gsneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____ w 90° Swing-Away Function. Mounts at top of back. Custom charge applies if ordered with short back. Kit Includes: One PERMOfix joint with Flat Plate, one PERMOfix joint with 45° swivel adjustment, two 6" rods, one 90° Swing-Away joint, one 9" gooseneck switch w 50" cable and mounting clamp.</i>	806.00	
<input type="checkbox"/> I105935-99-0	3.5mm Ext. Cable for Sensor/Switch, 0.3m Long	34.00	
<input type="checkbox"/> I105931-99-0	3.5mm Ext. Cable for Sensor/Switch, 1m Long	37.00	
<input type="checkbox"/> I105932-99-0	3.5mm Ext. Cable for Sensor/Switch, 2m Long	44.00	
<input type="checkbox"/> I105933-99-0	3.5mm Ext. Cable for Sensor/Switch, 3m Long	48.00	
<input type="checkbox"/> I105934-99-0	3.5mm Y-splitter Cable for Sensor/Switch <i>Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.</i>	37.00	

Additional Order Instructions (for Permobil):

Notes & Comments (not for Permobil):

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, Inc. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. *The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.*