

Chairman Koala Mini-Flex



*** Indicates a Required Field**

Dealer Information

Order ID: _____

*ATP: _____

Dealer Code: _____

*Dealer Name: _____

Address: _____

*City: _____

* State/Zip: _____

* Phone#: _____

PO#: _____

Email Address: _____

Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

* First Name: _____

* Last Name: _____

Diagnosis: _____

Funding Source: _____

Client Age: _____

Client Measurements

* Weight: _____

*Height: _____

A) Top of Shoulders: _____

B) Chest Depth: _____

C) Chest Width: _____

D) Seat Depth: _____

E) Top of Head: _____

F) Elbow to Hand: _____

G) Seat Pan to Elbow: _____

H) Hip Width: _____

*I) Knee to Foot: _____



Three column pricing on order/quote:

Display HCPCS Codes on order/quote:

Please send order/quote to fax#: **(800) 231-3256**

Email to: sales@permobil.com

Permobil Inc.

300 Duke Dr.

Lebanon, TN 37090

Tel: (800) 736-0925

Fax: (800) 231-3256

www.permobil.com

Prices effective **January 1, 2022**

Chairman Koala Mini-Flex

Chair Model and Colors

CHOOSE BASE:

Part Number	Description	Price	HCPCS Code
<input checked="" type="checkbox"/> I106680-99-0	Koala Mini-Flex w/ R-net <i>Price Includes: FWD Koala Power Base setup for one Power Seat Function, R-net Electronics, Tie-Down Hardware for Strap Systems, Charger and Standard Positioning Belt. Must choose Flat Free or Pneumatic Tires below. Max Speed = 4MPH, Weight Capacity = 80lbs, STFH = 15".</i>	10,504.00	K0890

CHOOSE DRIVE TIRES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103889-99-0	Tri Spoke Split Rim w/ FF Drive Tires - 2.5 x 8"	No Charge	
<input type="checkbox"/> I103890-99-0	Tri Spoke Split Rim w/ Pneumatic Tire - 2.5 x 8"	No Charge	

CHOOSE SHROUD COLOR:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105756-99-0	Shroud for Koala - Onyx Black	No Charge	
<input type="checkbox"/> I105757-99-0	Shroud for Koala - Pop Star Pink	No Charge	
<input type="checkbox"/> I105758-99-0	Shroud for Koala - Platinum Silver	No Charge	
<input type="checkbox"/> I106369-99-0	Shroud for Koala - Volcano Red	No Charge	
<input type="checkbox"/> I106370-99-0	Shroud for Koala - Cobalt Blue	No Charge	
<input type="checkbox"/> I106371-99-0	Shroud for Koala - Purple Rain	No Charge	
<input type="checkbox"/> I106372-99-0	Shroud for Koala - Extreme Yellow	No Charge	
<input type="checkbox"/> I106945-99-0	Shroud for Koala - Green Apple	No Charge	

CHOOSE HUBCAP COLOR:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10871	Match Hub Cap Color to Shroud: _____ <i>Note: Hub caps will be matched to shroud color if no selection is made.</i> Solid Color Hub Caps - Pair	No Charge	
	<i>Platinum Silver: _____ Onyx Black: _____ Cobalt Blue: _____ Phantom Grey: _____ Pop Star Pink: _____ Volcano Red: _____ Purple Rain: _____ Titan Blue: _____ Sublime Green: _____ Crystal Blue: _____ Vintage Wine: _____ Extreme Yellow: _____ Chili Pepper Red: _____</i>		

Base Options

MUST CHOOSE BATTERIES:

Part Number	Description	Price	HCPCS Code
<input checked="" type="checkbox"/> I101979-99-0	Batteries, Koala (40AH) Sealed Gel, x2, Installed	914.00	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I106860-99-0	Light Kit - Koala	1,067.00	
<input type="checkbox"/> I100865-99-0	Front Bumper Koala	623.00	
<input type="checkbox"/> I106919-99-0	Enhanced Steering Perform Unit (ESP), R-net <i>Recommended item with switched drive inputs for improved tracking.</i>	1,580.00	
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	547.00	

Joystick Options

CHOOSE JOYSTICK TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I108319-99-0	PLEASE NOTE: <i>For additional input devices, please refer to the Alternative Drive Systems Section.</i> Permobil Joystick Module	1,125.00	
<input type="checkbox"/> I105246-99-0	R-net LED Remote Joystick - No Display <i>This Joystick does NOT have a Display Screen, does NOT contain mono jack ports, and does NOT include Bluetooth or Assignable Buttons. PC Programming required and profiles limited to five.</i>	No Charge	

Chairman Koala Mini-Flex

Joystick Options

CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard	No Charge	
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	132.00	E2324
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	186.00	E2323
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	132.00	E2323
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	174.00	E2323
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	174.00	E2323
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	174.00	E2323
<input type="checkbox"/> I105442-99-0	1.75" Black Foam Ball For Joystick Knob	132.00	E2323
<input type="checkbox"/> I105846-99-0	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	174.00	E2323
<input type="checkbox"/> I105847-99-0	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	174.00	E2323
<input type="checkbox"/> I105848-99-0	Bodypoint J/S Handle Dome Shaped, Rubber	174.00	E2323

CHOOSE JOYSTICK MOUNT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10127	Fixed Joystick Mount - Right	No Charge	
<input type="checkbox"/> I10136	Fixed Joystick Mount - Left	No Charge	
<input type="checkbox"/> I100849-99-0	Joystick Mount-Right Swing Away	387.00	E1028
<input type="checkbox"/> I100848-99-0	Joystick Mount-Left Swing Away	387.00	E1028
<input type="checkbox"/> I101260-99-0	Midline Control Mount	356.00	E1028

Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I11068	Expandable Controller - R-net <i>This item is required when ordering PN ID51153, PN ID51149, or when ordering an Alternative Drive Control.</i>	717.00	E2377
<input type="checkbox"/> I10513	Harness for Expandable Controller <i>Must be ordered with any Expandable Controller.</i>	645.00	E2313
<input type="checkbox"/> I105408-99-0	Attendant Control R-net	1,427.00	E2331
<input type="checkbox"/> I10582	Output Module (IOM) for R-net Control System <i>Enables a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). Each module can provide up to 4 switch closures. May require specialty cable or wiring. Up to three output modules can be installed on one wheelchair.</i>	1,344.00	

Seating System

STANDARD SEATING SYSTEM:

Part Number	Description	Price	HCPCS Code
<input checked="" type="checkbox"/> IKOALA-MFII	Mini-Flex Seat w/ Power Tilt <i>Single actuator elevates into tilt (up to 25°) to allow seating system to clear base. Seating system also includes an adjustable seat frame, adjustable backrest and legrest angle, positioning belt and one-piece footplate.</i>	8,627.00	E1002

Seating System Options

CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10418	Seat Function Control Switch for R-net <i>Does not include thru drive control operation of power seat functions.</i>	No Charge	
<input type="checkbox"/> I10419	Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control. Order with R-net electronics or remote joystick options.</i>	2,009.00	E2310

CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104135-99-0	Push Buttons on ICS Switchbox (3m)	No Charge	
<input type="checkbox"/> I104136-99-0	Toggle Switches on ICS Switchbox (3m)	340.00	

Chairman Koala Mini-Flex

Seating System Options

CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105409-99-0	Do Not Install ICS Switchbox, Send as Part	No Charge	

CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge	
<input type="checkbox"/> I103471-99-0	Mount ICS Switchbox onto Armrest - VS, Koala <i>ICS Switchbox will be mounted to the end of the armrest pad. If this item is selected, the ICS Switchbox will NOT swing away with the Joystick Module.</i>	No Charge	
<input type="checkbox"/> I106256-99-0	Mount ICS Switchbox on Backrest <i>Mount on Left: _____ Mount on Right: _____</i>	163.00	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Seat function interface box that accommodates up to 8 momentary switches. Provides an additional method for accessing the seat functions. 1/8" mini-mono switches not included.</i>	1,329.00	
<input type="checkbox"/> I106743-99-0	Remote Stop System - All Electronics <i>Consists of a 315MHz Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Custom charges may be required when ordering with a PS/PS Jr seat, a low backrest height or a pediatric wheelchair. If ordering on a K450, push handles must also be selected.</i>	1,315.00	E2321

Backrest Options

CHOOSE BACKREST TYPE:

Part Number	Description	Price	HCPCS Code
BACKREST NOTE: <i>Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Backrest cushions are provided with the MF Backrests only.</i>			
<input type="checkbox"/> I103185-99-0	MF Backrest 10"W x 14"H	No Charge	
<input type="checkbox"/> I103186-99-0	MF Backrest 12"W x 16"H	No Charge	
<input type="checkbox"/> I103187-99-0	MF Backrest 14"W x 20"H	No Charge	
<input type="checkbox"/> I10113	Tubular Backrest Frame Mini-Flex, 12" Wide <i>Backrest canes to mount aftermarket backrest systems.</i>	No Charge	
<input type="checkbox"/> I10114	Tubular Backrest Frame Mini-Flex, 14" Wide <i>Backrest canes to mount aftermarket backrest systems.</i>	No Charge	

CHOOSE BACKREST CUSHION:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104309-99-0	MX/MF Backrest Cushion, 10"W x 14"T - Leatherette	No Charge	
<input type="checkbox"/> I104311-99-0	MX/MF Backrest Cushion, 12"W x 16"T - Leatherette	No Charge	
<input type="checkbox"/> I104313-99-0	MX/MF Backrest Cushion, 14"W x 20"T - Leatherette	No Charge	
<input type="checkbox"/> I104308-99-0	MX/MF Backrest Cushion, 10"W x 14"T - Black Mesh	No Charge	
<input type="checkbox"/> I104310-99-0	MX/MF Backrest Cushion, 12"W x 16"T - Black Mesh	No Charge	
<input type="checkbox"/> I104312-99-0	MX/MF Backrest Cushion, 14"W x 20"T - Black Mesh	No Charge	

Seat Sizing

CHOOSE SEAT TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103188-99-0	MF Seat 10"W x 10"D	No Charge	
<input type="checkbox"/> I103189-99-0	MF Seat 12"W x 12"D	No Charge	
<input type="checkbox"/> I103190-99-0	MF Seat 14"W x 14"D	No Charge	
<input type="checkbox"/> I10115	Tubular Seat Frame Mini-Flex, 12" Wide <i>Seat frame tubes to mount aftermarket seating systems.</i>	No Charge	

Chairman Koala Mini-Flex

Seat Sizing

CHOOSE SEAT TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10116	Tubular Seat Frame Mini-Flex, 14" Wide <i>Seat frame tubes to mount aftermarket seating systems.</i>	No Charge	

CHOOSE SEAT CUSHION:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104145-99-0	Leatherette Seat Cushion MX/MF - 10"W x 10"D	621.00	
<input type="checkbox"/> I104146-99-0	Leatherette Seat Cushion MX/MF - 12"W x 12"D	621.00	
<input type="checkbox"/> I104147-99-0	Leatherette Seat Cushion MX/MF - 14"W x 14"D	621.00	
<input type="checkbox"/> I104305-99-0	Black Mesh Seat Cushion MX/MF - 10"W x 10"D	625.00	
<input type="checkbox"/> I104306-99-0	Black Mesh Seat Cushion MX/MF - 12"W x 12"D	625.00	
<input type="checkbox"/> I104307-99-0	Black Mesh Seat Cushion MX/MF - 14"W x 14"D	625.00	

Armrest Options

CHOOSE LEFT ARMREST PAD:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I100770-99-0	Armrest Left, 8" Leatherette	No Charge	
<input type="checkbox"/> I100774-99-0	Armrest Left, 10" Leatherette	No Charge	

CHOOSE RIGHT ARMREST PAD:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I100769-99-0	Armrest Right, 8" Leatherette	No Charge	
<input type="checkbox"/> I100773-99-0	Armrest Right, 10" Leatherette	No Charge	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10137	Height Adjustable Armrest Assembly	360.00	E0973
<input type="checkbox"/> I10016	Height Adjustable Armrest Assembly - No Charge <i>No charge with power tilt and/or recline.</i>	No Charge	

Legrest Options

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I101088-99-0	Footplate, Small w/ 5-12" Pan to Ftplate Adj. <i>Footplate is 10"W x 6.75"D.</i>	No Charge	
<input type="checkbox"/> I101089-99-0	Footplate, Large w/ 5-12" Pan to Ftplate Adj. <i>Footplate is 12"W x 7.875"D.</i>	No Charge	
<input type="checkbox"/> I100599-99-0	Foot Strap Kit	102.00	E0951,E0952
<input type="checkbox"/> I107503-99-0	BP Ankle Huggers, Med Side Release, Adj (Pair) <i>This item will be sent as a part with the wheelchair. It will not be installed.</i>	300.00	

Accessories

<input type="checkbox"/> I100518-99-0	Positioning Belt	No Charge	
<input type="checkbox"/> I100519-99-0	Positioning Belt Retractable	316.00	

HEADREST OPTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I100749-99-0	MF Headrest for Koala/PS Jr - Leatherette <i>Must also select PN I10552.</i>	356.00	E0955
<input type="checkbox"/> I10552	Adj/Removable Headrest Hardware	291.00	E1028
<input type="checkbox"/> I107499-99-0	Universal Headrest Adapter <i>Please Note: Recline is limited to 160° or less (depending on backrest height) when ordering a Universal Headrest Adapter on a 23" or 24" tall Corpus backrest.</i>	287.00	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I100753-99-0	Thigh Support Kit Leatherette <i>Must also select PN I10550.</i>	328.00	E0953

Chairman Koala Mini-Flex

Accessories

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I100751-99-0	Lateral Supports Leatherette <i>Must also select PN I10551.</i>	311.00	E0956
<input type="checkbox"/> I10550	Adj Removable - Thigh Support Hardware (Pair)	573.00	E1028
<input type="checkbox"/> I10551	Adj Removable - Trunk Support Hardware (Pair)	576.00	E1028
<input type="checkbox"/> I100867-99-0	Knee Block Leatherette	583.00	K0108
<input type="checkbox"/> I107100-99-0	RAM X-Grip 10" Tablet Holder w/Swing Away Mount <i>Compatible with most 10" screen tablets with or without a protective case. Includes adjustable swing away mount (I103476-99-0). Customization may be required if ordered with an Omni.</i>	444.00	
<input type="checkbox"/> I107222-99-0	USB Charger <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets. The charger features an isolated circuit design for improved safety; and by pulling power from the ICS, prevents battery drainage when the chair is powered off. Includes mounting hardware and ICS connector. No charge with 4 power seat function configurations. Mount on left armrest: _____ Mount on right armrest: _____</i>	220.00	
<input type="checkbox"/> I105805-99-0	Medical Necessities Bag Hooks <i>Please Note: Recline is limited to 165° or less (depending on backrest height) when ordered with 23" or 24" tall Ergo backrest. Custom required if ordered with a ROHO AGILITY backrest. May ship as a part at shorter backrest heights.</i>	105.00	
<input type="checkbox"/> I108276-99-0	Bodypoint PivotFit Shoulder Harness - Medium	275.00	E0960
<input type="checkbox"/> I100578-99-0	Small Angle Adjustable Tray <i>Must also select PN I10556.</i>	322.00	E0950
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each)	291.00	

Alternative Drive Systems

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10751	Mini-Joystick for Omni - (9-pin connector) <i>Operating Force = 0.5 oz (13 g) Neutral to Full Deflection = 0.2" (5mm)</i>	3,447.00	E2312
<input type="checkbox"/> I103893-99-0	Compact Joystick for Omni - (9-pin connector) <i>Operating Force = 7.8 oz (220 g) Neutral to Full Deflection = 1.1" (28mm) Joystick includes two 1/8" mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,686.00	E2373
<input type="checkbox"/> I103894-99-0	Compact Joystick LITE for Omni - (9-pin connector) <i>Operating Force = 1.8 oz (50 g) Neutral to Full Deflection = 0.7" (17mm) Joystick includes two 1/8" mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,907.00	E2373

ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104690-99-0	Alt JS Swing Away Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	295.00	E1028
<input type="checkbox"/> I107135-99-0	Alt JS / Omni Combo Swing Away Mount Kit - Armrest <i>Used to mount Omni and Alternative Proportional Joystick to the same armrest. Includes Swing Away Omni/Tablet Mount, with a longer mounting rod, two PERMOfix modules, one 3" rod and one 4" rod. Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	558.00	E1028
<input type="checkbox"/> I104143-99-0	Alt JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	920.00	E1028
<input type="checkbox"/> I10752	Haims Harness for Chin Control Mount, Large <i>For neck width up to 7½".</i>	654.00	E1028
<input type="checkbox"/> I10753	Haims Harness for Chin Control Mount, Small <i>For neck width up to 5½".</i>	654.00	E1028

ALT PROPORTIONAL JOYSTICK ACCESSORIES:

Part Number	Description	Price	HCPCS Code
-------------	-------------	-------	------------

Chairman Koala Mini-Flex

Alternative Drive Systems

ALT PROPORTIONAL JOYSTICK ACCESSORIES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104904-99-0	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid arms. Switches set-up for Mode & On/Off function. Switch operation force = 4.5oz (130 g)</i> Mount at front of Joystick: _____ Mount at rear of Joystick: _____ Mount at left of Joystick: _____ Mount at right of Joystick: _____	274.00	
<input type="checkbox"/> I104641-99-0	Dual Switches for Compact Joystick -Gooseneck Mount <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches set-up for Mode & On/Off function. Switch operation force = 4.5oz (130 g)</i>	493.00	

SIP & PUFF SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10573	Sip & Puff System for Omni <i>Requires Omni Control Module. Includes breath tube kit.</i>	1,844.00	E2325,E2326
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Includes breath tube kit.</i>	2,454.00	E2325,E2326
<input type="checkbox"/> I106807-99-0	Sip & Puff Kit for Total Control Head Array <i>For use with Total Control Head Array. Requires Omni Control Module. Includes breath tube kit.</i>	951.00	E2325,E2326
<input type="checkbox"/> I10762	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573, I10574 or I106807-99-0.</i>	295.00	

Head Array Options

TOTAL CONTROL HEAD ARRAY:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I106805-99-0	Permobil Head Array - 3-Switch System <i>Features an occipital pad sensor and two lateral arms configured with proximity sensors. The arms are mounted to the upper attachment points on the head array. Must select hardware below to complete assembly. ESP recommended. Choose Omni Control Module separately.</i>	4,589.00	E2330
<input type="checkbox"/> I11575	Swing Away Hardware for Side Sensor - Left	250.00	E1028
<input type="checkbox"/> I11576	Swing Away Hardware for Side Sensor - Right	250.00	E1028
<input type="checkbox"/> I11072	Program Head Array as a 5-Switch System <i>Select additional arms/sensors below and indicate the intended use of each sensor in the Notes & Comments section. Choose Omni Control Module separately.</i>	No Charge	
CHOOSE HEAD ARRAY MOUNT HARDWARE (REQUIRED):			
<input type="checkbox"/> I105927-99-0	Stealth TWB480 Backrest Mounting Hdwr with Links	305.00	E1028
<input type="checkbox"/> I105941-99-0	Stealth TWB480 Backrest Mnt Hdwr with Flipback	492.00	E1028
SELECT ADDITIONAL LEFT ARM PACKAGE (OPTIONAL):			
<input type="checkbox"/> I10976	Add Left Lateral Arm w Proximity Sensor	655.00	E1028
<input type="checkbox"/> I10977	Add Left Lateral Arm w Egg Switch	536.00	E1028
<input type="checkbox"/> I10978	Add Left Lateral Arm w Pad (No Sensor or Switch)	437.00	E1028
SELECT ADDITIONAL RIGHT ARM PACKAGE (OPTIONAL):			
<input type="checkbox"/> I10979	Add Right Lateral Arm w Proximity Sensor	655.00	E1028
<input type="checkbox"/> I10980	Add Right Lateral Arm w Egg Switch	536.00	E1028
<input type="checkbox"/> I10981	Add Right Lateral Arm w Pad (No Sensor or Switch)	437.00	E1028

Alt Drive Systems Setup

CHOOSE OMNI (REQUIRED):

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I108609-99-0	Omni2 Control Module for R-net - Koala <i>Two-section design features 3.5" VGA high resolution display, USB charger port, built in IR control and integrated Bluetooth technology for control of up to 4 Bluetooth devices. Must also select mounting hardware below. Please specify which type of alternative drive input device you will be using. Input Device: _____</i>	1,962.00	
<input type="checkbox"/> I103476-99-0	Adjustable Swing-Away Hardware - Omni/Tablet Mount in front of Left Arm: _____ Mount in front of Right Arm: _____	295.00	

Chairman Koala Mini-Flex

Alt Drive Systems Setup

CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10080	Mount Seat Control Switchbox on Backrest	No Charge	

Alt Drive Accessories

PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105907-99-0	Proximity Sensor + Swivel Assembly <i>Only available with the Total Control Head Array. Swivel assembly allows sensor to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array. Since Proximity Sensors require power, they CANNOT be used as On/Off switches.</i>	378.00	E1028
<input type="checkbox"/> I105919-99-0	Egg Switch + Swivel Mount Assembly <i>Swivel assembly allows switch to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	255.00	
<input type="checkbox"/> I105912-99-0	Solid Oval Pad + Swivel Assembly <i>Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	154.00	
<input type="checkbox"/> I106675-99-0	Oval Plate + Swivel Assembly - No Switch <i>Blank plate provides mounting option for Micro Lite, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	135.00	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I101769-99-0	Microlight Switch Qty: _____	206.00	
<input type="checkbox"/> I101767-99-0	Buddy Button Qty: _____	174.00	
<input type="checkbox"/> I106293-99-0	Egg Switch Black Qty: _____	174.00	
<input type="checkbox"/> I101771-99-0	Mini Cup Switch Qty: _____	174.00	
<input type="checkbox"/> I107500-99-0	Proximity Switch Qty: _____	688.00	
<input type="checkbox"/> I107495-99-0	Mini Joystick 5 Switch OC for Omni - (9-pin)	747.00	
<input type="checkbox"/> I10772	Dome Switch, Backrest Mnt w Gooseneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____ Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 6" rod, one 8" rod and one 9" gooseneck switch with 50" cable and mounting clamp.</i>	545.00	
<input type="checkbox"/> I10846	Dome Switch, Backrest Swg Awy Mnt w Gsneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____ w 90° Swing-Away Function. Mounts at top of back. Custom charge applies if ordered with short back. Kit Includes: One PERMOfix joint with Flat Plate, one PERMOfix joint with 45° swivel adjustment, two 6" rods, one 90° Swing-Away joint, one 9" gooseneck switch w 50" cable and mounting clamp.</i>	806.00	
<input type="checkbox"/> I105935-99-0	3.5mm Ext. Cable for Sensor/Switch, 0.3m Long	34.00	
<input type="checkbox"/> I105931-99-0	3.5mm Ext. Cable for Sensor/Switch, 1m Long	37.00	
<input type="checkbox"/> I105932-99-0	3.5mm Ext. Cable for Sensor/Switch, 2m Long	44.00	
<input type="checkbox"/> I105933-99-0	3.5mm Ext. Cable for Sensor/Switch, 3m Long	48.00	
<input type="checkbox"/> I105934-99-0	3.5mm Y-splitter Cable for Sensor/Switch <i>Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.</i>	37.00	

Additional Order Instructions (for Permobil):

Notes & Comments (not for Permobil):

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, Inc. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. *The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.*