

M300 Corpus HD

permobil



Dealer Information

* Contact: _____

Dealer Code: _____

* Dealer Name: _____

Address: _____

* City: _____

* Province/Postal Code: _____

* Phone#: _____

* Fax#: _____

PO#: _____

Email Address: _____

Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

* First Name: _____

* Last Name: _____

Diagnosis: _____

Funding Source: _____

Client Age: _____

Client Measurements

* Weight: _____

* Height: _____

A) Top of Shoulders: _____

B) Chest Depth: _____

C) Chest Width: _____

D) Seat Depth: _____


E) Top of Head: _____

F) Elbow to Hand: _____

G) Seat Pan to Elbow: _____

H) Hip Width: _____

*I) Knee to Foot: _____



Three column pricing on order/quote:

Please send order/quote to fax#: (800) 231-3256

Email to: sales@permobil.com

Permobil Ltd.
12 - 75 Mary St.
Aurora, Ontario
L4G 1G3
Tel: (800) 736-0925
Fax: (800) 231-3256
www.permobil.ca

Prices effective **January 1, 2022** and payable in Canadian dollars.

Chair Model and Colors

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10908	M300 Base Corpus HD - SPO <i>Price Includes: M300 Power HD Base Set Up For One Power Seat Function, Corpus HD Seat Frame, Shock Absorbing Suspension System, 8" Flat-Free Casters, Tie-Down Hardware for Strap Systems, Charger, Flat Free or Pneumatic Tires and a Standard Positioning Belt. R-net 120amp Electronics. Max Speed = 5.0mph, Weight Capacity = 450lbs.</i>	12,452.00
<input type="checkbox"/> I105799-99-0	M300 Base Corpus HD - MPO <i>Price Includes: M300 HD Power Base Set Up For Multiple Power Seat Functions, Corpus HD Seat Frame, Shock Absorbing Suspension System, 8" Flat-Free Casters, Tie-Down Hardware for Strap Systems, Charger, Flat Free or Pneumatic Tires and a Standard Positioning Belt. R-net 120amp Electronics. Max Speed = 5.0mph, Weight Capacity = 450lbs.</i>	12,452.00

CHOOSE WHEEL PACKAGE:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I11081	Grey Wheel Package - MWD <i>Includes grey tires on 8" silver casters and drive wheel rims. Hub cap color is selected below. Please select drive wheel type: Flat free drive tires: _____ Pneumatic drive tires: _____</i>	No Charge
<input type="checkbox"/> I11082	Black Wheel Package - MWD <i>Includes black tires on 8" silver casters and drive wheel rims. Hub cap color is selected below. Please select drive wheel type: Flat free drive tires: _____ Pneumatic drive tires: _____</i>	No Charge
<input type="checkbox"/> I11083	Blackout Wheel Package - MWD <i>Custom Option. Includes black tires and custom matte black 8" casters, drive wheels and hub caps. Please select drive wheel type: Flat free drive tires: _____ Pneumatic drive tires: _____</i>	437.00

OPTIONAL DRIVE TIRES:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10782	Tri-Spoke Split Rim w/ Snow Tire Pneumatic 3x8" <i>Snow tires are designed for use in cold weather conditions. The nylon compound in these tires remains flexible at low temperatures, but wears more easily. Therefore, the snow tires (sold as wheel assemblies) are not recommended for year-long use and must be ordered along with standard drive wheels. Not available with black tires or black rims.</i>	567.00

CHOOSE SHROUD COLOR:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105885-99-0	M300 HD Shroud - Onyx Black	No Charge
<input type="checkbox"/> I105881-99-0	M300 HD Shroud - Platinum Silver	No Charge
<input type="checkbox"/> I105887-99-0	M300 HD Shroud - Sublime Green	No Charge
<input type="checkbox"/> I106169-99-0	M300 HD Shroud - Cobalt Blue	No Charge
<input type="checkbox"/> I106171-99-0	M300 HD Shroud - Phantom Grey	No Charge
<input type="checkbox"/> I106159-99-0	M300 HD Shroud - Pop Star Pink	No Charge
<input type="checkbox"/> I106170-99-0	M300 HD Shroud - Purple Rain	No Charge
<input type="checkbox"/> I106160-99-0	M300 HD Shroud - Volcano Red	No Charge

CHOOSE HUBCAP COLOR:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10871	Match Hubcap Color to Shroud: _____ <i>Note: Hub caps will be matched to shroud color if no selection is made.</i> Solid Color Hub Cap - Pair <i>Platinum Silver: _____ Onyx Black: _____ Cobalt Blue: _____ Phantom Gray: _____ Pop Star Pink: _____ Volcano Red: _____ Purple Rain: _____ Titan Blue: _____ Sublime Green: _____ Crystal Blue: _____ Vintage Wine: _____ Extreme Yellow: _____ Chili Pepper Red: _____</i>	No Charge

Base Options

MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

Part Number	Description	MSRP - CA
-------------	-------------	-----------

Base Options

MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105875-99-0	HD Power Adjustable Seat Height - 8" Travel <i>Please Note: Drive is restricted when seat is in an elevated position.</i>	4,031.00
<input type="checkbox"/> I105880-99-0	HD Seat Tube Fixed	No Charge

MUST CHOOSE BATTERIES:

Part Number	Description	MSRP - CA
BATTERY OPTIONS:		
<input type="checkbox"/> I108999-99-0	Permobil Batteries, Grp 24 AGM GEL, Pair <i>Includes Permobil VoltPro charger optimized for Permobil batteries.</i>	1,013.00
<input type="checkbox"/> I107032-99-0	Batteries, Grp 34 Gel, Pair <i>Includes Permobil VoltPro charger for non-Permobil gel batteries.</i>	1,001.00

CHOOSE SEAT PAN TO FLOOR HEIGHT:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105878-99-0	Set Seat to Floor Height at 18"	No Charge
<input type="checkbox"/> I105877-99-0	Set Seat to Floor Height at 19"	No Charge
<input type="checkbox"/> I105876-99-0	Set Seat to Floor Height at 20"	No Charge

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105428-99-0	Enhanced Steering Perform Unit (ESP), R-net <i>Recommended item with switched drive inputs for improved tracking.</i>	1,284.00
<input type="checkbox"/> I105297-99-0	LED Light Kit for M300 / M400 w/ R-net Electronics	1,247.00
<input type="checkbox"/> I10715	Compact Laptop Ventilator Tray - Corpus Seat <i>Mounts to rear of Corpus backrest hinge. Accommodates Pulmonetic Systems / CareFusion LTV Series and similar sized Laptop Ventilators. The tray is 330mm Tall, adjustable Width between 200mm - 330mm, and adjustable Depth between 50mm - 165mm. Overall backrest angle limited to 150° when mounting this ventilator tray. Custom Option. Additional drive restrictions may apply.</i>	1,819.00
<input type="checkbox"/> I10787	Battery Tray for Laptop Ventilator Tray <i>Mounts to rear of Laptop Ventilator Tray. Measures 290mm W x 170mm H x 42mm D. Accommodates Pulmonetic Systems Lithium-Ion Battery Pack, CareFusion SprintPack Battery System, or Richardson PowerTech Vent Power Center. Custom Option.</i>	291.00
<input type="checkbox"/> I10788	O2 Holder for Laptop Ventilator Tray <i>Mounts to Left or Right Side of Laptop Vent Tray. O2 Holder is 175mm Tall. Width infinitely adjustable 90mm - 170mm. Depth infinitely adjustable 95mm - 125mm. Custom Option.</i>	411.00
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	567.00

Joystick Options

CHOOSE JOYSTICK TYPE:

Part Number	Description	MSRP - CA
PLEASE NOTE: <i>For all remote joystick modules shown below, you must also select an expandable controller. For additional input devices, refer to Alternative Drive Systems Section.</i>		
<input type="checkbox"/> I107779-99-0	Color Joystick Module, (CJSM) <i>This R-net color joystick features: A Large High-Resolution 3.5" VGA color display, control of up to 4 Bluetooth devices, Infrared control, programmable buttons and a programmable switch jack; all in a compact, durable Aluminum housing.</i>	1,444.00
<input type="checkbox"/> I11186	Permobil Joystick Module, NC <i>Designed exclusively for Permobil, this R-net color joystick features: a 2.8" color display.</i>	No Charge
<input type="checkbox"/> I107829-99-0	3.5mm Y-Splitter Cable for PJSM/CJSM <i>For use with the PJSM or CJSM Mode/Profile jack. Allows 2 mechanical switches to control 4 separate functions. Requires PC Programming to assign the short hold and long hold commands to each switch. Switches sold separately.</i>	45.00

Joystick Options

CHOOSE JOYSTICK TYPE:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105246-99-0	R-net LED Remote Joystick - No Display <i>This Joystick does NOT have a Display Screen, does NOT contain mono jack ports, and does NOT include Bluetooth or Assignable Buttons. PC Programming required and profiles limited to five.</i>	No Charge
<input type="checkbox"/> I105700-99-0	R-net Remote Colour JS w/ Lights, Heavy Duty Gimbal <i>Heavy-Duty modification enables gimbal to hold up to extreme operating forces. Reinforced housing is attached to a standard Colour LCD R-net Joystick. 3" long joystick knob has 3/4" diameter. Not compatible with alternative joystick handles.</i>	2,208.00

CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard	No Charge
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	129.00
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	190.00
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	129.00
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	173.00
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	173.00
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	173.00
<input type="checkbox"/> I105442-99-0	1.75" Black Foam Ball For Joystick Knob	129.00
<input type="checkbox"/> I105846-99-0	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	166.00
<input type="checkbox"/> I105847-99-0	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	166.00
<input type="checkbox"/> I105848-99-0	Bodypoint J/S Handle Dome Shaped, Rubber	166.00

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10664	Joystick Mount Right - Fixed	No Charge
<input type="checkbox"/> I10665	Joystick Mount Left - Fixed	No Charge
<input type="checkbox"/> I108037-99-0	Slimline Retractable Joystick Mount - Right	489.00
<input type="checkbox"/> I108038-99-0	Slimline Retractable Joystick Mount - Left	489.00
<input type="checkbox"/> I11232	Adjustable Height Panel Bracket <i>Provides 1.77" of vertical adjustment below the armrest and 25° of angle adjustment. Compatible with all Permobil joystick mounts except for the M1 Fixed Joystick Mount.</i>	379.00
<input type="checkbox"/> I105210-99-0	Mount Standard Joystick on Corpus Back Hinge <i>Standard R-net Joystick is mounted to backrest hinge. Can be chosen instead of PN I105408-99-0. Mount on Left: _____ Mount on Right: _____</i>	339.00
<input type="checkbox"/> I108268-99-0	Bodypoint Tri-Lock Midline Mount Assembly <i>Includes L-shaped joystick arm and height adjustable mounting clamp. Joystick mount sold separately (select from the mounts listed below). Please indicate if you want midline mounted on left or right armrest: Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	747.00
<input type="checkbox"/> I108280-99-0	Bodypoint Tri-Lock Mount for R-net JS	306.00
<input type="checkbox"/> I108279-99-0	Bodypoint Tri-Lock Mount for Compact Joystick	306.00
<input type="checkbox"/> I10842	PERMOFix Mount for Mini/Compact Joystick <i>Includes a shim that allows the PERMOFix joint to be mounted on the joystick arm. Provides an adjustable mount for a Mini Joystick or Compact Joystick.</i>	306.00

Part Number	Description	MSRP - CA
<input checked="" type="checkbox"/> I10411	Expandable Controller - R-net <i>Must be ordered with all expandable/remote joysticks.</i>	No Charge
<input checked="" type="checkbox"/> I10513	Harness for Expandable Controller <i>Must be ordered with any Expandable Controller.</i>	No Charge
<input type="checkbox"/> I105225-99-0	Co-Pilot Attendant Control R-net - Corpus	5,444.00
<input type="checkbox"/> I105408-99-0	Attendant Control R-net	1,413.00

Joystick Options

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10582	Output Module (IOM) for R-net Control System <i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). Each module can provide up to 4 switch closures. May require specialty cable or wiring. Up to three output modules can be installed on one wheelchair.</i>	1,891.00

Seating System

CHOOSE SEATING SYSTEM:

Part Number	Description	MSRP - CA
	PLEASE NOTE: <i>Weight Limit 450lbs. For the Corpus HD Seat on the M300, there is a half speed drive limit when the Total Backrest Angle (TBA) > 135° and a drive inhibit at 150° of TBA. There is also a drive inhibit at 20° of tilt.</i>	
<input type="checkbox"/> I10911	HD Seat w/ No Tilt & 150° Power Recline	10,372.00
<input type="checkbox"/> I10912	HD Seat w/ 45° DuraCore Tilt & Manual Recline	10,017.00
<input type="checkbox"/> I10913	HD Seat w/ 45° DuraCore Tilt & 150° Power Recline	15,726.00

Seating System Options

CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10418	Single Seat Function Control Switch for R-net <i>Does not include thru drive control operation of power seat functions.</i>	No Charge
<input type="checkbox"/> I10419	Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control. Order with R-net electronics or remote joystick options.</i>	1,797.00
<input type="checkbox"/> I10585	Multiple Seat Function Control Kit for R-net <i>For operation of two or more seat functions through the wheelchair electronics via the drive control. Order with R-net electronics or remote joystick options.</i>	2,994.00

CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	MSRP - CA
	PLEASE NOTE: <i>ICS Switchbox is only available with chairs that include an R-Net LED Remote Joystick.</i>	
<input type="checkbox"/> I102798-99-0	Push Buttons on Seat Function Control Kit	No Charge
<input type="checkbox"/> I102799-99-0	Toggle Switches on Seat Function Control Kit	320.00
<input type="checkbox"/> I105409-99-0	Do Not Install ICS Switchbox, Send as Part	No Charge

CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge
<input type="checkbox"/> I104651-99-0	Mount ICS Switchbox onto Armrest <i>Mount on Left: _____ Mount on Right: _____</i>	No Charge
<input type="checkbox"/> I106256-99-0	Mount ICS Switchbox on Backrest <i>Mount on Left: _____ Mount on Right: _____</i>	180.00

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Seat function interface box that accommodates up to 8 momentary switches. Provides an additional method for accessing the seat functions. 1/8" mini-mono switches not included.</i>	1,413.00
<input type="checkbox"/> I106743-99-0	Remote Stop System - All Electronics <i>Consists of a 315MHz Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Custom charges may be required when ordering with a PS/PS Jr seat, a low backrest height or a pediatric wheelchair. If ordering on a K450, push handles must also be selected.</i>	1,348.00

Seating System Options

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10750	Memory Seat Program, Corpus w/ R-net	131.00
<input type="checkbox"/> I106209-99-0	Legrest Movement Linked to Recline Electronically <i>Power Legrest will raise when Power Recline moves backward and will lower when back moves forward. If Power Legrest is operated, Power Recline will NOT automatically reposition. Feature works through ICS Switchbox and Driver Control.</i>	53.00
<input type="checkbox"/> I10898	Enable Independent Repositioning Mode <i>IRM allows for sequential movement of power seat functions into a posterior position (35° tilt, 115° recline, 135° legrest) and returns to upright position (5° tilt, 100° recline, 110° legrest). This can be accessed using a standard drive control or a single switch. On the Permobil Joystick Module, IRM is accessible through the mode/profile jack (long switch hold). Please see the Alt Drive Accessories section for switch options.</i>	No Charge

Backrest Options

CHOOSE FIXED OR SLIDING BACKREST KIT:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105871-99-0	Sliding Backrest Kit - 5" Mechanical Shear	No Charge
<input type="checkbox"/> I105869-99-0	Fixed Backrest Kit - No Mechanical Shear	42.00

CHOOSE BACKREST TYPE:

Part Number	Description	MSRP - CA
BACKREST HEIGHT NOTE: <i>Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°.</i>		
<input type="checkbox"/> I106016-99-0	Corpus HD Ergo Backrest - 20"W x 21"T (C3GERGBCA)	1,251.00
<input type="checkbox"/> I106017-99-0	Corpus HD Ergo Backrest - 22"W x 21"T	1,251.00
<input type="checkbox"/> I106018-99-0	Corpus HD Ergo Backrest - 24"W x 21"T	1,251.00
PLEASE NOTE: <i>Taller backrests listed below include: foam lumbar pads*, a UT rail mounted at the top for lateral adjustment of the headrest, and split backrest feature with height adjustment in 1" increments. Please indicate below where you want the height set at so that a properly sized backrest cushion can be provided. * Foam lumbar pads not included with Solid Backs.</i>		
<input type="checkbox"/> I106019-99-0	Corpus HD Ergo Back -20"W x 24 - 29"T (C3GERGBCA)	1,251.00
<input type="checkbox"/> I106020-99-0	Corpus HD Ergo Back -22"W x 24 - 29"T	1,251.00
<input type="checkbox"/> I106021-99-0	Corpus HD Ergo Back -24"W x 24 - 29"T	1,251.00
	Set Backrest at 24" _____	
	Set Backrest at 25" _____	
	Set Backrest at 26" _____	
	Set Backrest at 27" _____	
	Set Backrest at 28" _____	
	Set Backrest at 29" _____	

CHOOSE BACKREST CUSHION:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10659	Leatherette Cushion - Corpus Ergo Back	No Charge
<input type="checkbox"/> I10660	Black Mesh Cushion - Corpus Ergo Back	60.00

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10725	Split Backrest Cushion - Corpus <i>Custom Option.</i>	384.00

Seat Sizing

CHOOSE SEAT WIDTH SETTING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105859-99-0	Set Seat Width at 19"	No Charge

Seat Sizing

CHOOSE SEAT WIDTH SETTING:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105860-99-0	Set Seat Width at 21"	No Charge
<input type="checkbox"/> I105861-99-0	Set Seat Width at 23"	No Charge
<input type="checkbox"/> I105862-99-0	Set Seat Width at 25"	No Charge

CHOOSE SEAT DEPTH:

Part Number	Description	MSRP - CA
PLEASE NOTE: <i>Odd numbered seat depths can be made 1" shorter without any additional parts. For example, a 21"D seat can also be set at 20"D.</i>		
<input type="checkbox"/> I105863-99-0	Corpus HD Seat Depth 19"	No Charge
<input type="checkbox"/> I105864-99-0	Corpus HD Seat Depth 20"	No Charge
<input type="checkbox"/> I105865-99-0	Corpus HD Seat Depth 21"	No Charge
<input type="checkbox"/> I105866-99-0	Corpus HD Seat Depth 22"	No Charge
<input type="checkbox"/> I105868-99-0	Corpus HD Seat Depth 23"	No Charge
<input type="checkbox"/> I105856-99-0	Corpus HD Seat Depth 24"	No Charge
<input type="checkbox"/> I105857-99-0	Corpus HD Seat Depth 25"	No Charge

CHOOSE SEAT CUSHION:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10651	Corpus Ergo Cushion - Leatherette (C3GERGSCA)	654.00
<input type="checkbox"/> I10652	Corpus Ergo Cushion - Black Fabric	654.00
<input type="checkbox"/> I11113	Incontinence Cover for Ergo Cushion <i>Made from a stretch Dartex fabric and suitable for all Ergo cushions including the Corpus Ergo Air. The cover is cut to the contours of the Ergo cushion and is designed to be used as an outer layer of protection. Cover size will be matched to cushion size. Sent as a part. Part Number will be provided on the quote / order confirmation based on the Seat Width and Seat Depth selected.</i>	208.00
<input type="checkbox"/> I11166	ROHO QUADTRO SELECT HIGH PROFILE	777.00
<input type="checkbox"/> I11168	ROHO HIGH PROFILE Single Compartment	636.00
<input type="checkbox"/> I11173	Incontinence Cover for ROHO Cushion <i>Cover size will be matched to ROHO cushion size. Sent as a part.</i>	145.00
<input type="checkbox"/> I105118-99-0	Omit Seat Cushion	No Charge

Armrest Options

CHOOSE ARMBAR MOUNTING POSITION:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I106022-99-0	Set Armbar at 20"W Between Armrests: _____	No Charge
<input type="checkbox"/> I106023-99-0	Set Armbar at 21"W Between Armrests: _____	No Charge
<input type="checkbox"/> I106024-99-0	Set Armbar at 22"W Between Armrests: _____	No Charge
<input type="checkbox"/> I106025-99-0	Set Armbar at 23"W Between Armrests: _____	No Charge
<input type="checkbox"/> I106026-99-0	Set Armbar at 24"W Between Armrests: _____	No Charge
<input type="checkbox"/> I106027-99-0	Set Armbar at 25"W Between Armrests: _____	No Charge

CHOOSE ARMREST HEIGHT MOUNTING POSITION:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10923	Height Adjustable Armrests	439.00
<input type="checkbox"/> I10924	Height Adjustable Armrests - NC	No Charge

SET ARMREST HEIGHT AT: _____

Corpus HD Armrests are infinitely adjustable from 8" - 16", measured from the seat pan to the top of the armrest pad. No additional parts are needed to get this full range of adjustment. In one inch increments (ex: 8", 9", etc.), please indicate where you would like the armrests set.

CHOOSE ARMREST PAD RIGHT:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105104-99-0	4"x13" Corpus Armrest Pad Right - Leatherette	No Charge

Armrest Options

CHOOSE ARMREST PAD RIGHT:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105997-99-0	4"x16" Corpus HD Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I105999-99-0	4"x18" Corpus HD Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I106047-99-0	4"x13" Corpus Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/> I105993-99-0	4"x16" Corpus HD Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/> I105995-99-0	4"x18" Corpus HD Armrest Pad Right - Black Fabric	No Charge
<input type="checkbox"/> I107593-99-0	4"x13" Corpus Arm w/ Permobil Gel Pad - Right	213.00
<input type="checkbox"/> I108246-99-0	4"x16" Corpus HD Arm w/ Permobil Gel Pad - Right	218.00
<input type="checkbox"/> I107597-99-0	4"x18" Corpus HD Arm w/ Permobil Gel Pad - Right	227.00
<input type="checkbox"/> I108861-99-0	Moldable Forearm Support - Right <i>11.3"L x 5"W. Not compatible with armrest pouch, trays (except Small Tray) or the VS Chest Support.</i>	360.00
<input type="checkbox"/> I108863-99-0	Contoured Hand Pad Medium - Right <i>6.4"L x 5.5" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.</i>	179.00
<input type="checkbox"/> I108865-99-0	Contoured Hand Pad Large - Right <i>7.7"L x 6.1" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.</i>	189.00

CHOOSE ARMREST PAD LEFT:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105102-99-0	4"x13" Corpus Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I105996-99-0	4"x16" Corpus HD Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I105998-99-0	4"x18" Corpus HD Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I106046-99-0	4"x13" Corpus Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/> I105992-99-0	4"x16" Corpus HD Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/> I105994-99-0	4"x18" Corpus HD Armrest Pad Left - Black Fabric	No Charge
<input type="checkbox"/> I107592-99-0	4"x13" Corpus Arm w/ Permobil Gel Pad - Left	213.00
<input type="checkbox"/> I108247-99-0	4"x16" Corpus HD Arm w/ Permobil Gel Pad - Left	218.00
<input type="checkbox"/> I107596-99-0	4"x18" Corpus HD Arm w/ Permobil Gel Pad - Left	227.00
<input type="checkbox"/> I108860-99-0	Moldable Forearm Support - Left <i>11.3"L x 5"W. Not compatible with armrest pouch, trays (except Small Tray) or the VS Chest Support.</i>	360.00
<input type="checkbox"/> I108862-99-0	Contoured Hand Pad Medium - Left <i>6.4"L x 5.5" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.</i>	179.00
<input type="checkbox"/> I108864-99-0	Contoured Hand Pad Large - Left <i>7.7"L x 6.1" W. For use with Moldable Forearm Supports only. Not compatible with Slimline or Fixed joystick mounts or Adjustable Height Panel Bracket.</i>	189.00

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I105463-99-0	Corpus Arm Trough Adapter 1" - Pair <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked when selecting this option.</i>	375.00
<input type="checkbox"/> I105462-99-0	Corpus Arm Trough Adpater 7/8" - Pair <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts. Please Note: Medial Swivel must be must be locked when selecting this option.</i> Install Arm Trough Adapters: _____ Send Arm Trough Adapters as Parts: _____	375.00
<input type="checkbox"/> I107339-99-0	Armrest Pouch - Right	273.00
<input type="checkbox"/> I107338-99-0	Armrest Pouch - Left	273.00
<input type="checkbox"/> I107840-99-0	UT Armrest Side Support - Right Outer <i>Includes 2"T x 6.5"W tapered pad.</i>	391.00

Armrest Options

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107863-99-0	UT Armrest Side Support - Left Outer <i>Includes 2"T x 6.5"W tapered pad.</i>	391.00
<input type="checkbox"/> I107864-99-0	UT Armrest Side Support - Right Inner <i>Includes 2"T x 6.5"W tapered pad. May need to increase width between the armrests to ensure armrests will flip-up.</i>	391.00
<input type="checkbox"/> I107865-99-0	UT Armrest Side Support - Left Inner <i>Includes 2"T x 6.5"W tapered pad. May need to increase width between the armrests to ensure armrests will flip-up.</i>	391.00
<input type="checkbox"/> I107839-99-0	UT Elbow Support - Right <i>Must also select PN I10881.</i>	210.00
<input type="checkbox"/> I107838-99-0	UT Elbow Support - Left <i>Must also select PN I10882.</i>	210.00
<input type="checkbox"/> I10881	Adjustable Elbow Support Hardware - Right	306.00
<input type="checkbox"/> I10882	Adjustable Elbow Support Hardware - Left	306.00

Legrest Options

CHOOSE LEGREST TYPE:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I106053-99-0	Corpus HD Manual Adj Legrest Elevation	No Charge
<input type="checkbox"/> I106052-99-0	Power Elevating Center Mount Legrests - 175° <i>Includes mechanical articulation.</i>	4,587.00

CHOOSE FOOTPLATE STYLE:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I106008-99-0	2pc Wide Corpus Footplates - 9"D x 7.5"W	No Charge
<input type="checkbox"/> I106009-99-0	2pc Extra Wide Corpus Footplates - 10"D x 9"W	308.00

FOOTPLATE HEIGHT POSITION:

Part Number	Description	MSRP - CA
	SET FOOTPLATE HEIGHT AT: _____ <i>This is measured from the SEAT PAN to the top of the FOOTPLATE, when the legrest assembly is at 90° DO NOT INCLUDE CUSHION THICKNESS IN YOUR MEASUREMENT. In one inch increments, please indicate desired footplate height. NOTE: Standard footplate height range is 6" - 18"; however, height may need to be reduced for shipping. Calf Pads are not available at heights less than 10".</i>	

CHOOSE LATERAL FOOTPLATE ROTATION OPTIONS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I106006-99-0	Lateral Footplate Rotation - Fixed	No Charge
<input type="checkbox"/> I106007-99-0	Lateral Footplate Rotation - Enabled	41.00

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I106190-99-0	VS Calf Support Kit (A) w Mesh Cover & Adj Hardware <i>5.5"W x 6"T Pad. Compatible with footplate heights of 10" - 18".</i>	No Charge
<input type="checkbox"/> I106188-99-0	VS Calf Support Kit (H) w Mesh Cover & Adj Hardware <i>7"W x 7"T Curved Pad. Compatible with footplate heights of 11" - 18".</i>	455.00
<input type="checkbox"/> I106186-99-0	VS Calf Support Kit (D) w Mesh Cover & Adj Hardware <i>8"W x 8"T Pad. Compatible with footplate heights of 12" - 18".</i>	455.00
<input type="checkbox"/> I106225-99-0	Stealth Gel Pads - Wide Footplates (Pair)	444.00
<input type="checkbox"/> I106224-99-0	Stealth Gel Pads - Extra Wide Footplates (Pair)	458.00

Accessories

<input type="checkbox"/> I104518-99-0	Bodypoint Non-Padded Hip Belt - Medium	No Charge
<input type="checkbox"/> I104520-99-0	Bodypoint Padded Hip Belt - Medium	155.00
<input type="checkbox"/> I10940	Bodypoint Non-Padded Hip Belt, UniTrack - Large	No Charge

Accessories

MUST CHOOSE POSITIONING BELT:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I104516-99-0	Bodypoint Padded Hip Belt - Large	213.00
<input type="checkbox"/> I106531-99-0	Bodypoint Padded Hip Belt - Extra Large	289.00
<input type="checkbox"/> I100519-99-0	Positioning Belt Retractable	322.00
<input type="checkbox"/> I100600-99-0	Positioning Belt Retractable Long	350.00
<input type="checkbox"/> I108267-99-0	Bodypoint Padded Hip Belt "Airplane" Latch - Med	273.00
<input type="checkbox"/> I108266-99-0	Bodypoint Padded Hip Belt "Airplane" Latch - Large	289.00

HEADREST OPTIONS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I107470-99-0	Permobil 7"W Headrest <i>Must also select headrest mounting hardware.</i>	295.00
<input type="checkbox"/> I107469-99-0	Permobil 10"W Headrest <i>Must also select headrest mounting hardware.</i>	295.00
<input type="checkbox"/> I107468-99-0	Permobil 14"W Headrest <i>Must also select headrest mounting hardware.</i>	295.00
<input type="checkbox"/> I107480-99-0	Adj/Removable Headrest Link Hardware - Standard	315.00
<input type="checkbox"/> I109178-99-0	Extra Link for Link Headrest Hardware <i>Used to extend range of standard link hardware. Sent as part.</i>	63.00
<input type="checkbox"/> I108721-99-0	Adj/Removable BodiLink Head Support Hdwr - MAX 8 <i>Must be selected when ordering more than 2 swing-away headrest accessories or chin boom assemblies.</i>	325.00
<input type="checkbox"/> I108722-99-0	BodiLink Headrest Bracket - Corpus 23-28 <i>For 23" - 28" Corpus Ergo backrest heights.</i>	148.00
<input type="checkbox"/> I108723-99-0	Bodilink Headrest Bracket Corpus Short/AGILITY <i>For 20" Corpus Ergo backrests and AGILITY Max Contour Back Systems.</i>	148.00
<input type="checkbox"/> I11305	Swing Away Egg Switch Assembly <i>Mount on Upper Left: _____ Mount on Upper Right: _____</i>	645.00
<input type="checkbox"/> I11306	Swing Away Microlight Assembly <i>Mount on Upper Left: _____ Mount on Upper Right: _____</i>	740.00
<input type="checkbox"/> I11307	Swing Away Facial Pad Assembly 4"L x 2"T <i>Mount on Upper Left: _____ Mount on Upper Right: _____</i>	525.00
<input type="checkbox"/> I11308	Swing Away Facial Pad Assembly 6"L x 2"T <i>Mount on Upper Left: _____ Mount on Upper Right: _____</i>	525.00
<input type="checkbox"/> I106331-99-0	Permobil 10"W Headrest, Leatherette <i>Not available with ROHO AGILITY Back Systems. Must also select PN I10552.</i>	672.00
<input type="checkbox"/> I10552	Adj/Removable Headrest Hardware	144.00
<input type="checkbox"/> I107499-99-0	Universal Headrest Adapter <i>Please Note: Recline is limited to 160° or less (depending on backrest height) when ordering a Universal Headrest Adapter on a 23" or 24" tall Corpus backrest.</i>	271.00

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I103070-99-0	UT Amp. Support Pad (H) w/ Mesh Cover - (1 Each) <i>7"W x 7"T Curved Pad. Must also select PN I10553.</i>	403.00
<input type="checkbox"/> I103071-99-0	UT Amp. Support Pad (D) w/ Mesh Cover - (1 Each) <i>8"W x 8"T Flat Pad. Must also select PN I10553.</i>	403.00
<input type="checkbox"/> I10553	Adj Removable - Amp Support Hardware (1 Each) <i>Please indicate if you would like the amputation support mounted on the right or left side of the seat frame: Mount Amp Suppt on Left: _____ Mount Amp Suppt on Right: _____</i>	144.00
<input type="checkbox"/> I109040-99-0	BodiLink Premium LPTS Pad (B), Stretch-Air - Pair <i>4.5"W x 3.5"T Pad w/additional foam insert and Stretch-Air cover. Must also select BodiLink LPTS Mounting Hardware.</i>	284.00

Accessories

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I109043-99-0	BodiLink Premium LPTS Pad (A), Stretch-Air - Pair <i>6"W x 5.5"T Pad w/additional foam insert and Stretch-Air cover. Must also select BodiLink LPTS Mounting Hardware.</i>	284.00
<input type="checkbox"/> I109041-99-0	BodiLink Premium LPTS Pad (C), Stretch-Air - Pair <i>8"W x 3.5"T Pad w/additional foam insert and Stretch-Air cover. Must also select BodiLink LPTS Mounting Hardware.</i>	284.00
<input type="checkbox"/> I109042-99-0	BodiLink Premium LPTS Pad (G), Stretch-Air - Pair <i>12"W x 3.5"T Pad w/additional foam insert and Stretch-Air cover. Must also select BodiLink LPTS Mounting Hardware.</i>	284.00
<input type="checkbox"/> I109044-99-0	Adj/Removable BodiLink LPTS Hdwr Med Corpus Pair <i>Medium TT Power Mount hardware can be extended upward to clear cushions of 3-5.5" thick.</i>	336.00
<input type="checkbox"/> I109194-99-0	Adj/Removable Lateral Thigh/Hip Hardware - Pair	144.00
<input type="checkbox"/> I105260-99-0	Transfer Handles, Corpus Seat	554.00
<input type="checkbox"/> I105292-99-0	Push Handles, Corpus Seat	No Charge
<input type="checkbox"/> I106210-99-0	RAM Self-Leveling Cup Holder	96.00
<input type="checkbox"/> I106806-99-0	RAM X-Grip Phone Holder	96.00
<input type="checkbox"/> I107100-99-0	RAM X-Grip 10" Tablet Holder w/Swing Away Mount <i>Compatible with most 10" screen tablets with or without a protective case. Includes adjustable swing away mount. Note: Customization may be required if ordered along with an Omni.</i>	529.00
<input type="checkbox"/> I107189-99-0	USB Charger <i>Provides up to 1.5A (at 5V) to charge most smartphones and tablets. The charger features an isolated circuit design for improved safety; and by pulling power from the ICS, prevents battery drainage when the chair is powered off. Includes mounting hardware and ICS connector. Mount on left armrest: _____ Mount on right armrest: _____</i>	218.00
<input type="checkbox"/> I105805-99-0	Medical Necessities Bag Hooks <i>Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with 23" or 24" tall Ergo backrest. Custom required if ordered with a ROHO AGILITY backrest. May ship as a part at shorter backrest heights.</i>	104.00
<input type="checkbox"/> I102059-99-0	Medical Necessities Bag <i>Includes hardware to attach bag to backrest. Please Note: Recline is limited to 165° or less (depending on backrest height) when ordering with a 23" or 24" tall Ergo backrest. Custom required if ordered with a ROHO AGILITY backrest.</i>	294.00
<input type="checkbox"/> I105840-99-0	Bodypoint Monoflex Medium Chest Support Belt Kit <i>18.25"L x 4"W Pad Size. Item will be shipped as a part.</i>	203.00
<input type="checkbox"/> I105839-99-0	Bodypoint Monoflex Large Chest Support Belt Kit <i>21.5"L x 4.5"W Pad Size. Item will be shipped as a part.</i>	203.00
<input type="checkbox"/> I105850-99-0	Bodypoint Monoflex X-Lrge Chest Support Belt Kit <i>24.5"L x 5"W Pad Size. Item will be shipped as a part.</i>	203.00
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each)	144.00
<input type="checkbox"/> I108227-99-0	Permobil Adjustable Tray - 17"W	362.00
<input type="checkbox"/> I108229-99-0	Permobil Adjustable Tray - 19"W	362.00
<input type="checkbox"/> I108232-99-0	Permobil Tray Hardware Right Side Adjustable	263.00
<input type="checkbox"/> I108233-99-0	Permobil Tray Hardware Left Side Adjustable	263.00
<input type="checkbox"/> I105224-99-0	Oxygen holder	358.00
<input type="checkbox"/> I105223-99-0	Crutch holder	372.00

Stealth Accessories

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I106221-99-0	Stealth 4" x 6" Lateral Sprrts, Original Hdwr, HD <i>Must also select I10723. Kit contains thin, Stealth hardware with swing away brackets and two height adjustment rails. The rails are bolted to the backrest to provide mounts for the lateral hardware. PLEASE NOTE: The spacing between the armrests must be at least 1" wider than the backrest width when ordering these laterals.</i>	487.00
<input type="checkbox"/> I10723	Adj Removable Stealth Lateral Hardware (Pair)	635.00

Alternative Drive Systems

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10751	Mini-Joystick for Omni - (9-pin connector) <i>Operating Force = 13 g Neutral to Full Deflection = 5mm</i>	3,278.00
<input type="checkbox"/> I103893-99-0	Compact Joystick for Omni - (9-pin connector) <i>Operating Force = 220 g Neutral to Full Deflection = 28mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,599.00
<input type="checkbox"/> I103894-99-0	Compact Joystick LITE for Omni - (9-pin connector) <i>Operating Force = 50 g Neutral to Full Deflection = 17mm Joystick includes two 3.5mm mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.75" (42mm) Foam Ball Knob.</i>	1,830.00

ALT PROPORTIONAL JOYSTICK MOUNTING HARDWARE:

Part Number	Description	MSRP - CA
	NOTE: <i>Please see the PERMOfix Mounting Solutions, within the Alt Drive Accessories section, for additional mounting options. The Swing Away Chin Boom Assemblies are Compatible with Permobil Headrests I107470-99-0, I107469-99-0 or I107468-99-0 and must be selected with I108721-99-0.</i>	
<input type="checkbox"/> I107946-99-0	Swing Away Chin Boom Assm for Compact Lite - Left	593.00
<input type="checkbox"/> I107948-99-0	Swing Away Chin Boom Assm for Compact Lite - Right	593.00
<input type="checkbox"/> I107857-99-0	Swing Away Chin Boom Assm for Mini Joystick - Left	593.00
<input type="checkbox"/> I107947-99-0	Swing Away Chin Boom Assm for Mini Joystick - Right	593.00
<input type="checkbox"/> I104690-99-0	Alt. JS Swing Away Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	311.00
<input type="checkbox"/> I107135-99-0	Alt JS/Omni Combo Swing Away Mount Kit - Armrest <i>Used to mount Omni and Alternative Proportional Joystick to the same armrest. Includes Swing Away Omni/Tablet Mount, with a longer mounting rod, two PERMOfix modules, one 3" rod and one 4" rod. Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	629.00
<input type="checkbox"/> I104143-99-0	Alt. JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	1,007.00
<input type="checkbox"/> I10752	Haims Harness for Chin Control Mount, Large <i>For neck width up to 19cm.</i>	719.00
<input type="checkbox"/> I10753	Haims Harness for Chin Control Mount, Small <i>For neck width up to 14cm.</i>	719.00

ALT PROPORTIONAL JOYSTICK ACCESSORIES:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I104904-99-0	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid arms. Switches set-up for Mode & On/Off function. Switch operation force = 130 g Mount at front of Joystick: _____ Mount at rear of Joystick: _____ Mount at left of Joystick: _____ Mount at right of Joystick: _____</i>	298.00
<input type="checkbox"/> I104641-99-0	Dual Switches for Compact Joystick - Gooseneck Mnt <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches set-up for Mode & On/Off function. Switch operation force = 130 g</i>	543.00

SIP & PUFF SYSTEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10573	Sip & Puff System for Omni <i>Requires Omni Control Module. Includes breath tube kit.</i>	1,923.00
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Includes breath tube kit.</i>	2,559.00
<input type="checkbox"/> I106807-99-0	Sip & Puff Kit for Total Control Head Array <i>For use with Total Control Head Array. Requires Omni Control Module. Includes breath tube kit.</i>	1,062.00

Alternative Drive Systems

SIP & PUFF SYSTEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10762	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573, I10574 or I106807-99-0.</i>	311.00

PROPORTIONAL FOOT CONTROL:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10767	Proportional Foot Control <i>Not available with Active Reach Packages. Mount on Left: _____ Mount on Right: _____</i>	3,376.00
<input type="checkbox"/> I10766	Adjustable Mounting Hrdw - Foot Control <i>Required when selecting PN I10767.</i>	311.00

Total Control Head Array

TOTAL CONTROL HEAD ARRAY:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I106805-99-0	Permobil Head Array - 3-Switch System <i>Features an occipital pad sensor and two lateral arms configured with proximity sensors. The arms are mounted to the upper attachment points on the head array. Must select hardware below to complete assembly. ESP recommended. Choose Omni Control Module separately.</i>	5,933.00
<input type="checkbox"/> I11072	Program Head Array as a 5-Switch System <i>Select additional arms/sensors below and indicate the intended use of each sensor in the Notes & Comments section. Choose Omni Control Module separately.</i>	No Charge
CHOOSE HEAD ARRAY MOUNT HARDWARE (REQUIRED):		
<input type="checkbox"/> I105927-99-0	Stealth TWB480 Backrest Mounting Hdwr with Links	330.00
<input type="checkbox"/> I105941-99-0	Stealth TWB480 Backrest Mnt Hdwr with Flipback	545.00
SELECT LEFT ARM PACKAGE (OPTIONAL)		
<input type="checkbox"/> I10976	Add Left Lateral Arm w Proximity Sensor	449.00
<input type="checkbox"/> I10977	Add Left Lateral Arm w Egg Switch	299.00
<input type="checkbox"/> I10978	Add Left Lateral Arm w Pad (No Sensor or Switch)	183.00
SELECT RIGHT ARM PACKAGE (OPTIONAL):		
<input type="checkbox"/> I10979	Add Right Lateral Arm w Proximity Sensor	449.00
<input type="checkbox"/> I10980	Add Right Lateral Arm w Egg Switch	299.00
<input type="checkbox"/> I10981	Add Right Lateral Arm w Pad (No Sensor or Switch)	183.00

ASL FUSION HEAD ARRAY:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I109582-99-0	ASL Fusion Head Array (available from 4/1/2022) <i>Proportional and digital head array includes occipital (ASL 935) and lateral (ASL 933) pads, and pulse mechanical switch (ASL 320BK) for reset/mode changes. Must select hardware listed below to complete assembly.</i>	8,515.00
<input type="checkbox"/> I11580	Swing Away Hardware for Lateral Sensor - Left	425.00
<input type="checkbox"/> I11581	Swing Away Hardware for Lateral Sensor - Right	425.00
<input type="checkbox"/> I11582	Adj/Removable Headrest Mounting Hardware	425.00
<input type="checkbox"/> I109583-99-0	ASL Fusion Display and Programmer with Mount <i>Includes ASL 165 and mounting bracket.</i>	905.00

Alt Drive Systems Setup

CHOOSE OMNI (REQUIRED):

Part Number	Description	MSRP - CA
<input type="checkbox"/> I108500-99-0	Omni2 Control Module for R-net - Corpus <i>Two-section design features 3.5" VGA high resolution display, USB charger port, built in IR control and integrated Bluetooth technology for control of up to 4 Bluetooth devices. Must also select mounting hardware below. Please specify which type of alternative drive input device you will be using. Input Device: _____</i>	2,095.00
<input type="checkbox"/> I103476-99-0	Adjustable Swing-Away Hardware - Omni/Tablet <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	305.00

Alt Drive Systems Setup

CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I10080	Mount Seat Control Switchbox on Backrest	No Charge

Alt Drive Accessories

PERMOFIX MOUNTING SOLUTIONS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I106655-99-0	PERMOfix UniTrack Fixed Mount Kit <i>Includes two base joints, two 4" rods and one 8" rod.</i>	334.00
<input type="checkbox"/> I106656-99-0	PERMOfix UniTrack Removable Mount Kit <i>Includes quick release mount and thigh support hardware, two base joints and three different length rods.</i>	587.00
<input type="checkbox"/> I106657-99-0	PERMOfix UniTrack Swing Away Mount Kit <i>Includes PERMOfix Swing Away module, two base joints, one 2" rod and one 4" rod.</i>	846.00
<input type="checkbox"/> I105907-99-0	Proximity Sensor + Swivel Assembly <i>Only available with the Total Control Head Array. Swivel assembly allows sensor to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array. Since Proximity Sensors require power, they CANNOT be used as On/Off switches.</i>	437.00
<input type="checkbox"/> I105919-99-0	Egg Switch + Swivel Mount Assembly <i>Swivel assembly allows switch to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	291.00
<input type="checkbox"/> I105912-99-0	Solid Oval Pad + Swivel Assembly <i>Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	177.00
<input type="checkbox"/> I106675-99-0	Oval Plate + Swivel Assembly - No Switch <i>Blank plate provides mounting option for Micro Lite, Buddy Button or other switch options. Swivel assembly allows pad to be mounted with PERMOfix or directly to lateral arm of Total Control Head Array.</i>	153.00

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA
<input type="checkbox"/> I101769-99-0	Microlight Switch Qty: _____	192.00
<input type="checkbox"/> I101767-99-0	Buddy Button Qty: _____	145.00
<input type="checkbox"/> I106293-99-0	Egg Switch Black Qty: _____	160.00
<input type="checkbox"/> I101771-99-0	Mini Cup Switch Qty: _____	142.00
<input type="checkbox"/> I107500-99-0	Proximity Switch Qty: _____	672.00
<input type="checkbox"/> I107495-99-0	Mini Joystick 5 Switch OC for Omni - (9-pin)	1,141.00
<input type="checkbox"/> I10772	Dome Switch, Backrest Mnt w Gooseneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____ Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights. Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 150mm rod, one 200mm rod and one 225mm gooseneck switch with 1.2m cable and mounting clamp.</i>	599.00
<input type="checkbox"/> I10846	Dome Switch, Backrest Swg Awy Mnt w Gsneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____ w 90° Swing-Away Function. Mounts at top of back. Custom charge applies if ordered with short back. Kit Includes: One PERMOfix joint with Flat Plate, one PERMOfix joint with 45° swivel adjustment, two 6" rods, one 90° Swing-Away joint, one 9" gooseneck switch w 50" cable and mounting clamp.</i>	898.00
<input type="checkbox"/> I105935-99-0	3.5mm Ext. Cable for Sensor/Switch, 0.3m Long	41.00
<input type="checkbox"/> I105931-99-0	3.5mm Ext. Cable for Sensor/Switch, 1m Long	45.00
<input type="checkbox"/> I105932-99-0	3.5mm Ext. Cable for Sensor/Switch, 2m Long	49.00
<input type="checkbox"/> I105933-99-0	3.5mm Ext. Cable for Sensor/Switch, 3m Long	55.00
<input type="checkbox"/> I105934-99-0	3.5mm Y-splitter Cable for Sensor/Switch <i>Allows two Proximity Sensors/Mechanical Switches to connect to the same direction jack.</i>	45.00

Notes & Additional Instructions:

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, LTD. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgment:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty.